Table of Contents by Language

ArthroCare Atlas System User's Manual – English	1
Manuel de l'opérateur du système Atlas ArthroCare – Français	31
Gebrauchsanweisung für ArthroCare System Atlas – Deutsch	61

ArthroCare® Atlas™ System

User's Manual

This equipment has been tested and found to comply with the limits for medical devices to IEC 60601-1-2:1994. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. If this equipment does cause harmful interference to other devices, which can be verified by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device,
- Increase the separation between the affected equipment and the Controller,
- Connect the affected equipment to an outlet or circuit different from that to which the Controller is connected, or
- Consult the manufacturer or field service technician for help.

Table of Contents

Description, Indications for Use, and Contraindications	7
Description	7
Indications for Use	7
Contraindications	7
System Overview	8
Connection Diagram	8
Principle of Operation	9
Warnings, Precautions, and Adverse Events	10
Warnings	10
Precautions	10
Adverse events	11
Controls, Indicators, and Alarms	12
Controls & Indicators	12
Alarms	13
Diagram of Controls, Indicators, and Alarms	14
Unpacking, Assembly, and System Check	15
Unpacking	15
Assembly and System Check	15
Instructions for Use	16
Operator Training Requirements	16
General System Operation	16
Voltage Outputs	17
System Preparation and Care	18
System Preparation	18
Wand Selection	19
System Shut Down	19
System Storage and Transport	19
Equipment Disposal	19

Cleaning and Sterilization20
Controller20
Foot Control20
Wand21
Patient Cable21
 Maintenance and Troubleshooting22
Maintenance22
Fuse Replacement22
Troubleshooting Guide22
 Product Specifications24
Technical Specifications24
Controller Output Graphs25
 Controller Classification and Safety Verification28
Classification28
Safety Verification28
 Customer Service29
Warranty Information29
Product Complaints29
 Symbols Key30

Description, Indications for Use, and Contraindications

Description

The ArthroCare® Atlas™ System is a bipolar, radiofrequency (RF) electrosurgical system designed for use in arthroscopic and orthopedic procedures. The System consists of the following components:

- 1) a bipolar radiofrequency controller;
- 2) a reusable, non-sterile Power Cord;
- 3) a reusable, non-sterile Foot Control;
- 4) a reusable, non-sterile Patient Cable; and
- 5) a disposable, sterile ArthroWand or Microdebrider Wand.

An optional reusable non-sterile Hand Control is also available for use with the System.

The Controller is the voltage source that delivers RF energy to the treatment site via a reusable Patient Cable and sterile Wand, or an Integrated Cable Wand (ICW), a sterile Wand/Cable combination.

The reusable Patient Cable is supplied non-sterile and is designed for sterilization prior to use.

The sterile disposable Wand is available in various single or multi-electrode configurations and is supplied separately. It may have a Cable component incorporated into its design to replace the Patient Cable.

The Controller is activated by either a reusable Foot Control or an optional reusable Hand Control.

Indications for Use


The ArthroCare Atlas System is indicated for resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels in arthroscopic and orthopedic procedures.

Contraindications

The ArthroCare Atlas System is contraindicated in any procedures where a conductive solution is not used. The System is also contraindicated for patients who have cardiac pacemakers or other electronic implants without specific instructions from the manufacturer of the cardiac pacemaker or implant. Please refer to the Wand Instructions for Use for a more comprehensive list of contraindications regarding specific procedures.

System Overview

Connection Diagram

1. Controller
2. Power Cord
3. Foot Control
4. Ablation Pedal
5. Coagulation Pedal
6. Set Point Adjustment Button
7. Patient Cable
8. Controller Connector
- 9a. Patient Cable Connector
- 9c. Cable Connector with Black Plug
10. Wand
11. Shaft
12. Return Electrode
13. Active Electrode Tip
14. Suction Tube (optional on Wand Style)
15. Irrigant Tube (optional on Wand Style)
16. Handle
17. Integrated Cable Wand (ICW)

Principle of Operation

The ArthroCare Atlas Controller is designed to deliver RF energy to the electrode elements located at the distal end of the sterile single-use Wands. Current flows between the active electrodes elements and the return electrode element, both elements or poles located on the Wand itself in a bipolar configuration, providing a localized energy field. The result of this arrangement is controlled energy delivery with minimal collateral tissue damage.

The ArthroCare Atlas System works by passing RF energy through a conductive solution (such as normal saline or Ringer's lactate) in close proximity to or in contact with the tissue to be treated. The conductive solution forms a thin layer between the active and return electrode elements. In the ablation mode, when sufficient energy is applied, the conductive solution is converted into a vapor layer (plasma) containing energized charged particles. When the high-energy charged particles come in contact with tissue, they cause its disintegration through molecular dissociation.

This mode of operation results in relatively low treatment site temperatures when compared to conventional electro-surgical and monopolar RF systems, thus yielding limited collateral thermal damage to the surrounding untreated tissue.

The system can also function when a lower voltage is applied between the active and return electrode(s). In this case, the electrical field is below the threshold required to create a plasma layer and resistive tissue heating occurs. This mode is useful when a greater thermal effect is needed, i.e. for coagulation of blood vessel. The appropriate voltage setting will depend on the design of Wand used, tissue type, and desired tissue effect.

Warnings, Precautions, and Adverse Events

The following is a list of Warnings and Precautions that apply to the general operation of the ArthroCare Atlas System. For specific warnings and precautions, please refer to the Wand and the Patient Cable Instructions for Use.

WARNINGS

- Failure to follow all applicable instructions may result in serious surgical consequences.
- Fire Hazard: **DO NOT** place active accessories near or in contact with flammable materials (such as gauze or surgical drapes).
- Electrosurgical accessories, which are activated or hot from use, can cause a fire.
- Accessory tips may remain hot enough to cause burns after the electrosurgical current is deactivated.
- Inadvertent activation or movement of Wands outside the field of vision may result in injury to the patient.
- Localized burns to the patient or physician may result from electrosurgical current carried through other instruments and conductive objects.
- Electrosurgical current may be generated in conductive objects by direct contact with the active electrode or by the active or return electrode being in close proximity to a conductive object.
- If excessive heating or physical forces cause damage to the Wand tip, foreign body fragments may result, possibly requiring extended surgery for removal.
- **DO NOT** use the ArthroCare Atlas System with non-conductive media (e.g. sterile water, dextrose, air, gas, glycine, etc.). Use only conductive media such as normal saline or Ringer's lactate.
- Electric Shock Hazard: **DO NOT** connect wet accessories to the Controller.
- Controller failure could result in an unintended increase in output power.

PRECAUTIONS

- Prior to initial use, ensure that all package inserts, warnings, precautions, and Instructions for Use are read and understood.
- Safe and effective electrosurgery is dependent not only on equipment design, but also, to a large extent, on factors under the user's control. Only persons having adequate training and familiarity with orthopedic surgery should perform procedures with the ArthroCare Atlas System.
- Consult medical literature relative to techniques, complications, and hazards prior to performance of any procedure.
- Evaluate patients for predisposing medical problems that may be aggravated by the stress of surgery.
- A thorough understanding of the principles and techniques involved in electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device and other medical instruments. Ensure that insulation or Controller grounding is not compromised.
- When instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- When not in use, remove the Wand from the surgical site and place away from metallic objects. Wands should remain separated from other electrosurgical equipment to avoid inadvertent electrical coupling between devices. Inadvertent activation may cause injury to patient and/or user or equipment damage.
- **DO NOT** wrap Patient Cable around metal objects. Wrapping cords around metal objects may induce currents that could lead to shocks, fires, or injury to the patient or surgical personnel.

- Use caution when using Wand tips to probe or manipulate tissue. Forceful contact between Wand tips and tissue or other instruments may result in damage to the instrument.
- **DO NOT** use the Wand as a lever to enlarge surgical site or gain access to tissue.
- **DO NOT** allow fluid to contact any electrical connectors on the Wands, Controller, or Cables during use.
- Maintain the lowest Controller Set Point necessary to achieve the desired tissue effect.
- Confirm proper activation of the Wand if a Controller Set Point is chosen outside of the selected, default settings.
- **DO NOT** allow patient contact with grounded conductive objects, such as a surgical table frame or an instrument table, to avoid risk of shock. Grounding pads should not be used.
- **DO NOT** contact metal objects with an activated Wand.
- Observe fire precautions at all times. Sparking and heating associated with electrosurgery may be an ignition source.
- **DO NOT** use flammable agents for cleaning and disinfection of the Controller or Cables.
- As with other electrosurgical units, electrodes and Cables can provide paths for high frequency current. Position the cables to avoid contact with the patient or other electrical leads.
- High frequency (HF) electrosurgical equipment such as the ArthroCare Atlas System may adversely affect the operation of other electronic equipment.
- Electrodes should remain separated from other electrosurgical equipment to avoid inadvertent electrical coupling between devices.
- Monitoring electrodes should be positioned as far as possible from the surgical electrodes when HF surgical equipment and physiological monitoring equipment are used simultaneously on a patient. Monitoring electrodes are not recommended.
- Monitoring equipment incorporating high frequency current-limiting devices is recommended.
- **DO NOT** remove the cover of the Controller. Refer servicing to qualified personnel.
- **DO NOT** obstruct the exhaust fan (located at rear of Controller).
- **DO NOT** touch the Controller's fan and/or speaker while touching the patient.
- Before each use, check that all Controller indicator lights and audio signals are functional. Make sure that the power cable plug is properly connected to the Controller receptacle.
- To avoid risk of fire, only replace the Controller fuses with the same type and rating.
- The ArthroCare Atlas System is designed to be operated exclusively as a unit. Only use accessories provided by ArthroCare.

ADVERSE EVENTS

As a consequence of electrosurgery, damage to surrounding tissue through iatrogenic injury could occur.

Controls, Indicators, and Alarms

Controls & Indicators

The Atlas Controller incorporates the following controls and indicators:

FRONT PANEL

1. On/Off Switch

This switch turns the AC power on and off. When the switch is on, green numbers will be visible in the Set Point Display Window. This display may remain active for up to 5 seconds following power turn-off.

2. Warning Indicator

A red indicator illuminates and an audible alarm when a Controller-specific failure or malfunction occurs.

3. Foot Control / Hand Control Connected Indicator

A green Foot Control / Hand Control Connector Indicator will illuminate when a Foot Control or Hand Control is properly connected.

4. Coagulation Activation Indicator

A blue indicator will illuminate when the Foot Control/Hand Control Coagulation function is depressed with a Wand attached.

5. Wand Connected Indicator

A green Wand Connector Indicator will illuminate when the Patient Cable and the Wand are properly connected. For an ICW, the Indicator will illuminate when the Wand cable connector end is properly connected to the Controller.

6. Set Point Adjustment

Increment and Decrement Arrow buttons control the Ablation Set Point. This output level can also be adjusted from the Ablation Set Point Adjustment function on the Foot Control or the Hand Control. Refer to the Instructions for Use section of this manual for the corresponding voltage levels for each setting.

7. Set Point Display Window

This display indicates the output level. A nominal setting will automatically be displayed when the System is powered up with a Wand attached or when a new Wand is connected. When the Controller is turned on, and prior to connection of a Wand, this display will show 0.

8. Ablation Activation Indicator

A yellow indicator will illuminate when the Foot Control/Hand Control Ablation function is depressed with a Wand attached.

9. Foot Control / Hand Control Receptacle

The Foot Control / Hand Control plugs into the tan receptacle located on the front of the Controller.

10. Symbol for Defibrillator-Proof Type BF Equipment

This equipment provides a degree of protection against electric shock to TYPE BF applied parts as defined in IEC 601-1, it also has an F type applied part capable of withstanding the effects of defibrillator discharge.

11. Black Cable Receptacle

The Cable Receptacle with the black ring will accept ICWs with black mating end.

12. Tan Cable Receptacle

The Cable Receptacle with the tan ring will accept the reusable Patient Cable and ICWs with gray mating end.

13. Sliding Door

This door slides left and right and allows access to either the Black Cable Receptacle or the Tan Cable Receptacle.

BACK PANEL**14. Tone Volume Control**

The Tone Volume Control regulates tone volume. To increase volume, turn the knob clockwise. To reduce volume, turn the knob counterclockwise.

15. Non-Ionizing Radiation Symbol

This symbol indicates that this equipment intentionally emits RF energy during activation.

16. Fuse Rating Symbol

This symbol indicates that only fuses with the appropriate rating should be used. Check the Controller back label for the appropriate fuse rating. See Maintenance and Troubleshooting section for fuse replacement instructions.

17. Power Cord Receptacle / Fuse Holder

The Controller Power Cord plugs into this receptacle. The fuse holder is behind the receptacle.

18. Equipotential Ground Symbol

This symbol identifies the conductor that is used to bond the equipment to earth ground.

19. CE Mark - European Certification Symbol

This indicates compliance with the European Commission Medical Device Directive (93/42/EEC).

20. Attention Symbol

This symbol alerts the user to read and understand this manual and accompanying instructions before operating the equipment.



21. Auxiliary Input

For manufacturer use only.

Alarms

1. Intermittent dual tone: current limit alarm or overload
2. Intermittent monotone: connection fault alarm (Patient Cable, Wand, Integrated Cable Wand, Over Temperature Condition, Over Power Condition)

Diagram of Controls, Indicators, and Alarms

Unpacking, Assembly, and System Check

Unpacking

Verify that all items have been received and are not damaged. Damage should be reported at once to the Customer Service Department. Save all containers and packaging material; they will be required if it is necessary to return the equipment.

Assembly and System Check

1. Connect the Power Cord to the receptacle on the rear panel of the Controller. Connect the other end of the Power Cord to the electrical outlet. If it is necessary to use a Power Cord other than the one supplied with the System, the additional Power Cord should comply with appropriate electrical standards and be suitable for hospital use.
2. Turn the On/Off Switch on the front panel of the Controller to the 'On' position. The Set Point Display should indicate a setting of '0'. All other indicator lights should be off.
3. Attach the Foot Control or Hand Control to the tan receptacle on the front of the Controller. The Foot Control / Hand Control Connector Indicator on the front panel of the Controller should illuminate.
4. Depress the Ablation function on the Foot Control or the Hand Control. The red Warning Indicator on the front panel of the Controller should illuminate and the Controller should emit an intermittent monotone alarm. The yellow Ablation Activation Indicator on the front panel should be off.
5. Attach the Patient Cable to the left cable receptacle of the two cable receptacles on the lower left of the Controller. The Patient Cable and Wand indicators should not illuminate. No alarms should sound.
6. Attach an ArthroWand to the Patient Cable. The green Wand Indicator should illuminate.

CAUTION: DO NOT contact metal objects with an activated Wand.

CAUTION: DO NOT place active accessories near or in contact with flammable materials (such as gauze or surgical drapes).

CAUTION: Electrosurgical accessories, which are activated or hot from use, can cause a fire.

CAUTION: Accessory tips may remain hot enough to cause burns after the electrosurgical current is deactivated.

CAUTION: Localized burns to the patient or physician may result from electrosurgical current carried through other instruments and conductive objects.

7. Taking care not to touch the Wand end, depress the Ablation (yellow) pedal of the Foot Control. The yellow Ablation Activation Indicator should illuminate. If the Controller does not function as described above, please contact Customer Service immediately.

Instructions for Use

Operator Training Requirements

The operator should be experienced in electrosurgical techniques. It is recommended that the user remain current with advances in orthopedic procedures. Additional training on the use of the ArthroCare Atlas System from an ArthroCare representative is recommended.

General System Operation

Wands are activated using either the Foot Control or the Hand Control. The Foot Control and the Hand Control have three functions as described below:

Ablation Activation

Pressing the Ablation Activation function operates the Controller's normal ablation mode and will activate the Wand.

Coagulation Activation

Pressing the coagulation function will activate the Controller's preset coagulation mode. And will activate the Wand for hemostasis of blood vessels.

Ablation Set Point Adjustment

Pressing the Ablation Set Point Adjustment function will adjust the Ablation voltage level on the Controller. Each time the Ablation Set Point Adjustment function is pressed, the Ablation voltage level increases by one level up to the maximum set point for each Wand style. Once the maximum level for the connected Wand has been reached, the System will cycle back to Set Point 1.

NOTE: The Foot Control and Hand Control cannot be used simultaneously.

Voltage Outputs

The voltage outputs correspond to the Set Point Displays as listed in the Table below. Depending the particular Patient Cable in use, the display corresponds to the output voltage in open circuit as follows:

Tan Patient Cable Receptacle

Display	Output Voltage (Vrms $\pm 10\%$)
1	100
2	125
3	150
4	175
5	200
6	225
7	250
8	275
9	300
C	70

Black Patient Cable Receptacle

Display	Output Voltage (Vrms $\pm 10\%$)
1	100
2	129
3	158
4	186
5	215
6	244
7	273
8	301
9	330
C	70

When the Controller is first powered on with no Wand attached, the Set Point Display digit '0' appears on the front panel display as a default setting. If a Wand is attached at or after power-up, the Controller will adjust the Display and voltage output to nominal settings. These settings will usually provide the best effect in most situations.

NOTE: If a Set Point is selected outside of the default range (between the initial Set Point and maximum Set Point), proper activation of the Wand should be confirmed.

Once a Wand has been properly connected, the voltage level can be increased or decreased by pressing the voltage level adjustment button located on the front panel. Ablation level can also be adjusted by pressing the Ablation set point adjustment function on the Foot Control or the Hand Control. The ablation setting levels may be adjusted to levels throughout the appropriate range, depending on the maximum voltage permitted by the particular Wand in use.

System Preparation and Care

System Preparation

1. Prior to each use, inspect the ArthroCare Atlas System for possible damage to the Controller Casing and cables.
 2. Insert the receptacle end of the Power Cord into the Power Cord Receptacle located at the rear of the Controller. Insert the plug end of the Power Cord into an appropriately grounded electrical outlet. Position the Controller so that the exhaust fan located in the rear of the Controller is unobstructed and directed away from the patient.
 3. Press the On/Off Switch on the front panel of the Controller to the 'On' (I) position. The Voltage Output Level Display should indicate a setting of '0'. All other indicator lights should be off.
 4. Connect the Foot Control or the Hand Control to the gray Foot Control / Hand Control Receptacle on the front of the Controller. The Foot Control / Hand Control Connected Indicator (green) on the Controller front panel will illuminate.
 5.
 - a. If using the reusable Patient Cable, ensure that the Patient Cable has been cleaned and sterilized and that the connectors are dry before use. Slide the sliding door to the right and connect the Patient Cable to the tan Receptacle on the front of the Controller. Connect the appropriate Wand to the other end of the Patient Cable.
NOTE: If sliding door has been removed, contact Customer Service.
 - b. If using an ICW with gray mating end, slide the sliding door to the right and connect the ICW to the tan Receptacle.
 - c. If using an ICW with black mating end, slide the door to the left and connect the ICW to the black Cable Receptacle on the front of the Controller.
- NOTE:** Refer to each Wand Instructions for Use for specific instructions regarding surgical preparation and procedures.
6. Set appropriate Set Point Display required to obtain desired effects.

Wand Selection

Select the Wand type most appropriate for the procedure. The Controller will preset nominal and maximum ablation set points for each Wand style to as a guide for safe and effective operation.

NOTE: Initial and maximum Set Points are suggested settings. Proper activation of the Wand should always be confirmed.

In order to adjust the ablation set points, the Wand and Patient Cable must be connected to the Controller.

System Shut Down

1. Turn the power switch to the 'OFF' position. After a brief (less than 4 seconds) delay, all lights on the Controller and the Set Point Display will go off.
2. Disconnect the suction tubing if appropriate.
3. If using the reusable Patient Cable, remove the Wand from the Patient Cable and disconnect the Patient Cable from the Controller. Dispose of the Wand and prepare the Patient Cable for sterilization and further use.
4. If using the ArthroWand with Integrated Cable, disconnect the Wand cable connector from the Controller. Do not attempt to separate the ArthroWand from the Cable component. Discard the Wand with Integrated Cable.

System Storage and Transport

All ArthroCare Atlas System components may be safely stored and transported at an ambient temperature range of -40°C (-40° F) to +70°C (158° F). The relative humidity should be between 10 and 85%.

Equipment Disposal

The ArthroCare Atlas System contains electronic printed circuit board assemblies. It should be disposed of in accordance with any applicable national or institutional policies relating to obsolete electronic equipment. Dispose of the Wand, Wand with Integrated Cable, and the Patient Cable according to normal institutional practices for potentially contaminated items.

Cleaning and Sterilization

Controller

DO NOT sterilize or immerse in liquid. Wipe clean with a soft cloth and mild detergent as needed.

Foot Control

DO NOT sterilize. Clean with detergents and disinfectants according to standard practices. Disinfect with liquid chemical disinfectants such as chlorine solutions, iodophors, glutaraldehydes and hydrogen peroxides. Follow manufacturer guidelines for concentration and length of exposure.

Wand

The Wand is supplied sterile. The Wand is intended **for single use only**. **DO NOT clean, resterilize, or reuse the Wand as this may result in product malfunction, failure, or patient injury, which may also expose the patient to the risk of transmitting infectious diseases.** Please refer to the Instructions for Use associated with each Wand type for specific information concerning Wand use.

Patient Cable

The Patient Cable is supplied NON-STERILE. It is reusable if cleaned and sterilized prior to each use. Ensure that a validated and approved steam sterilization process is used during the sterilization process.

Cleaning

- Before cleaning, secure the protective caps over the Patient Cable connectors to protect them from possible damage. The protective caps should fit snugly over the Patient Cable connector. If the protective caps are loose, check to see if the silicone o-ring is positioned within the groove on the cable connector. Replacement o-rings are available through ArthroCare. The Patient Cable should then be wiped down with a mild detergent solution using a soft brush or sponge to remove any gross contaminants from the Patient Cable. The Patient Cable can be rinsed under running water to remove any cleaning residue. The Patient Cable should not be immersed in water or any other solution. Remove the protective caps from the Patient Cable connectors and inspect the Patient Cable for any damage such as cuts or nicks. Damaged Patient Cable should not be reused. Sterilize the clean Patient Cable by following one of the recommended sterilization methods:

Pre-Vacuum Steam Sterilization

If the pre-vacuum steam sterilizer has a prefixed cycle, use the "Hard Goods" cycle. If it does not have a prefixed cycle, the following parameters are recommended:

- Set temperature at 270-272°F (132-133°C);
- Set exposure time for 10 minutes for wrapped Patient Cables;
- Set exposure time for 4 minutes for unwrapped Patient Cables;
- Set drying time for 5 minutes minimum.

Cleaning and Sterilization (Cont.)

Gravity Displacement Steam Sterilization

The following parameters are recommended for gravity displacement steam sterilization:

- Set temperature at 270-272°F (132-133°C);
- Set exposure time to 15 minutes for wrapped Cables;
- Set exposure time to 10 minutes for unwrapped cables;
- Set drying time for 8 minutes minimum.

Or

- Set temperature at 250-254°F (121-123°C);
- Set exposure time to 35 minutes for wrapped or unwrapped cables;
- Set drying time for 8 minutes.

CAUTION: Failure to properly clean the Patient Cable may lead to inadequate sterilization.

CAUTION: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and calibrated.

CAUTION: To insure adequate sterilization, make sure that the protective caps do not cover the cable connectors during sterilization.

CAUTION: Make sure Cable is thoroughly dry before use. Wet Cable may damage the Cable and Controller.

Maintenance and Troubleshooting

Maintenance

Other than fuse replacement, the Controller has no user-serviceable parts. It is designed to provide consistent output levels and is calibrated by clock crystals, voltage references, and fixed resistors. There are NO internal adjustments in the instrument and, due to the integrated calibration methods, no annual maintenance check is required. There is NO software incorporated in the ArthroCare Atlas Controller. If any component malfunctions, call Customer Service for a return authorization.

Fuse Replacement

The fuse holder is located on the back of the Controller. To replace a fuse, turn off the power to the Controller and unplug the Power Cord from the power outlet at the rear of the Controller. After waiting at least 10 seconds for internal circuitry to discharge, use a screwdriver or similar tool to remove the fuse holder by depressing the locking tabs. Replace both fuses with the same type and rating as specified on the rear panel of the Controller. Reinsert the fuse holder until the locking tabs snap into place. Reconnect the Power Cord and restore power to the Controller. If a fuse fails again, disconnect all power to the Controller and contact Customer Service.

Troubleshooting Guide

If you are experiencing problems with the ArthroCare Atlas System, you may want to use the following troubleshooting guide to help identify or eliminate the problem before contacting Customer Service:

- **System does not power up after the power switch is pressed**

Check that the Power Cord is properly connected to the Controller and plugged into an appropriately grounded outlet. If the unit is plugged in properly, check if the fuses have blown. To change the fuses, follow the instructions for Fuse Replacement.

- **Green Foot Control / Hand Control Connected Indicator light does not illuminate**

Check that the Foot Control / Hand Control is properly connected to the Controller and that the cord is not nicked, cut or frayed. Do not use the Foot Control / Hand Control if the cord has been damaged. If the problem persists, change the Foot or Hand Control. If the Foot Control/Hand Control Connected Indicator Light is still not illuminated, return the System for Service.

- **Green Wand Connected Indicator light does not illuminate when a Wand is connected to the reusable Patient Cable.**

Make sure that the Wand is securely seated in the Patient Cable, and that the Patient Cable is properly connected to the Controller. If the problem persists, first change the Patient Cable and then the Wand. If the Wand indicator light is still not illuminated, return the System for service.

Maintenance and Troubleshooting (Cont.)

- **Nothing happens when one of the device activation functions on the Foot Control or the Hand Control is depressed**

Verify that both the Foot Control / Hand Control Connected Indicator light and the Wand Connected Indicator light are illuminated when the Foot Control or Hand Control is depressed. Check if the voltage level has been adjusted to a level appropriate for operation (usually 1 or greater). Confirm that the Wand tip and shaft are covered by a conductive irrigant. Make sure that the Wand is securely seated in the Patient Cable, and that the Patient Cable is properly connected to the Controller. If the problem persists, first replace the Patient Cable and then the Wand. If the System still fails to operate, return for service.

- **The Wand does not activate; an intermittent monotone alarm sounds and a red warning light illuminates**

This generally indicates a connection problem. Either the Wand is not fully seated in the Patient Cable, or the Patient Cable is not properly connected to the Controller. Check all connections. If the alarm continues to sound when the Foot Control or the Hand Control button is depressed, first replace the Patient Cable and then the Wand. Return the System for service if the problem continues.

- **A dual tone alarm sounds and a red warning light illuminates when the Wand is activated.**

This is a safety feature of the ArthroCare Atlas System and may occur if the Wand is activated for an extended period of time without contacting tissue. To reset the unit, step off the foot pedal if the System is connected to the Foot Control. If the System is connected to the Hand Control, release the button from the Hand Control. Make sure that the Wand is in good contact with the target tissue, and depress the Foot Control or Hand Control again. If the alarm continues to sound, first replace the Patient Cable and then the Wand. If the problem persists, return the System for service.

Product Specifications

Technical Specifications

Patient Cable

Overall Length	3.3 m (10 feet)
Sterilization Method	Steam

Controller

Input Power Requirements

Voltage	100-240 V~
Frequency50/60 Hz
RMS Current8 Amps Max
Fuse Rating	T8 A 250V for 110 V~ T4 A 250V for 220 V~

Output Power

Fundamental Frequency100 kHz
Voltage Range0-330 Vrms @100 kHz
Max. Output Power400 W @ 250 Ω
Operating Temperature10° C to 40° C

Controller Dimensions


Weight (max)	<12 kg (<24 lbs)
Height14.2 cm (5.65 inches.)
Width34.9 cm (13.75 inches.)
Length40.9 cm (16.1 inches.)

Foot Control and Hand Switch


Cable Length	4.7 m (15 ft.)
Sterilization (Hand Control only)	Steam

Controller Output Graphs

Output power at each set point with specified load resistance (per IEC 60601-2-2, subclause 6.8.3) is given in the graphs below:


Output Power vs. Set Point at 300 Ohms Resistive Load, Tan Patient Cable Receptacle


Output Power vs. Set Point at 250 Ohms Resistive Load, Black Patient Cable Receptacle

The output power (at full and half settings) versus load resistance (per IEC 60601-2-2, subclause 6.8.3) is given in the graphs belows:


Output Power vs. Load Resistance, Tan Patient Cable Receptacle


Output Power vs. Load Resistance, Black Patient Cable Receptacle

The open circuit voltage for each set point is given in the graphs below:

Output Voltage vs. Set Point, Tan Patient Cable Receptacle

Output Voltage vs. Set point , Black Patient Cable Receptacle

Controller Classification and Safety Verification

Classification

According to IEC 60601-2-2, Specification for High Frequency Surgical Equipment, the Controller is classified as follows:

- Type of protection against electrical shock:
Class I equipment.
- Degree of protection against electrical shock:
Defibrillation proof, type BF (Isolating/floating).
- Degree of protection against harmful ingress of water:
 - Controller meets requirements of IEC 60601-2-2, subclause 44.6.
 - Foot Control meets requirements of IEC 60601-2-2, subclause 44.6, watertight construction (IPX8).
- Equipment not suitable for use in the presence of a flammable anesthetic mixture.
- Mode of operation: capable of continuous operation.

Safety Verification

The ArthroCare Atlas System meets the requirements of IEC 601-1, IEC 60601-2-2. It is recommended that the biomedical engineering department test the System to ensure that it meets the following leakage levels:

- Leakage current $\leq 20\mu\text{A}$ at 100-240 V~, 50/60 Hz
isolated patient connections
- Leakage current $\leq 200\mu\text{A}$ at 100-240 V~, 50/60 Hz
non patient applied parts

If the System fails to meet the following specifications, please contact ArthroCare Customer Service for a return merchandise authorization.

Customer Service

Warranty Information

The ArthroCare Atlas Controller and Foot Control are warranted for one year and the warranty for the reusable Patient Cable extends for a period of 90 days from the date of shipment to the original purchaser. Any component of the System, which develops defects resulting from defective material or workmanship during these time periods will be replaced or repaired without charge.

Product Complaints

All questions or concerns related to the quality, reliability and/or durability of this product should be directed to Customer Service or an authorized ArthroCare representative. Please contact Customer Service for an authorized ArthroCare representative for a return authorization.

Manufacturer

ArthroCare Corporation
Sunnyvale, CA U.S.A.
Toll Free: (800) 797-6520 (Customer Service)
www.artrocare.com

Manufactured By

ArthroCare Corporation
Sunnyvale, CA U.S.A

ArthroCare Costa Rica
La Aurora Heredia, Costa Rica

Authorized European Representative

ArthroCare Europe AB
Stockholm, Sweden
+46 8 442 72 30

Symbols Key

	See Instructions for Use
	Date of Manufacture
	Coagulation
	Ablation
	Wand Connected
	Foot Control/Hand Control Connected
	Defibrillator- Proof Type BF Equipment
	Fuse Rating
	Non-Ionizing Radiation
	Equipotential Ground
	Tone Volume Control
	Fragile, Handle with Care
	Temperature Limitations
	Keep Dry
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

The ArthroCare Atlas System is covered by the following U.S. Patents: 5,697,909; 5,679,281; 5,697,536; 5,697,882; 5,683,366; 5,681,282; 5,766,153; 5,810,764; 5,843,019; 5,871,469. Additional patents issued and pending.

ArthroCare and Atlas are registered trademarks or trademarks of ArthroCare Corporation.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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ArthroCare und Atlas sind eingetragene Warenzeichen oder Warenzeichen der ArthroCare Corporation.