ENGLISH..... 1  
FRANÇAIS .. 24  
ITALIANO .. 48  
POLSKI..... 72  
TÜRKÇE..... 96

# NIOX MINO®

User Manual  
Manuel de l'utilisateur  
Manuale d'uso  
Instrukcja obsługi  
Kullanıcı Kılavuzu

CE

Aerocrine

## Remember

- It is of vital importance to adhere to the specified operating conditions:  
 Ambient temperature: +16 to 30°C  
 Humidity: 20 to 60% RH (non-condensing)
- Mobile phones and cordless phones might interfere with the instrument and should therefore be kept away from the instrument. Interference could make it impossible to perform a measurement.
- When within the specified operating conditions it is recommended to keep NIOX MINO® with installed Sensor connected to a power outlet whenever possible.
- It is recommended to allow a new Sensor to be connected for two hours before the first measurement is done.
- It is recommended not to use NIOX MINO in the proximity of areas where volatile substances such as organic fluids or disinfectants are being used. Special attention should be paid to aerosols and disinfection baths (either open vessels or ultrasonic baths).
- Always use a NIOX MINO Bag (closed) for transportation and storage of NIOX MINO when it is disconnected from the power outlet. (The NIOX MINO Bag is sold separately.)
- Shelf-life - NIOX MINO instrument: Minimum 3 years at time of delivery, or 3,000 measurements.  
 NIOX MINO Sensor: Maximum 12 months after opening package and mounted in NIOX MINO, whichever comes first.

**Note!**

This manual is intended for medical care personnel and for use when teaching patients how to use the NIOX MINO instrument. To aid this, a demonstration mode is available in the instrument. Please refer to the "Demonstration mode" section, page 8.

Contact for Middle-East:

Stallergenes Middle-East Regional Office  
 Al Madina Al Munawar St., Al Baddad Bldg, No.219  
 3rd Floor, Office No.:307  
 Amman-Jordan  
 Telefax:00962-6-5532008/5542008

Country	Telephone	Telefax
Jordan	96265353761	96265354593
Lebanon	9614403410/11	9614520164
Saudi Arabia	96626720888	96614771374
Bahrain	96626720888	96614771374
Kuwait	96524312240/41	96524312239
Qatar	97444410270	97444430094
United Arab Emirates	97143474800	97143477909
Oman	97143474800	97143477909
Iraq	96265353761	96265354593

# Table of contents

<b>Warning!</b> .....	<b>2</b>	Install Visual Incentive .....	<b>20</b>
<b>Intended use (Europe)</b> .....	<b>2</b>	Start Visual Incentive.....	<b>20</b>
<b>Training requirements</b> .....	<b>2</b>	Download instrument data .....	<b>20</b>
<b>Presentation</b> .....	<b>3</b>	Shut down Visual Incentive.....	<b>20</b>
<b>Installation and set-up</b> .....	<b>4</b>	<b>Research applications</b> .....	<b>21</b>
<b>Measurement</b> .....	<b>6</b>	Special 6s FeNO measurement .....	<b>21</b>
Enter patient ID (optional).....	6	Nasal measurement .....	<b>21</b>
Perform FeNO measurement.....	7	<b>Cautions</b> .....	<b>22</b>
<b>View stored results</b> .....	<b>8</b>	<b>Clinical documentation</b> .....	<b>22</b>
<b>Ambient NO measurement</b> .....	<b>8</b>	<b>Vigilance</b> .....	<b>23</b>
<b>Demonstration mode</b> .....	<b>8</b>	<b>Guidance and manufacturer's declaration</b> .....	<b>23</b>
<b>Change settings</b> .....	<b>9</b>	Electromagnetic emissions .....	<b>23</b>
Time and date.....	9		
Sound volume .....	9		
Top light intensity .....	9		
Bluetooth on/off.....	9		
QC on/off .....	9		
<b>Information menu</b> .....	<b>10</b>		
<b>Turn off NIOX MINO®</b> .....	<b>10</b>		
<b>General care</b> .....	<b>11</b>		
Preventive inspections .....	11		
Change Sensor.....	11		
Change NO scrubber.....	11		
Disposal of used/expired products .....	12		
Return shipments.....	12		
Support.....	12		
Limited warranty.....	12		
<b>Troubleshooting</b> .....	<b>12</b>		
Alert codes .....	13		
<b>Technical data</b> .....	<b>14</b>		
Display buttons and symbols .....	14	User manual EPM-000200, version 01, October	
Symbols explanation.....	16	2013, <b>for instruments with software version 20XX</b>	
Backplate.....	16	<b>and 21XX.</b> X can be any number between 0 and 9.	
Base label .....	16	The version number for your instrument can be seen	
Dimensions and weight .....	16	in the Information menu, see page 10.	
Electrical data .....	16		
Noise level (stand by) .....	16	Information in this document is subject to change.	
Shelf-life .....	16	Amendments will be made available by Aerocrine	
Operating conditions .....	16	AB as they occur.	
Exhaled NO - performance data .....	17		
Linearity .....	17		
Precision.....	17		
Accuracy .....	17		
Method comparison.....	17		
Inhalation parameters .....	17		
Exhalation parameters.....	17		
Memory capacity.....	17		
Transport and storage.....	17		
Patient filter (mouthpiece) .....	17		
Responsible manufacturer.....	17		
<b>NIOX MINO parts and accessories</b> .....	<b>18</b>		
<b>Connect NIOX MINO to a PC using USB</b> .....	<b>19</b>		
<b>Connect NIOX MINO to a PC using Bluetooth</b> .....	<b>19</b>		
<b>NIOX MINO Visual Incentive</b> .....	<b>20</b>		

User manual EPM-000200, version 01, October 2013, **for instruments with software version 20XX and 21XX.** X can be any number between 0 and 9. The version number for your instrument can be seen in the Information menu, see page 10.

Information in this document is subject to change. Amendments will be made available by Aerocrine AB as they occur.

- NIOX MINO is CE-marked according to In Vitro Diagnostic Device Directive 98/79/EC and approved for clinical use in EEC countries.
- NIOX MINO is 510(k) cleared, k101034, by FDA.
- NIOX MINO is RoHS compliant.
- Copyright © 2013 Aerocrine AB, Solna, Sweden.
- Aerocrine, NIOX MINO and NIOX are registered trademarks of Aerocrine AB.

## Warning!

- Handle the NIOX MINO® instrument as stated in this manual. Aerocrine accepts no responsibility for damaged equipment or faulty results if the equipment is not handled according to this manual.
- Do not use a damaged NIOX MINO instrument or damaged components.
- Use only the power supply unit provided.
- Keep the instrument out of water. Ensure that no liquid is spilled or dripped on the instrument.
- Do not heat or dispose of the instrument or Sensor in fire. Please refer to the "Disposal of used/expired products" section.
- Take care not to drop the instrument or subject it to strong impact.
- It is recommended not to use NIOX MINO in the proximity of areas where volatile substances such as organic fluids or disinfectants are being used. Special attention should be paid to aerosols and disinfection baths (either open vessels or ultrasonic baths).
- NIOX MINO should not be used adjacent to or stacked with other equipment.
- The NIOX MINO Sensor contains chemicals that could be harmful if swallowed.
- Touch only the grey cap when exchanging the Sensor.
- Do not clean the Sensor. Cleaning of the Sensor with ethanol or similar disinfectant might destabilize it for a non-predictable time period.
- Keep the Sensor out of water. Ensure that no liquid is spilled or dripped on the Sensor.
- The NO scrubber contains potassium permanganate and should be disposed of as hazardous waste in accordance with the local waste disposal regulations.
- When selecting an accessory for your NIOX MINO please keep in mind that an accessory not recommended by Aerocrine may result in loss of performance, damage to your NIOX MINO, fire, electric shock, injury or damage to other property. The product warranty does not cover product failure or damage resulting from use with non-approved accessories. Aerocrine takes no responsibility for health and safety problems or other problems caused by the use of accessories not approved by Aerocrine.
- No modification of the NIOX MINO instrument or the Sensor is allowed.

Also see Cautions page 22.

## Intended use (Europe)

NIOX MINO® Airway Inflammation Monitor measures Fractional Exhaled Nitric Oxide (FeNO) in human breath. The level of Nitric Oxide is elevated when eosinophilic airway inflammation is present and decreases as a response to anti-inflammatory treatment. The FeNO concentration in expired breath can be measured by NIOX MINO with assurance that such measurements are repeatable and according to guidelines for NO measurement established jointly by the European Respiratory Society and the American Thoracic Society<sup>1</sup>.

Measurement of FeNO by NIOX MINO is a quantitative, non-invasive, simple and safe method to measure the FeNO concentration in patients with asthma-like airway symptoms in order to establish and monitor the presence of eosinophilic airway inflammation.

NIOX MINO is suitable for children, approximately 4 - 17 years, and adults 18 years and older.


FeNO measurements should be used as part of the regular assessment and monitoring of asthma patients to improve the patient's asthma management. NIOX MINO should be used as directed in the NIOX MINO User Manual. NIOX MINO may not be used with infants or by children approximately under the age of 4, as measurement requires patient cooperation.

<sup>1</sup> Am J Respir Crit Care Med 2005; 171: 912-930.

## Training requirements

NIOX MINO should only be used as directed in this manual, by trained physicians, nurses, respiratory therapists and laboratory technicians. "Trained" status is achieved only after careful reading of this manual.

# Presentation


## Installation and set-up


### **Caution!**

Be careful when opening the Sensor can. The inside of the opening has sharp edges. Do not touch or clean the white Sensor membrane


1. Insert the Sensor into the compartment at the base of the instrument


2. Turn the orange swivel until locked

3. Spread the latches apart, one at a time, (Optionally: Use the red tool delivered with the instrument to spread the latches apart)


... and carefully pull out the cover


4. Insert the new NO scrubber into the cover


Make sure to use a new NO scrubber, with an unbroken hole


5. Replace the cover and make sure it snaps in place

6. Attach the power supply unit to the instrument and to a power outlet

7. Wait for the instrument to start up and press the Forward button


8. Reminder to insert a new NO scrubber. Press the Forward button

9. Select OTHER on the display

### **Note!**

This screen is only displayed at first start-up of the instrument. The Japan icon is only available for SW 21XX.


10. Select Settings

11. Select Clock

12. The time and date settings are shown

13. Click the time format icons to change between US and ISO time and date

14. Use the plus and minus buttons to set the time and date

15. Select the forward button to change between hour, minute, year, month and day


16. Select OK to accept the changes


17. Select NO scrubber

18. Again select NO scrubber

19. Input the passcode 0000 using the number buttons to confirm that a new NO scrubber is installed

20. Select OK to accept the changes

21. Optional: Activate Bluetooth for wireless communication with PC or connect a USB cable, see page 19

22. Select Return to go back to the Main screen


23. Allow the instrument to stabilize.

**Note!**

Stabilization of the instrument normally may take up to 30 minutes following connection of the power supply unit to the power outlet. However, it is recommended to allow a new NIOX MINO® Sensor to be connected for two hours before the first measurement is done.

## Measurement


Depending on what is shown on the display, proceed as follows:

Stabilization in progress  
Wait until ready

Instrument in sleep mode  
Touch the display

Ready for use

## Enter patient ID (optional)

### Note!

If Patient ID is used, it has to be entered before each measurement (even if it is the same patient).

Input a patient specific ID number, up to 10 digits.


1. Select the ID button on the main screen

2. Input the patient specific ID using the number buttons

Use the shift buttons to change between digits 0 to 4 and 5 to 9

3. Select OK to accept the patient ID


## Perform FeNO measurement

Make sure that the instrument is ready for measurement (activate by pressing the display if in sleep mode)

### Important!

Always attach a new patient filter for each new patient

### Choose one of three options in order to follow the measurement

Listen to the sound signals and view the top light.


Look at the display using a mirror.

If the optional PC based NIOX MINO® Visual Incentive is used, look at the PC screen. (see page 20)

#### 1. Empty lungs

#### 2. Inhale deeply through the filter to total lung capacity


The cloud on the display is inflated and the top light is turned off while you inhale

The sun/arrow rises  
Inhalation time is displayed if meter is selected

#### 3. Exhale slowly through the filter until the top light and sound turns off

Continuous sound and steady light = OK


Intermittent high frequency sound and flickering light = exhalation too hard

Intermittent low frequency sound and flickering light = exhalation too weak

#### 4. Wait for result

FeNO value  
Patient ID (if entered)  
Sequence number  
(each measurement result is given a unique identification number)

## View stored results

A measurement result is automatically displayed at the end of a measurement. All previous results are stored in the instrument, and can be viewed at any time.


1. Select Mode

2. Select Measurement results

3. The latest stored measurement is displayed, showing:

A. FeNO value  
B. Measurement mode  
C. Patient ID (if entered)  
D. Sequence number (each measurement result is given a unique identification number)  
E. Time and date of measurement

4. Use the previous and next buttons to step through the stored measurements

5. Select Return to go back to the Mode screen


## Ambient NO measurement

### Note!

A completed ambient NO measurement is counted as one measurement for the Sensor and NIOX MINO®.

1. Attach a new patient filter


2. Select Mode

3. Select Ambient measurement

4. Wait for result (approximately 4 minutes)

## Demonstration mode

Use this mode for a new patient in order to demonstrate the different steps during a measurement and the corresponding illustrations on the display as well as light and audio feed-back. (Sound volume has to be activated.)

1. Select Mode

2. Select Demo

3. The inhalation screen is shown and the top light is turned off

4. Use the Forward button to advance to the next screen

5. The normal exhalation screen is shown

Correct exhalation pressure: The cloud in a steady centered position, the top light is lit, and the audio emits a constant sound

6. Use the Forward button to advance to the next screen

7. The exhalation too weak pressure screen is shown

8. Use the Forward button to advance to the next screen


9. The exhalation too strong pressure screen is shown

10. Use the Forward button to advance to the next screen

11. The wait for result screen is shown (static view)

12. Select Return to go back to the Mode screen


# Change settings

## Time and date

1. Select Settings
2. Select Clock
3. The time and date settings are shown
4. Select the time format icons to change between US and ISO time and date
5. Use the plus and minus buttons to set the time and date
6. Select the forward button to change between hour, minute, year, month and day
7. Select OK to accept the changes

4. The top light setting is shown


5. Use the plus and minus buttons to set the light intensity
6. Select OK to accept the changes


## Bluetooth on/off

1. Select Settings
2. Select Bluetooth
3. The Bluetooth setting is shown
4. Use the buttons to deactivate or activate Bluetooth
5. Select OK to accept the changes

## Sound volume

1. Select Settings
2. Select Sound and Light
3. The sound volume setting is shown
4. Use the plus and minus buttons to set the sound volume
5. Select OK to accept the changes

## QC on/off

### Note:

An optional External Quality Control (QC) procedure is available for NIOX MINO®.

Refer to the NIOX MINO External Quality Control user manual. Contact your local Aerocrine representative or Aerocrine Customer Support for further information.


## Top light intensity

1. Select Settings
2. Select Sound and Light
3. Select top light

## Information menu

1. Select Settings
2. Select Information menu
3. The information screen shows:
  - A. Number of remaining sensor measurements
  - B. Sensor expiration date
  - C. Sensor serial No.
  - D. Number of remaining instrument measurements
  - E. Instrument expiration date
  - F. Instrument serial No.
  - G. Instrument software version
  - H. Ambient RH (relative humidity), measured inside the instrument (the value may differ from RH measured outside the instrument)


4. Select Return to go back to the Settings screen

... alternatively select Alert to view the alert codes

5. The alert code information screen is shown, showing the 16 latest alerts

### Note!

The alert codes in the list are only to be used in contact with Customer support

- Date of the alert
- Alert code (for customer support purpose only)
6. Select Return to go back to the Information screen

## Turn off NIOX MINO®

### Note!

When within the specified operating conditions (see page 16) it is recommended to keep NIOX MINO with installed Sensor connected to a power outlet whenever possible.

1. Pull out the power supply unit connector from NIOX MINO

## Moving NIOX MINO

The instrument can be moved to another power outlet with optimized initiating startup.

### Note!

The instrument must be connected to the new power outlet within a few minutes.


1. Select Settings


2. Select Controlled power off


3. Select OK to confirm

4. Move the instrument and connect it to a power outlet

5. Wait until the instrument is ready to use

## General care

- Use a new NIOX® Filter (patient filter) for each new patient.
- Clean the instrument with a cloth dampened with 70% ethanol, isopropanol or similar disinfectant. Be careful when using disinfectants, as excess of alcohol (ethanol) might permanently damage the Sensor. Do not use spray detergents.

**Note!**

Never attempt to open or service the NIOX MINO® instrument or Sensor.

## Preventive inspections

Before each measurement, visually inspect:

- that NIOX MINO is not damaged
- that the Sensor is in place
- that a NIOX Filter is attached
- that the Power supply cord is undamaged and correctly connected to the instrument and to the power outlet.


If any item is missing or damaged, contact your local Aerocrine representative or Aerocrine AB.

## Change Sensor


**Caution!**

Be careful when opening the Sensor can. The inside of the opening has sharp edges. Do not touch or clean the white Sensor membrane


1. Press and hold the blue button...


...while turning the orange swivel


2. Remove the Sensor

3. Insert the new Sensor


4. Turn back the orange swivel until locked


## Change NO scrubber


1. Spread the latches apart, one at a time, (Optionally: Use the red tool delivered with the instrument to spread the latches apart) ... and carefully pull out the cover

2. Remove the used NO scrubber from the cover

3. Insert the new NO scrubber into the cover


Make sure to use a new NO scrubber, with an unbroken hole

Unbroken      Broken

4. Replace the lid and make sure it snaps in place


5. Select Settings


6. Select NO scrubber

7. Again select NO scrubber

8. Input the passcode 0000 using the number buttons to confirm that a new NO scrubber is installed

9. Select OK

### **Disposal of used/expired products**

- Used or expired NIOX MINO® Sensors should be recycled according to the local recycling program for electronic equipment. Note that there is a LiMnO<sub>2</sub> battery inside the Sensor unit\*.
- An expired NIOX MINO should be recycled, according to the local program for electronic equipment. Note that there is a LiMnO<sub>2</sub> backup battery inside the instrument\*.
- The NO scrubber contains potassium permanganate and should be disposed of as hazardous waste in accordance with the local waste disposal regulations.

NIOX MINO is RoHS compliant.

\* The batteries are not user-replaceable parts.

### **Return shipments**

Contact your local Aerocrine representative or Aerocrine AB.

### **Support**

Please contact your local Aerocrine representative or Aerocrine Customer Support if you encounter problems, which you can not solve with the actions stated in this manual.

For contact details, see back-page, and provide the following information:

- Your name, address and telephone number.
- Serial No. (both instrument and Sensor).
- Problem description (as thorough as possible).
- Alert codes or lists.

### **Limited warranty**

Aerocrine AB provides a Limited Warranty for this instrument and original accessories delivered with this instrument. Conditions are defined when the items are purchased.

Do NOT try to repair the instrument. It is NOT permitted to open the instrument. Any attempt will make the warranty invalid and performance according to the specifications can not be guaranteed.

## **Troubleshooting**

<b>Warning</b>	<b>Action</b>
	The inhalation was too weak to initiate a measurement or an exhalation into the instrument was performed prior to an inhalation - Stop the procedure immediately when this warning appears. Wait until the main menu screen is displayed and repeat the inhalation with a stronger inhalation force.
	NO scrubber reminder. The symbol is shown at first start-up of the instrument as a reminder to insert and set the software for a new NO scrubber. See <i>Installation and set-up</i> section page 4.
	NO scrubber almost expired - Order a new NO scrubber. The symbol is shown when 10% of the measurements are remaining or 2 weeks before expiration date and continue until the NO scrubber has expired. A NO scrubber can be used for 1000 measurements or 1 year. Refer to the Change NO scrubber section page 11.
	No Sensor connected - Connect a Sensor.
	Sensor almost expired – Order a new Sensor. The symbol is shown when 10% of the measurements are left or 2 weeks before expiration date and continue until the Sensor has expired. Refer to the Change Sensor section page 11.

Warning	Action	Code	Action
	<p>Instrument almost expired – Order a new instrument.</p> <p>The symbol is shown 4 months before the instrument expires or when 10% of the measurements are remaining. The instrument will not work after the indicated date, or after the indicated number of measurements.</p> <p>It is still possible to view measurements stored in the instrument memory and download data to a PC.</p> <p>Make sure that the ambient temperature is between +16 and +30°C.</p> <p>Wait for the Sensor to stabilize.</p>	A13	Select return and repeat the measurement and do not breathe through the patient filter during analysis.
		A14	Wrong passcode for NO scrubber exchange.
	Remove any sources of disturbance (such as cordless or mobile telephones, or gas emitting appliances). Wait for the Sensor to stabilize.	A20	Check that the ambient temperature is within specification. If necessary, shut the instrument down, move it to another location and restart the instrument.
	Wait for the Sensor to stabilize. < 4 minutes (countdown started).	A21	Remove any sources of disturbance (such as cordless/mobile telephones, or gas emitting appliances). When the instrument is ready try to repeat the measurement. If the alert persists, unplug and reconnect the power supply unit to restart the instrument.
		A22	Unplug and reconnect the power supply unit to restart the instrument.
		A23	Remove any sources of disturbance (such as cordless/mobile telephones, or gas emitting appliances). When the instrument is ready try to repeat the measurement. If the alert persists, unplug the power supply unit, remove and reinsert the Sensor, reconnect the power supply unit and restart the instrument.
		A24	Check that the supply voltage is within specification. If necessary replace the power supply unit.
<b>Alert codes</b>		<b>Connection alerts</b>	
		A30	Check the Bluetooth connection to the PC
		A31	Check the USB connection to the PC
<b>Instrument and Sensor expiration alerts</b>		<b>Instrument and Sensor expiration alerts</b>	
		A90	Instrument expiration date has passed or all instrument measurements have been used. It is still possible to view measurements stored in the instrument memory and download data to a PC.
		A91	Sensor expiration date has passed or all measurements on the Sensor have been used. Replace the Sensor.
Code	Action		
<b>User alerts</b>			
A10	Exhalation too strong. Select return and repeat the measurement with less force.		
A11	Exhalation too weak. Select return and repeat the measurement with greater exhalation force and exhale until the signal for completed exhalation is heard.		
A12	No exhalation detected. Select return and repeat the measurement and exhale into the instrument directly after inhalation.		

## Technical data

### Display buttons and symbols

Button	Description	Button	Description
	Enter patient ID		Shift to keypad 0-4
	Measurement modes		Shift to keypad 5-9
	Settings		Bluetooth on
	QC settings		Bluetooth off
	Standard 10 second exhalation FeNO measurement		US time and date settings
	Special 6 second exhalation FeNO measurement (research appl.)		ISO time and date settings
	Ambient NO measurement		Alert info
	Nasal measurement (research application)		Sound and top light
	Stored measurement results		Sound
	Nasal flow 2 ml/s (research application)		Top light
	Nasal flow 5 ml/s (research application)		NO scrubber setting
	Demonstration mode		Controlled power off
	Clock		
	Information		
	Forward		
	OK		
	Return		
	Previous		
	Next		
	Decrease		
	Increase		
to	Number		
	Backspace		

Symbol	Description	Symbol	Description
	Bluetooth activated		Wait
	USB cable connected		Wait, time
	Sleep mode		Wait for temperature stabilization
	Standard FeNO measurement		Wait for Sensor stabilization
	Special 6s FeNO measurement (research application)		
	Nasal measurement (research application)		
	Ambient NO measurement		
	Attach nasal olive to nostril (research application)		
	Insert Sensor		
	NO scrubber		
	NO scrubber		
	General warning		
	Instrument expiration warning		
	Sensor warning		
	Temperature warning		
	NO scrubber warning		

## Symbols explanation

Responsible manufacturer

The product meets the requirements of applicable European directives

ETL-Listed. Conforms to ANSI/AAMI ES60601-1 and certified to CAN/CSA C22.2 No. 60601-1.

Electrical safety Type B applied part 100-240 V ~50-60 Hz

The product should be recycled according to the local program for electronic equipment.

Consult instructions for use

Refer to instruction manual/booklet

Expiration date

Temperature limitation

For single use only

In vitro diagnostic medical device

Durability - period after opening


Hazardous material

Class II equipment

## Backplate

1 2 3 4 5 6 7

1. Responsible manufacturer
2. Refer to instruction manual/booklet
3. In vitro diagnostic medical device
4. Type B applied part
5. The product should be recycled according to the local program for electronic equipment
6. The product meets the requirements of applicable European directives
7. ETL-Listed. Conforms to ANSI/AAMI ES60601-1 and certified to CAN/CSA C22.2 No. 60601-1

## Base label

1. Instrument expiration date (year:month)
2. Instrument article number
3. Instrument serial number

## Dimensions and weight

Height: 230 mm

Width: 128 mm

Depth: 96 mm

Weight of instrument including Sensor: 0.8 kg

## Electrical data

Electrical safety classification: The equipment complies with the requirements according to IEC 61010-1 and 60601-1 Class II Type B.

### NIOX MINO instrument

Secondary voltage (external power adapter): 6 V ===

Power consumption: < 15 VA

### Power supply unit

Mains Voltage: 100-240 V ~50-60 Hz

## Noise level (stand by)

< 35 dBA, at a distance of 1 m

## Shelf-life

NIOX MINO® instrument: Minimum 3 years at time of delivery, or 3,000 measurements.

NIOX MINO Sensor: Maximum 12 months when installed in NIOX MINO or expiration date as stated on the Sensor, whichever comes first.

Always use the NIOX MINO Bag (closed) for transportation and storage of NIOX MINO when it is disconnected from the power outlet. (The NIOX MINO Bag is sold separately).

## Operating conditions

Ambient temperature: +16 to +30°C

Humidity: 20 to 60% RH (non-condensing)\*

Atmospheric pressure range: 700 to 1060 hPa

Ambient NO concentration: <300 ppb

Exhaled flow rate during measurements at 50 ml/s ± 10% BTPS (Body Temperature Pressure Saturated, 37°C, 1013 hPa, 100% RH)

\* The ambient RH in the instrument can be viewed in the information menu, see page 10 and may differ from RH measured outside the instrument.

## **Exhaled NO - performance data**

Performance data typical for temperature range +16 to +30°C, Humidity range 20 - 60% RH, Pressure 1013 hPa

Measurement range:

FeNO: 5 to 300 ppb

Lowest Detection Limit: 5 ppb

Determination by analyzing gas concentrations around and below the detection limit. 5 ppb was the lowest detectable level.

Measurement time: < 2 min

Temperature stabilization time: ≤ 30 min

Ambient NO measurement time: < 4 min

## **Linearity**

Squared correlation coefficient  $r^2 \geq 0.998$ , slope 0.95 -1.05, intercept  $\pm 3$  ppb.

Determination based on the pooled regression analysis from 10 instruments using standard gas reference samples at 7 different concentration levels covering the operating measurement range.

## **Precision**

< 3 ppb of measured value < 30 ppb,  
< 10% of measured value ≥ 30 ppb.

Expressed as one standard deviation for replicate measurements with the same instrument, using a certified gas concentration of Nitric Oxide reference standard.

## **Accuracy**

±5 ppb or max 10%.

Expressed as the upper 95% confidence limit, based on absolute differences for concentrations ≤ 50 ppb and relative differences for concentrations > 50 ppb, from certified gas concentration of Nitric Oxide reference standard.

## **Method comparison**

< 10 ppb for values ≤ 50 ppb,  
< 20 % for values > 50 ppb

Expressed as the difference between a NIOX MINO® FeNO value and the corresponding FeNO value measured with NIOX® instrument from Aerocrine.

## **Inhalation parameters**

Inhalation to TLC (Total Lung Capacity) before start of exhalation, to secure proper activation:

Minimum -3 cm H<sub>2</sub>O (trigger only)

## **Exhalation parameters**

Exhalation time:

Standard mode: 10 s (clinical use)

Special 6 s. mode: 6 s (research use)

All exhalations are to be performed at an exhalation pressure of 10 - 20 cm H<sub>2</sub>O, to maintain a fixed flow rate of 50 ±5 ml/s.

## **Memory capacity**

3000 measurements

## **Transport and storage**

### **NIOX MINO instrument (transport and storage in its unopened original package)**

Relative humidity range: 10 to 100%, including condensation.

Temperature range: -20 to +60°C

Atmospheric pressure range: 500 to 1060 hPa

### **NIOX MINO Sensor (transport and storage in original package)**

Relative humidity range: 10 to 100%, including condensation.

Temperature range: -20 to +60°C

Recommended temperature range: +2 to +25°C

Atmospheric pressure range: 500 to 1060 hPa

### **NIOX MINO, including sensor (transportation in NIOX MINO Bag)**

Relative humidity range: 10 to 95%, non-condensing

Temperature range: +2 to +25°C

Atmospheric pressure range: 700 to 1060 hPa

When transporting the instrument from one location to another, a prolonged stabilization period might be required before measurements can be performed.

## **Patient filter (mouthpiece)**

Disposable filter (NIOX Filter) to be changed for each new measurement session and for each new patient.

Bacterial, viral filter according to Medical Device Directive Class IIa.

## **Responsible manufacturer**

Aerocrine AB,  
Råsundavägen 18,  
SE-169 67 Solna  
Sweden

[www.aerocrine.com](http://www.aerocrine.com)  
[www.niox.com](http://www.niox.com)

## NIOX MINO® parts and accessories

### **Warning!**

When selecting an accessory for your NIOX MINO product please keep in mind that an accessory not recommended by Aerocrine AB may result in loss of performance, damage to your NIOX MINO product, fire, electric shock, injury or damage to other property. The product warranty does not cover product failure or damage resulting from use with non-approved accessories. Aerocrine AB takes no responsibility for health and safety problems or other problems caused by the use of accessories not approved by Aerocrine.

Bacterial, viral filter according to Medical Device Directive Class IIa.

### **Parts included in NIOX MINO package (Item No. 09-1000)**

NIOX MINO instrument (09-1100)  
NIOX MINO Power Supply Unit (incl. 4 supply adapters) (09-1005)  
NIOX MINO USB cable (09-1010)  
NIOX MINO QC plug (09-1015)  
NIOX MINO Visual Incentive (on USB stick) (09-1025)  
NIOX MINO NO scrubber (09-1300)  
NIOX MINO User Manual (EPM-000200)

### **Accessories**

NIOX MINO Test Kit 100 (03-4000)  
(contains: 1 NIOX MINO Sensor\* for 100 measurements and 100 NIOX® Filters\*\*)  
NIOX MINO Test Kit 300, (03-4001)  
(contains: 1 NIOX MINO Sensor\* for 300 measurements and 300 NIOX Filters\*\*)  
NIOX MINO Test Kit 500, (03-4007)  
(contains: 1 NIOX MINO Sensor\* for 500 measurements and 550 NIOX Filters\*\*)  
NIOX MINO Test Kit 1000, (03-4010)  
(contains: 1 NIOX MINO Sensor\* for 1000 measurements and 1050 NIOX Filters\*\*)  
NIOX MINO Data Management Program (09-1026)  
NIOX MINO Bag (09-1020)

\* NIOX MINO Sensor

Pre-calibrated disposable sensor for 100, 300, 500 or 1000 measurements. Shelf-life: Maximum 12 months when installed in NIOX MINO or expiration date as stated on the Sensor, whichever comes first. (Use only NIOX MINO Sensors with Aerocrine item number 03-2005, 03-2010, 03-2030, 03-2050 or 03-2100.)

\*\* NIOX Filter


Disposable filter to be changed for every patient.

## Connect NIOX MINO® to a PC using USB

**Caution!**

- A PC connected to the USB connector has to be certified for one of the standards IEC 60601-1, IEC 61010-1, IEC 60950 or comparable with safety extra low voltage on the USB ports.
- The connected PC should be placed out of reach from the patient. Do not, simultaneously, touch the connected PC and the patient.

1. Connect the USB cable from the instrument to the PC

2. The USB icon  is shown on the NIOX MINO display to indicate that a USB connection is established

## Connect NIOX MINO to a PC using Bluetooth


1. Select Settings

2. Select Bluetooth

3. The Bluetooth setting is shown

4. Use the button to activate Bluetooth

5. Select OK to accept the changes

6. Follow the instructions delivered with the PC or from the Bluetooth device manufacturer for how to set up a Bluetooth connection:

- Create a new Bluetooth connection
- Select NIOX MINO
- Enter passcode **0000**

## NIOX MINO® Visual Incentive

A USB memory stick is delivered with NIOX MINO containing the Visual Incentive program.

### Operating System Requirements:

- Windows® XP\*, Vista, Windows® 7 or 8
- Windows® 32 or 64 bit
- .NET Framework 3.5 or higher

\*Windows® XP SP2 (Service Pack two). Installation files included on USB memory stick.

### Install Visual Incentive


1. Connect the USB memory stick to a USB connector on the PC
2. Open the **NIOX MINO Visual Incentive x.x.x.x** folder
3. Doubleclick the **setup.exe** file
4. Follow the instructions on the screen to install the software on the PC

### Start Visual Incentive

#### Note!

Make sure to have a connection between NIOX MINO and the PC. Refer to the Connect NIOX MINO to PC sections.

1. Doubleclick the **NIOX MINO Visual Incentive** file
2. Select the **Measurement** tab

3. Select display mode
4. Optional: Enter a Patient ID (max 10 characters)  
The ID is displayed at the top of the NIOX MINO

Main screen and saved with the measurement in NIOX MINO

5. Click **OK**
6. Perform measurements using Visual Incentive, see Measurement section page 7


### Download instrument data

#### Note!

Only to be used after contact with Customer support.

Make sure to have a connection between NIOX MINO and the PC. Refer to the Connect NIOX MINO to PC sections.

1. Doubleclick the **NIOX MINO Visual Incentive** file
2. Select the **Instrument Log** tab

3. Check the Instrument log box
4. Click **Download** to download data from NIOX MINO
5. Select the directory where to save the Instrument Log file and click **Save**

#### Note!

The Instrument Log file is named with serial number, date and time when it was downloaded.

### Shut down Visual Incentive

1. Click the **X** icon to close the program

## Research applications

The research applications are not validated for clinical use.

### ***Special 6s FeNO measurement***

This is a research application with an exhalation time of 6 seconds.

1. Select Mode

2. Select Special mode

3. Make a measurement according to the instructions on page 7, with an exhalation time of 6 seconds.

### ***Nasal measurement***

Contact your local Aerocrine representative or Aerocrine Customer Support for further information.

## Cautions

- **Mobile phones and cordless phones might interfere with the instrument and should therefore be kept away from the instrument. Interference could make it impossible to perform a measurement.**
- **Make sure that you use the required measurement mode, otherwise incorrect exhaled NO results might be obtained.**
- Use of substances containing alcohol close to NIOX MINO® may cause erroneous measurement results.
- Do not reuse the NIOX® Filters.
- The NIOX MINO instrument might produce some heat during normal operation, the temperature could increase up to 5°C above the ambient temperature. Make sure that the base ventilation slots are not blocked. Do not place the instrument on a bed, sofa, carpet, or other soft surface.
- A maximum of 10 measurements/hour can be performed during continuous use. However, it will be possible to perform 20 measurements per hour if the instrument is allowed to rest for a minimum of 30 minutes after that.
- The Special 6 s. exhalation and nasal modes are intended for research use only.
- The NIOX MINO Sensor shall be kept in its original unopened package before installation. For transportation and storage conditions, refer to the corresponding section in this manual.
- The NIOX MINO Sensor is sensitive to changes in ambient temperature and humidity.
- Be careful when opening the Sensor can. The inside of the opening has sharp edges. Do not touch or clean the white Sensor membrane.
- Best performance is achieved if the ambient conditions are stable. Refer to the "Operating conditions" section. Keep the instrument away from windows, sunshine, radiators, stove or open fire in order to avoid unstable conditions. When transporting the instrument from one location to another, a prolonged stabilization period before measurement might be required. Refer to the recommended transportation conditions in the "Transport and Storage" section. Always use a NIOX MINO Bag for transportation. Keep the NIOX MINO with mounted Sensor connected to the power outlet whenever possible.
- NIOX MINO has been tested and found to comply with the limits for medical devices according EN/IEC 60601-1 Safety requirements for medical electrical systems and electromagnetic compatibility. These limits are designed to provide protection against harmful interference in a typical medical installation. However, because of the increased use of radio-frequency transmitting equipment and other sources of electrical noise emitters in the health-care and

home environments, such as base stations for radio, cellular/cordless telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast, it is possible that high levels of such interferences due to close proximity or strength of a source, may result in disruption of performance of the device. If abnormal performance is observed, it may be necessary to reorient or relocate NIOX MINO.

- A PC connected to the USB connector has to be certified for one of the standards IEC 60601-1, IEC 61010-1, IEC 60950 or comparable with safety extra low voltage on the USB ports.
- The connected PC should be placed out of reach from the patient. Do not, simultaneously, touch the connected PC and the patient.

## Clinical documentation

NIOX MINO follows the American Thoracic Society (ATS) and European Respiratory Society (ERS) 2005 equipment recommendations for measurement of exhaled NO.<sup>1</sup>

The recommendations are based on analysis of NO with the chemiluminescence method. However, NIOX MINO is using a different method (electrochemistry). One consequence of this is that the NIOX MINO Sensor does not need any field calibration\*. Agreement has been demonstrated for the mean of two valid exhaled NO measurements in NIOX and the first valid NO measurement in NIOX MINO.<sup>2-5</sup>

Supported by these data, one valid NO measurement should be considered sufficient, instead of two as recommended in the guideline.

<sup>1</sup> Am J Respir Crit Care Med 2005; 171: 912-930.

<sup>2</sup> Alving K, Jansson C, Nordvall L: Performance of a new hand-held device for exhaled nitric oxide measurement in adults and children. Respiratory Research 2006;7:67-5.

<sup>3</sup> McGill C, Malik G, Turner S. Validation of a hand-held exhaled nitric oxide analyzer for use in children. Pediatr Pulmonol. 2006 Nov;41(11):1053-7.

<sup>4</sup> Khalili B, Boggs PB, Bahna SL. Reliability of a new hand-held device for the measurement of exhaled nitric oxide. Allergy 2007;62:1171-1174.

<sup>5</sup> Menzies D, Nair A, Lipworth BJ. Portable exhaled nitric oxide measurement: Comparison with the "gold standard" technique. Chest 2007;131(2):410-4.

\* The Sensor is already calibrated during manufacturing and does not need to be re-calibrated.

## Vigilance

Aerocrine, as a medical device manufacturer, must have a system in place to report to health authorities, any adverse incidents that have occurred with its medical products, a so-called Vigilance system.

The purpose of the Vigilance system is to ensure the health and safety of patients, users and others using medical products by reducing the likelihood of the same type of adverse incident being repeated. This is achieved by immediate notification of incidents to enable corrective and preventive actions.

An adverse incident is defined as: Any malfunction or deterioration in the characteristics and/or performance of a device, or any inadequacy in the labelling or instructions for use which might lead to or might have led to the death of a patient or user or to a serious deterioration in his/her state of health.

Manufacturers of medical devices are obliged to report adverse incidents to national health authorities within 10 or 30 days, dependent on the severity of the incident.

Any user of Aerocrine's products, who experience an adverse incident related to the product, must therefore immediately report this to Aerocrine or Aerocrine's local market representative. The report can be made by e-mail, fax, or telephone.

The report should contain the following information:

- When and where did the incident occur?
- What product / accessory was involved?
- Was the incident related to instructions for use of the product?
- Was the risk foreseeable and clinically acceptable in view of potential patient benefit?
- Was the outcome adversely affected by a patient's pre-existing condition?

The report should be forwarded to Aerocrine's local representative in the country where the incident occurred.

## Guidance and manufacturer's declaration

### *Electromagnetic emissions*

#### ***Caution!***

NIOX MINO® has been tested and found to comply with the limits for medical devices according EN/IEC 60601-1 Safety requirements for medical electrical systems and electromagnetic compatibility. These limits are designed to provide protection against harmful interference in a typical medical installation. However, because of the increased use of radio-frequency transmitting equipment and other sources of electrical noise emitters in the health-care and home environments, such as base stations for radio, cellular/cordless telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast, it is possible that high levels of such interferences due to close proximity or strength of a source, may result in disruption of performance of the device. If abnormal performance is observed, it may be necessary to reorient or relocate NIOX MINO.

#### ***Warning!***

NIOX MINO should not be used adjacent to or stacked with other equipment.

Adjacent or stacked use together with NIOX® Nitric Oxide Monitoring System has been tested by Aerocrine and proven as acceptable.

The technical description on the electromagnetic immunity test levels and electromagnetic environment guidance can be found on [www.niox.com](http://www.niox.com).

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Information in this document is subject to change.

Amendments will be made available by Aerocrine AB as they occur.

© Copyright 2013: Aerocrine AB, SOLNA, Sweden.

Based on the company's intellectual property, Aerocrine develops and commercializes product for the monitoring of nitric oxide (NO) as a marker of inflammation, to improve the management and care of patients with inflammatory disease in the airways.

Patents:

For information about patents relating to Aerocrine products see [www.aerocrine.com](http://www.aerocrine.com)

## A NEW DIMENSION IN ASTHMA CARE

**Aerocrine**

Aerocrine AB, Råsundavägen 18, SE-169 67 Solna, Sweden  
Phone: +46 8 629 07 80, Fax: +46 8 629 07 81, E-mail: [info@aerocrine.com](mailto:info@aerocrine.com)  
[www.aerocrine.com](http://www.aerocrine.com)

NIOX MINO® is CE marked according to European In Vitro Diagnostic Device Directive 98/79/EEC and is cleared 510(k) by FDA, 101034, (US).

Copyright © 2013 Aerocrine AB, Solna, Sweden.

Aerocrine, NIOX MINO® and NIOX® are registered trademarks of Aerocrine AB.

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- Copyright © 2013 Aerocrine AB, Solna, Svezia.