

## Thermo Scientific Accela PDA Detector

The Thermo Scientific Accela High Speed PDA Detector uses Thermo Scientific LightPipe technology to ensure fast photo diode array detection with the highest sensitivity.

- Proprietary LightPipe™ technology for highest sensitivity
- Specialized 1 cm path length, 2  $\mu$ L LightPipe flowcell for U-HPLC separations
- 5 cm path length, 10  $\mu$ L LightPipe flow cell provides the highest sensitivity available for HPLC applications
- Fiber-optic beam shaper provides maximum peak resolution
- Easy maintenance with integrated wavelength validation
- Reduced flow cell volume improves chromatographic efficiency

### High Speed Sample Detection

The Thermo Scientific Accela PDA is optimized for the detection of high speed chromatographic separations. The short 1 cm path length, combined with the incredibly minimized cell volume of 2  $\mu$ L, renders this LightPipe flowcell as the premium choice for fast separations. The extremely low level of dispersion in the LightPipe enables the flowcell to retain excellent peak shape and chromatographic resolution.

The innovative fiber-optic beam shaper converts the circular beam of light exiting the flowcell into a vertical beam, emulating a slit while providing nearly 100% light throughput to the diodes.

The combination of the Thermo Scientific LightPipe technology with the fiber-optic beam shaper provides the high sensitivity and resolution needed for even the shortest chromatographic separations.

### Complementary to LC-MS

Compounds that are poorly ionized may be difficult to detect and analyze solely with a mass spectrometer. The Accela™ PDA detector is an ideal LC-MS complement for many applications.

### Easy Maintenance

The Accela PDA has an integrated holmium oxide solution-filled cuvette for validating wavelength accuracy. This holmium oxide in perchloric acid solution is immensely advantageous to the solid holmium oxide filters used by most manufacturers, because it provides wavelength verification across the entire spectral range.

### Software Control

The Accela PDA can be controlled as a component of the Accela High Speed LC system by the Thermo Scientific ChromQuest Chromatography Data System, the Thermo Scientific Xcalibur Mass Spectrometry Data System, and other multi vendor software solutions. Thermo Scientific software solutions provide computer control of the PDA to ensure automated data handling of the high speed system.



# Product Specifications

## Thermo Scientific Xcalibur Data System

Xcalibur™ is an extremely versatile, easy-to-use mass spectrometry data system with tools for 21 CFR Part 11 compliance. Whether the user is a lab manager or a technician, our Xcalibur software offers easy navigation through the process of instrument setup, sequence setup, and data acquisition. Data analysis is simplified using any of the three browser environments: Qual Browser, Quan Browser, and Library Browser. Reporting is straightforward and effortless with the intuitive XReport software. Thermo Scientific LCQUAN quantitation software provides a structured and secure workflow. When operating the Accela system as part of an LC-MS solution, use Xcalibur for complete instrument control and mass spectrometer data handling.

## Thermo Scientific ChromQuest Chromatography Data System

ChromQuest™ offers maximum flexibility and efficiency in setting up chromatographic methods and optimizing the performance of the Accela LC system. ChromQuest protects your results by ensuring data traceability and integrity with enhanced security features, electronic documentation, and features for 21 CFR Part 11 compliance. Sophisticated analytical tools facilitate high sample throughput and powerful data processing for liquid chromatography.

ChromQuest maximizes laboratory efficiency with project management for information storage and transfer. ChromQuest uses robust integration, flexible calibration routines, and advanced spectral analysis tools to accelerate data analysis and results generation. When operating the Accela LC system as a stand alone LC solution, use ChromQuest for complete instrument control and data handling.

## Certification and Qualification Services

Every PDA is tested at the factory prior to shipment to ensure that it meets or exceeds performance specifications. A Declaration of Conformity (DC) accompanies each PDA and documents test results. On-site field testing (optional) qualifies system performance.

### Ordering Information

Part Number	Description
60057-60010	Accela Pump
60057-60110	Accela 600 Pump
60057-60020	Accela Autosampler
60057-60140	Accela PDA with 1 cm LightPipe flow cell
60057-60141	Accela PDA with 5 cm LightPipe flow cell
60057-60060	Accela System Kit
60057-60070	Thermo Scientific Hypersil GOLD, 50 × 2.1 mm, 1.9 µm column

## Specifications

### Wavelength Range

190–800 nm at 1 nm increment

### Wavelength Accuracy

± 1 nm at 254 nm and 640 nm

### Wavelength Calibration

Using Holmium oxide solution

### Wavelength Resolution

1.2 nm (512 pixel array)

### Absorbance Non-Linearity

< 5% at 2.0 AU at 257 nm

### Absorbance Range

-2.0 to +4.0 AU, 20-bit resolution

### Drift

< 1 mAU/hr after warm-up at 254 nm at a stable temperature (± 1 °C)

### Light Source

Pre-aligned Deuterium and Tungsten lamps

### Rise Time

User selectable: 0, 0.02, 0.05, 0.1, 0.2, 0.5, 1, 2, 5, or 10 sec

### Scan Rate

User selectable: 0.5, 1, 2, 4, 5, 10, 20, 40, or 80 Hz

### Discrete Channels

Three wavelength selectable channels

### Short-Term Noise (50 mm flow cell)

< 6 µAU/cm at 254 nm at 5 Hz

(20 Hz scan rate)

< 1.5 uAU/cm at 254 nm at 5 Hz

(40 Hz scan rate); typical

< 0.3 uAU/cm at 254 nm at 5 Hz

(80 Hz scan rate); typical

### Warm-up Time

90 min required to meet noise and drift specifications

### Cell Dimensions

10 mm, 2 µL LightPipe flow cell

50 mm, 10 µL LightPipe flow cell

### Cell Pressure Range

0–1000 psi

### Analog Outputs

20-bit digital/analog conversion,

unattenuated at 10 mV/AU,

100 mV/AU, or 1.0 V/AU

### Diodes

512

### Diode Spacing

1.2 nm

### Remote Controls

Ethernet interface for PC-based software control

### Operating Temperature

10–30 °C

### Ambient Environment

10–40 °C, 5–95% relative humidity (non-condensing)

### Dimensions

18 × 36 × 47 cm (H × W × D)

### Weight

19.5 kg

### Power Requirements

100/115 or 230 VAC, 50/60 Hz, 225 VA max

### Product Certification

EMC

CE

TUV CUS

FCC (EMI)

[www.thermoscientific.com](http://www.thermoscientific.com)

©2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Africa-Other** +27 11 570 1840  
**Australia** +61 3 9757 4300  
**Austria** +43 1 333 50 34 0  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** +86 10 8419 3588  
**Denmark** +45 70 23 62 60

**Europe-Other** +43 1 333 50 34 0  
**Finland/Norway/Sweden** +46 8 556 468 00  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Russia/CIS** +43 1 333 50 34 0  
**South Africa** +27 11 570 1840

**Spain** +34 914 845 965  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752



Thermo Fisher Scientific,  
San Jose, CA USA is ISO Certified.

PS62125\_E 08/11M

**Thermo**  
SCIENTIFIC