REF CW-USLS

# CyberWand™ Dual Ultrasonic Lithotripsy System

## Instructions for Use

EN

(FR)	Système de lithotritie ultrasonique à double sonde CyberWand	Mode d'emploi
(DE)	CyberWand Dual-Ultraschall-Lithotripsiesystem	Gebrauchsanweisung
(IT)	Systema ultrasonico a doppia sonda CyberWand per litotripsia	Istruzioni per l'uso
(ES)	Sistema dual de litotricia ultrasónica CyberWand	Instrucciones de uso
(SV)	CyberWands dubbla ultraljudslithotripsi-system	Bruksanvisning
(NL)	CyberWand dubbel ultrasoon lithotripsiesysteem	Gebruiksaanwijzing

Products and information available from

**Gyrus ACMI, Inc. • 136 Turnpike Road • Southborough, MA 01772-2104 • USA**  
**1-888-524-7266 or 1-763-416-3000 • www.gyrusacmi.com**

CyberWand is manufactured by

**Cyberonics Inc. • Erie, PA 16510 • USA**

**0123 EC REP**

MediTech Strategic Consultants B.V. • Maastrichterlaan 127-129 • 6291 EN Vaals • The Netherlands • Tel: +31.43.306.3320 • Fax: +31.43.306.3338

**CAUTION: Federal (US) Law restricts this device to sale by or on order of a physician.**

## Table of Contents

<b>Return/Repair/Warranty .....</b>	<b>3</b>
<b>WARNINGS AND CAUTIONS.....</b>	<b>3</b>
<b>Section 1.0 Unpacking and Initial Inspection .....</b>	<b>5</b>
<b>Section 2.0 General Description and Intended Use .....</b>	<b>6</b>
<b>2.1 Specifications .....</b>	<b>6</b>
<b>2.2 Safety Specifications .....</b>	<b>7</b>
<b>2.3 Front Panel.....</b>	<b>8</b>
<b>2.4 Back Panel .....</b>	<b>8</b>
<b>2.5 Footswitch (REF CW-USLFS) .....</b>	<b>8</b>
<b>2.6 Transducer (REF CW-USLT) .....</b>	<b>9</b>
<b>2.7 Probe Set (REF CW-USLRBP and CW-RPBX).....</b>	<b>9</b>
<b>2.8 Cleaning Stylet (REF CW-USLCS) and Wrench (REF CW-USLW) .....</b>	<b>9</b>
<b>2.9 Symbol Identifications on labeling and unit .....</b>	<b>10</b>
<b>Section 3.0 Assembling the Unit .....</b>	<b>10</b>
<b>3.1 Probe/System Tuning .....</b>	<b>11</b>
<b>Section 4.0 Operating Instructions .....</b>	<b>11</b>
<b>Section 5.0 Cleaning, Disinfection and Sterilization.....</b>	<b>12</b>
<b>5.1 Generator .....</b>	<b>12</b>
<b>5.2 Probes .....</b>	<b>12</b>
<b>5.3 Transducers .....</b>	<b>12</b>
<b>5.4 Steam Sterilization.....</b>	<b>13</b>
<b>5.5 Ethylene Oxide Sterilization.....</b>	<b>13</b>
<b>5.6 STERRAD™ Sterilization .....</b>	<b>13</b>
<b>5.7 Caution.....</b>	<b>13</b>
<b>Section 6.0 Operating and storage environment .....</b>	<b>14</b>
<b>Section 7.0 Preventive Maintenance.....</b>	<b>14</b>
<b>7.1 General Checks .....</b>	<b>14</b>
<b>7.2 Periodic Safety Checks .....</b>	<b>14</b>
<b>Section 8.0 Troubleshooting .....</b>	<b>14</b>
<b>8.1 Troubleshooting .....</b>	<b>14</b>
<b>8.2 Fuse Replacement.....</b>	<b>14</b>
<b>Section 9.0 Disposal of Waste Electrical Equipment.....</b>	<b>15</b>
<b>Section 10.0 Repair .....</b>	<b>15</b>

## Return/Repair/Warranty

### Return/Repair

**All returns must have prior authorization.** To comply with OSHA Bloodborne Pathogen Standard 29 CFR 1910.1030 regulations and U.S. Postal and Transportation law, all used medical devices returned for repair or replacement must be properly cleaned and decontaminated with a chemical germicide that has been cleared for use as a "Hospital Disinfectant." To ensure that the product has been properly decontaminated, a signed Decontamination Certificate should be enclosed in the package.

All medical devices returned to Gyrus ACMI for any reason must be shipped in accordance with Gyrus ACMI's return procedures (available upon request) and all applicable regulations. To obtain a return authorization number (RMA), return address, and instructions, please call Gyrus ACMI Customer Service at 1-888-524-7266 or 1-763-416-3000.

Please provide the product part number and lot/serial number. Product should be returned in its original packaging (when possible) and marked with the RMA number on the exterior of the package.

### Limited Express Warranty

SHOULD THE PRODUCT BECOME INOPERABLE DURING NORMAL AND PROPER USE IN ACCORDANCE WITH APPLICABLE INSTRUCTIONS AND WITHIN THE TIME FRAME SPECIFIED BELOW FROM THE DATE OF SHIPMENT, GYRUS ACMI WILL REPAIR OR REPLACE THE PRODUCT, AT ITS SOLE OPTION, AT NO CHARGE. GYRUS ACMI MAKES NO OTHER WARRANTIES WITH RESPECT TO THE PRODUCT AND EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AS TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER MATTER. IN NO EVENT SHALL GYRUS ACMI BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES. IN NO EVENT SHALL GYRUS ACMI BE LIABLE FOR ANY BREACH OF WARRANTY IN ANY AMOUNT EXCEEDING THE PURCHASE PRICE OF THE PRODUCT.

This warranty applies only to the original purchaser and will be voided if the product(s) are serviced or repaired by anyone other than Gyrus ACMI or an organization duly authorized by Gyrus ACMI for such purpose.

Ultrasonic Lithotripters .....	1 year
Accessories .....	90 days

## WARNINGS AND CAUTIONS

The following warnings and cautions apply to the use, care, and/or maintenance of the Gyrus ACMI® CyberWand Ultrasonic Lithotripsy System.

### Warnings

*(Indicate that danger to life or health can result from misusing the equipment)*

1. This device should be operated only by or under the direct supervision of a physician experienced in ultrasonic lithotripter procedures. The user should be thoroughly familiar with the operation of this device prior to use.
2. This device is designed to function with single-use ultrasonic lithotripsy probes designed specifically for the CW-USLS generator. **Probes are single-use should not be reprocessed for reuse.**
3. Perform the prescribed inspections and operational tests prior to first use and regularly thereafter to ensure continued satisfactory performance. Do not use a CyberWand Ultrasonic Lithotripter that fails to meet the criteria stated in the labeling or that has been damaged. Otherwise, injury to patient, personnel and/or an adverse effect on the procedure could result.
4. When used in an environment containing flammable gasses and high oxygen concentrations, be aware of the potential for a flammable event. Take routine cautions to prevent this kind of event. The product is not explosion proof.
5. A potential shock hazard exists when the instrument cabinet is opened. Refer all repairs to authorized service personnel.
6. Do not adjust electronic circuitry. Contact your local Gyrus ACMI representative or call Gyrus ACMI customer service at 1-888-524-7266 or 1-763-416-3000.

EN

7. A continuous outflow of irrigation from the body cavity (operative site) is necessary to prevent excessive pressure buildup.
8. Never attempt to service the device when it is connected to a power source. Hazardous voltages inside the device may cause severe electrical shock. Disconnect the power cord before servicing.
9. As a Type BF applied part, this instrument must never be applied directly to the heart and must not be used in procedures involving cardiac observation or surgery. Do not use electrically active devices, such as electrosurgical units, on the patient during lithotripsy.
10. The housing of the generator must be properly grounded to ensure safe operation. To reduce the risk of electric shock, be sure to connect the power cord to a properly grounded 3-pin receptacle. Do not use a 3-pin /2-pin adapter.
11. The receptacle to which this product is connected should have enough power capacity to fulfill the total power requirements of all equipment connected to it. If the capacity is insufficient, the circuit breaker of the medical facility may be tripped, cutting off the power supply to all equipment connected to the same power source.
12. To minimize the risk of electric shock and generator damage, keep all liquids away from the generator. If liquid is spilled into this product, immediately stop the procedure and contact Gyrus ACMI. If the power cord gets wet, fully dry it before use. Otherwise the user may receive an electric shock.
13. This product is provided with air vents for heat release on the rear panel and bottom. Do not block the air vents.
14. This product is a precision device; handle it with care. Avoid rough or violent handling, which may cause equipment damage.
15. If you suspect any abnormality while using this product, immediately stop the procedure, and if the trouble persists, contact Gyrus ACMI.
16. Keep the tip of probe in the field of view at all times when the CyberWand system is being used.
17. This equipment has been tested and found to comply with the EMC limits for the Medical Device (EN55011 Class A and EN 60601-1-2 2004). The limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with another device, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following methods:
  - a. Reorient or relocate the receiving device.
  - b. Increase the separation between instruments.
  - c. Connect the equipment into an outlet on a circuit different from that to which the other devices are connected.
  - d. Consult the manufacturer or field service technician for help.
18. As with any lithotripter device, adverse events may occur, such as, bleeding and/or damage in the kidney, and other potential events associated with endoscopy.
19. If a probe breaks during a procedure, use a grasping device to remove it. Always maintain a spare probe set to finish the procedure.

## Cautions

*(Indicate that equipment or other property may be damaged or may malfunction by misuse)*

1. Damage (mechanical and electrical) may result if the unit is dropped or struck against another object.
2. Thoroughly inspect all electrical cables and probes before each use. Do not use if there is any evidence of deterioration. Replace if any damage or excessive wear is observed.
3. Using the wrench provided, ensure that the probe is tightened on to the transducer before operation. Before attaching, be sure probe and transducer threads are clean.
4. Do not allow fluids to enter the cabinet. Damage can result.
5. Continuous irrigation and suction must be applied to prevent overheating of the transducer and probe. Clogs within the probe or transducer must be cleared immediately.
6. Do not allow the probe to press against the scope during the procedure. Pressure on the probe may stall the transducer and heat the probe within the scope's working channel.

## Section 1.0 Unpacking and Initial Inspection

Proper care and maintenance are critical for safe and effective operation of sophisticated medical equipment. We recommend careful inspection of all equipment upon receipt and prior to each use as a safeguard against possible injury to patient or operator.

1. To avoid inadvertent damage, study this manual thoroughly before handling assembling, testing, using or cleaning the CyberWand unit.
2. Examine the shipping carton and instrument for signs of damage. Any breakage or other apparent damage should be noted, the evidence retained, and the carrier or shipping agency notified.
3. Verify that the shipping carton contains the items listed below:
  - a. CW-USLG Generator
  - b. CW-USLT Transducer
  - c. CW-USLPC Power cord
  - d. CW-USLFS Foot switch
  - e. CW-USLRBP Renal/bladder probe sets (2) (non-sterile)\*  
\* CW-RBPBX 3-pack, Sterile Renal/Bladder Probes sold separately
  - f. CW-USLW Wrench
  - g. CW-USLCS Cleaning Stylet
  - h. 80834 Instructions for Use
4. Notify Gyrus ACMI Customer Service immediately if any damage or discrepancies are noted.

## Section 2.0 General Description and Intended Use

The Gyrus ACMI CyberWand Ultrasonic Lithotripsy System is an electromechanical device capable of fragmenting calculi and aspirating stone debris. It is classified as Class 1 equipment with a Type BF applied part, per IEC 60601-1. This product conforms to EU Council Directive 93/142/EEC for medical devices.

The CyberWand Ultrasonic Lithotripsy System is intended to be used to fragment urinary tract calculi in the kidney, ureter, and bladder.

The following items are included in the CyberWand System:

REF	Description	
CW-USLS	CyberWand Dual Ultrasonic Lithotripsy System: Includes Generator, Transducer, Footswitch, Power Cord, 2 Renal/Bladder probe sets, Wrench and Cleaning Stylet (Cleaning rod)	
	CW-USLG	Generator (100 ~ 240 VAC)
	CW-USLT	Transducer (Hand piece)
	CW-USLRBP	Renal/Bladder Probe Sets (2) 3.8 mm
	CW-USLFS	Footswitch
	CW-USLPC	Power Cord
	CW-USLW	Wrench
	CW-USLCS	Cleaning Stylet

NOTE: The above items can also be ordered separately from Gyrus ACMI.

### 2.1 Specifications

#### Generator

Power supply	Voltage	100 ~ 240 VAC
	Frequency	50/60 Hz
Physical	Dimensions	31x 12x 28 cm W x H x D
	Weight	6 kg
Classification (Electro-Medical)	Protection against electric shock	Class I
	Type of applied part	Type BF
Ultrasonic Energy	Frequency	21 kHz±1 kHz
	Maximum output	100 W

#### Footswitch

Physical Specifications	Dimension	20 x 4 x 15 cm W x H x D
	Weight	3.18 kg, 7lbs
	Length of Cord	3 m
Classification	Degree of water resistance	IPX8 watertight type (except plug)

#### Transducer

Housing Dimensions	Diameter	4 cm
	Length	17 cm
	Weight	300 g
Cable	Length	3 m
Sterilization	Steam or Ethylene Oxide	

### Probe Set (Single Use Only)

<b>Physical Dimensions</b>	Outside Diameter	3.75 mm
	Working Length	40 cm
	Tip Excursion	80 ± 20 µm
<b>Material</b>	Probe (Fixed & floating)	300 Series Stainless steel
	Free mass ring	Stainless steel
	Spring	Stainless steel
<b>Sterilization</b>	Steam or Ethylene Oxide	

## 2.2 Safety Specifications

Meets IEC 60601-1, UL60601-1, CSA 22.2 No 601-1.

EN60601-1-2, EN55011.

### Electromagnetic Emission Declaration

The CyberWand system is intended for use in the electromagnetic environment specified below. The customer or user of the CyberWand system should assure that it is used in such an environment

Emission Test	Compliance	Electromagnetic Environment
RF emission CISPR 11	Group 1	The CyberWand system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electric equipment.
RF emission CISPR 11	Class A	
Harmonic emission IEC 61000-3-2 and 61000-3-2	Complies	The CyberWand system is suitable for use in all establishments, other than domestic establishments and those directly connected to a low voltage power supply network that supplies building used for domestic purposes.


### Electromagnetic Immunity Declaration

The CyberWand system is intended for use in the electromagnetic environment specified below. The customer or user of CyberWand system should ensure that it is used in such an environment



Immunity Test	IEC 60601 test level		
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	C	Floors should be wood, concrete or tiled (ceramic). If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply line	A	Power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV differential ± 2 kV common mode	A	Power quality should be that of a typical commercial or hospital environment
Voltage dip, short interruptions and voltage variations on power sup- ply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycle 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle <5% $U_T$ (95% dip in $U_T$ ) for 5 sec	A	Power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during power mains interruptions. It is recommended that the CyberWand system be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	A	Magnetic fields in the network frequency should correspond to the typical values found in business and hospital environments.

Note  $U_T$  is the A.C mains voltage prior to application of test level

## General Description of CyberWand Unit and accessories

EN

### 2.3 Front Panel

#### 1. Power Indicator

Power LED is illuminated when power is applied to the system.  
I indicates power is "ON," O indicates power is "OFF."

#### 2. Tune Indicator

Tune LED, flashing when power is applied, indicates that the CyberWand system is not tuned and is not ready for use.

#### 3. Ready Indicator

When Ready LED is steady on, the system is tuned and ready for use.

#### 4. Fault Indicator

Fault LED is illuminated when either the footswitch or the transducer is disconnected.

#### 5. Footswitch Receptacle

Align the white dots to connect the footswitch to the generator

#### 6. Transducer Receptacle

Align the red dots to connect the transducer to the generator.

### 2.4 Back Panel

#### 1. Cooling Fans (2) with finger guard

#### 2. Power Receptacle

UL Listed 3-prong hospital grade plug

#### 3. Fuse Holder

**CAUTION:** Do not block fan during operation; unit will overheat.

### 2.5 Footswitch (REF CW-USLFS)

#### 1. LARGE Pedal Mode (gray)


Depressing LARGE pedal causes the CyberWand system to produce maximum pulse energy to fragment large stones into smaller stone pieces.

#### 2. SMALL Pedal mode (green)

Depressing SMALL pedal causes the CyberWand system to produce a pulse energy level that keeps the probe in contact with the smaller stone fragments allowing greater fragmentation. Small stones are repulsed less when the SMALL pedal is used.

#### 3. Connector plug

#### 4. Cable


## 2.6 Transducer (REF CW-USLT)

1. Transducer cover (Nose Cone)
2. Transducer housing contains crystals that produce the piezo-electric ultrasonic energy.
3. Suction Outlet
4. Cable
5. Connector plug

EN

## 2.7 Probe Set (REF CW-USLRBP and CW-RPBX)

1. Fixed (Inner) Probe
2. Free Mass Ring
3. Floating (outer) Probe
4. Return Spring

## 2.8 Cleaning Stylet (REF CW-USLCS) and Wrench (REF CW-USLW)

NOTE: these are very important for each case.


	I	O	~
1	2	3	4
	6	7	<b>SN</b>
5			8
	10		
9		11	
<b>Rx Only</b>	<b>Distributed by</b>		
13	14	15	

## 2.9 Symbol Identifications on labeling and unit

1. Type BF applied part
2. Power ON
3. Power OFF
4. Alternating current
5. Fuse
6. Footswitch
7. Transducer output
8. Serial Number
9. Instructions for Use
10. Caution, Consult Accompanying Documents.
11. Tested by TUV Product Service to UL60601-1: 2003, CSA 22.2 No. 601.1-M90
12. Quantity symbol: Numeral in symbol (in place of #) indicates the quantity of units in package.
13. Caution: Federal (US) Law restricts this device to sale by or on the order of a physician.
14. Distributed by
15. Reference Directive 2002/96/EC.

## Section 3.0 Assembling the Unit

1. Disinfect and sterilize transducer, probe, cleaning stylet, and wrench according to the instructions in Sections 5.2, 5.3, and 5.4. Note: The stylet must be removed from within the probe prior to disinfection and sterilization.
2. Disconnect the transducer cover (nose cone) from transducer body to expose the probe attachment screw threads
3. Firmly wipe the threads in the transducer and the mating screw end of the probe with a sterile 4" x 4" gauze pad to remove any moisture and to permit a firm seating.
4. Check to ensure that the Probe set parts are in the correct order before attaching probe to transducer.
  1. Floating (outer) Probe
  2. Return Spring
  3. Fixed (Inner) Probe
  4. Free Mass Ring


First, finger-tighten the Fixed (inner) Probe to the Transducer. Then use the wrench to be sure that the Fixed Probe is screwed tightly onto the Transducer.

**CAUTION:** Probe must be secured tightly in transducer to assure good coupling. Loosening of the probe during operation will decrease or halt lithotripsy action of the probe and may cause damage.

**Caution:** If the probe is not properly assembled, seated or tightened into the transducer, it will not perform properly. Check to ensure that the probe is assembled and seated correctly prior to use. If the probe does not appear to work effectively, use a cleaning stylet to make sure it is not obstructed by debris.

Do not use a damaged probe. Probes (REF CW-USLRBP) are supplied non-sterile\* and must be sterilized prior to use. CyberWand probes are designed as single use devices. DO NOT REUSE. BREAKAGE MAY OCCUR.

\* Sterile 3-pack Probe Sets (REF CW-RBPBX) are available, sold separately.

5. Remove the connector cap covering pins on the transducer cable and connect the transducer to the generator by aligning the red dots on the transducer connector with those on the transducer receptacle on the front panel (see photo at left). Push straight in. Check condition of cord and transducer body. Should repair be necessary, contact Gyrus ACMI Customer Service.
  6. Connect aspiration tubing (OD 9mm, ID 6mm) to the port at the rear of the transducer and connect the other end of the tubing to an appropriate suction source.
  7. Connect the footswitch to the generator by aligning the white indicators on the footswitch connector with those on the footswitch receptacle on the front panel (see photo at left). Push straight in.
- NOTE: It is recommended that the footswitch be protected from contamination by containing it in a plastic bag.
8. Turn the CyberWand system power ON and verify that the green POWER and TUNE LEDs on the front panel are illuminated. The TUNE LED should be flashing, indicating that the system is not tuned.
  9. Verify that both fans in the back of the generator are running and that nothing is blocking/interfering with air intake.

### 3.1 Probe/System Tuning

Depress either the LARGE or SMALL pedal on the footswitch for approximately 2 seconds or until the tuning cycle begins, then take foot off the pedal and allow the tuning cycle to run until complete. The probe vibrates in response to the generator computer tuning program. After pressing the footswitch, the flashing LED will become solid for the duration of the tuning cycle. Once tuned, the Ready LED will illuminate, the CyberWand system is then tuned and ready for use.

**NOTE:** During auto tuning, the tuning sound keeps changing. At the end of tuning, the sound should be strong. If the sound is weak, it means it is not tuned. Retighten or change the probe set if needed and retune again. In order to re-tune, it is necessary to turn the generator power off and then back on. If the problem persists, change the transducer and re-tune.

## Section 4.0 Operating Instructions

This section is intended to provide the operator with instructions on the use of the CyberWand system, not on how to perform ultrasonic lithotripsy. **Be sure a backup probe and backup transducer are sterilized and available before beginning any ultrasonic lithotripsy procedure.**

1. The tip of the probe should be visible at all times while in use.
2. The probe tip should be in contact with the stone while activated.  
Note: Position the probe gently against the stone for best fragmentation.
3. Do not contact tissue with the tip or sides of the probe; tissue damage could result.
4. Depending on the stone size, select LARGE stone mode or SMALL stone mode by depressing right or left foot pedal, respectively.

#### Notes:

LARGE Stone Mode: Depressing LARGE pedal (gray) causes the CyberWand system to produce maximum pulse energy to fragment large stones into smaller stone pieces.

SMALL Stone Mode: Depressing SMALL pedal (green) causes the CyberWand system to produce a pulse energy level that keeps the probe in contact with the smaller stone fragments allowing greater fragmentation. Small stones are repulsed less when the SMALL pedal is used.

5. Cleaning stylets are provided to clear the probe should obstructions occur. Should the ultrasonic action lessen during use, remove the suction tubing from the transducer. Pass the cleaning stylet through the transducer to the probe tip to dislodge any obstruction.
6. Continuous irrigation and suction should be used to avoid overheating the probe and transducer.  
**Note:** Normal physiological saline is the recommended irrigation for ultrasonic lithotripsy.

7. Continue fragmentation and suctioning of stone debris until remaining pieces are deemed small enough to pass, or patient is stone free.

**Note:** if the performance is weak, or after long periods of running, the system generator should be turned off and then back on and re-tuned as described in Section 3.1 above.

**Warning:**

- The tip of the probe must be kept visible at all times when the transducer is activated.
- Do not contact mucosal tissue with the tip of the probe as tissue damage may occur.

**Caution**

- A back-up transducer and probe set should be sterilized and available prior to beginning a procedure.
- Continuous irrigation and suction is necessary for proper operation of the CyberWand system.

EN

## **Section 5.0 Cleaning, Disinfection and Sterilization**

The following cleaning, disinfection, and sterilization procedures are recommended by Gyrus ACMI. Use of any procedures not expressly recommended by Gyrus ACMI may adversely affect or damage Gyrus ACMI devices and may void Gyrus ACMI's Limited Express Warranty.

Follow all applicable bloodborne pathogen procedures as indicated by OSHA and/or your hospital requirements when cleaning, disinfecting, and/or sterilizing instruments and accessories.

### **5.1 Generator**

1. Turn off power to the unit and remove the power cord before cleaning.
2. Wipe generator cabinet, front panel, power cord, and footswitch with a sponge (or towel) lightly moistened with a mild soap solution. Follow with a sponge (or towel) lightly moistened with clear tap water. Dry completely.
3. Avoid using chemical or abrasive cleaning agents, which can damage the finish.
4. Wipe the transducer cable with a gauze pad (or towel) dampened with a mild soap solution. Remove the soap solution with a gauze pad (or towel) lightly moistened with clear tap water. Dry completely.

### **5.2 Probes**

CyberWand probes are intended for single use. Discard after each procedure. Reuse may result in breakage of the probe.

### **5.3 Transducers**

After each use or prior to placing in a disinfectant or sterilizer, carefully inspect transducer and cable for tears, cracks or other signs of damage. Do not use damaged equipment.

1. Remove the transducer connector from the generator and place the connector cap over pins.
2. Remove the cover (nose cone) of the transducer (REF CW-USLT) by unscrewing counterclockwise.
3. Fully immerse both the transducer and its cover (nose cone) in an enzymatic cleaning solution that is formulated to dissolve proteinaceous material. Follow the manufacturer's written instructions for dilution, temperature, and time. Warm water (100–120°F or 38–49°C) should be used; water above these temperatures can make proteins difficult to remove from the instruments. Do not use dishwashing detergents, skin cleansers, or any other cleaning agent not designed to dissolve protein. Avoid using chemical or abrasive cleaners that could damage the finish.
4. Scrub instrument with a soft brush.
5. Rinse all components thoroughly in warm tap water to remove all traces of cleaning solution; then towel dry.

## 5.4 Steam Sterilization

Thoroughly clean, rinse, and dry transducers, probes, cleaning stylet, and wrench prior to sterilization. CyberWand components may be sterilized in a steam sterilizer controlled within the following parameters:

		Probe Sets Prevacuum	Transducer Prevacuum (polypropylene double wrapped)
<b>Temp:</b>	°F	270	270
	°C	132	132
<b>Vacuum</b>		20-28 in Hg	20-28 in Hg
<b>Pressure</b>		25-28 psig	25-28 psig
<b>Exposure (minimum)</b>		4 minutes	60 minutes

### 5.4.1 Caution

Never rinse the CyberWand components with cooled water for speedy cooling after sterilization. Sudden changes in temperature may damage the components.

NOTE:

- Both sterile (REF CW-RBPBX) and non-sterile (REF CW-USLRBP) probes are identical.
- All sterile and non-sterile probes are single-use.
- Sterile probes can be steam sterilized with the same method as non-sterile probes if opened and unused.
- See the Instructions for Use for the sterile probes for more details.

## 5.5 Ethylene Oxide Sterilization

The transducer, probe, cleaning stylet, and wrench can be sterilized in an ethylene oxide (EO) sterilizer.

1. Thoroughly clean, rinse, and dry prior to sterilization.
2. The devices can be placed in an ethylene oxide (EO) sterilizer controlled within the following parameters:

**EO (100%)**

**Temp °F (°C):** 131 (55)

**Rel. Humidity:** > 35%

**Vacuum:** 2–3 psia (0.14 – 0.20 bar)

**Preconditioning Time:** ≥ 30 minutes

**Gas Conc.:** 725–750 mg/l

**Exposure:** ≥ 60 minutes

3. After EO exposure, aerate for a minimum of 12 hours at 131°F (55°C) in a mechanical aerator.

## 5.6 STERRAD™ Sterilization

**DO NOT USE STERRAD TO STERILIZE TRANSDUCER.** STERRAD sterilization has not been tested for, submitted to, or cleared by the FDA for transducer, probe, or component sterilization.

### 5.7 Caution

CyberWand transducer components can be processed until they exhibit evidence of deterioration, such as cracking on the transducer or broken insulation on the cable. At such time they must be removed from use and discarded. Please refer your facility's procedure for medical waste disposal. Check carefully for cracking and broken cable after each sterilization.

## Section 6.0 Operating and storage environment

<b>Operating environment</b>	Ambient temperature	10 – 40°C
	Relative humidity	30% – 70%
	Air Pressure	70 – 106 kPa
<b>Storage environment</b>	Ambient temperature	-25 – +55°C
	Relative humidity	10% – 85%
	Air Pressure	70 – 106 kPa

## Section 7.0 Preventive Maintenance

### 7.1 General Checks

The operator should ensure the generator and all cables and all components on transducer are undamaged prior to beginning any procedure.

### 7.2 Periodic Safety Checks

**The safety checks should be performed at least every 24 months.** The inspection can only be performed by a person who has adequate training, knowledge, and practical experience. All the equipment used in safety check must be calibrated. Gyrus ACMI recommends returning the unit for safety inspection or contact your local representative for information about maintenance contracts.

## Section 8.0 Troubleshooting

### 8.1 Troubleshooting

Irregular Description	Possible cause	Solution
Power indicator LED is not illuminated.	The power cord is not securely plugged into a receptacle	Connect power cord into receptacle. Note: Ensure the power switch is at OFF position before plugging in.
	The power cord is not securely plugged into generator	Connect power cord into generator. Note: Ensure the power switch is at OFF position before plugging in.
	Power switch is not turned ON	Turn the switch ON.
	The fuse has blown	Replace fuse as described in 8.2 Fuse Replacement below.
Fault Indicator light is illuminated.	Footswitch is not plugged into generator	Plug footswitch into generator. Note: Ensure the power switch is at OFF position before plugging in.
	Transducer is not plugged into generator	Plug transducer into generator. Note: Ensure the power switch is at OFF position before plugging in.
Poor performance or low power	<ul style="list-style-type: none"> <li>• Mis-tuned probe</li> <li>• Loose Probe</li> </ul>	<ul style="list-style-type: none"> <li>• Re-tune system (Section 3.1).</li> <li>• Retighten probe and re-tune system.</li> </ul>

### 8.2 Fuse Replacement

#### Warning:

**Ensure that the power cord is disconnected from the instrument before opening the fuse holder.**

1. Remove the fuse holder with a flat screwdriver, to expose fuse
2. Replace the blown fuse with a new fuse (TYPE T 250V 5A).
3. Replace the fuse holder securely into the generator.

## Section 9.0 Disposal of Waste Electrical Equipment

It is the sole responsibility of the user/user facility to properly dispose of equipment in a safe manner consistent with all compliance guidelines, applicable policies and procedures, statutes, and regulations that govern waste electrical equipment discarded by their facility in conjunction with all applicable local, city, county, town, state, or country specific codes for waste management.

The safe disposal of waste electrical equipment should be undertaken only at the end of the equipment's useful life (i.e., the device is nonfunctional & nonrepairable), and should be undertaken in an environmentally conscious manner minimizing environmental impact.

It is the sole responsibility of the user/user facility to determine whether waste electrical equipment should be collected and disposed separately from or together with unsorted municipal waste. Public or private solid waste collection and disposal facilities that offer various forms of recovery (dismantling, reuse, and recycling, where possible) should be selected reducing the amount of waste being discarded and providing the highest level of protection for human health and the environment.

It is the sole responsibility of the user/user facility to determine the conditions under which the waste electrical equipment is to be collected, noting that items contaminated by the user may require confirmed decontamination by means of the best available treatment options.

## Section 10.0 Repair

This product does not contain any user serviceable parts. Do not disassemble, modify or attempt to repair it; patient or user injury and/or equipment damage can result.

Some problems that appear to be malfunctions may be correctable by referring to Section 8.0 Troubleshooting.

A technical manual is available upon request for Biomedical Technicians.

NOTE: Do not open transducer or generator systems. A technical manual is to help in troubleshooting the problems or to help in the diagnosis if the part or system is returned for repairs.

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