PEAK PlasmaBlade™ TnA

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
PEAK PlasmaBlade™ TnA

Device Description

The PEAK PlasmaBlade™ TnA is a single-use, monopolar RF device. It consists of a single handpiece and two interchangeable electrode tips; one for tonsils and one for adenoids. The handpiece connects to the PULSAR® Generator and may be operated with the integrated handswitch or the PULSAR Footswitch. The tonsil tip features a bendable midsection and a curved and tapered blade with an opening in the center to allow for the evacuation of smoke and fluids. The adenoid tip consists of a wire electrode housed in a plastic tip with a bendable suction lumen that allows for the evacuation of tissue, fluids, and smoke. Both tips connect to the suction shaft of the handpiece.

Indications

The PEAK PlasmaBlade TnA device is only indicated for cutting and coagulation of soft tissue during otolaryngology (ENT) surgery including adenoidectomy and tonsillectomy (Pharyngeal, Tubal, Palatine).

Adverse Events

As a consequence of electrosurgery, damage to surrounding tissue through iatrogenic injury could occur.

Warnings

- Observe fire precautions at all times. An electrosurgical device may provide an ignition source due to sparking and heating.
- Do not use in the presence of flammable anesthetics or oxidizing gases such as nitrous oxide and oxygen. Do not activate the device until vapors from alcohol-based skin prepping agents have dissipated. Naturally occurring gases that accumulate in body cavities can also be an ignition source.
- Proper suction connection is important to prevent a build-up of flammable or oxidizing gases at the surgical site. During use, monitor the device for inadequate suction performance, which might be apparent by excessive /

visible smoke or burnt odor. Use the Cleaning Brush provided to remove excess eschar build up and maintain a clear channel for suction.

- Do not contact metal objects and instruments with the PEAK PlasmaBlade while power is being applied as unintended tissue damage and electrode tip damage could occur.
- Activation of the PEAK PlasmaBlade outside the field of view could cause patient injury.
- Activation of the handpiece when not in contact with target tissue may cause capacitive coupling.
- Do not allow patient contact with grounded metal objects, as such contact may result in patient or user injury.
- Activation of the handpiece simultaneously while aspirating fluid may alter the path of electrical energy away from target tissue.
- Inadvertent patient contact may result in burns. When not in use, place the device in a dry and nonconductive area away from the patient.
- Ensure that only the active tip of the device is in contact with the patient during use.
- Monopolar devices require a Patient Return Electrode. The unit must detect proper Patient Return Electrode impedance before generator output can be active. The impedance is continuously monitored and displayed while in Monopolar Mode. The unit presents audible and visible alarms if the unit detects improper impedance with the Patient Return Electrode in Monopolar Mode and will disable generator output. Refer to the manufacturer's instructions for application site and placement procedures when applying the Patient Return Electrode. Do not rely entirely on the impedance sensing feature as it can be affected by a damaged (shorted) Patient Return Electrode. It is recommended that the operator verify appropriate placement and contact of the Patient Return Electrode. Inadequate contact of the Patient Return Electrode may result in patient alternate site burns or injury.
- Do not reuse, resterilize or reprocess the PEAK PlasmaBlade as it is supplied sterile and intended for Single Use Only. A device that has been resterilized or reprocessed may not perform properly and may result in patient or user injury.
- Avoid fluid contact with the handle and its interfaces and connections as this may result in patient or user injury, or device failure.
- Position the cable to avoid patient contact to protect against high frequency current paths to the patient as this may result in patient or user injury.
- Position the cable to avoid creating a tripping hazard.
- The use of electrosurgery in the presence of internal or external cardiac pacemakers is potentially hazardous. Interference from the electrical current can cause device malfunction. Consult the cardiac pacemaker manufacturer for further information before proceeding with the surgery.

- The device should not be used near electrocardiograph electrodes as it can cause interference.
- Ensure that the tip is fully seated into the suction shaft as any gaps may cause unintended tissue damage, and injury could occur.
- Excessive bending may compromise performance and cause device failure, resulting in patient or user injury.
- Postoperative bleeding is inherent to surgical procedures and may occur regardless of the electrosurgical hemostatic device used. Appropriate and thorough post-op surgical care should be administered to the patient.

Precautions

- The PEAK PlasmaBlade should only be used by qualified medical personnel possessing training in the surgical procedures to be performed.
- Take care when handling the electrode tip to prevent possible damage to the tip and to prevent user injury.
- Take care when cleaning around the adenoid wire electrode, as excessive force may cause damage or breakage.
- Use the lowest power setting and the shortest activation time possible to achieve the desired end effect.
- Unless the product is being used to spot coagulate a vessel, it is recommended to keep the electrode tip in motion while activated to avoid excessive eschar buildup. Excessive eschar buildup can compromise device performance, including reduced or clogged suction.
- Do not use sharp or abrasive instruments or materials to clean eschar buildup on the electrode tip as this may damage the tip.
- When not in use, the device should be placed in a secure sterile location to prevent accidental damage to the tip.
- Prior to initial use, ensure that all package inserts, including warnings, cautions, instructions for use, as well as the *PULSAR Generator Operator's Manual*, are read and understood.
- Prior to use, carefully inspect the package before opening. Do not use the device if the package appears damaged in any way. Return any damaged packages to Medtronic Advanced Energy LLC.
- Prior to use, inspect the PEAK PlasmaBlade for any defects. Do not use if insulation or connectors are damaged.
- The PEAK PlasmaBlade is intended for use only with the PULSAR Generator.

Equipment List

- PULSAR Generator
- Optional PULSAR Footswitch
- Any standard Patient Return Electrode

Instructions For Use

- Read all instructions carefully. Failure to comply may lead to electrical or thermal injury, or cause device malfunction.
- Please refer to the *PULSAR Generator Operator's Manual* for step-by-step instructions regarding initial setup and

preparation for surgery. The PEAK PlasmaBlade may only be used with the PULSAR Generator.

Before Surgery

1. Inspect the PEAK PlasmaBlade package to ensure that it is unopened and undamaged. If it is damaged or the sterile seal is broken, do not use the device.
2. Remove the PEAK PlasmaBlade from the packaging using aseptic technique.
3. Inspect the device for any signs of damage. If the device is damaged it must not be used.
4. Connect suction to the suction connector on the handpiece.
5. For detailed instructions on the PULSAR Generator use, please refer to the operator's manual. Do not use the PEAK PlasmaBlade until reading and understanding the directions for the PULSAR Generator.
6. Switch on the PULSAR Generator to perform the initial setup and self-test.
7. Plug the PEAK PlasmaBlade into the Monopolar Connector of the PULSAR Generator. Upon connection, the unit will alarm if the device is not recognized as a PEAK PlasmaBlade device or if it fails a unit test. In event of an alarm, replace the device.

The PULSAR Generator will display default power level settings for the PEAK PlasmaBlade upon connection. These settings must be verified and adjusted by the user to ensure the appropriate settings are selected. The lowest power setting possible should be used to achieve the desired end effect. Table 1 indicates the PEAK PlasmaBlade settings available.

Table 1: PEAK PlasmaBlade Settings

Tissue Effect	Level (LCD display)	Mode	Power [W-RMS] @ Load (Ohms)	Maximum Output Voltage (Vp-p)
CUT low hemostasis, low collateral damage ↓ high hemostasis, higher collateral damage	1	Low Cut	6 @ 40	1660
	2		10 @ 40	1750
	3		20 @ 40	1830
	4	Medium Cut	30 @ 600	780
	5		40 @ 600	780
	6		50 @ 600	780
	7	High Cut (Blend 1)	30* @ 600	1560
	8	High Cut (Blend 2)	40 @ 600	1560
	9		40 @ 600	1930
	10		50 @ 600	2110
COAG low hemostasis, low collateral damage ↓ high hemostasis, higher collateral damage	1	Low Coag	20 @ 600	1550
	2		30 @ 600	1880
	3		35 @ 600	2020
	4		42 @ 600	2200
	5		50 @ 600	2370
	6	High Coag	25** @ 1000	3630
	7		30 @ 1000	3960
	8		35 @ 1000	4260
	9		40 @ 1000	4520
	10		50 @ 1000	5000


* When Cut level is changed from 6 to 7, Cut mode changes from Medium Cut to High Cut and there is a decrease in power. Hence, output power will increase when the Cut level is changed from 7 to 6.

** When Coag level is changed from 5 to 6, Coag mode changes from Low Coag to High Coag and there is a decrease in power. Hence, output power will increase when the Coag level is changed from 6 to 5.

- If using the optional footswitch, plug the wireless receiver into the Footswitch Connector of the PULSAR Generator. The footswitch can only be used to activate output on a monopolar device if a bipolar device is not connected.
- Position the Patient Return Electrode on the patient and then connect the Patient Return Electrode to the Patient Connector receptacle of the PULSAR Generator.
- Monopolar devices require a Patient Return Electrode. The PULSAR Generator must detect proper Patient Return Electrode impedance before it can be active. The operator should confirm the presence of a properly connected Patient Return Electrode.
- The PULSAR Generator has an Impedance indicator display (Single Foil = , Split Foil = ) that indicates the type and impedance level of the Patient Return Electrode in monopolar mode. The indicator will remain green with proper return electrode impedance. The indicator will turn red when the unit does not detect proper return electrode impedance. Ensure proper return electrode impedance before starting surgery. Refer to the *PULSAR Generator Operator's Manual* for more information.

Using the Adenoid Tip

- Ensure power is *not* applied while handling the tip, as this may result in user injury.
- Remove the Adenoid Tip from the package by peeling open the ends of the pouch.
- Inspect the tip for damage. If the tip is damaged, do not use it.


Push to white line on shaft.

- Connect the tip by holding the finger grip and pushing the tip straight into the suction shaft until the proximal end of the finger grip is aligned with the white line on the shaft.
- Ensure that the tip is fully seated and the finger grip overlaps the end of the suction shaft. If the tip is not fully seated unintended tissue damage and/or injury may occur.
- The Adenoid Tip can be rotated in any direction. Use finger force to bend the tip to a maximum of 60° from the suction shaft. Bend the tip only in the direction perpendicular to the suction opening. Excessive bending may compromise performance and cause device failure, resulting in patient or user injury.
- The Adenoid Tip can be cleaned using the Cleaning Brush by gently brushing the eschar off the electrode from the front. Take care when cleaning around the wire electrode as excessive force may cause damage or breakage.
- To remove the Adenoid Tip, grasp the finger grip and pull the tip straight off.

Using the Tonsil Tip

- Ensure power is *not* applied while handling the tip, as this may result in user injury.
- If applicable, remove the Tonsil Tip from the package by peeling open the ends of the pouch.
- Inspect the tip for damage. If the tip is damaged, do not use it.

Push to white line on shaft.

4. Connect the tip by squeezing the finger grip and pushing the tip straight onto the suction shaft until the proximal end of the finger grip is aligned with the white line on the shaft. Take care not to damage the cutting blade.
5. Ensure that the tip is fully seated and the finger grip overlaps the end of the suction shaft.
If the tip is not fully seated, unintended tissue damage and/or injury may occur.
6. The Tonsil Tip can be rotated in any direction by squeezing the finger grip and rotating. The midsection can be bent using finger force to a maximum of 60° from the suction shaft. Bend the tip only in the direction perpendicular to the suction opening. Excessive bending may compromise performance and cause device failure, resulting in patient or user injury.
7. To remove the Tonsil Tip, squeeze the finger grip and pull the tip straight off.

During Surgery

Place the device on a sterile tray when not in use.

1. To cut, press the yellow button/pedal on the handswitch or footswitch.
2. To coagulate, press the blue button/pedal on the handswitch or footswitch.
3. Eschar buildup on the tip can be removed manually with gloved fingers or gauze pads.

After Surgery

1. Turn off the PULSAR Generator.
2. Disconnect the PEAK PlasmaBlade, footswitch, and Patient Return Electrode from the PULSAR Generator.
3. Discard the PEAK PlasmaBlade after use. It is intended for Single Use Only. Do not reuse, resterilize, or reprocess. A device that has been resterilized or reprocessed may not perform properly and may result in patient or user injury.

How Supplied

The PEAK PlasmaBlade is supplied sterile and is intended for Single Use Only. Contents are sterile unless the package is opened or damaged. Do not reuse, resterilize, or reprocess the PEAK PlasmaBlade. The product is sterilized using Ethylene Oxide (EO).

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Trademarks: PEAK, PEAK PlasmaBlade, and PULSAR are registered trademarks or trademarks of Medtronic Advanced Energy LLC.
Marchi: PEAK, PEAK PlasmaBlade e PULSAR sono marchi depositati o marchi di fabbrica di Medtronic Advanced Energy LLC.
Marken: PEAK, PEAK PlasmaBlade und PULSAR sind eingetragene Marken oder Marken von Medtronic Advanced Energy LLC.
Védjegyek: A PEAK, a PEAK PlasmaBlade és a PULSAR a Medtronic Advanced Energy LLC bejegyzett védjegyei vagy védjegyei.
Medtronic B.V.