

# MicroBeta Trilux – Your total solution

**wallac**



**Radiolabels**

**Luminescence**

**PerkinElmer™**  
life sciences.

**PerkinElmer™**  
life sciences.

# MicroBeta® TriLux is th

**MicroBeta® TriLux is a truly universal counting instrument. It can count beta, gamma or luminescent labels on filters, in microplates or in tubes. There are 12, 6, 3, 2 or 1 detector versions and the samples can be in 384, 96 or 24 format. In terms of detector design, sample handling, counter operation and running costs – all pivotal factors when deciding on a new counter – MicroBeta wins every time.**



#### **For basic research**

MicroBeta is an "OPEN" system, designed for compatibility with any counting vessel manufactured by any other company. Detector construction, cassette sample handling, crosstalk correction and software compatibility are all designed to put the user in control of the instrument rather than have the limitations of the counter dictate the assay method.

#### **For high throughput screening**

As both a high sample capacity stand-alone and as part of a fully automated system, MicroBeta has gained acceptance in high throughput screening laboratories around the world. With up to 12 detectors and only 9 readings per 96-well plate, it is the fastest plate counter available.








#### **The Wallace logo - hallmark of quality and innovation**


MicroBeta® TriLux comes from the company which pioneered plate counting. In 1987 the Wallace Betaplate™ was the world's first multi-detector liquid scintillation counter for filtermats. Later came MicroBeta, the very first microplate counter.

While leading the field, PerkinElmer Life Sciences has secured patents on numerous innovative techniques that help MicroBeta TriLux users to achieve reliable results while saving work, time, money and waste.


# e best possible choice



*Today's leading multidetector liquid scintillation and luminescence counter comes in five different detector versions.*



*A selection of the filters, plates and other counting vessels commonly used with MicroBeta® TriLux*




## Easy GLP

Compliance with Good Laboratory Practice (GLP) is an essential feature of any analytical instrument, and like all other Wallac counters MicroBeta includes customized protocols for monitoring instrument performance over an extended period.

# Built around patent technology solution

At the very heart of the instrument,  
the detector



An important feature of MicroBeta TriLux is the twin photomultiplier tube detector design. This is the basis for coincidence counting, a robust technique which assures unsurpassed counting performance for beta and gamma labels. For luminescence counting too, separate high speed electronic circuitry provides a high dynamic range for samples in 384, 96 or 24-formats.

For counting opaque white plates, it is possible to count with just the upper PMT.


*The upper detectors of MicroBeta® TriLux move to accommodate plates, filters, or tubes, which can thus be counted in the vertical position.*



# ed and proven ns

## Cassettes, for reliability and freedom of choice in sample handling

In the Wallac MicroBeta, samples are supported by a cassette-based sample changing mechanism. Samples are never stacked on top of each other and so the risks of contamination or sticking and jamming are minimized. Cassettes provide the only practical method of counting flexible plates, filtermats and tubes.



Counting protocol identification and other MicroBeta® commands are read from a reusable bar code sticker that is fixed to the cassette. This provides the most convenient and economical identification method.

For high throughput screening, when it is often imperative that each plate is positively identified, a second bar code reader is focused to read directly from the plate. This optional accessory is compatible with Code 39 and several other bar code languages.

# A truly versatile inst and luminescence

## Liquid samples – from 25 $\mu$ L to 4 mL

The cassette handling system provides support for flexible plates, rigid plates and tubes and so every possible combination of requirements for sample volume and solvent resistance can be catered for. All types of liquid cocktail are equally suited to the detector system.

For isotopes that may give rise to isotopic crosstalk, like  $^{32}\text{P}$  or  $^{125}\text{I}$ , our exclusive crosstalk correction program can be enabled to improve sensitivity. Results can be reported in CPM or DPM, single or dual label. Supporting 384, 96 or 24 well formats, only MicroBeta<sup>®</sup> can provide this level of flexibility.

	Cocktail	CPM	DPM	Eff.%	Bkd	x-talk%
H-3 *	Betaplate Scint	194666	334222	58	8	0.00
C-14 *	Betaplate Scint	278168	307972	90	9	0.00
P-32	SuperMix	378437	382760	99	7	0.25
I-125	SuperMix	138247	227047	61	9	0.30
Cr-51	SuperMix	69223	248314	28	3	0.03


\* sample in organic solvent

*MicroBeta<sup>®</sup> TriLux crosstalk, background and efficiency figures for commonly measured nuclides. All measurements using the Wallac Isoplate<sup>™</sup>, 25  $\mu$ L sample with 150  $\mu$ L cocktail.*

## A popular choice for Chromium Release

Chromium release is a typical assay that is widely counted on the MicroBeta. Counting efficiency is five to six times more efficient than a gamma counter whilst the 8 x 12 format is maintained. The procedure is very straightforward. Supernatant is transferred to a new plate and after the addition of cocktail, sealing and shaking the samples are ready for counting.

Percent specific marker release



# Instrument for LSC

Restrictions on the use of radioactive materials and increased disposal costs are adding pressure to the search for non-radioactive alternatives. The MicroBeta TriLux provides an opportunity to develop the increasingly popular "glow type" luminescence assays while routine use of radiolabels continues to be necessary. Separate, ultra fast electronic circuitry is employed in luminescence counting. Linear responses of up to 24 million counts per second are obtainable. All types of plates, filters and tubes can be counted. Up to 12 samples are counted at a time in 384, 96 or 24-well formats.

## Temperature Control Option for Luminescence


With MicroBeta, temperature control is not necessary for liquid scintillation counting. For some luminescence assays, however, a low and constant background may be desirable. The temperature control unit is a conveniently sized, solid state Peltier device that maintains a constant temperature of the upper PM tubes. The active temperature can be included in the printout.

## Measure flash-type luminescence with MicroBeta® JET



*The MicroBeta® JET, a sister product to TriLux, has built in reagent dispensers to allow flash-type luminescence.*

*Please contact your local representative for further details.*




*Measurement of ALP-reporter gene activity AMPPD® substrate measured using MicroBeta®*

## Wallac luminescence kits

New luminescence kits for use with MicroBeta TriLux include the Wallac GeneLux™ enhanced luciferase assay and the Wallac CytoLux™ assay. The latter represents a genuine high performance non-radioactive method for cytotoxicity and cell viability assay.

# Consumables, cassettes and


Plate	How to use it	Applications	Plate or vial	How to use it		
	<b>1450-401</b> 96-well sample plate, 250 mL, round bottom wells, flexible, made of PET-G.	During sample preparation the plate may be used with support (1450-481). It is counted in 96-well cassette (1450-101) with either permanent seal (1450-461) or removable seal (1450-462). SuperMix or HiSafe scintillator is recommended.	Microvolume radiolabelled assays, Chromium Release, SPA, Cerenkov <sup>32</sup> P, Coated plate assays		<b>1450-402</b> 24-well sample plate, 1 mL, flat bottom wells, flexible, made of PET-A.	During sample preparation the plate may be used with support (1450-481). It is counted in 24-well cassette (1450-102) with either permanent seal (1450-461) or removable seal (1450-462). SuperMix or HiSafe scintillator is recommended.
	<b>1450-514 or 1450-515</b> Isoplate high performance 96-well plate, polystyrene clear 350 mL flat bottom wells, rigid white exterior.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462). SuperMix scintillator is recommended.	Microvolume radiolabelled assays, SPA, Cerenkov <sup>32</sup> P		<b>1450-408</b> Heatsealable 24-well sample plate, 1 mL, flat bottom wells, flexible, made of PET.	During sample preparation the plate may be used with support (1450-481). It is counted in 24-well cassette (1450-102) with either permanent seal (1450-461) or removable seal (1450-462). SuperMix or HiSafe scintillator is recommended.
	<b>1450-571 or 1450-572</b> Black Isoplate, high performance 96-well plate, polystyrene clear 350 mL flat bottom wells, rigid black exterior.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence assays		<b>1450-601 or 1450-602</b> Visiplate 24-well sample plate, 2 mL, polystyrene, flat bottom wells	Used with rigid 96-well cassette (1450-105) with SuperMix scintillator.
	<b>1450-581 or 1450-582</b> B & W Isoplate, high performance 96-well plate, polystyrene white 350 mL flat bottom wells, rigid black exterior.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence assays		<b>1450-603 or 1450-604</b> Visiplate TC, tissue culture.	Used with rigid 96-well cassette (1450-105) with SuperMix scintillator and permanent seal.
	<b>1450-518 or 1450-519</b> Isoplate HB, high binding.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462). SuperMix scintillator is recommended.	Radiolabelled coated plate assays		<b>1450-605 or 1450-606</b> Black Visiplate TC, tissue culture.	Used with rigid 96-well cassette (1450-105) with SuperMix scintillator and permanent seal.
	<b>1450-575 or 1450-576</b> Black Isoplate HB, high binding.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence coated plate assays		<b>FlashPlate® and ScintiPlate</b> scintillant coated 96-well microplate, or 1450-501 or 1450-502 ScintiPlate 96-well microplate manufactured with scintillating solid phase.	Counted in 96-well cassette (1450-101) with permanent seal.
	<b>1450-585 or 1450-586</b> B & W Isoplate HB, high binding.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence coated plate assays		<b>Flashplate</b> 384-well 384-well is used.	Flashplate 384-well 384-well is used.
	<b>1450-516 or 1450-517</b> Isoplate TC, tissue culture.	Used with rigid 96-well cassette,(1450-105) with SuperMix scintillator and either permanent seal (1450-461) or removable seal (1450-462).	Radiolabelled adherent cell assays		Standard 96-well plates	May be counted in 96-well cassette (1450-101) and either permanent seal (1450-461) or removable seal (1450-462). SuperMix scintillator is recommended.
	<b>1450-573 or 1450-574</b> Black Isoplate TC, tissue culture.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence cell studies		24-well cell culture plates	Counted in 24-well cassette (1450-102) with insert tubes (1450-108) or permanent seal (1450-461). SuperMix scintillator is recommended.
	<b>1450-583 or 1450-584</b> B & W Isoplate TC, tissue culture.	Used with rigid 96-well cassette,(1450-105) with either permanent seal (1450-461) or removable seal (1450-462).	Luminescence cell studies		Eppendorf tubes	Counted in 1450-101 SuperMix scintillator is recommended.



Cassette 96-well  
1450-101



Cassette 24-well  
1450-102



Rigid 96-well  
cassette  
1450-105



Corning cassette  
1450-110 used  
with insert tubes  
1450-109




24-well cassette  
1450-102 with  
1450-108 adapters



4 mL tube cassette  
1450-117



# and scintillators for MicroBeta®

Use it	Applications	Filter or tube	How to use it	Applications
Sample preparation the plate sealed with support (1450-481). Counted in 24-well cassette (1450-102) ②. For permanent seal (1450-461) useable seal (1450-462). For x or HiSafe scintillator recommended.	Microvolume radiolabelled assays, Chromium Release, SPA, Cerenkov <sup>32</sup> P	 4 mL tubes, e.g. 1200-421	Counted using the 4mL tube cassette (1450-117) ⑧. SuperMix scintillator is recommended.	Microvolume LSC
Sample preparation the plate sealed with support (1450-481). Counted in 24-well cassette (1450-102) ②. For permanent seal (1450-461) useable seal (1450-462). For x or HiSafe scintillator recommended.	Microvolume radiolabelled assays	 1450-421 Filtermat A, 8 x 12 format glassfibre.	Counted in filter cassette (1450-104) ⑦ using Betaplate Scint or MultiLex solid scintillator.	Labelled cell assays
For rigid 24-well cassette (1450-110) ④. For SuperMix scintillator and permanent seal (1450-461)	Radiolabelled assays requiring large sample sizes	 1450-422 Filtermat A, 4 x 6 format glassfibre.	Counted in filter cassette (1450-116) ③ using Betaplate Scint or MultiLex solid scintillator.	Labelled cell assays
For rigid 24-well cassette (1450-110) ④. For SuperMix scintillator and permanent seal (1450-461)	Radiolabelled adherent cell assays	 1450-521 Filtermat B, 8 x 12 format glassfibre, double thickness.	Counted in filter cassette (1450-104) ⑦ using Betaplate Scint or MultiLex B/HS solid scintillator.	Receptor ligand binding assays
For rigid 24-well cassette (1450-110) ④. For permanent seal (1450-461)	Luminescence cell studies	 1450-424 Filtermat B, 4 x 6 format glassfibre, double thickness.	Counted in filter cassette (1450-116) ③ using Betaplate Scint or MultiLex B/HS solid scintillator.	Receptor ligand binding assays
For using rigid 96-well cassette (1450-105) ⑨ with permanent seal (1450-461).	See Wallac Application Note 1450-1000	 Betaplate™-format filtermat, e.g. 1205-401, with 6 x 16 format	Cut into 4 parts and counted in filter cassette (1450-107) ⑨ using Betaplate Scint or MultiLex solid scintillator.	All types of filter application
Plate is also available in 8 x 12 format. With this the cassette (1450-130) ⑫		 1450-522 DEAE filtermat for negatively charged compounds, 8 x 12 format glassfibre.	Counted in filter cassette (1450-104) ⑦ using Betaplate Scint or MultiLex solid scintillator.	Enzyme assays, e.g. reverse transcriptase
Counted using cassette (1450-105) ⑨. For permanent seal (1450-461) or removable seal (1450-462) with cocktail.	General applications	 1450-523 P30 filtermat for positively charged compounds, 8 x 12 format glassfibre.	Counted in filter cassette (1450-104) ⑦ using Betaplate Scint or MultiLex solid scintillator.	Binding positively charged compounds, e.g. protein kinases.
For using rigid 24-well cassette (1450-110) ④ with permanent seal (1450-461) and with cocktail. For x scintillator is recommended.	See Wallac Application Note 1450-0003	 1450-423 Nylon membrane, 8 x 12 format.	Counted in filter cassette (1450-104) ⑦ using Betaplate Scint or MultiLex solid scintillator, or in <sup>32</sup> P cassette ⑩ without scintillator.	RNA/DNA dot blots, see Wallac Application Note 1450-1010
For using 24-well cassette (1450-102) with adapters (1450-108) ⑤. For x scintillator is recommended.	Microvolume LSC	 Millipore® Multiscreen® plate	Counted in Millipore cassette (1450-106) ⑪ with cassette liner (1450-433) using SuperMix or MultiLex solid scintillator.	See Wallac Application Note 1450-984

⑦



Filter cassette  
1450-104

⑧



Filter cassette  
1450-116

⑨



Filter cassette  
1450-107

⑩



<sup>32</sup>P cassette  
1450-118

⑪



Millipore cassette  
1450-106  
Cassette Liner  
1450-433

⑫



384-well cassette  
1450-130

# MicroBeta® TriLux with gives you a dedicated

## MicroBeta® and SPA

Amersham's Scintillation Proximity Assay (SPA) represents a unique sample type for which a radical new counting method, ParaLux™ Count Mode, has been developed.




*Reliable counting of 384-well plates places great demands on any counting system. Only MicroBeta® and ParaLux™ Count Mode can maintain counting efficiency and accurate quench correction of coloured samples.*


In SPA and CytoStar-T™ samples the activity is much closer to the lower PM tube than the upper one and so, especially in the case of coloured samples, MicroBeta® detector design provides the ultimate counting geometry.

# h ParaLux™ d SPA counter

ParaLux Count Mode combines twin PM tube detection with dual multi-channel analyser (MCA) data analysis. Counting efficiency is up to 500% better than that achieved with any other counting method. A completely new quench parameter, AQP(I), is a genuine advance in counting technique. Superior DPM calculations at considerably higher levels of colour quench mean that time consuming and costly dilution and re-counting can be avoided. Improved DPM precision can also lead to overall shorter counting times.



H-3 Efficiency / %



Comparison of tritium SPA colour quench series counted with conventional coincidence circuitry and ParaLux™ high efficiency mode.


## An optimum solution for Cytostar-T™ plates, too

MicroBeta counters are ideal for work with monolayer cells on Cytostar-T plates. The scintillating base of these plates is directly over the MicroBeta's lower PM tube, to allow maximum capture of signal.


$$\text{AQP(I)} = \text{MCA-1} / \text{MCA-2}$$

$$\text{High Efficiency CPM} = \text{MCA-1} + \text{MCA-2}$$

$$\text{Low Background CPM} = \text{MCA-1}$$




ParaLux™ low background mode. Tritium SPA quench curves. AQP(I) provides overlapping colour quench curves. This permits the use of a single quench series for DPM calculations.



(<sup>14</sup>C) thymidine uptake into synchronous hamster lung V79 cells on 96 well Cytostar-T scintillating microplates.

# MicroBeta<sup>®</sup>, the highest lowest cost option

There are **two methods** of preparation  
for filter samples.



The MicroBeta<sup>®</sup> TriLux is the ideal solution for labelled cell studies, ligand receptor binding assays, enzyme assays and DNA/RNA hybridization assays. There are two important reasons why. Firstly, the superior sample handling afforded by the cassette support mechanism offers true versatility when considering samples harvested or spotted on to filtermats. Secondly, the twin PM tube detectors provide superior counting geometry. Compared to all other methods, tremendous savings in consumables costs, sample handling, counting time and waste volume can be made.

## With liquid scintillator

Using liquid scintillator, the sample is put in a bag, a few millilitres of cocktail are spread over the filtermat, the bag is sealed and then placed in a cassette.

## With solid scintillator

MeltiLex<sup>®</sup> is an easy-to-use alternative if the sample needs to be permanently fixed to the filter or if solid waste is preferred. You simply place the filtermat or filter plate on a sheet of MeltiLex and apply gentle heat with, for example, the Wallac MicroSealer or a hotplate.



# Best efficiency, for filter counting


## The method of choice for $^{32}\text{P}$ counting

For user convenience and superior counting efficiency, the  $^{32}\text{P}$  cassette is used for  $^{32}\text{P}$  labelled filter samples. Absolutely no sample preparation is required because solid scintillator is permanently fixed in the cassette. Filters can be counted wet or dry and because this is a non-destructive method, membranes used in DNA hybridization assays can be re-probed


Count **any type** of  
filter plate - **any way**  
you want

MicroBeta detector geometry is ideal for counting filter plates.


As when counting from filters, samples can be prepared with either liquid cocktail or, if long term stability and solid waste is preferred, MeltiLex solid scintillator may be used.



$^{32}\text{P}$  cassette and  $^{32}\text{P}$  labelled filter.



In addition to a complete range of Wallac ready cut, printed filtermats, you can also use most standard filtermats including Skatron-type formats.



MicroBeta® is a particularly strong option for the MultiScreen Assay System. For fully automated systems or simply for user convenience, the MicroBeta offers the only facility to count with the underdrain intact.



# Adherent cell studies

## “Freeplate” counting with MicroBeta®

Specially treated Wallac Tissue culture grade microplates are available in 96- and 24-formats for both radiolabel and luminescence assays.

A variety of plates with opaque and clear wells are supplied individually wrapped with a lid.



For adherent cell assays that require traditional 24-well TC plates counting geometry can be optimised by the addition of white tubes prior to counting.

*Specifications for all Wallac plates and most other commonly used plates are supplied with your MicroBeta®. For other plates you can easily input the plate dimensions.*


A screenshot of the 'Plate/Filter' software window. The window has a title bar 'Plate/Filter' and a close button. It is divided into two main sections: 'Identification' and 'Dimensions'.  
**Identification:**  
- A text box contains '24-format TC plate'.  
- Below it, a text box contains 'Use with 1450-109 cassette, 1450-110 inserts and 100ul SuperMix'.  
**Size:**  
- Three radio buttons: '24 (4x6)' (selected), '96 (8x12)', and '384 (16x24)'.  
**Plate/Filter:**  
- Two radio buttons: 'Plate' (selected) and 'Filter'.  
**Dimensions:**  
- Two radio buttons: 'Wallac' and 'Custom' (selected).  
- Four input fields: 'A1 X' (17.3), 'A1 Y' (13.4), 'A6 X' (113.6), and 'D1 Y' (71.5).  
- A checkbox 'Use all detectors' is unchecked.  
At the bottom are 'OK' and 'Cancel' buttons.

# Scintillating plates for all your binding assays


ScintiPlate® and FlashPlate® are designed for homogeneous radiometric assays. Employing the principle of scintillation proximity, only radiolabel that is bound to the scintillating surface of the well produces light.

FlashPlates are available in 96- and 384-formats and with a wide variety of surface coatings for immunoassay, enzyme, functional and binding assays.

ScintiPlates, in 96-format, are either plain or streptavidin coated. A tissue culture grade ScintiPlate is recommended for adherent cell assays.



*Measurement of cyclic AMP on FlashPlate with MicroBeta.*




# Whatever your application, whatever your lab type

The Wallac MicroBeta® comes as part of a comprehensive package including instrumentation, software, reagents and application support. Its users are found in academic research, and in industrial laboratories carrying out, for example, high throughput screening.

Whatever your specialization, as a MicroBeta user, the wide-ranging expertise of the PerkinElmer Life Sciences team in all aspects of plate counting is at your disposal.

*For high throughput screening, MicroBeta forms a complete and seamless system integrated with the Beckman Biomek 2000 BioRobotics System.*



AMPPD is a registered trademark of Tropix Corp.  
Biomek is a registered trademark of Beckman Instruments Inc.  
CytoStar is a trademark of Nycomed Amersham plc.  
Falcon is a trademark of Becton Dickinson, Inc.  
FlashPlate is a registered trademark of Packard Instrument Company, exclusively licensed to PerkinElmer Life Sciences.  
Millipore is a registered trademark and MultiScreen is a trademark of Millipore Corp.  
Products for scintillation proximity assay (SPA) are available from Nycomed Amersham plc.

**Worldwide Headquarters:** PerkinElmer Life Sciences, 549 Albany Street, Boston, MA 02118-2512 USA (800) 551-2121

**European Headquarters:** PerkinElmer Life Sciences, Imperiastraat 8, BE-1930 Zaventem Belgium +32 2 717 7911

**Technical Support:** in Europe: techsupport.europe@perkinelmer.com in US and Rest of World: techsupport@perkinelmer.com

**Belgium:** Tel: 0800 94 540 • **France:** Tel: 0800 90 77 62 • **Netherlands:** Tel: 0800 02 23 042 • **Germany:** Tel: 0800 1 81 00 32 • **United Kingdom:** Tel: 0800 89 60 46  
**Switzerland:** Tel: 0800 55 50 27 • **Italy:** Tel: 800 79 03 10 • **Sweden:** Tel: 020 79 07 35 • **Norway:** Tel: 800 11 947 • **Denmark:** Tel: 80 88 3477 • **Spain:** Tel: 900 973 255

[www.perkinelmer.com/lifesciences](http://www.perkinelmer.com/lifesciences)