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Fig. Front cover: BIOSTAT® RM 20 front view

Sartorius Stedim Systems GmbH

Technical Documentation
Schwarzenberger Weg 73–79
34212 Melsungen
Germany

www.sartorius-stedim.com
E-mail: info.fermentation@sartorius-stedim.com

Printed in Germany on chlorine-free paper.
If you require mandatory information regarding specific properties of your BIOSTAT® RM that goes beyond that which is described here, please contact us. No part of this publication may be reprinted, reproduced or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH.

This operating manual is valid for RM devices.

BIOSTAT® is a registered trademark of Sartorius Stedim Biotech GmbH.

Content

Introduction	6	4. Menu Operation	21
Symbols	7	4.1 The Main Menu	21
Safety Instructions	8	4.2 Preparing and Inflating the CultiBag RM	22
Setup Location Requirements	9	4.2.1 Unpacking the CultiBag RM	22
Other Operating Conditions	9	4.2.2 Securing the CultiBag RM to the Bag Holder RM	22
Term Definitions	10	4.2.3 Installing the Filter Heater	23
1. Design and Function	11	4.2.4 Connecting the Air Supply	23
1.1 Intended Use of the BIOSTAT® RM 20 50 basic	11	4.2.5 Inflating the CultiBag RM	23
1.2 CultiBag RM	11	4.3 Cultivation Start	24
1.3 System Design	12	4.3.1 Filling the CultiBag RM with the Culture Medium	24
1.4 Rocker 20 50	13	4.3.2 Setting the Cultivation Parameters and Starting Cultivation	25
1.4.1 Operation	13	4.4 Operating the Control Elements	25
1.4.2 Aeration Control	14	4.4.1 Aeration Menu	25
1.4.3 Air and CO ₂ Supply	14	4.4.1.1 Aeration with Air	26
1.4.4 Aeration Rate	14	4.4.1.2 Aeration with CO ₂	26
1.4.5 CO ₂ Aeration	14	4.4.2 Speed and Angle Menu	27
1.4.6 Rocker Speed	14	4.4.3 Temperature Menu	28
1.4.7 Rocker Angle	14	4.4.4 Position Menu	29
1.4.8 Temperature Control	14	4.4.5 User Level Menu	30
1.4.9 Temperature Range	15	4.4.6 Alarm Menu	32
1.4.10 Serial Interface	15	4.4.7 Settings Menu	32
1.4.11 Housing Material	15	4.4.7.1 Setting the Date and Time	32
1.4.12 Dimensions Weight	15	4.4.7.2 Service Information	33
1.4.13 Power Supply Power Disconnect	15	4.4.7.3 Setting the Bag Size	34
1.4.14 Room Temperature Range	15	4.4.7.4 Selecting a Control Mode	34
1.5 Bag Holder	15	4.4.7.5 Sampling Time	36
1.5.1 Properties, Specifications	15	4.4.7.6 Batch Time	36
1.6 Heating Elements	16	4.5 Cultivation Parameters	37
1.7 Hood for Bag Holder 20 or 50	16	4.6 Sampling and Inoculation	37
1.7.1 Dimensions	16	4.6.1 Preparing Sampling or Inoculation	37
1.8 Surface Temperature Sensor	16	4.6.2 Inoculation	38
1.9 Filter Heater	16	4.6.3 Sampling	38
2. Delivery and Installation Instructions	17	4.7 Completing the Process, Cell Harvesting and Shutting Down the Device under Safe Conditions	39
2.1 Delivery Check	17	5. Cleaning and Maintenance	41
2.2 Space Requirement	17	5.1 Intermediate Cleaning	41
2.3 Connecting to Supply Lines in the Laboratory	17	5.2 Servicing the Device by Operating Personnel	41
2.3.1 Work Area Preparation	17	5.3 Servicing and Calibration by Authorized Service Personnel	42
2.3.2 Power Connection	17	6. Transportation and Storage	43
2.4 BIOSTAT® RM 20 50 basic	18	7. Appendix	45
2.4.1 Connections, Interfaces	18	7.1 Schematic Representation of Aeration	45
3. Initial Startup	19	7.2 Overpressure Control	45
3.1 Overview	19	7.3 Alarm Messages	46
3.2 Installation Kit	19	7.4 Malfunctions and Countermeasures	52
3.3 Rocker Design	19	7.5 Agreements	54
3.3.1 Installing the Bag Holder	19	7.5.1 Agreements and Intended Use	54
3.3.2 Installing the Surface Temperature Sensor	19	7.5.2 Service	54
3.3.3 Power Supply Connection	20	7.5.3 Device Disposal	55
3.3.4 Connecting the Compressed Air and CO ₂ Supply	20	7.6 Decontamination Declaration	55
3.3.5 Device Startup	20		

Introduction

These operating instructions describe the installation and operation of the BIOSTAT® RM 20 | 50 basic bioreactor.

Detailed information about the measurement and control element as well as other optional upgrades can be found in the corresponding additional documentation.

About this Documentation

The BIOSTAT® RM 20 | 50 basic was designed for the development and optimization of manufacturing processes for the cultivation of cells and microorganisms in single-use systems with low shearing forces. The operation of the bioreactor requires special knowledge in the aseptic area to ensure safe handling of cells and microorganisms. Special handling of cell and microorganism cultures may be required to prevent biological hazards. This also applies to the disposal of cultures and single-use components.

These operating instructions contain safety information regarding possible hazards as well as corresponding countermeasures, which only apply to the devices being described and which supplement other workplace regulations of the operator for the respective process. For applications associated with special hazards, additional safety equipment or measures for the protection of personnel and the work environment may be required. This documentation does not go into detail regarding such circumstances as well as legal or other binding regulations.

The laboratory model of the single-use bioreactor system BIOSTAT® RM 20 | 50 basic is part of the product line of Sartorius Stedim Systems GmbH.

For additional questions regarding the BIOSTAT® RM 20 | 50 basic or other Sartorius Stedim Systems GmbH products, please contact:

Sartorius Stedim Systems GmbH
Schwarzenberger Weg 73–79
D-34212 Melsungen, Germany
Tel.: +49.0.5661.713400
Fax: +49.0.5661.929945
info@sartorius-stedim-systems.com
ww.sartorius.com

Version Information

© Sartorius Stedim Systems GmbH. No guarantee is assumed for the information provided. Sartorius reserves the right to make technical changes to this document without any prior notice. No part of this document may be changed, or reproduced, nor is any other use permitted without obtaining prior written permission.

Biological hazard

This symbol refers to process-dependent risks to personnel. The possible hazards – damage to devices or health effects – that can result from this depends on the respective process and the microorganisms and cells being used.

This symbol indicates electrical hazards – damage to devices or health effects – due to voltage or current.

The surface in this area can have very high temperatures.

Skin contact can result in burns | injury.

Avoid direct contact with skin in this area and wear corresponding protective clothing | equipment.

Warning of hand injuries!

This symbol indicates risks to operating personal from moving parts.

Keep hands and other body parts away from this area.

Caution!

This indicates that non-compliance can cause damage to the device.

Notes with the "IMPORTANT!" warning indicate risks that can lead to material damage.

This symbol marks work steps that must be carried out with special care or it indicates important aspects that must be observed in order to avoid improper operation or malfunctions.

▷ This symbol indicates other noteworthy aspects or information.

1., 2., ... Numbered paragraphs indicate operating steps that should be performed in the order given.

This symbol references additional information in these instructions or in other documentation.

The operator (or responsible persons) must ensure that suitable conditions are in place for the respective process:

1. Only operating personnel qualified to run the process may be employed for working with the bioreactor; they must be aware of potential hazards, and be thoroughly familiar with the handling of all devices associated with the bioreactor.
2. The device must be operated and used according to the operating instructions. This is essential in order to protect the device and personnel.
3. Unauthorized persons should not be allowed access to the bioreactor. If the respective process is a potential hazard source (e.g. from cultures, media), the work area must be identified with suitable hazard symbols (such as "Biological hazard") and be marked off or blocked off. ► "Biological Safety Regulations."
4. The work area must be suitable for the process. It must be able to withstand, e.g. acids, alkalines or media and, in case of contamination, e.g. from cultures, it must be easy to decontaminate, disinfect or clean.
5. The operator (or responsible persons) must instruct personnel about hazards to persons and the work environment that may occur during the process, including mandatory safety requirements as well as provide safety equipment.
6. In general, personnel must wear suitable work clothing and appropriate personnel protective equipment (gloves, safety goggles and a facemask if necessary).
7. The device may only be opened by Sartorius service personnel.
8. Please note that the device may be operated under IP 23 protection conditions.
9. Please note that interference with the device is prohibited while it is operation.
10. Make sure that no unauthorized objects are attached to the rocker and the hood.
11. Make sure that the device is not connected to the mains supply during setup, installation and disassembly.
12. Make sure that setup, installation and disassembly work is performed by trained personnel.
13. Note the details on the ID label before installation and start up.
14. Make sure that the mains connection is equipped with a ground cable connection.
15. Only use the device for the purpose specified in the operating instructions. All safety measures on the device may become ineffective if used in any other way.

Setup Location Requirements

The work area and laboratory equipment must be suitable for the requirements of the bioreactor according to its specifications.

Please note the information under ► "2.2 Space Requirement" when selecting the setup location

1. Please note the weight and dimensions of the device and its accessories. Suitable tools and equipment must be used to transport the bioreactor to its setup location.
 2. Check and ensure that there is a supply of air to the underside of the device to keep the electronics cool.
- The work area must be robust enough to hold the bioreactor with all of its peripheral devices.

Chemical hazard | Gas leak

If the device is used improperly or is damaged, this can lead to a gas leak during processing.

For this situation, make sure that the BIOSTAT® RM 20 | 50 basic is only installed and operated in rooms with a suitable and proper ventilation system.

- You must make sure that the laboratory power lines and gas supply lines meet the specifications of the bioreactor.
- You can only use systems and accessories approved by Sartorius Stedim Systems for use with the BIOSTAT® RM 20 | 50 basic.
- Make sure that all connections between laboratory lines, the CultiBag RM and accessories are installed properly. Check for leaks and unintentional escapes.
- Make sure that all parts of the system and accessories are in full working order. This applies especially to the CultiBag RM. Damaged parts may not be used (this also applies if you only have reasonable suspicion that there is damage).

Other Operating Conditions

- You must make sure that the system is disconnected from the electrical power supply during all setup, assembly and removal activities.
- You must make sure that all setup, assembly and removal activities to the system are only carried out by trained personnel authorized by Sartorius.

Please observe all information on the manufacturer's ID label during installation.

Term Definitions

Spare part	Spare parts are mechanical components and parts which can only be replaced by SSB Service (e.g housing components).
Wearing part	Wearing parts are mechanical components and parts which undergo mechanical wear and can be replaced by the user (e.g. seals).
Disposables	These are components that are required used for the process and can be replaced by the user (e.g. bags).

1. Design and Function


1.1 Intended Use of the BIOSTAT® RM 20 | 50 basic

Fig. 1:
Design of the BIOSTAT® RM 20 basic

1.2 CultiBag RM

1. Design and Function

The BIOSTAT® RM 20 | 50 basic is a rocker platform for single-use bioreactors (CultiBag RM). The shaft mixing process of this platform uses mechanical energy to ensure homogenous mixing of cells. The energy supply required for this is achieved via the low-shear, rocking motion of the rocker to which the CultiBag RM is attached. This type of cell-friendly motion ensures homogeneous cell distribution in the cultivation medium whereby the surface of the medium is regenerated without causing a buildup of bubbles.

The BIOSTAT® RM 20 | 50 basic is a medium-scale platform with replaceable bag holders (size 20 or 50) for bags with an operational or cultivation volume of max. 10 or max. 25 liters. The system can be operated standalone on a laboratory table or in the incubator. An optional internal aeration unit is available that can produce an air | CO₂ mixture. Both compressed air and ambient air can be used for this. An IR sensor is used to measure the CO₂ concentration. The BIOSTAT® RM 20 | 50 basic can also be combined with the BIOSTAT® RM Tower control unit, which is operated via an intuitive touch screen and controls the process via single-use sensors. It also includes the software for automatic measured value recording (MFCS/DA).

If the rocker is connected to the BIOSTAT® RM Tower it can be controlled by it. "BIOSTAT® RM serial port active" appears on the rocker display.


All target values and control functions can be displayed and controlled via the BIOSTAT® RM Tower.

For more information about connecting the rocker to the BIOSTAT® RM Tower, please see the relevant operating instructions.

Cell cultivation takes place in the chamber of the single-use CultiBag RM bioreactor, which was specially developed for use with the BIOSTAT® RM 20 | 50 basic. CultiBags RM come gamma-sterilized and ready for use. The film that the CultiBag RM is made of corresponds to the USP Class VI requirements. The volume data for the CultiBag RM always refers to the gross volume. The maximum operating volume is 50% of the gross volume because the remaining space acts as the headspace. The CultiBag RM is designed for single-use. This does away with the cleaning and sterilization required for a conventional bioreactor, saving considerable effort and labor. Thanks to the single-use design of the CultiBag RM, cross-contamination can be prevented and process reliability increased.

The BIOSTAT® RM 20 | 50 basic is designed to only use CultiBags RM.


CultiBag RM 20L, view from the top


- 1 Dip tube for cell removal
- 2 Luer septum for sampling or inoculation
- 3 Luer connection
- 4 C-Flex tubing with a Luer connector (8" × 4"),
- 5 C-Flex tubing with a MPC connector (3/8" × 5/8"), with dip tube
- 6 Sterile filter supply air
- 7 Sterile filter exhaust air
- 8 Valve

1.3 System Design

The next figure shows the schematic design of the BIOSTAT® RM 20 | 50 basic:

A BIOSTAT® RM 20 | 50 basic package consists of the following components:


- ▷ Rocker 20 | 50
Power Connection: 100-240 VAC | 600 W, fuse: 6.3A slow-acting, frequency: 50-60 Hz
Noise level according to IEC 61672:2003:
 - System switched on, idle, aeration off: 39 dB
 - System switched on 42 RPM, aeration off: 62dB
 - System switched on, idle, aeration on: 54 dB
 - System switched on, 42 RPM, aeration on: 63dB
 Material: housing ABS, structural elements chromium steel and aluminum
Protection class: IP23


- ▷ Bag Holder 20 or
▷ Bag Holder 50
Material: ABS
Clamping bars: chromium steel

DZ020L-R2BH
DS050L-RMBH

- ▷ Aeration module
with air | CO₂ mixture
Mixture regulation: 0.8 - 15 % CO₂
Gas flow rate: 1 × 1000 ml/min oder 2 × 500 ml/min

DZ050L-R2CO

- ▷ Heating Panel 20 | 50
Power supply: 48 VDC, Device: H1 | H2
Output: 140 W | Heating Panel
Overheating protection: bimetal trip 70°C +/- 5°C

DZ020L-R2HP


- ▷ Hood:
 - White hood for Holder 20 and | or
 - White hood for Holder 50 and | or
 - Material: ABS

DS020L-R2WL
DS050L-R2WL

- ▷ Surface temperature sensor
 - Type Pt100/3 conductor | Class A
 - Protection class: IP67
 - Connection: Device | Temp1 | Temp2
 - Material: POM white, cable PFA, sensor surface chromium steel 1.4404

DZ-----R2PT

- ▷ Filter heater for Rocker 20 | 50
 - Electric resistance heating surface, 24 VDC | 6 W,
 - with automatic temperature control for approx. 48 °C

DZ-----R2FH

Caution!

Only use in combination with Rocker 20 | 50. Device connection FH1 | FH2.

- ▷ Tubing set for aeration
- ▷ Country-specific connection cable
- ▷ Quick reference guide
- ▷ CD-ROM with operating instructions and quick reference guide

CultiBags RM are not included in the scope of supply.

Bags must be ordered separately. When ordering, please specify the type (Basic, Optical or Perfusion).

For more information, please contact your Sartorius Stedim Biotech representative.

1.4 Rocker 20 | 50

1.4.1 Operation

The rocker is a fully-automatic, PLC-controlled system which is operated via a touch screen. The touch screen is on the front of the device and uses a hierarchical menu structure with a backlight. The integrated functions include:

- Setting the rocking rate and angle
- Individual temperature control for 2 bags (2, 10L) or 1 bag (20, 50L) depending on room temperature
- Aeration of 2 bags (0–500 ml/min) or 1 bag (0–1000 ml/min)
- Setting the bag configuration
- Integrated air | CO₂ mixture (optional)
- Platform positioning for harvesting and sampling
- Alarm display
- 3 different user levels
- Trend display for the display of data
- Time and date display
- Selection of control mode: local or DCU
- Service interval display

1.4.2 Aeration control (only with optional aeration module)

The target value for the flow rate is freely selectable. A suitable tubing set with connectors for the CultiBag RM are included in the package. An integrated pressure sensor controls the overpressure in the CultiBag RM and sets the aeration to achieve a max. of 30 mbar of overpressure.

1.4.3 Air and CO₂ Supply (only with optional aeration module)

Only available with the optional aeration module.

For the air intake, you can select either an integrated air pump that supplies the bag with room air or a compressed air supply on the rear side of the rocker housing for external air intake. The overpressure in the CultiBag RM is monitored via an integrated pressure gauge (for details, see ► "7.2 Overpressure Control" on page 47).

You must observe the gas pressure values listed in the operating instructions. This is the only way to ensure proper functioning of the system.

- ▷ Room air supply: via a membrane pump which compresses the air at the rear "AIR IN" connection
- ▷ Compressed air: integrated pressure regulator for connection to an external compressed air source, + 1.5 bar
- ▷ Quick release coupling: for the direct connection of a rigid tube (AD 4.0 mm)
- ▷ CO₂ integrated pressure regulator for connecting to an external CO₂ source.
1.5 bar quick release coupling for the direct connection of a rigid tube (AD 4.0 mm)
- ▷ The gasses used must be dry as well as dust and grease free.
Use a filter if necessary.

1.4.4 Aeration Rate (only with optional aeration module)

Mass flow controller for flow speeds of 50–1000 ml/min;
Reading accuracy ± 5 %

1.4.5 CO₂ Aeration (only with optional aeration module)

The CO₂ concentration controller has a target value range of 0.8 % – 15%, the display shows the % share of CO₂ in the total gas flow. Reading accuracy ± 5%.

1.4.6 Rocker Speed

8–42 ± 1 rocker movements/min, can be set on the display

1.4.7 Rocker Angle

4–10°, ± 0.3°, can be set on the display.

1.4.8 Temperature Control

PLC-programmed PID controller. The heating plates can be used to heat an individual 2L or 10L bag on one side or two 2L or 10L bags or one 20L or 50L bag on both sides. The temperature of the heating panel can be controlled via the touch screen on the front of the device. Only room air is used for cooling. No additional device is required. A heating | cooling function with water cooling is available as an option for the BIOSTAT® RM 20 | 50 optical and perfusion models

1.4.9 Temperature Range Temperature input|Converter + 15 ... + 40 °C, ±0.2 °C
Minimal temperature depends on the ambient temperature

1.4.10 Serial Interface 1 × RS-232
Optional: Profibus DP, Ethernet

1.4.11 Housing Material Frame: Stainless steel
Outside housing: ABS (UL compliant)

1.4.12 Dimensions | Weight
Incl. Bag Holder 20 and hood for Bag Holder 20:
W × D × H 765 × 600 × 400 mm | 30.0 kg
Incl. Bag Holder 50 and hood for Bag Holder 50:
W × D × H 1085 × 600 × 450 mm | 31.3 kg

1.4.13 Power supply | Power Disconnect Power switch on the back of the device. Automatic voltage selection.
100 – 240 VAC | Fuse 6.3 AT

Caution!
Make sure that the power supply is equipped with a ground lead (IEC 61001-1).
Never use the device without a ground lead

1.4.14 Room Temperature Range 0°C to +40 °C

1.5 Bag Holder

The bag holder is a platform that holds the Cultibag during operation. The bag holder includes both heating elements and the bag is attached to both sides using the lever. Bags can be replaced at any time without the need for tools. The bag holder is made from ABS and comes in a 20 and 50 version.

1.5.1 Properties, Specifications

Bag Holder Dimensions in mm Weight	CultiBag RM basic, optical or perfusion	Number of mountable bags	Min. operating volume in the CultiBag RM [Liter]**	Max. operating volume in the CultiBag RM [Liter]**
Holder 20 765 × 600 × 60 5.5 kg	CultiBag RM 2L	1 (l*) or 2 (r*)	0.15	1
	CultiBag RM 10L	1 (l*) or 2 (r*)	0.5	5
	CultiBag RM 20L	1 (c*)	1	10
Holder 50 1085 × 600 × 60 7.8 kg	CultiBag RM 50L	1 (c*)	5	25


* l = left side, r = right side, c = center

** "Optical" and "Perfusion" CultiBags RM require a higher min. volume due to the integrated single-use sensors.

The BIOSTAT® RM 20 | 50 basic comes standard with a dual temperature measurement and control system. If the rocker is used with a BIOSTAT® RM DCU Tower "Single" version, only the left temperature measurement function can be used. A Twin (single rocker) version tower is required for dual bag control.

1.6 Heating Elements

The heating elements H1 and H2 are integrated directly into the bag holders. They use electric resistance heating that has an integrated temperature safety switch.

The reactor has 2 heating circuits that can be controlled separately (depending on the bag configuration – see ► "4.4.7 Settings Menu" on page 34).

The target value is specified via the "Set" button. Once set, the controller can be activated. If the Pt100 is not connected, this will trigger an alarm.

Control values can range from 15.0°C to 40.0°C.

Heating output: two heating elements, each with 140 W.

Power supply: 48 VDC, connection to device: H1 and H2

Overheating protection: bimetal trip 70°C +/- 5°C

Only use heating elements in combination with Rocker 20 | 50.


Please note that in an uncontrolled state the surfaces of the heating elements can reach a temperature of up to 75°C!

1.7 Hood for Bag Holder 20 or 50

The hood is a safety cover with an opening on the front for bag handling (e.g. sampling). The hood protects the attached bag from mechanical influences during operation and reduces heat loss.


The hood is made from ABS and is also available in white. The hood should always be used during operation.

1.7.1 Dimensions

Hood for Bag Holder 20:

Dimensions: 765 × 600 × 255 mm,
Weight: 2.5 kg

Hood for Bag Holder 50:

Dimensions: 1085 × 600 × 255 mm,
Weight: 3.7 kg.

1.8 Surface Temperature Sensor


The three-stage, Class A Pt100 resistance thermometer (EN60751) measures the CultiBag RM temperature. The thermometer including the cable is run through the opening into the Bag Holder 20 | 50 and connected to the Rocker 20 | 50 via a plug connection. The head is secured in the opening via an O-ring. Connection to device: Temp 1 and Temp 2, protection class: IP67

Make sure that the surface of the Pt100 sensor makes full contact with the underside of the installed bag in order to ensure correct and safe system operation.

1.9 Filter Heater

The heating of the CultiBag RM produces slight condensation on the exhaust filter. The filter heater keeps the filter dry and prevents it from being blocked. A green LED indicates when the filter heater is in operation. When the LED flashes, this means it is heating up. Rapid flashing indicates: overheating or defective electronics.

The rocker PLC checks to ensure that the filter heater is inserted as soon as aeration starts. If this is not the case, an alarm message is triggered.

Temperature range: 40–50 °C, ± 5 °C

Material: housing, polycarbonate, heating output: 6 Watt

Power supply to device: connection FH1 and FH2

24 VDC | 6 W. The filter heater is only suitable for use with the rocker.

Make sure that the filter heating on the exhaust filter is correctly installed and connected to the rocker, before the bag is aerated!

This is the only way to ensure problem-free operation of the system.

2. Delivery and Installation Instructions

2. Delivery and Installation Instructions

2.1 Delivery Check

The BIOSTAT® RM 20 | 50 basic is delivered after conducting extensive function tests. Please proceed as follows if there is transport damage or any device defects appear after installation:

1. Check the order confirmation and shipping documents for correctness.
2. Check all system components for damage.
3. Inform Sartorius Stedim Biotech Customer Service regarding which parts are defective or missing.
4. Check the suitability of all units and individual parts before the first cultivation process.
5. Document the defects and malfunctions carefully and provide this information to Sartorius Stedim Biotech Customer Service.

2.2 Space Requirement

The following table contains information about the space requirement for the BIOSTAT® RM 20 | 50 basic:

Bioreactor configuration:	Space requirement W × H × D (mm):
BIOSTAT® RM 20 basic	765 × 600 × 580
BIOSTAT® RM 50 basic	1085 × 600 × 580

Make sure that the setup location has the required dimensions and is level. This is the only way to ensure problem-free operation of the system.

2.3 Connecting to Supply Lines in the Laboratory

The connections required are included in the scope of delivery. Additional connections are available for special requirements. Please contact Sartorius Stedim Biotech Customer Service for further information.

2.3.1 Work Area Preparation

1. Make sure that there is enough work area available for the setup and installation of the bioreactor:
 - ▷ Placement of the BIOSTAT® RM 20 | 50 basic
 - ▷ Connection to laboratory lines
 - ▷ Placement of the CultiBag RM
 - ▷ Additional peripheral devices
2. Check the power connections:
 - ▷ The main connections in the laboratory | work area must conform to the specifications of the bioreactor.

2.3.2 Power Connection

Danger of electrical damage to device!

The power supply must correspond to the requirements of the device. Voltage fluctuations may not occur. Make sure that the power connection is equipped with a ground lead. Make sure that a correct and disruption-free power supply is guaranteed. Only the named components may be connected to each other. Never connect to the power supply when the device is switched ON.

Check the labels on the devices to ensure that the correct power connections were delivered and that the power cables required for the supply connections in the laboratory are equipped with the proper connections. Please contact Sartorius Stedim Biotech Customer Service, if the devices are equipped with incorrect plugs and cables.

2.4 BIOSTAT® RM 20|50 basic

2.4.1 Connections, Interfaces

Fig. 1: Schematic representation of the back of the BIOSTAT RM

Fig. 2: Schematic representation of the left side of the BIOSTAT RM

Fig. 3: Schematic representation of the right side of the BIOSTAT RM

Key:

- SIG 1: Reserve (LED connection or Loadcell)
- D-link: Optional Ethernet
- SIG 2: Alarm output, potential-free, max. 0.5 A
- SIG 3: Reserve (weight measurement)
- D-LINK 2: DCU/ MFCS connection D-SUB
- D-LINK 3: Optional Profibus DP D-SUB
- AIR: Compressed air connection max.1.5Bar | air pump intake
- CO2: Gas connection max. 1.5 (possibly O2)
- T6.3A: Main switch | Power connection, 100-240VDC broad band. Mains adapter, fuse 6.3A slow-acting
- GAS1: Gas mixture outlet to RM bag
- TEMP 1: PT100 temperature sensor connection
- H 1: Heating panel 48V | 3.15A
- FH 1: Filter heater 24V | 800mA
- GAS 2: Gas mixture outlet to RM bag
- TEMP 2: PT100 temperature sensor connection
- H 2: Heating panel 48V | 3.15A
- FH 2: Filter heater 24V | 800mA

3. Initial Startup

3.1 Overview

3. Initial Start-up

The following steps are included in the initial start-up and operation of the BIOSTAT® RM 20 | 50 basic:

1. Rocker setup
2. Installation and inflation of the CultiBag RM
3. Cultivation start
4. Inoculation, sampling and cultivation process
5. Replacement of the culture medium (optional)
6. Completing the process, cell harvesting and shutting down the device under safe conditions
7. Cleaning and maintenance

3.2 Installation Kit

The scope of delivery includes the required connections and connecting pieces.

1. Only use connections and accessory parts for the BIOSTAT® RM 20 | 50 basic that have been approved by Sartorius Stedim Biotech.
2. Only replace connections and accessory parts with products that have been supplied and/or approved by Sartorius Stedim Biotech.

Malfunctions or defects that are caused by the use of unapproved parts for the BIOSTAT® RM 20 | 50 basic are not covered by the warranty.

3.3 Rocker Design

3.3.1 Installing the Bag Holder

1. Place the bag holder on the rocker and make sure that it is positioned between both holding points on the front and back of the basic unit in the appropriate openings to prevent the holder from shifting.
2. Place the heating mats in the openings of the bag holder. Guide the cable through the appropriate openings. Connect the right and left cables to the sockets labeled H1 and H2.

3.3.2 Installing the Surface Temperature Sensor

1. Attach the left sensor to the bag holder. Do this by guiding the sensor cable through the round opening in the center of the heating element. Press down on the sensor from above until it is attached in the opening via the o-ring. The sensor is positioned correctly when its silver side is pointing up so that it will later lie directly on the CultiBag RM as shown in the figure.

Note

If the sensor is positioned incorrectly, this can lead to overheating of the heating panel. In this case, the heating will stop automatically via a bimetal trip when the temperature reaches 70 +/-5°C.

2. Now connect the sensor plug to the socket reserved for it on the left side of the basic unit.
3. Procedure as above for the right PT100 sensor.

3.3.3 Power Supply Connection

1. Connect the power cable to the rocker and to the power source as shown in the figure. Please also observe the note in section ► "2.3.2 Power Connection" on page 19.


3.3.4 Connecting the Compressed Air and CO₂ Supply

1. If you are using a compressed air supply, connect it to the "Process Air Inlet" connection on the back of the basic unit.
2. Connect the CO₂ source to the corresponding input of the rocker.
3. Turn on the pressure for the gasses being used according to the info in section ► "1.4.3 Air and CO₂ Supply (only with optional aeration module)" 1.4.3 Air and CO₂ Supply (only with optional aeration module)" on page 14.

Please note that the device begins operating again automatically using the set parameters after a power outage. This function can be activated/deactivated in the System menu under "Automatic System Restart."

3.3.5 Device Startup


- ▷ Make sure that the device has been set up in accordance with the steps listed above and that all components have been properly installed.
- ▷ The device is now operational.
- ▷ The "0-1" switch on the back left of the device must be used to turn the device on/off and/or disconnected it from the power supply.

Warning to avoid hand injuries!
Note that inserting your hands between the bag holder and the rocker is prohibited during operation. Failure to comply with this may result in material damage and personal injury.

4. Menu Operation


4.1 The Main Menu

4. Menu Operation

Welcome Screen


This is displayed for approx. 5 seconds when the device is turned on. It includes the firmware version of the CPU, the configuration version of the display and the bag configuration. The latter includes information about which bags can be used with the current bioreactor settings (control parameters).

Main Screen


This is the main operating display. It can be used to access each submenu. For example, when the Speed button is pressed, it takes you immediately to the Speed menu where you can start or stop the controller. Then you return to the Main screen. The Main screen always displays the nominal values of the various parameters. Buttons that are greyed indicate that the function is deactivated or that the user does not have authorization to this function.

Note: The symbol "*" on the screen illustrations stands for the actual value.

Overview page:

When you press the "Overview" button, it takes you to this page.

It displays the most important parameters as actual values. A green LED next to the nominal values indicates whether or not the respective controller has been started.

Pressing the blue Home button returns you to the Main page.

The rocker display is only used in conjunction with the BIostat® RM 20 | 50 basic for data entry. If the rocker is attached to the BIostat® RM Tower, the Tower touch screen can be used for operation. However, the control mode of the bioreactor must be set to "DCU Tower" in the menu item ► "4.4.7 Settings Menu" on page 34.

4.2 Preparing and Inflating the CultiBag RM

4.2.1 Unpacking the CultiBag RM

Remove the sterile single-use CultiBag RM bioreactor from the plastic protective film.


Caution!
Do not damage the sterile cultivation chamber when removing it. Carefully open the outside packaging. Do not use a sharp instrument as this could damage the CultiBag RM.

4.2.2 Securing the CultiBag RM to the Bag Holder RM

1. Open the guide rails on the bag holder by releasing the lever on both sides of the holder as shown in the figure.
2. Now place one or two CultiBags RM in the bag holder and press the plastic bars on both sides of the bags into the now open guide rails of the bag holder. If you are only placing one CultiBag RM in the holder, place it on the left side of the holder to ensure error-free temperature measurement. For information regarding the number and size of CultiBags RM that can be attached to the rocker, please see section ► "1.5.1 Properties, Specifications" on page 15.
3. Now press the lever down again to secure the CultiBag(s) RM.

4.2.3 Installing the Filter Heater

1. Insert the exhaust filter between both parts of the exhaust filter heater. Both parts are held together by magnets. (The exhaust filter is equipped with an overpressure valve.)

2. Two magnets close at the opening of the filter heater to ensure a secure hold.

3. Now connect the filter heater to the corresponding socket on the left and/or right side of the rocker (see figure). The exhaust filter heater of the left bag is connected to the left side. The exhaust filter heater of the right bag is connected to the right side if it is being used.

A green LED indicates when the heater is in operation. If the LED flashes, this indicates the warm-up phase. If the LED lights up continually, it indicates that it has reached operating temperature.

4. The rocker controller queries whether or not the filter heater is connected as soon as aeration starts. If this is not the case, an alarm is triggered.

The exhaust filter is now heated to >40 °C. This helps to ensure that the filter does not become blocked.

4.2.4 Connecting the Air Supply

1. Connect the "Gas 1 and/or Gas 2" connection of the BIOSTAT® RM 20 | 50 basic to the sterile exhaust filter on the CultiBag RM using the supplied tubing. (The exhaust filter is the filter without the overpressure valve).
If two CultiBags RM are being used, use both rocker outputs.

Caution!
Connecting the bag (CultiBag RM) directly to a gas source or using another air pump can create overpressure. This overpressure could burst the bag because the integrated pressure sensor now cannot send a signal to the external gas source to stop the aeration for overpressure in the CultiBag RM. This situation would not be covered by the warranty.

4.2.5 Inflating the CultiBag RM

Make sure that the tubing clamps on the intake air filter and exhaust filter are open and close all additional tubing clamps on the bag.

1. Turn the device on via the main switch on the back. The Start menu appears on the display for a few seconds showing the respective program version of the device before switching to the Main menu.

2. Start the aeration as described in the following chapter ► "4.4.1 Aeration Menu" on page 27.

Caution!

The internal pump should not be selected if another compressed air source is connected to the device. If the internal pump is activated and compressed air is connected to the device, this can cause damage to the pump.

3. Wait until the CultiBag RM has inflated completely and make sure that it is secured to the bag holder via the lever.

The CultiBag RM should be fully inflated, i.e. there should be no folds or kinks. Make sure that the exhaust filter is functioning by pressing down gently on the bag chamber and making sure that air escapes. Once the CultiBag RM is inflated, the air intake can be shut down when the desired process values have been reached.

4.3 Cultivation Start

4.3.1 Filling the CultiBag RM with the Culture Medium

The CultiBag RM should be inflated before being filled with the culture medium (see section ► "4.2.5 Inflating the CultiBag RM" on page 25).

This reduces the foaming of the medium. Regulate the air intake so that the CultiBag RM remains inflated during filling and does not collapse. The rocker function should be turned off during filling.

1. Make sure that there is an aseptic connection between the chamber of the CultiBag RM and the media container by using a BioWelder, Luer connections or MPC connections to connect both tubes. If you are using Luer or MPC connections, you must work at a sterile workbench.

To ensure aseptic conditions, the CultiBag RM must be filled at a sterile workbench if necessary. The Luer connection for sampling can also be used for this.

2. Pump the desired amount of the culture medium into the bioreactor and make sure that the media container is also separated from the bioreactor under aseptic conditions.

Make sure that the tubing clamps of the air intake filter and exhaust filter are open when filling the medium.

3. Place the hood on the bioreactor.

4.3.2 Setting the Cultivation Parameters and Starting Cultivation

1. Turn the rocker on at the main switch. The main switch is located on the back of the rocker.
2. Select the size of the CultiBag RM you are using in the CultiBag menu as described in chapter ► "4.4.7 Settings Menu" on page 34.
3. Set the rocker speed as described in the menu ► "4.4.2 Speed and Angle Menu" on page 29 to ensure visible wave movement. Turn on the rocker. Bei starker Schaumbildung muss die Wippgeschwindigkeit reduziert werden. Der CultiBag muss komplett aufgeblasen sein, da sich die Schaumbildung durch Falten oder Unebenheiten verstärkt.
4. Set the cultivation temperature in the Temperature (T) menu as described in ► "4.4.3 Temperature Menu" on page 30.

Start the heating process after you have made sure that the surface temperature sensor is properly placed under the CultiBag RM. Otherwise, there is a danger of overheating. Uniform heating of the culture medium is ensured when the rocker movements are started before activating the heating panel.

4.4 Operating the Control Elements


4.4.1 Aeration Menu

4.4.1.1 Aeration with Air

Flow 1 and 2 menu:


There are 2 flow cycles on the reactor that can be controlled independently. The target value is specified via the "Set" button. The target value can range from 50 ml/min. to 1000 ml/min. If both flow cycles are active, the max. target value per control cycle is limited to 500 ml/min. Once set, the "Gasflow control" can then be activated. The air supply is defined using "supply over." "Process air" means that the reactor is supplied by external process air (1.0 ... 1.5 bar). Pressing the button switches to the internal air pump. If this pump is to be used, the user must make sure that no external compressed air source is connected because this input is used as the intake line for the pump.

The "Bag Inflation" function can be used to automatically inflate an unfilled bag. To do this, the correct bag configuration must be specified in the Settings menu. A timer function is then used to pump the proper amount of air into the bag. The aeration will also stop when the pressure in the bag reaches 30 mbar.

The red bar indicates the present activity of the valve.


The trend for Gasflow1 and Gasflow2 is set up with a time window of 30 min.


4.4.1.2 Aeration with CO₂

CO₂ menu:

The target value is specified via the "Set." The target value range is 0.8% to 15.0%. The % display of CO₂ is the amount in the entire gas flow.

The trend for CO₂ is set up with a time window of 30 min.

4.4.2 Speed and Angle Menu

Speed and Angle menu:


The "rocking process" can be started and stopped in both menus. The controller regulates the rocking rate as well as the angle for the desired target values. Press "Set" to access the keypad page and enter target values.


Speed control range: 8–42 rocks

Angle control range: 4.0–10.0 °deg


1 rock corresponds to one movement from the center to the front position – then to the back position – and then back to the center.

The "Sample Position" function automatically moves to the sampling position. When this button is pressed, the platform moves to a suitable 10° forward position. At the same time, the heating is turned off to avoid inhomogeneous temperature distribution. When the button is pressed again, the process continues. The Sample function is also equipped with an automatic safety timer. Once the time set in the Sampling Time menu ► "4.4.7.5 Sampling Time" on page 38 expires, the rocker turns back on automatically.


The trend for speed and angle is displayed together. A max. of 30 min. is displayed for this. This scaling is fixed and cannot be changed by the user.


4.4.3 Temperature Menu


Temperature 1 and 2 menu:

The reactor has 2 heating circuits that can be controlled separately (depending on the bag configuration – see ► “4.4.7 Settings Menu” on page 34). The target value is specified via the “Set” button. Once set, the controller can be activated. If the Pt100 is not connected, this will trigger an alarm. Target values can range from 15.0 °C to 40.0 °C.

The trend for Temperature1 and Temperature2 is set up with a time window of 30 min.


4.4.4 Position Menu

Position menu:

This is accessed via the Main page. The Position process is used to move and hold the pan to a specific tilted position. The Rocking process must be deactivated in order to start positioning.

Then "Start pos process" can be activated. The target value is specified via the "Set" button. Then the direction is selected. If the positioning operation must be stopped, "Start pos process" must be deactivated. This returns the pan to the center position.

Warning to avoid hand injuries!

Note that inserting your hands between the bag holder and the rocker is prohibited during operation. Failure to comply with this may result in material damage and personal injury.

4.4.5 User Level Menu

User level page


The User level is displayed on each screen. If the system is "locked," no changes can be made to the system. You must press the "Enter password" button to switch to the User level. This accesses the keypad page. You can enter up to 8 characters. After entering the password, press the "OK" button to confirm.

The system then switches automatically to the User level associated with the password.

The "Lock user" button is used to re-lock the system. You can change standard passwords with supervisor rights. Use "Create new passwords" to switch to this page.

User level

Keypad page

Press "Enter" to enter the password or "ESC" to cancel.

Standard passwords factory setting:

User: SART02

Supervisor: SART01

There are 3 user levels.

1. Locked ... no system changes can be made.
Most menus are also locked.
2. User ... all parameters can be changed and all menus can be accessed except for the PID parameter settings.
3. Supervisor ... all functions are released incl. the PID parameter settings.
The supervisor can also create new passwords for users and supervisors.
Calibration functions can only be access by service technicians.

Create new passwords:

This is used to assign a new password for the User level.


To enter a new password, press "Enter new password"

to open the entry page. You must re-enter the password in "Re-Enter new password" to confirm.

Then press "Save new password" to apply it.

Pressing this button switches from "User" to "Supervisor."

4.4.6 Alarm Menu

Alarm page:


This is accessed from the Main page via the "Alarm" button. You can view the last alarms to occur by using the arrow keys to move through the list.

This screen is displayed automatically by the system when an alarm or message occurs.

The alarm or message is confirmed by pressing the Home button.

For example, if an unfavorable status is detected for the reactor, i.e. overpressure in the bag, the flow regulator is shut down and the Alarm page is displayed with the current alarm along with its entry.


After leaving the Alarm page, the reason for the overpressure must be found and the flow regulator restarted after correcting the problem.

4.4.7 Settings Menu

Settings:

The "Settings" can be accessed from the Main page.


4.4.7.1 Setting the Date and Time

The date and time can be set here.

To do so, press the entry fields. This opens the keypad page for entry. Once all entries have been made, press the "Set date and time" button to confirm.


This transfers the data to the CPU.

Automatically Summer- | Wintertime change:
 This setting is used to automatically adjust the time for daylight savings.
 Please note that the device must be in operation at the time for the automatic change to occur.

4.4.7.2 Service Information

The Service page provides information about operating hours, number of movements and when the next service is due.


These displays are reset by the service technician once service is complete. (The reset takes place at a password-protected technician level).

- ▷ If service does not take place on time, "Service required" scrolls across the screen every 3 minutes.
 This does not deactivate any functions. It is just a notice regarding the need for service.

4.4.7.3 Setting the Bag Size

The bag configuration is defined here and displayed on each screen. The following configurations are available:

- | | |
|----------|---|
| 1 × 1L: | Only Heater1 (left heating panel) in operation, Pt100_1 and Temp_1 controller in operation |
| 1 × 2L: | Only Heater1 (left heating panel) in operation, Pt100_1 and Temp_1 controller in operation |
| 1 × 10L: | Only Heater1 (left heating panel) in operation, Pt100_1 and Temp_1 controller in operation |
| 1 × 20L: | Heater1 and 2 (left and right heating panel) in operation, Pt100_1 and Temp_1 controller in operation |
| 1 × 50L: | Heater1 and 2 (left and right heating panel) in operation, Pt100_1 and Temp_1 controller in operation |
| 2 × 1L: | Pt100_1 and Temp_1 controller controls Heater1 (left heating panel), Pt100_2 and Temp_2 controller controls Heater2 (right heating panel) |
| 2 × 2L: | Pt100_1 and Temp_1 controller controls Heater1 (left heating panel), Pt100_2 and Temp_2 controller controls Heater2 (right heating panel) |
| 2 × 10L: | Pt100_1 and Temp_1 controller controls Heater1 (left heating panel), Pt100_2 and Temp_2 controller controls Heater2 (right heating panel) |


4.4.7.4 Selecting a Control Mode

The Control mode via the CPU can be selected here.


In this case, the display is selected as the master, i.e. all commands and parameters are entered via the display as well as all alarms and warning messages.

When you switch the "ON" button, the DCU Tower is activated as the master.

After 3 seconds, it automatically switches to the next page.


This displays the versions, bag configuration and the connection status. A successful connection is indicated by a green OK checkmark. The device automatically switches to this page after a device restart. All commands, target values and alarms are now managed via the DCU Tower.


When service is performed or calibration work is carried out on the rocker, this must be carried out via the rocker display. To do this, press the "Choose controlmode" button to return to Control selection (not possible via the DCU Tower due to the amount of data involved).

This defines whether or not controllers for Speed, Temperature 1 and 2, Gasflow 1 and 2 and CO₂ are restarted automatically after the device is restarted.


4.4.7.5 Sampling Time

This specifies the time when the rocker is restarted automatically when using the Sampling function ► "4.6.1 Preparing Sampling or Inoculation" on page 39.

4.4.7.6 Batch Time

The device has a timer that can be started in this menu and displays the batch time that has expired. The batch time is also displayed in the Overview menu.

4.5 Cultivation Parameters

Each cell line requires optimization of the cultivation conditions.

The following table contains recommendations for the aeration rate and the rocker angle to start the cultivation process (room air with an oxygen content of 20 %).

Liquid volume (L)/CultiBag RM	Rocker speed [r/min]	Aeration rate [l/min]	Rocker angle
1/2	15–20	0,1	6–8
5/10	15–20	0,1–0,2	6–8
10/20	15–20	0,2	6–8
25/50	15–20	0,3–0,5	6–8

Air intake (aeration rate) and CO₂ concentration

An oversized air intake only has a minimal effect on the oxygen intake. The aeration rate should be kept to a minimum to reduce the evaporation of the medium. The recommended aeration rate is 0.1–0.2 l/min. We recommend a CO₂ concentration of 5% because this is used for most mammalian cells. However, this must be adjusted to each cell line and medium.

Rocker speed

The rocker speed must be set to the minimum speed that ensures thorough mixing and oxygen intake without leading to excessive foam buildup and the associated shearing forces. At the start, the rocker speed is normally 15–20 rocker movements/min. The rocker speed required for the oxygen intake is designed for the respective cell population and its metabolism. Normally, 15–25 rocker movements/min. are enough to ensure the oxygen requirement of cell thicknesses up to 5×10^6 cells/ml.

Rocker angle

A smaller rocker angle can minimize foaming of the medium. The rocker angle can be reduced 2–3° for high foaming media. The rocker angle can be increased up to 8° for cells with a high oxygen requirement, e.g. insect cells. The maximum recommended rocker angle is 10°.


4.6 Sampling and Inoculation

4.6.1 Preparing Sampling or Inoculation

This procedure requires a standard laboratory syringe with a Luer septum without a needle. This procedure must not be carried out on a sterile workbench.

1. Stop the rocker movement. The rocker has an easy-to-use Sampling function for this, which is explained in more detail in section 4.4.2.
2. Do not interrupt the surface aeration.
3. Remove the protective cover from the Luer septum of the sampling tube.

4. Disinfect the Luer septum: Wipe or spray the upper part of the Luer septum with 70% alcohol (or similar medium).

5. Connect the syringe: connect a sterile, single-use syringe without a needle to the Luer connection under aseptic conditions.

6. Open the tubing clamp of the Luer septum.

4.6.2 Inoculation

1. Fill the CultiBag RM with the syringe contents via the Luer septum. Then remove the syringe from the Luer septum.
2. Disinfect the Luer septum again and re-close it using the protective cover.
3. Let the remaining liquid in the tubing run back into the bioreactor. Press the tubing together a few times to ensure that any remaining liquid runs back into the bioreactor. Then close the tubing clamp.
4. Now restart the rocker movement. Set the speed so that a wave motion is visible on the surface of the liquid. Decrease the speed if too much foam builds up.

4.6.3 Sampling

Sampling is carried out in the same way as inoculation (syringe without a needle is also used here). However, you must make sure that the cells which are inside the tubing are removed before taking your representative sample. For this reason, take two samples as described in this procedure. The first one you will throw away and the second is your representative sample.

The Luer septum is designed for 50 sterile sampling and/or inoculation procedures. It contains a silicon septum that seals the tubing automatically when the connection is severed. The sample is also protected completely because no aerosols can get in from the outside. At lower rocker speeds (< 15 Rocker movements/min), the cells can settle so that the sampling may not be representative. Therefore, the rocker speed should be increased to 15 rocker movements/min. approx. 5-10 minutes before the sampling. Then it can be reduced again.

If you are working with large micro carriers or large cell aggregates, you should not use the sampling equipment because it has a relatively small opening. You should use the tubing with the MPC connection (outside diameter 11.1 mm).

4.7 Completing the Process, Cell Harvesting and Shutting Down the Device under Safe Conditions

1. Stop the rocker movement as described in section ► "4.4.2 Speed and Angle Menu" on page 29. The bag holder stops automatically at a slightly forward angle so that the dip tube remains in the medium. For bags with the dip tube at the back, the rocker can be moved to another stop position as described in section ► "4.4.4 Position Menu" on page 31.
2. Turn the temperature control off as described in section ► "4.4.3 Temperature Menu" on page 30.
3. Turn the aeration off as described in section ► "4.4.1 Aeration Menu" on page 27.
4. Remove the hood from the rocker platform.
5. Clamp off the intake air and exhaust filter.
6. Remove the air intake tubing from the sterile intake air filter.
7. There are several options available for emptying the CultiBag depending on the process conditions:
 - ▷ To remove the entire contents via the harvest tubing (with a dip tube), use a suitable peristaltic pump and hold the CultiBag RM in an upright position until it is completely empty.
 - ▷ To release the overpressure, remove the holder with the attached CultiBag RM from the rocker platform and wait 10–15 minutes until the cells have settled. Pump out the excess and add fresh medium or move the cell culture to a larger BIOSTAT® RM bioreactor.
8. Store the CultiBag RM under the prescribed process conditions before you begin with downstream processing.
9. Turn off the aeration with the compressed air (if connected).
10. Turn off the main switch of the rocker.

For your notes

5. Cleaning and Maintenance

Dangerous voltage when opening the device or due to contact with the device during servicing!

The "0-1" switch must be turned off and the power cable disconnected from the device for maintenance and servicing.

5.1 Intermediate Cleaning

The BIOSTAT® RM 20 | 50 basic is a single-use bioreactor system that does not require cleaning or sterilization contrary to conventional, reusable bioreactors. However, the device should be cleaned on a regular basis because soiling from daily operation, e.g. splashes or spilled liquids are unavoidable:

1. The bag holder can be cleaned with alcohol or other common cleaning agents. In addition, the holder can also be carefully rinsed down with water. However, you must disconnect the temperature sensor. Make sure that the temperature sensor always reassembled in its holder after cleaning.
2. Wipe down the hood of the bag holder using a soft moistened cloth. Do not use solvents | liquids that can corrode or damage the material (see also section ► "1.7 Hood for Bag Holder 20 or 50" on page 17).

3. **First disconnect the device from the power supply. Wipe down the rocker housing using a moist cloth. For cleaning, you can use either water or a common disinfecting agent such as 70 % ethanol. Do not use solvents | liquids that can corrode or damage the material (see also section ► "1.4.10 Serial Interface" on page 15). Do not immerse the device in water nor rinse it down with water. Only use a moistened cloth for cleaning. Make sure that water does not get into the housing because this can cause serious damage to the electrical components.**

5.2 Servicing the Device by Operating Personnel

The maintenance activities to be carried out by operating personnel are limited to:

- ▷ Replacing the aeration tubing between the "AIR OUT" connection and the sterile air intake filter of the CultiBag
- ▷ The immediate replacement of the following defective parts:
 - Filter heater for the Rocker 20 | 50
 - Bag Holder 20 and 50
 - Hood for Bag Holder 20 and 50
 - Surface temperature sensor
 - Heating panel 20 | 50
- ▷ Regular lubrication is not required.

- ▷ **Service and calibration of the BIOSTAT® RM 20 | 50 basic are reserved for qualified, trained and authorized service personnel.**

5.3 Servicing and Calibration by Authorized Service Personnel

Several components such as the bearing require regular maintenance. We recommend maintenance every 8000 operating hours, including the required calibration of all parameters.

The battery for operating the PLC integrated into the BIOSTAT® RM 20 | 50 basic has a lifespan of more than five years. However, we recommend that it be replaced within this 5-year-period. This service is included within the scope of maintenance if necessary.

Service and calibration of the BIOSTAT® RM 20 | 50 basic are reserved for qualified, trained and authorized service personnel only (as per ► "5.2 Servicing the Device by Operating Personnel").

For more information, please contact your Sartorius Stedim Biotech service technician.

6. Transportation and Storage

The following points must be observed with transporting and storing the BIOSTAT® RM 20 | 50 basic:

**Heavy weight of the BIOSTAT® RM 20 | 50 basic (> 20 kg)!
Improper lifting and transporting of the BIOSTAT® RM 20 | 50 basic can cause injury to the user and damage to the device. Please make sure that the BIOSTAT® RM basic is only lifted using suitable lifting gear. All applicable and country-specific safety regulations regarding the lifting and transport of heavy loads must be observed.**

- ▷ Please note the dimensions and weight of the BIOSTAT® RM 20 | 50 basic for transport purposes (as per ► "1.4.12 Dimensions|Weight").
 - ▷ The BIOSTAT® RM 20 | 50 basic should only be stored in air-conditioned rooms suitable for this purpose.
-

Use the original package for long-term storage.


For your notes

7. Appendix

7.1 Schematic Representation of Aeration

The BIOSTAT® RM 20 | 50 basic can be used to aerate the cell culture with either external compressed air or compressed room air.


A pressure sensor monitors the pressure and is equipped with an alarm signal (see section ► "7.2 Overpressure Control" on page 47).

7.2 Overpressure Control

A pressure sensor records the pressure value of the CultiBag RM (see section ► "7.1 Schematic Representation of Aeration" on page 47).

If the pressure rises to $> +30$ mbar, the control value closes automatically until the pressure decreases to $+20$ mbar. This protects the cell culture from overpressure and protects the CultiBag RM from damage.

7.3 Alarm Messages

The BIOSTAT® RM 20 | 50 basic displays alarm messages on the screen to inform the operator of various malfunctions. When in doubt, please contact your Sartorius Stedim Biotech field employee.

Alarm message	Possible causes	Measures
Battery low	This is displayed when the PLC backup battery is too weak. In this case, you should immediately call a service technician to replace the battery. You can acknowledge this alarm by turning the system off and on.	The battery must be replaced within two weeks to prevent a program loss.
Error Pt100_1	This is displayed when an error occurs in the Pt100 sensor during heating operation or when the sensor is disconnected from the system during heating operation. Temp controller 1 is deactivated for safety reasons in case of malfunction. The controller must be restarted manually after acknowledging the alarm. Press the Home button to acknowledge this alarm.	Connect the surface temperature sensor. If required, you should also test the other Pt100 sensor.
Error Pt100_2	This is displayed when an error occurs in the Pt100 sensor during heating operation or when the sensor is disconnected from the system during heating operation. Temp controller 2 is deceived for safety reasons in case of malfunction. The controller must be restarted manually after acknowledging the alarm. Press the Home button to acknowledge this alarm.	Connect the surface temperature sensor. If required, you should also test the other Pt100 sensor.

Alarm message	Possible causes	Measures
Error Heater1	<p>This is displayed when no power consumption for heating plate 1 is being measured during heating operation. Either the bimetal (integrated into the heating panel) has triggered due to over temperature or the heater was disconnected from the system during heating operation.</p> <p>Temp controller 1 is deactivated for safety reasons in case of malfunction. The controller must be restarted manually after acknowledging the alarm.</p> <p>Press the Home button to acknowledge this alarm.</p>	<p>Check the contact between the surface temperature sensor and the Rocking bag.</p> <p>Wait until the heating surface has cooled before repeating the test.</p>
Error Heater2	<p>This is displayed when no power consumption for heating plate 2 is being measured during heating operation. Either the bimetal (integrated into the heating panel) has triggered due to over temperature or the heater was disconnected from the system during heating operation.</p> <p>Temp controller 2 is deactivated for safety reasons in case of malfunction. The controller must be restarted manually after acknowledging the alarm.</p> <p>Press the Home button to acknowledge this alarm.</p>	<p>Check the contact between the surface temperature sensor and the Rocking Bag.</p> <p>Wait until the heating surface has cooled before repeating the test.</p>

Alarm message	Possible causes	Measures
No CO₂	<p>If the CO₂ controller is active (requires that the Gasflow1 controller is also active), this alarm message is displayed after 5 min when the CO₂ actual value is less than 0.6% (target value min. 0.8%).</p> <p>This does not deactivate the CO₂ controller.</p> <p>You can acknowledge this alarm message via the "Home" button.</p>	<p>Please check whether or not the connection tubing to the device has been crushed or kinked.</p>
Drive System Err	<p>If no movement is detected in the motor within a period of 15 s or an error is detected in the motor amplifier (e.g. overload), this alarm is triggered.</p> <p>This alarm is acknowledged by turning the device off and on (min. 5 s off). This message is also displayed on the screen as a ticker when there is an error.</p>	<p>If the error continues to be displayed even after a device restart, you should contact a service technician.</p>
Error Filter-H1	<p>The filter heater is always in operation when connected to the system. However, it is only monitored during active gas aeration. In this case, it is important that the filter heater functions and this is measured via the power consumption. If filter heater 1 is disconnected from the system (or defective) while the Gasflow1 controller is active, this alarm is triggered after a few seconds. This does not deactivate the controller. This alarm message is only used 1x per gas flow controller start. You can acknowledge this alarm message via the "Home" button.</p>	<p>Make sure that the filter heater is connected to the correct output. If this error should continue to occur, test the function using another filter heater.</p>

Alarm message	Possible causes	Measures
Error Filter-H2	The filter heater is always in operation when connected to the system. However, it is only monitored during active gas aeration. In this case, it is important that the filter heater functions and this is measured via the power consumption. If filter heater 2 is disconnected from the system (or defective) while the Gasflow2 controller is active, this alarm is triggered after a few seconds. This does not deactivate the controller. This alarm message is only used 1x per gas flow controller start. You can acknowledge this alarm message via the "Home" button.	Make sure that the filter heater is connected to the correct output. If this error should continue to occur, test the function using another filter heater.

Alarm message	Possible causes	Measures
Overpr Flow1	<p>If the bag pressure 1 should exceed the upper limit set in the Technician menu (normally 30 mbar), the Gasflow1 valve will close and the alarm message is triggered. This does not deactivate the Gasflow1 controller. If the bag pressure decrease below the lower limit set in the Technician menu (normally 20 mbar), the Gasflow1 valve is reopened and the control process continues.</p> <p>You can acknowledge this alarm message via the "Home" button.</p>	<p>Please check whether or not the connection tubing between the device and bag has been crushed or kinked.</p> <p>Is the filter heater of the respective bag functional and in operation?</p> <p>This can be checked by touching the filter heater – it should feel warm.</p> <p>Did you set a gas flow rate that is too high for a bag that is too small?</p>
Overpr Flow2	<p>If bag pressure 2 should exceed the upper limit set in the Technician menu (normally 30 mbar), the Gasflow1 valve will close and the alarm message is triggered. This does not deactivate the Gasflow2 controller. If the bag pressure decrease below the lower limit set in the Technician menu (normally 20 mbar), the Gasflow2 valve is reopened and the control process continues.</p> <p>You can acknowledge this alarm message via the "Home" button.</p>	<p>Please check whether or not the connection tubing between the device and bag has been crushed or kinked.</p> <p>Is the filter heater of the respective bag functional and in operation?</p> <p>This can be checked by touching the filter heater – it should feel warm.</p> <p>Did you set a gas flow rate that is too high for a bag that is too small?</p>

Alarm message	Possible causes	Measures
No Airflow1	If the Gasflow1 controller is active and no airflow higher than 40 ml/m is measured for 2 min (target value is at least 50 ml/m), this alarm message is displayed. This does not deactivate the controller. You can acknowledge this alarm message via the "Home" button.	When operating with process air – please check the input pressure range. Is it between 1.0 and 1.5 bar? Is the supply tubing to the device crushed or kinked? When operating with an internal air pump – please contact Sartorius Service.
No Airflow2	If the Gasflow2 controller is active and no airflow higher than 40ml/m is measured for 2 min (target value is at least 50ml/m), this alarm message is displayed. This does not deactivate the controller. You can acknowledge this alarm message via the "Home" button.	When operating with process air – please check the input pressure range. Is it between 1.0 and 1.5 bar? Is the supply tubing to the device crushed or kinked? When operating with an internal air pump – please contact Sartorius Service.

7.4 Malfunctions and Countermeasures

Problem	Possible causes	Measures
The CultiBag seems to inflate too tightly.	The CultiBag should inflate so that there are no folds. It should also not be too tight that overpressure occurs and kinks form at the hold points.	<p>If the bag inflates too tightly, check that the sterile exhaust filter is not blocked. Do this by connecting a tube to the filter and placing the other end in a container with water. The filter is blocked if no bubbles appear. Replace the filter on a sterile workbench.</p> <p>In some cases, the return valve to the sterile exhaust filter may be blocked. Unscrew it and replace it with a new one.</p>
The CultiBag seems to be underinflated.	If the CultiBag RM is underinflated, it creates excessive foaming and poor mixing. "No airflow" may appear on the display.	<p>Check the "AIR OUT" output on the rocker to ensure that aeration is functioning. The sterile exhaust filter could be blocked, the closing clamp could be shut or the aeration tubing not properly secured. Replace the exhaust filter on a sterile workbench or fill the culture medium into a new CultiBag.</p>
Extreme foam buildup in the CultiBag RM	A certain amount of foam is normal for the BIOSTAT® RM 20 50 basic. However, the surface of the medium should not be covered with more than 50% foam.	<p>If there is too much foam buildup, check to make sure that CultiBag aeration is functioning properly. Under inflation can lead to strong foam buildup. Normally, the foam will recede after a few hours. If the foam remains unchanged after this time, reduce the rocker speed. However, make sure that it is fast enough to ensure a sufficient oxygen supply for the medium. If there is still too much foam after this, reduce the rocker angle. A smaller rocker angle produces less foam and the rocker speed can again be increased slightly. For more information, please contact your Sartorius Stedim Biotech representative.</p>

Problem	Possible causes	Measures
The bioreactor heats up too slowly.	The culture medium does not heat up or only very slowly.	Check that the surface temperature sensor is positioned correctly. Make sure that the desired temperature is set correctly. See chapter ► "4.4.3 Temperature Menu" on page 30. Make sure that the correct bag configuration is selected. See chapter ► "4.4.7 Settings Menu" on page 34. The LEDs must light up continuously or flash. Make sure that the CultiBag RM is in motion. The heating process for the medium takes place slowly and continually. The device requires over one hour to heat up a 25 liter culture medium from room temperature to 37 °C.
The temperature controller is not functioning properly.	The temperature of the culture medium is not being controlled correctly.	Make sure that the surface temperature sensor is positioned correctly under the bag. Make sure that the CultiBag RM is in motion. If the temperature is still not working properly, recalibration of the temperature measurement may be necessary. Please contact your Sartorius Stedim Biotech representative.
The temperature seems to be incorrect.	The temperature displayed in the Main menu does not match that of the room temperature or an independent measurement.	Check that the surface temperature sensor is positioned correctly. Make sure that the CultiBag RM is in motion. If the temperature is still not working properly, recalibration of the temperature measurement may be necessary. Please contact your Sartorius Stedim Biotech representative.

For any other questions, please contact your Sartorius Stedim Biotech representative.

7.5 Agreements

7.5.1 Agreements and Intended Use

All components of the BIOSTAT® RM 20 | 50 basic are designed for biotechnological processes in a corresponding laboratory setting. If, for example, there are doubts concerning the use of aggressive (corrosive) media together with the CultiBag or in the respective laboratory setting, all components must undergo corrosion resistance testing. Warranty agreements for all Sartorius Stedim Biotech products are geared towards the General Terms and Conditions. Any deviations from this must be declared in writing.

- ▷ The warranty covers design, manufacturing or material errors as well as their resulting malfunctions.
 - Defective material will be repaired or replaced.
- ▷ The agreements do not apply to disposables and other parts that undergo normal wear (e.g. electrodes, o-rings, seals, filters).

The warranty becomes void:

- ▷ if there is any improper handling of the bioreactor and its device parts, if they are used outside of their specifications or the accompanying information is ignored
- ▷ if the bioreactor is set up under unsuitable environmental conditions, e.g. in the presence of corrosive media
- ▷ if third-party components are used that are not approved for use by Sartorius Stedim Biotech.

7.5.2 Service

Defective devices cannot be repaired in the workshop of the customer himself. Defective devices can be returned to Sartorius Stedim Biotech. Repairs are carried out according to the maintenance conditions of Sartorius Stedim Biotech. Please contact your Sartorius Stedim Biotech representative regarding device returns.

The sender must prove compliance with the safety regulations.

A corresponding "Declaration of Decontamination and Cleaning of the Device" (located in the appendix) must be filled out before return shipment (see also section ► "7.6 Decontamination Declaration" on page 57).

7.5.3 Device Disposal

CultiBags are single-use items that must be disposed of after each cultivation. Sterilization may be required before you can dispose of the CultiBag. This procedure must correspond to the respective national regulations for the disposal of biological harmful single-use devices.

1. Turn off the device and disconnect it from the power supply.
2. Disconnect the compressed air supply.
3. Disinfect or clean any contaminated parts if necessary.
4. Break down the BIOSTAT® RM 20 | 50 basic into its individual parts:
 - ▷ Stainless steel components
 - ▷ Plastic components
 - ▷ Electrical components
 - ▷ Electronic components
 - ▷ Battery
 - ▷ Please observe any applicable national regulations.

7.6 Decontamination Declaration

Defective devices and accessory parts can be sent to Sartorius Stedim Systems GmbH for warranty claims and service. This shipment must include a decontamination declaration noting which media, microorganisms and/or cells the device or components have been in contact with and what type of disinfection and/or decontamination it has undergone.

- ▷ The recipient (i. e. the Service Department of Sartorius Stedim Systems GmbH) must be able to view the decontamination declaration before unpacking the goods.
- ▷ A copy of the decontamination declaration for sending devices to Sartorius Stedim Systems GmbH is included in this appendix or available as a separate document upon request. If you required additional copies of this document, please contact Sartorius Stedim Systems GmbH.

Declaration about decontamination and cleaning of equipment and components

To protect our personnel, we require all equipment or components be free of biological, chemical, or radioisotopic contaminants. We will only accept such equipment or components when:

- the equipment or components have been adequately CLEANED and DECONTAMINATED.
- this declaring document has been completed, signed and returned by an authorized person.

Please help us in assuring a safe, hazard-free work environment.

Description of the Equipment or Component(s)

Description / Cat. No.:	
Serial no.:	
No. of invoice/delivery note:	
Date of delivery:	

Contamination / Cleaning

Attention: Please specify exactly the biological, chemical, or radioisotopic contaminant.	Attention: Please describe the cleaning and decontamination procedure/method.
The equipment was contaminated with:	and it has been cleaned and decontaminated by:

Legally binding declaration

I (we) certify that all information given in this form is correct and complete. The equipment and components have been adequately decontaminated and cleaned according to the legal requirements. No chemical or biological or radioisotopic risks remain that can endanger exposed persons' safety or health.

Company / Institute:			
Address / Country:			
Tel.:		Fax (with area code):	
Name of the authorized person:			
Position:			
Signature / Date:			

Please pack the equipment properly and send it to your local service representative or to Sartorius Stedim Systems GmbH, Germany (carriage paid to receiver).

Sartorius Stedim Systems GmbH Technical Service Schwarzenberger Weg 73-79 34212 Melsungen Germany

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The screenshot shows a menu titled 'Userlevel-menu' with a yellow header. Below the header, it displays 'User-level: Locked' and 'TT.MM. 20JJ/HH:mm Bag-config:'. The main area contains a 'Create new password for:' section with a button labeled 'User'. Below this are three input fields: 'Enter new password' (masked with asterisks), 'Re-Enter new password' (masked with asterisks), and 'Save new password'. At the bottom, there is a 'Back' button and a blue house icon.

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The screenshot shows a numeric keypad interface with a yellow header. The keypad consists of a grid of buttons labeled with numbers 1-9, letters A-Z, and special keys like ESC, CLR, BS, and ENT. Below the keypad, there is a row of asterisks (*****).

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The screenshot shows the 'System settings' menu. At the top, it displays 'System settings' and 'User-level: Locked'. On the right, it shows the date and time 'TT.MM. 20JJ/HH:mm' and 'Bag-config:'. The main content area has a grey background and contains the following elements:

- A question: 'Automatically controller-start after system power-restart?' with an 'OFF' button.
- A label: 'Time period to hold the sampling position:'
- A numeric input field with four asterisks and a decimal point, followed by the text 'seconds'.
- A 'Back' button at the bottom left.
- A home icon at the bottom right.

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The first screenshot shows the 'Settings-menu'. At the top, it displays 'Settings-menu' and 'User-level: Locked'. On the right, it shows the date and time 'TT.MM. 20JJ/HH:mm' and 'Bag-config:'. The main content area has a grey background and contains several buttons:

- 'Date and time'
- 'Service'
- 'Bag configuration'
- 'Technician level'
- 'Control mode'
- 'System'
- 'Batch-timer'

At the bottom, there is a 'Back' button and a home icon.

The second screenshot shows the 'Batch timer' configuration screen. At the top, it displays 'Batch timer' and 'User-level: Locked'. On the right, it shows the date and time 'TT.MM. 20JJ/HH:mm' and 'Bag-config:'. The main content area has a grey background and contains the following elements:

- Labels for time units: 'hh:', 'mm:', and 'ss:'.
- Three numeric input fields, each with two asterisks, separated by colons.
- A 'Start batch-timer' button.
- A 'Back' button at the bottom left.
- A home icon at the bottom right.

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Sartorius Stedim Biotech GmbH
August-Spindler-Str. 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or by any
means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make changes
to the technology, features,
specifications and design of the
equipment without notice.

Status:
November 2011,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany