🖲 FUJIFILM

FLUORESCENT IMAGE ANALYZER FLA-5000

Operation Guide

Science System, Industrial Materials Dept.

Fuji Photo Film Co., Ltd.

Version 1.0 October 2001

Thank you for purchasing our fluorescent image analyzer Model FLA-5000.

This Operator's Guide is intended to provide a brief description of the operating procedures for the image analyzer unit incorporated into your system together with the functional options available from the unit. For a detailed description of the operation of the unit, please refer to the Operating Manual provided.

Index

1	Features of the FLA-5000	3
2	Turning ON/OFF the FLA-5000	3
3	Startup/Shutdown the FLA-5000 Image Reader	4
4	Exposing IPs and Precautions	5
5	Reading IP	6
6	Reading Fluorescent Samples	9
7	Reading Chemiluminescent Samples	16
8	Reading digitize Samples	19

Features of the FLA-5000

* The Fluoro Image Analyzer FLA-5000 employs solid-state lasers of wavelengths suitable for the excitation of fluorescent pigments originally developed by FUJI PHOTO FILM Co., Ltd. It combines several lasers of different wavelengths with four types of filters and may read gel, membranes, etc. dyed with various fluorescent pigments. In addition, it may read imaging plates (IPs) originally developed as radiation energy sensors by FUJI PHOTO FILM Co., Ltd.

<Main features>

- * Makes high-sensitivity, high-resolution images of gel, membranes, etc. dyed with fluorescent pigments.
- * Produces ultra high-resolution images equivalent in the quality to X-ray film;
- * Accommodates the IP featuring ultra-high sensitivity, wide dynamic range, reliable linearity, high resolution and capability of reuse;
- * Ensures high resolution, high sharpness and linearity not attainable by existing filmless systems; and
- * Obviates the need for a darkroom and automated processing lab for configuration of the system.

2 Turning ON/OFF the FLA-5000

- 1 Turning on procedures
 - 1. Press the POWER switch on the right side of the unit to the "I" position.



2. The POWER indicator lamp should light up as shown below to indicate that the unit is correctly turned ON.



Note: If the alarm buzzer sounds during the power-ON sequence (with the ERROR lamp blinking), check if the door of the stage setting block is securely closed.

2 Turning off procedures

Check that scanning is not in progress in the unit, then press the POWER switch on the right side of the unit to the "0" position.

3 Startup/Shutdown the FLA-5000 Image Reader

- 1 Startup procedures
 - 1. Turn on the FLA-5000 and peripheral devices.
 - 2. Turn on the computer (DOS/V PC or Macintosh).
 - Make sure that the FLA-5000 has warmed up. (Only the power lamp on the upper left panel on the front of the FLA-5000 is lit when warming-up is completed.) Start the FLA-5000 Image Reader from the startup menu or using the shortcut key. (On the Macintosh, double-click the alias to start the software.)





4. The main window of the FLA-5000 Image Reader is displayed.



2 Shutdown procedures

Click on Exit from the File menu. (On the Macintosh, click on Quit from the File menu.)

4 Exposing IPs and Precautions

<Procedures for exposure>

- 1. Immediately before the start of exposure, complete a thorough erasing process.
- 2. Wrap the sample in Saran Wrap® or other appropriate wrapping film.
- Put the sample into the cassette with the surface to be scanned facing up.
 For proper position alignment between the sample and the IP, use the grids (inscribed at inter vals of 25 mm) as a guide.



4. Referring to the diagram below, load the IP into the cassette making sure that the surface of the IP for exposure faces the sample.

Use the notch on the IP as a guide to identify the correct orientation of the sample, the IP and the loading of the sample-containing IP into the carrier.



5. Close the cassette lid fully until a click is heard.

<Precautions for handling the IP>

- 1. Keep the IP away from any moisture
 - Although the IP is designed to be water resistant, moisture makes it susceptible to a lowering of sensitivity, causing incomplete linearity. Allow the sample to dry out completely before expo sure.

Tightly wrap any sample that is not dried completely in Saran Wrap to ensure proper moisture tightness.

2. Keep the IP away from any volatile solvents

Overnight or longer exposure of any sample containing dichloromethane, chloroform, acetone, acetic acid, or other solvent can cause shrinkage of the protective layer on the surface of the IP, leading to deformation of the plate. A deformed IP can cause jamming during the scanning. Wrap such samples with a double layer of Saran Wrap to ensure proper air tightness before exposure.

3. Wear cotton gloves

When handling the IP, wear cotton gloves to protect the surface of the IP against any contami nation. When removing the IP from the cassette, use a suction disc. Using the bare hands (fingernails) can cause the edges to gradually peel off, eventually resulting in an unusable IP (due to the failure of scanning caused by jamming in the unit).

- 4. Use a dry cotton cloth for cleaning the IP Clean the surface of the IP with a dry cotton cloth or other appropriate material. Never use a water-dampened cloth. For any heavy contamination, use a cloth with a small amount of anhydrous ethanol (Extra Pure or Guaranteed Reagent). In actuality, however, it is necessary to ensure that such appropriate anhydrous ethanol to be used has been kept in a brown re agent bottle or stored in the proper environment recommended by the maker, since any anhy drous ethanol stored under incorrect storage conditions can cause deterioration of the IP.
- Keep the IP away from any light source after completion of exposure After exposure, keep the IP strictly hidden from any direct light source. If the IP is exposed to light, the image information in the IP may be lost.



Click on this button to return to the Main window.





For the setting of any additional conditions other than those described below, refer to the instructions for Setting Conditions in the IP S Mode (page 6).

(1) CH1: 800 V

This voltage is applied to the first PMT.

Input a voltage value directly in this box. (You may input a value between 250 and 1000.) The greater the value is, the higher the voltage to be applied to the PMT and the reading sensitivity are.



- 6. Place the IP on the IP stage.
 - * You may use only IPs with magnetic absorption layers on the back on the IP stage. Since an IP is stuck on the IP stage by the magnetic force, any IPs without magnetic absorption layers may not be attached to the IP stage. Such IPs are unusable.
 - 6-1 Put the IP cassette with an exposed IP on the side of the IP stage.



6-2 Turn down the room lighting (below 20 luxes).(Keep the room lighting turned down until the IP stage is set and the main cover is closed.)

6-3 Take the IP out of the IP cassette and set it on the IP stage immediately.



Set the IP on the back side (with a numeric value and alphabetic characters) of the IP stage so that the IP reading surface (white or blue side with fluorescent material applied) faces up.

- 7. Load the IP stage into the FLA-5000 body.
 - 7-1 Open the door of the stage setting block, and set the IP stage with the stuck IP facedown.



- 7-2 Push the IP stage to the very end (until it butts).
- 7-3 Close the door of the stage setting block.
- 8. Starting reading
 - 8-1 Click the Read... button, and reading starts.



8-2 When reading finishes normally, the following dialog box appears. Click the total button.



8-3 You may finish reading at any time before the whole reading area is read. Click the button when you want to finish reading.

6 Reading Fluorescent Samples

Set the reading 1. Turn on the FLA-5000 and peripheral devices. conditions. 2. Turn on the computer (DOS/V PC or Macintosh). 3. Make sure that the FLA-5000 has warmed up. (Only the power lamp on the upper left panel on the front of the FLA-5000 is lit when warming-up is Set a Fluorescent completed.) Start the FLA-5000 Image Reader from the startup menu Samples on the or using the shortcut key. (On the Macintosh, double-click the alias to start the software.) stage. 4. The main window of the FLA-5000 Image Reader is displayed. 5. Select a reading mode and set the reading conditions. Set the stage on the FLA-5000. The number of modes available for reading fluorescent samples 5-1 Start reading. depends on the number of photo-multiplier tubes (PMTs) built in the FLA-5000 as shown below. One PMT is built in the FLA-5000 Two PMTs are built in the FLA-5000 \longrightarrow Six reading modes are available. 1 laser 1 image ... One laser is used to excite fluorescent samples. Fluorescent light from a sample Available passes through a single filter when one PMT and forms a single image. is built in. 1 laser 1 image Cyclic ... 1-laser, 1-image reading is repeated (up to four times). 1 laser 2 image ... One laser is used to excite Available when two fluorescent samples. PMTs are built in. Fluorescent light from a sample is divided and filtrated through two filters to form two images. 1 laser 2 image Cyclic ... 1-laser, 2-image reading is repeated (up to two times). 2 laser 2 image ... Two lasers are used to excite fluorescent samples simultaneously. Fluorescent light from a sample is divided and filtrated through two filters to form two images. 2 laser 2 image Cyclic ...2-laser, 2-image reading is repeated (up to two times).



11 aser 2i mage

11 aser 2image Cyclic







2) aser 2) mage Cyclic

- (1) File Name Assign a file name to the read image data for saving.
- (2) Comments Enter any comments required for saving with the assigned file name.
- (3) Change Filter... Use this button when changing any filter in formation. (No need to use this option in the normal IP reading operation.)
- (4) Laser Select the laser to be used from the pulldown menu.
- (5) CHT: 400 V Enter the voltage value to be applied to the first PMT. (Voltage values must be an integer between 250 and 1000.)
- (6) FilterSelect the filter to be used for reading.(7) Root Folder Change
- Specify where to save the file for saving the image data.
- (8) Gradations Click to select the number of gradations of a read image.
- (9) Resolution Click to select the pixel size for reading.
- (10) Sample Area Click to select the method for setting the reading area.

(11) Reading Area

The actual reading area will be shown in the form of a red-colored box. As necessary, move the cursor onto the red frame or to any point in the red box and drag the box to change the size and/ or position of the reading area.

- (12) Read Click on this button to start the reading operation.
- (13) Save Template
 Use this button to save the reading conditions in a file. (For details, see the Operation Manual.)
- (14) Return Click on this button to return to the Main window.
- (15) Read Cycle (a)
 Specify how many times to repeat the reading cycle. It is possible to repeat reading up to four (4) times in the "1 laser 1 image Cyclic" mode.
- (16) CH1: 400 V CH2: 400 V

Woltage to be applied to the first PMT

CH2: 400 ▼ →Voltage to be applied to the second PMT

Enter the voltage values. (The voltage values must be an integer between 250 and 1000.)

- (17) Read Cycle (b) Specify how many times to repeat the reading cycle. It is possible to repeat reading up to two (2) times in the "1 laser 2 image Cyclic" or "2 laser 2 image Cyclic" mode.
- (18) Laser 1 : 2 473nm -Laser 2 : 2 532nm -

Select the laser to be used from the pulldown menu.

- 6. Set a fluorescent sample on the fluorescent sample stage or multi-stage. Setting a gel sample:
 - 6-1 Set a fluorescent gel sample on fluorescent sample stage. When reading a fluorescent gel sample, set the fluorescent sample stage on the FLA-5000 first and set the fluorescent gel sample on the fluorescent sample stage then as shown below. This is an easy way.
 - 6.1.1 Open the door of the stage setting block, and set the fluorescent sample stage with the printed frame side faceup on the FLA-5000.





Set the stage in the front position. Do not pull it to the very end.

6.1.2 Set the gel stopper as shown below.



6.1.3 Put a gel sample so that the side with charging holes faces up. Fix it with the gel stopper.



6.1.4 Press the fluorescent sample stage slowly to the very end (until it butts).



6.1.5 Close the door of the stage setting block.

Setting the tighter plates:

- 6.2 Set a tighter plate on the multi-stage.
 - 6.2.1 Set a tighter plate in any position of the tighter plate frame set on the multi-stage.



6.2.2 Open the door of the stage setting block, and set the multi-stage with the printed frame side faceup on the FLA-5000. Push the multi-stage to the very end.



6.2.3 Close the door of the stage setting block.

Setting gel with glass:

- 6.3 Set a gel sample with glass on the multi-stage.
 - 6.3.1 Detach the tighter plate plug-in as shown below, if it has been set.



6.3.2 Insert one or two gel samples with glass from the front side. Adjust the movable stays to the glass edges (by moving them from the back to the front).



6.3.3 Turn the levers on both sides horizontally to lock them.

6.3.4 Fix the glass with two glass holders in the front position and one glass holder in the rear position. (Turn the dial of each glass holder clockwise to fix the glass.)



There are spacers on the side of the glass holders. Each spacer has two levels. Slide the spacers to adjust the balance of the glass holders.

6.3.5 Open the door of the stage setting block, and set the multi-stage with the printed frame side faceup on the FLA-5000. Push the multi-stage to the very end (until it butts).



6.3.6 Close the door of the stage setting block.

- 7. Start reading.
 - 7,1 Click the Read... button, and reading starts. The condition of the preset reading area, where reading is completed, is displayed in the reading status real-time display window. Reading is executed from the front side of the stage to the back. On the monitor screen, read data is displayed from below to above.
- a) 1-laser, 1-image



7.2 When reading finishes normally, the following dialog box appears. Click the button.



7.3 You may finish reading at any time before the whole reading area is read.Click the ______ button when you want to finish reading.

Reading Chemiluminescent Samples 1. Turn on the FLA-5000 and peripheral devices. Set the reading 2. Turn on the computer (DOS/V PC or Macintosh). conditions. 3. Make sure that the FLA-5000 has warmed up. (Only the power lamp on the upper left panel on the front of the FLA-5000 is lit when warming-up is completed.) Start the FLA-5000 Image Reader from the startup menu Set the stage on the FLA-5000. or using the shortcut key. (On the Macintosh, double-click the alias to start the software.) 4. The main window of the FLA-5000 Image Reader is displayed. Set a Chemilumines- 5. Setting the reading conditions cent Sample on the stage. 5-1 Click the button. Start reading. The following window opens. 5-2 . 🗆 🗙 FLA-5000 Chemiluminescence Mode (6) (1)(7) (8) (9) (3)

Set the reading conditions by following the instructions shown below.

(1) File Name Assign a file name to the read image data for saving.

(4)

(5)

- (2) Comments Enter any comments required for saving with the assigned file name.
- (3) Change Filter... Use this button when changing any filter information. (No need to use this option in the normal chemiluminescent sample reading operation.)
- (4) OHL FOR Y
 Enter the voltage value to be applied to the first PMT. (Voltage values must be an integer between 250 and 1000.)
- (5) Filter Select the filter to be used for reading. In the chemiluminescent sample reading mode, no filter is used as a general rule. In the chemiluminscence mode, the filter tray must contain at least one throughposition where no filter module is installed. (For details, refer to the Operation Manual.)
- (6) Root Folder Change Specify where to save the file for saving the image data.

(7) Gradations Click to select the number of gradations of a read image.

Read.

Save As Te

(10)

(11)

(13)

(12)

(8) Resolution Click to select the pixel size for reading.
(9) Sample Area

Click to select the method for setting the reading area.

- (10) Reading Area The actual reading area will be shown in the form of a red-colored box. As necessary, move the cursor onto the red frame or to any point in the red box and drag the box to change the size and/ or position of the reading area.
 (11) Read
 - Click on this button to start the reading operation.
- (12) Save Template Use this button to save the reading conditions in a file. (For details, see the Operation Manual.)
- (13) Return Click on this button to return to the Main window.

- 6. Set the fluorescent sample stage on the FLA-5000.
 - 6.1 Open the door of the stage setting block, and set the fluorescent sample stage with the printed frame side faceup on the FLA-5000.



Set the stage in the front position. Do not pull it to the very end.



- 7. Set the chemiluminescent sample on the fluorescent sample stage. The following descriptions show an example of setting a membrane sample.
 - 7.1 Set a membrane sample with the chemiluminescent side facedown on the fluorescent sample stage.



Tape the corners of the membrane to reduce influences of physical vibrations upon the image.

7.2 Press the fluorescent sample stage slowly to the very end (until it butts).



7-3 Close the door of the stage setting block.

8. Starting reading

- 8.1 Click the Read... button, and reading starts.
- 8.2 When reading finishes normally, the following dialog box appears. Click the button.

mage Reader FLA-5000 🛛 🔀	
Reading completed successfully.	
<u> </u>	

8.3 You may finish reading at any time before the whole reading area is read. Click the button when you want to finish reading.

8 Reading Digitize Samples



- Click to select the filter to be used for reading. (In the digitizing mode, two options are available, Y510 and O575.)
 Root Folder Change
- (7) Root Folder Change Specify where to save the file for saving the image data.
- Click on this button to return to the Main window.

Operation Manual.)

Return

(14)

6. Set the fluorescent sample stage on the FLA-5000.

6.1

- sample stage with the printed frame side faceup on the FLA-5000.

Open the door of the stage setting block, and set the fluorescent

Set the stage in the front position. Do not pull it to the very end.

- 7. Set the fluorescent sample stage on the FLA-5000.
 - 7.1 Set a sample with the migration side facedown on the fluorescent sample stage.



- 7.2 Put the fluorescent plate for digitizing supplied with the fluorescent sample stage on the digitize sample.
 - * Put the fluorescent plate with the mat (lusterless) side facedown.
- 7.3 Press the fluorescent sample stage slowly to the very end (until it butts).



7-4 Close the door of the stage setting block.

- 8. Starting reading
 - 8.1 Click the Read... button, and reading starts.



8.2 When reading finishes normally, the following dialog box appears. Click the total button.

Inner Deeden ELA 5000				
٩	Reading completed successfully.			
	<u>OK</u>			

8.3 You may finish reading at any time before the whole reading area is read. Click the button when you want to finish reading.

Science System, Industrial Materials Dept.

Fuji Photo Film Co., Ltd.

Version 1.0 October 2001