# BiliChek

USER MANUAL | MANUEL DE L'UTILISATEUR | MANUALE DELL'UTENTE |  
BENUTZERHANDBUCH | MANUAL DEL USUARIO | MANUAL DO UTILIZADOR

Manufactured for  
Children's Medical Ventures, LLC  
1001 Murry Ridge Lane  
Murrysville, PA 15668 USA

**EC REP**  
Respironics Deutschland  
Gewerbestrasse 17  
82211 Herrsching, Germany

REF 1062275

\*+H90610622750A\*

1005175  
HM 9/7/09

**BiliChek**  
user manual

**English**

BiliChek, BiliCal and BilEclipse are trademarks of Respiromics, Inc. and its affiliates.

© 2009 Respiromics, Inc. and its affiliates.

# Table of Contents

1. Introduction .....	1
About This Manual .....	1
BiliChek Description.....	1
Intended Use.....	1
Jaundice .....	2
Theory of Operation .....	2
Protecting a Measurement Site.....	4
Components of the BiliChek system.....	5
2. Warnings, Cautions, Contraindications, and Symbols.....	6
Warnings.....	6
Cautions.....	7
Contraindications .....	7
Notes .....	8
Symbols .....	8
3. Set Up.....	10
Initial Set Up.....	10
Attaching the Hand-Held Unit Battery.....	10
Connecting the Charger Base Power Cord .....	11
Connecting the Charger Base Ethernet Cable .....	12
Charging the Hand-Held Unit Battery .....	12
4. Viewing and Changing Settings .....	13
BiliChek Hand-Held Unit (HHU) .....	13
BiliChek Keypad .....	13
Turning the Unit On.....	14
Navigating the Menu Screens .....	14
Accessing On-Screen Help .....	15
Returning to the Home Screen .....	15

Changing the Settings.....	15
Customizing the Settings.....	16
Selecting a Language.....	16
Selecting a Time and Date Format.....	17
Entering the Current Time and Date .....	19
Additional Settings .....	20
Calibrating and Taking a Reading .....	22
Using Other Features .....	24
History .....	24
Stand-by .....	25
Recommendations for Use During and After Phototherapy .....	25
During Phototherapy .....	25
BilEclipse .....	26
After Phototherapy .....	27
5. Cleaning, Maintenance, and Service .....	28
Cleaning .....	28
Maintenance .....	29
Changing the Measurement Light Bulb .....	28
Changing the Battery .....	31
Service .....	31
Customer Service Information.....	31
6. Specifications.....	32
Environmental.....	32
Physical .....	32
BiliChek System.....	32
BiliChek Hand-Held Unit.....	32
Charger Base.....	32
Light Bulb.....	32
Standards Compliance .....	33
Electrical Requirements.....	33
WEEE/RoHS Recycling Directives .....	33
Disposal .....	33
Appendix A. EMC Information .....	34
Appendix B.Warranty .....	38

# BiliChek

## user manual

# 1. Introduction

## About This Manual

This User manual has been prepared for both first-time and experienced users of the BiliChek System. This manual provides basic instructions for use and troubleshooting procedures. To ensure safe operation of this product the user should read the entire instruction manual before operating the unit. Save this manual for future reference.

English

## BiliChek Description

The BiliChek provides a numerical measurement of predicted bilirubin count in mg/dL or  $\mu\text{mol}/\text{L}$  within a clinically beneficial range that has been correlated with total serum bilirubin concentration measured by High Pressure Liquid Chromatography (HPLC). The device is intended for use in the hospital or institutional environment, including clinics, physician's offices, and home nursing agencies, to assist clinicians in monitoring the status of neonates for the development of hyperbilirubinemia.

## Intended Use

The BiliChek Non-Invasive Bilirubin Analyzer is a non-invasive transcutaneous bilirubinometer, which is intended as an index to predict serum bilirubin levels prior to phototherapy (pre-phototherapy), during, and post phototherapy in neonates, without regard to gender, gestational age, or body weight. Neonates whose BiliChek test results are indicative of hyperbilirubinemia are evaluated by their physician(s) for appropriate patient management.

## Jaundice

Jaundice is a condition characterized by yellowing of the skin and eyes that is caused by an excess of bilirubin in the blood. Bilirubin is a normal waste product resulting from the breakdown of red blood cells and is removed from the body by the liver. Prior to birth, bilirubin in an infant is processed by the mother's liver and excreted. After birth, an infant must eliminate bilirubin without the mother's help. It may take the infant's system several days to begin eliminating bilirubin from the blood faster than it is produced. Infants who are born prematurely, who are underfed, or who belong to certain ethnic groups are at increased risk of developing jaundice. The initial screening of jaundice is the observation of yellow skin. This is a subjective determination prone to errors due to differing skin colors and gestational ages. Bilirubin levels normally peak in full-term infants 4 to 5 days after birth and may peak later in premature infants. If the infant's serum bilirubin levels continue to rise the infant is at risk for kernicterus (the deposit of bilirubin in the brain at toxic levels which can cause permanent neurologic impairment).

## Theory of Operation


The BiliChek works by directing white light into the skin of the newborn and measuring the intensity of the specific wavelengths that are returned. By knowing the spectral properties of the components within the skin, one can subtract out the interfering components and determine the concentration of bilirubin. This is called a transcutaneous bilirubin measurement.

Photons of specific wavelengths are preferentially absorbed by certain molecules. By plotting the absorption against the wavelength, one can visualize characteristic absorption spectra of the particular molecules. For example, melanin has a near-linear absorption spectrum in the visible spectrum and, like the scattering phenomenon, there is greater absorption of photons with shorter wavelengths than in the red region of the spectrum. Conversely, hemoglobin is a much more complicated absorber which is compounded by the fact that oxyhemoglobin and deoxyhemoglobin have different profiles. The peak absorption of photons by bilirubin occurs at a wavelength of 460 nm. This is in the blue portion of the spectrum and is the reason why blue lights are sometimes preferred for phototherapy. It is also in a region of the spectrum at which hemoglobin absorption is relatively low.

The light reflected from the skin of neonates and collected by the BiliChek is analyzed with a highly sophisticated, proprietary algorithm to generate a serum bilirubin measurement. The major skin components (see figure) which impact the spectral reflectance in newborns are (1) dermal maturity, (2) melanin, (3) hemoglobin, and (4) bilirubin. The intensity of the reflected light is converted to absorbance units, Optical Density (OD), for analysis. The total OD at each measured wavelength is equal to the sum of the component ODs. This is represented mathematically as:

$$\text{OD}_{\text{total}} = \text{OD}_{\text{derm}} + \text{OD}_{\text{mel}} + \text{OD}_{\text{heme}} + \text{OD}_{\text{bili}}$$

By utilizing known spectral characteristics of each component, the contribution to the total OD can be sequentially subtracted until bilirubin is all that is remaining.

### Non-Invasive Skin Spectral Components<sup>1</sup>

<sup>1</sup>Jacques, S., Saidi, Il, Ladner, Al, and D. Oelberg. *Developing an optical fiber reflectance spectrometer to monitor bilirubinemia in neonates.*

Each of these components have also been considered in the development of the algorithm by which the BiliChek computes the bilirubin measurement. By measuring the intensity of the reflected light at more than 100 individual wavelength bands across the visible spectrum, a spectral analysis can be performed. By normalizing against a reference standard, the BiliCal, the variations over time of the bulb and the photo-detectors are eliminated. The peak absorption of the bilirubin at approximately 460 nm is measured by subtracting the absorption of the other components from the total skin absorption. This is represented by the equation below.


$$\text{OD}_{\text{bili}} = \text{OD}_{\text{total}} - \text{OD}_{\text{mel}} - \text{OD}_{\text{heme}} - \text{OD}_{\text{derm}}$$

The absorption of light by a molecule is proportional to its concentration. Therefore, the absorption of bilirubin, as calculated above, is proportional to the concentration of bilirubin in the subcutaneous capillary beds and subcutaneous tissue. By subtracting the contribution of the aforementioned components, the bilirubin absorbance can be quantified and is correlated to the laboratory total serum bilirubin levels.

<sup>1</sup>SPIE Proceedings 2975, Laser-Tissue Interactions, San Jose, CA Feb. 1997. 1-7.


## Protecting a Measurement Site

Exposure of a neonate to phototherapy creates changes in the skin due to the therapeutic effects of phototherapy lights. If a small portion of the newborn skin is protected from the therapeutic effects of phototherapy lights, the relationship between the skin bilirubin in the protected skin site and blood bilirubin is likely to be maintained. The accuracy of the BiliChek transcutaneous bilirubin measurement taken from a protected site during and after phototherapy has been determined in a controlled clinical study. The recommendations in the User Instructions are based on the results of this study.

## Components of the BiliChek System

The BiliChek system includes the following components. Some components and accessories will be packaged in separate boxes.

To reorder parts and accessories, visit the Respiromics web site at [www.respiromics.com](http://www.respiromics.com) or contact Respiromics Customer Service at 1-724-387-4000 or +49 8152 93060.

1. BiliChek Hand-Held Unit
2. Battery Door Cover
3. Battery
4. Protective Tip
5. Power Supply
6. BiEclipse Phototherapy Protective Patches (Bag of 10)
7. Charger Base
8. User Manual (not shown)
9. Quick Start Guide (not shown)

**English**

 **WARNING:**

The BiliChek Hand-Held Unit is to be used only with these and approved Respiromics replacement parts and accessories.

## 2. Warnings, Cautions, Contraindications, and Symbols

*Caution! US federal law restricts this device to sale by or on the order of a physician.*

### Warnings

- The BiliChek Hand-Held Unit is to be operated only in conjunction with Respiration approved charger base, battery packs, power supply, and BiliCal Individual Calibration Tips.
- Do not point at eyes while taking a measurement.
- Do not use on areas of the skin with visible bruising.
- If you suspect BiliChek is interfering with other electronic devices or vice versa, try the following measures to correct it:
  - Relocate the device until the interference stops.
  - Move the BiliChek farther away from the device.
  - Plug the BiliChek Power Supply into a different power outlet so that the BiliChek Power Supply and the device are on different breaker or fuse circuits.
- Do not use in the presence of flammable substances such as anesthetics, cleaning agents, or gases that support combustion (e.g. oxygen, nitrous oxide).
- Do not puncture or incinerate the battery pack or dispose of in a landfill. Batteries should be disposed of properly according to local regulations.
- Do not use outdoors or place in direct sunlight.

#### **WARNING**

A warning indicates the possibility for injury to the user or operator.

- The charger base unit (battery charger) should be plugged into a properly grounded outlet. Do not operate with damaged plug or cord.

## Cautions

- The BiliCal Individual Calibration Tip should not be exposed to light for extended periods of time. Keep the BiliCal in its package prior to use.
- If the product is dropped and frequent error messages occur, it should be returned to Respirationics.
- Do not immerse in water or other liquid.
- Do not let the power cord hang over the edge of counter or table as this may pose a tripping hazard.
- Pins of connectors identified with the ESD warning symbol should not be touched. Connections should not be made to these connectors unless ESD precautionary procedures are used. Precautionary procedures include methods to prevent build-up of electrostatic discharge (e.g., air conditioning, humidifications, conductive floor coverings, and non-synthetic clothing), discharging one's body to the frames of the equipment or system or to earth or a large metal object, and bonding oneself by means of a wrist strap to the equipment or system or to earth.

### CAUTION

*A caution indicates the possibility for damage to the device.*

English

## Contraindications

- The BiliChek should not be used in patients who have undergone an exchange transfusion.
- Avoid using measurement sites with the following conditions since they can produce erroneous results:
  - Bruising
  - Birthmarks
  - Hematomas
  - Excessive Hairiness

## Notes

**Note:** A note places emphasis on an operating characteristic.

Certain medical conditions and medications are known to disrupt bilirubin-albumin carrying capacity and binding. These factors may alter the transcutaneous bilirubin measurement by increasing or decreasing the bilirubin deposition in the skin especially at high total serum bilirubin concentrations (>15 mg/dL or 255 µmol/L). Therefore, additional caution should be used in interpreting the BiliChek measurement in patients with serious illnesses or complications such as septicemia, Patent Ductus Arteriosus (PDA), Intraventricular Hemorrhage (IVH), and Respiratory Distress Syndrome (RDS). In these patients, a separate serum bilirubin measurement should be considered if the BiliChek indicates elevated bilirubin levels.

There is limited data on use of the BiliChek on patients under intensive phototherapy (>30 µW/cm<sup>2</sup>/nm or 1.5 mW/cm<sup>2</sup>), double-sided phototherapy, or after exchange transfusions. Caution should be used in interpreting BiliChek results under these conditions.

## Symbols

The following symbols appear on the BiliChek unit or accessories.

Symbol	Description	Symbol	Description
	Attention, consult instructions for use		Indoor Use
	Use With		Type BF Applied Part
	Do Not Reuse		Enter Key

Symbol	Description	Symbol	Description
	Use By		Compliant with the Waste Electrical and Electronics Equipment/Restriction of the Use of Certain Hazardous Substances in Electrical and Electronics Equipment (WEEE/RoHS) recycling directives
	Catalog Number		Canadian/US Safety Certification. Conforms to UL Std. 60601-1. Certified to CAN/CSA Std. C22.2 No. 601.1
	Batch Code		European Representative
	Serial Number		Temperature Range
	European Declaration of Conformity		Relative Humidity Range
	AC Power		DC Power
	Help Key		Navigation Arrow Keys


# 3. Set Up

## Initial Set Up

Visually inspect all components for physical damage. If damage is noted, please contact your local dealer/distributor or Resironics.


### Attaching the Hand-Held Unit Battery

1. Clip the battery connector into the connector on the hand-held unit.

2. Fold the battery wires underneath the light bulb compartment. Slide the battery inside the hand-held unit.

Light Bulb Compartment

### **WARNING**

*The battery is designed to work only with the BiliChek. Do not attempt to use the battery pack with other products.*


### Connecting the Battery

### Inserting the Battery


3. Attach the battery compartment cover by sliding the cover up until it clicks securely in place.

#### Attaching the Cover

### Connecting the Charger Base Power Cord

1. Attach the appropriate plug style to the power supply, if necessary.

#### WARNING


Periodically inspect the power cord for damage or signs of wear. Discontinue use and replace if damaged.

English

#### Attaching the Plug

2. Plug the socket end of the power cord into the power inlet on the bottom of the charger base.


#### Connecting the Power Supply to the Charger Base

3. Plug the pronged end of the cord into a wall outlet. Ensure all connections are secure.

## Connecting the Charger Base Ethernet Cable

The BiliChek unit can send data to your laboratory information system if you have the necessary software. To allow this capability, you must connect the ethernet cable included with your BiliChek system.

1. Connect one end of the Ethernet cable into the connector on the bottom of the charger base.


**Connecting the Ethernet Cable**

2. Connect the other end of the cable to your network connection.

## Charging the Hand-Held Unit Battery

Before initial use, the battery should be fully charged.

1. Place the BiliChek onto the charger base. A lightning bolt icon will appear on the battery icon to indicate the battery is charging.
2. Charge for eight hours or until the battery status icon on the screen indicates a full charge.


# 4. Viewing and Changing Settings

## BiliChek Hand-Held Unit (HHU)

The BiliChek Hand-Held Unit consists of the following components:


- Display window
- Keypad
- Calibration Tip
- Barcode Scanner
- Battery compartment
- Charging Pins

BiliChek Hand-Held Unit - Components

## BiliChek Keypad


Use the keypad to navigate through the on-screen menus.

Keypad

## Turning the Unit On

1. Ensure the battery is connected to the unit and fully charged.
2. If the unit is in sleep mode, press the ✓ (Enter) key to begin.
3. The unit will power up and then will display the home screen.

**Note:** The BiliChek Hand-Held Unit does not have an off switch. The BiliChek will automatically turn off if it is idle for 2 minutes.

**Note:** Dashes will display for the current date and time until they are set.


Home Screen

## Navigating the Menu Screens

To navigate through all of the menu screens and settings:

- Use the Up and Down arrow keys to scroll through the menu options and settings.
- Use the Left and Right arrow keys to perform the actions specified on the screens.
- Use the ✓ (Enter) key to save the setting on the screen.


## Accessing On-Screen Help

### Help

From any screen, press the **?** (Help) key on the BiliChek keypad. A help message for that action will appear.

### Returning to the Home Screen

You can return to the home screen at any point by pressing the **House** (Home) key. If you return to the home screen in the middle of a reading, you must start the measurement process over.

**BiliChek Keypad**


**Note:** To exit the Setup mode at any time, press the **House** (Home) key on the bottom left of the keypad. Any unsaved data will be lost.

## Changing the Settings

Before you use the BiliChek for the first time, you may review or change certain factory settings. You can also change any of the settings at a later time. Choose Setup on the home screen to change settings. At any time during the setup process, you can press the **?** (Help) key on the bottom right of the keypad to access on-screen instructions.

### IMPORTANT!

Your institution will determine which information you are required to enter.

**Setup Selection**


## Customizing the Settings

Before changing any settings, it is recommended that you first change the BiliChek display to your native language. Next, you should select the appropriate Time Format and Date Format, and then enter the current Time and Date. The sections below will walk you through each of these settings.


### Selecting a Language

Before you begin to customize your BiliChek settings, you must change the display screen to your native language. The BiliChek display screen is defaulted to the English language. To change the language, you will need to do the following:

1. From the Home screen, use the down arrow key to select **Setup**.

2. Once in the Setup screen, use the down arrow key to select **Language**.

3. Scroll down until your language is highlighted. Press the ✓ (Enter) key to save the setting.

4. A screen will display to verify the language you have selected.

## Selecting a Time and Date Format

### Time Format

1. From the Setup menu, select **Time Format**.

**Note:** After changes are made to the sound level, date, or time selections, you will automatically be exited back to the Home screen. After changes are made to any other selection the unit will exit back to the main Setup screen.

2. Select 12 or 24 hour clock (12 hour is the default).

## Date Format

1. From the Setup menu, select **Date Format**.

2. Select either Month/Day/Year (mm/dd/yyyy) or Day/Month/Year (dd/mm/yyyy) (mm/dd/yyyy is the default).

## Entering the Current Time and Date

### Current Time

- Select **Current Time** from the Setup menu.

Current Time Selection

- Using the keypad, type in the Time.

Entering the Time

- Press the up arrow to place the cursor on the on-screen display.  
If you are using the 12-hour clock, go to step 4. Otherwise, go to step 5.
- Using the arrow keys, scroll to highlight **am** or **pm** and then press the ✓ (Enter) key to select am or pm.

Selecting am/pm

- Press the down arrow to highlight **Done** and then press the ✓ (Enter) key. The display will exit back to the Home screen and the time will now be displayed.

## Current Date

- Select **Current Date** from the Setup menu.

- Using the keypad, type in the Date.

- Press the ✓ (Enter) key to select **Done**. The display will exit back to the Home screen and the Date will now be displayed.

## Additional Settings

The following features can be customized from the Setup screen. To select a setting, highlight the selection and press the ✓ (Enter) key.

- Sound level:** Select 0 for no sound, select 1, 2, or 3 for sound with 3 being the loudest.
- Units of measure:** Select mg/dL or µmol/L.

**Note:** If you are going to connect to your laboratory information system, you must configure each of the following settings.

- **Nurse ID:** Select a method of entry. Enter the ID by scanning a barcode or typing information using the BiliChek keypad.
- **1st and 2nd Patient ID:** Select how you would like to enter the patient ID. You can enter the ID by scanning a barcode or typing the information using the on-screen keyboard. In addition, you can select None if you do not want to enter a patient ID.
- **1st and 2nd ID Type:** Select which type of patient ID you will be entering. The options are: medical record number, patient name, or date of birth. If you entered None for the 1st and or 2nd patient ID setting above, you will not need to enter the ID type.
- **Network alert:** If the unit loses network connectivity or is unable to connect to the network, this is the amount of time during which the unit tries to reconnect to the network before generating an on-screen alert (none, 1 minute, 5 minutes, or 30 minutes).

In addition, when the BiliChek unit is on the charger base, the following options are visible in the Setup menu and required for connectivity. Contact your Information Technology department to receive the proper settings for connectivity.


- IP address
- Subnet mask
- Gateway
- Host IP
- Host Port

**Important!** Network settings are stored in the charger base. If you use a different BiliChek unit with the charger base you will not need to configure the network settings again.

**Note:** To set the unit so it works similar to the original BiliChek unit, set the Nurse ID and 1st and 2nd Patient ID to None.

## Calibrating and Taking a Reading

- From the Home screen, select **TcB Reading**.

**Note:** If you set Nurse ID and Patient ID to None in the Setup menu, Steps 2 and 3 are not required.

- You can enter a Nurse ID and two of the following forms of Patient ID:
    - Patient Name
    - Patient Medical Record Number
    - Patient Date of Birth
  - To enter information, scan the appropriate barcode or type the information using the on-screen keyboard. To scan a barcode:
    - Highlight **Scan In** to activate the barcode scanner. Point the scanner on the hand-held unit towards the barcode and press the **✓ (Enter)** key.
    - The scanned information will appear on the screen. If it is correct, select **Confirm** by pressing the **✓ (Enter)** key.
- To type in the information:
- Using the on-screen keyboard, type in the appropriate Nurse and/or Patient ID. Use the arrow keys to move to the appropriate position on the keyboard and press the **✓ (Enter)** key to select each letter/number.
  - Highlight **Done** and press the **✓ (Enter)** key to save the settings.

**Note:** Press the **?** (Help) key on the BiliChek keypad for more information during any of these steps.

**Note:** This unit is equipped with a Critical Power-off feature. After completing the calibration, the unit will beep (unless the sound is set to level 0) if no measurement is performed within a preset period of time. The unit will power-off 2 minutes later if no key is pressed. If the unit automatically shuts off, it is necessary to perform another calibration after restarting.

**Note:** To exit during calibration or measurement, press the **Home** (Home) key to return to the home screen. Data will not be saved.

- Press the **?** (Help) key to see the on-screen instructions for additional information.
4. Apply a new BiliCal tip to the hand-held unit. Firmly press it on to the unit to ensure proper seating.
  5. Select **Done** to begin calibration. The unit will beep and the display will indicate Calibration successful.

**Note:** The unit will not beep if the sound is set to level 0.


6. Pull on the BiliCal tab and peel away the protective covering and calibration material from the tip and dispose, leaving only the clear film on the BiliCal tip. Press the **✓** (Enter) key to select **Done**.
7. Gently press the BiliCal tip against the patient's forehead or sternum. Lift the BiliCal tip from the patient's forehead or sternum and then repeat the process four more times.
8. A single tone sounds after each reading to indicate that it was successful. If an error occurs during the reading, a yellow error screen will appear with a description of the error. Follow the on-screen troubleshooting instructions to correct the problem.

**Note:** The unit will not beep if the sound is set to level 0.

10. When the readings are complete, the BiliChek displays and saves all entered information.
11. Remove and dispose of the BiliCal tip and return the hand-held unit to the charger base. If your unit is linked to a laboratory information system, all new data entered since the last transfer will be sent once the unit is placed back on the charger base.

## **WARNING**

*It is critical to calibrate the unit and apply a new BiliCal Individual Calibration Tip immediately prior to performing a test. Inaccurate test results will occur if this is not done. Removal of the tip after calibration and prior to a patient measurement will result in loss of calibration. The entire calibration sequence must be repeated with a new BiliCal Tip.*

**Note:** The BiliChek unit should always be stored on the charger base when it is not in use.


**TcB Records Sent Screen**

## Using Other Features

### History

Select History to access patient data that is stored in the hand-held unit. The unit stores up to 30 records. Once the record limit is reached, the oldest record will be deleted.

The top right-hand corner of each screen will contain a sent status, which indicates one of the following:

1. The record has been sent.

2. The record has not been sent.

3. The record contains an error and cannot be sent.

## Stand-by

Select Stand-by from the Home screen to turn the screen off between readings.

Stand-by Selection

## Recommendations for Use During and After Phototherapy

### During Phototherapy

- Select the site at which the BiliChek measurement will be taken. The preferred sites are the flat area of the forehead between the infant's eyebrows or the sternum.
- Be sure to hold the BiliCal Individual Calibration Tip flat against the patient's forehead or sternum until the test is complete.

Place Disposable Tip Flat on Forehead or Sternum

- A photo-opaque material (such as BilEclipse Phototherapy Protective Patch or adhesive skin patches) must be positioned over the measurement site prior to the start of phototherapy. Any necessary precautions should be taken to ensure that the material does not move and expose the measurement site to the phototherapy lights.

**Note:** Contact Respiration for recommended accessories to be used with phototherapy.

- All phototherapy lights should be turned off while a BiliChek measurement is taken. Ensure that the preferred measurement site is still covered. Open the BilEclipse protective patch light blocking material and take the BiliChek measurement from the preferred site. Close the BilEclipse flap and continue the phototherapy or remove the BilEclipse and discontinue the phototherapy as indicated. (You may choose to leave the BilEclipse in place after phototherapy is discontinued to identify the correct measurement site for future measurements.)

**Note:** Studies on the repeatability of the BiliChek measurement have established precision specifications. The source of variation in the measurement is the heterogeneous nature of skin and the bilirubin deposition. Small changes in successive transcutaneous bilirubin measurement of less than 2.0 mg/dL or 34 µmol/L (increases or decreases) should not be interpreted as definite trends.

## BilEclipse

1. Always clean the skin site before you place the BilEclipse on the infant. Respiromics makes no recommendations on the type of cleaning agent that is used. Any liquids, lotions, or creams will reduce the adhesiveness of the patch and increase the likelihood that it will fall off during phototherapy.
2. When you lift the LID of the patch to take a measurement, it may be helpful to place your finger on the edge of the BASE in order to make sure that the entire patch is not pulled off of the infant's skin.
3. When you close the LID, make sure that it is in good contact with the BASE around the entire patch. Although it is unlikely that light will penetrate the sides, any light that does reach the skin could affect the BiliChek measurement.

### BilEclipse

**Note:** Clinical studies indicate that as many as 48 hours may be required before the skin treated by phototherapy returns to the bilirubin level of an unexposed site.

**Note:** Newborns that are placed near windows with high exposure to sunlight may experience "natural phototherapy". This may be particularly true when the newborn has already been discharged home and is returning to a hospital, clinic, or office facility for testing.

4. After phototherapy has been completed, it is recommended that the BilEclipse be left on the skin until no future BiliChek measurements are likely to be made. Although no further changes in the skin will occur after the phototherapy lights are turned off, leaving the patch in place will ensure that future measurements are made from the same, protected location.

## After Phototherapy

Appropriate measures should be taken to ensure that any measurements taken from newborns after they have received phototherapy are taken from a site on the forehead or sternum that was not exposed to the phototherapy lights. The exact procedures that should be used will depend on the type of phototherapy protection used. Although leaving a BilEclipse Phototherapy Protective Patch in place is not necessary once the phototherapy lights are removed, it will ensure that the same site is used for subsequent BiliChek measurements.

# 5. Cleaning, Maintenance, and Service

## Cleaning

The BiliChek System is made of a smooth, easy to clean surface. Follow the instructions in this section any time the unit or charger base is dirty.

1. Use soapy water, a 10% bleach solution or full strength ammonia.
2. Use a damp, soft sponge or cloth to apply the cleaner.
3. Apply the cleaning solution to the sponge or cloth and wipe down the unit or charger base.
4. Allow the equipment to air dry.
5. Wipe the unit and charger base with a dry cloth.

### **WARNING**

When cleaning, DO NOT USE:

- Phenolic compound based germicide cleaner/disinfectant
- Gluteraldehyde disinfectant/sterilants
- Regular commercial cleaners or laundry detergents
- Iodine solutions, strong acids or strong alkali solutions

These solutions could leave residue on the surfaces, and /or be abrasive or harmful to the infant.

### **WARNING**

*Do not immerse the BiliChek system in water or other liquid. If liquid spills onto the unit, wipe with a damp cloth and let dry before using.*

### **WARNING**

*Do not attempt to clean and/or reuse BiliCal Individual Calibration Tips.*

### **CAUTION**


*Clean the optical measurement tip (only) with 90% or higher Isopropyl Alcohol using a soft optical surface cleaning wipe*

## Maintenance


### Changing the Measurement Light Bulb

If an error screen indicates the light bulb is not operating, replace it using a Replacement Bulb Kit (contact Resironics for ordering information).


1. Remove the battery compartment cover. Pry the cover open by inserting a small flat head screwdriver into the top of the cover and sliding the cover off.

2. Disconnect the battery.

3. Disconnect the lamp harness connector.

#### CAUTION

The BiliChek is to be used only with an approved Replacement Bulb Kit from the manufacturer. Do not over-tighten the lamp cap.

#### Removing the Battery Cover


#### Disconnecting the Battery

#### Disconnecting the lamp harness connector


**Note:** You may need to slightly lift the lamp housing to remove the lamp harness connector.

4. Remove the lamp cap using a flat tip or Phillips head screw driver.


Removing the lamp cap

5. Pull the light bulb assembly to remove.

Removing the bulb assembly

6. Insert the factory specified replacement lamp assembly.


Inserting the new lamp

7. Install the replacement lamp cap supplied with the bulb. Do not over-tighten.
8. Reconnect the lamp harness connector.
9. Reconnect the battery and close the battery compartment.
10. Apply a new BiliCal tip.
11. Turn the power on.
12. Perform multiple calibrations until the unit accepts calibration.

## Changing the Battery

1. Remove the battery compartment cover. Pry the cover open by inserting a small flat head screwdriver into the top of the cover and sliding the cover off.

**Removing the Battery Cover**

2. Disconnect the battery.
3. To install the new battery, follow the *Attaching the Battery* instructions earlier in this manual.

## Service

The measurement light bulb is the only field serviceable component which can be accessed from the battery pack opening. Do not attempt to separate the upper and lower bodies or remove the junction screw. Doing so could damage the fiber-optics and will invalidate the manufacturer's warranty.

## Customer Service Information

If you need to contact Respiromics directly, call the Respiromics Customer Service department at 1-800-345-6443, 1-724-387-4000, or +49 8152 93060.

You can also use one of the following addresses:

Respiromics	Respiromics
1001 Murry Ridge Lane	Deutschland
Murrysville, Pennsylvania	Gewerbestrasse 17
15668-8550 USA	82211 Herrsching Germany

# BiliChek

## user manual

## 6. Specifications

### Environmental

	Storage	Operating
Temperature	-20 to +50° C (-4 to 122° F)	5 to 35° C (41 to 95° F)
Relative Humidity	15 to 95% Non-condensing	15 to 95% Non-condensing

### Physical

#### BiliChek System

Total Weight: 390 g

#### BiliChek Hand-Held Unit

Weight  
(with battery pack): 250 g

Dimensions: 5.23 cm W x 20.45 cm L x 5.94 cm H (2.1 in W x 8.1 in L x 2.3 in H)

#### Charger Base

Weight: 140 g

Dimensions: 9.14 cm W x 9.01 cm L x 4.30 cm H (3.6 in W x 3.5 in L x 1.7 in H)

### Light Bulb

Bulb Type: Tungsten

Bulb Life: Minimum 30,000 measurements

## Standards Compliance

This device is designed to conform to the following standards:

- IEC 60601-1: Medical electrical equipment - Part 1: General requirements for safety.
- IEC 60601-1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests.

## Electrical Requirements

AC Power (Charger Base)	Input: 100-240 V AC, 50/60 Hz, 0.4 A Output: 5V DC, 2.0 A
Battery (Hand-Held Unit)	Type: Lithium Ion Voltage: 3.7V
Type of Protection Against Electrical Shock	Internally Powered
Degree of Protection Against Electrical Shock	Type BF Applied Part
Degree of Protection Against Ingress of Water	Ordinary Protection, IPX0

## WEEE/RoHS Recycling Directives

If you are subject to the WEEE/RoHS recycling directives, refer to [www.respirronics.com](http://www.respirronics.com) for the passport for recycling this product.

## Disposal

Dispose of device in accordance with local regulations.

# BiliChek

## user manual

# Appendix A. EMC Information

## Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for supply mains ±1 kV for input/output lines	Mains power quality should be that of a typical home or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV for common mode	Mains power quality should be that of a typical home or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical home or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment.
NOTE: $U_T$ is the a.c. mains voltage prior to application of the test level.			

## Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test	Compliance		Electromagnetic Environment - Guidance
		Level	Level	
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms		Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  <b>Recommended separation distance:</b> $d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m		$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <b>a</b> , should be less than the compliance level in each frequency range <b>b</b> .  Interference may occur in the vicinity of equipment marked with the following symbol: 

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

**a:** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

**b:** Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

## Recommended Separation Distances between Portable and Mobile RF Communications Equipment and this Device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Output of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter; where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

# **BiliChek**

## user manual

# **Appendix B. Warranty**

Respironics, Inc. warrants your BiliChek System against defects in material and workmanship for a period of one (1) year from the date of purchase. This warranty does not cover any damage caused by accident, misuse, tampering, or negligence such as failure to follow the instructions provided in this guide. In the event your unit fails to give satisfactory performance within the warranty period and conditions, Respironics, Inc. will repair or replace your device at no charge for parts or labor. The foregoing warranties are in lieu of all other warranties expressed or implied, including without limitation any implied warranty of merchantability or fitness for a particular purpose.

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