LEEP System 1000[®] Workstation, 220 VAC Model KH1000A

**Operating Manual (English)
Bedienungsanleitung (Deutsch / German)
Manual de funcionamiento (Español / Spanish)
Manuel d'utilisation (Français / French)
Manuale di funzionamento (Italiano / Italian)
Manual de Funcionamento (Português / Portuguese)**

English	Page 1
Deutsch / German	Seite 33
Español / Spanish	Página 65
Français / French	Page 97
Italiano / Italian	Pagina 129
Português / Portuguese	Página 161

CooperSurgical

LEEP System 1000[®] Workstation, 220 VAC

Model KH1000A

Operating Manual

(English)

This page left blank intentionally

**LEEP System 1000® Workstation, 220 VAC
Model KH1000A**

Table of Contents

Section	Content	Page
1.	Description.....	5
1.1	Introduction	5
1.2	Smoke Evacuator Description	5
1.3	LEEP System 1000 Electrosurgical Generator Description	5
1.4	LEEP Cart Description	5
2.	Unpacking and Assembly.....	6
2.1	Unpacking the LEEP Cart (box #1)	6
2.2	Unpacking and Installing the Smoke Evacuator (box #2)	7
2.3	Unpacking and Installing the LEEP System 1000® Electrosurgical Generator (box #3)	7
2.4	Installing the Filters and Tubing to the Smoke Evacuator	8
2.5	Connecting the Smoke Evacuator and LEEP System 1000® Electrosurgical Generator together and to the Wall Outlet.....	9
2.6	Installing the Foot Pedal Switch on the LEEP System 1000 Electrosurgical Generator	9
2.7	Installing the Electrodes and Dispersive Patient Plate	9
3.	LEEP System 1000® Electrosurgical Generator Features	11
4.	Front Panel of the LEEP System 1000® Electrosurgical Generator	12
5.	Professional Use Guide	13
5.1	LEEP Cart Use	13
5.2	Indications	14
5.3	Contraindications	14
5.4	LEEP Procedure and Technique	14
5.5	Safety Precautions	14
5.6	Electrosurgical Procedures	15
5.7	Turning on the Smoke Evacuator and LEEP System 1000® Electrosurgical Generator.....	15
6.	Operating the System	16
6.1	Setting the Controls and Output Mode	16
6.2	Guidelines for Power Settings	17
6.3	Thermal Effects on Tissue Treated with Loop Electrodes	18
7.	Electrosurgical Precautions.....	18
8.	KH1000A Accessories.....	21
9.	Periodic Safety Checks and Maintenance	21
9.1	LEEP System 1000® Electrosurgical Generator	21
9.2	Smoke Evacuator	21
10.	Safety Circuits	21
11.	Practical Suggestions	22
12.	Cleaning the LEEP System 1000® Electrosurgical Generator and Smoke Evacuator	22

Table of Contents (continued)

Section	Content	Page
13.	Troubleshooting	22
14.	Accessories	23
15.	Liability Statement	23
16.	Warranty	23
17.	Service and Repair	24
18.	LEEP System 1000 [®] Electrosurgical Generator EMC Compliance Information	25
19.	Specifications	27
19.1	Electrical – LEEP System 1000 [®] Workstation	27
19.2	General Specifications – LEEP System 1000 [®] Electrosurgical Generator	28
19.3	Smoke Evacuator	30
19.4	LEEP Cart	30
20.	Explanation of Symbols.....	31

CooperSurgical

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105
www.coopersurgical.com

International

Phone: (203) 601-9818
Fax: (203) 601-4747

Section 1 Description

1.1 Introduction

Congratulations on owning the LEEP System 1000[®] Workstation, which includes the following components:

LEEP System 1000 Workstation, 220 VAC (Model KH1000A)			Accessories (filters, tubing, etc.)
LEEP System 1000 [®] Electrosurgical Generator 220 VAC	Smoke Evacuator 220 VAC	LEEP Cart	

Before you can use this equipment, you must assemble the LEEP System 1000 Electrosurgical Generator and the Smoke Evacuator onto the LEEP Cart. Refer to Section 2 for Unpacking and Assembly instructions.

1.2 Smoke Evacuator Description

The CooperSurgical Smoke Evacuator System three-stage air filtration system is used to remove airborne particulate plume produced during office and surgical procedures and has the following features:

- Low noise level
- Triple filtration of air provides efficiency level for 0.014 microns rated at 99.999%. This includes a pre-filter, a charcoal filter for odor removal and a final safety filter placed after the charcoal filter
- Adjustable high air flow for effective collection of plume
- Virtually maintenance-free
- Conveniently attaches to the CooperSurgical LEEP System 1000 Workstation

1.3 LEEP System 1000 Electrosurgical Generator Description

The Electrosurgical Generator has the following features:

- Isolated power output and LED display located in the front for precise power selection, delivery and ease of use
- Flush faceplate membrane facilitates operation and cleaning
- Microprocessor-controlled for increased precision, accuracy, reproducibility and safety
- Provides a choice of CUT, BLEND and COAG waveforms to accommodate subtle differences in technique and electrode performance
- Pneumatic Foot Pedal for maximum safety
- Audible safety features include distinct tones for each operating setting
- Automatic self-test mechanism ensures accurate system operation
- Integrated Smoke Evacuator controls

1.4 LEEP Cart Description

Advantages of this style cart:

- Added mobility and functionality for transporting the Electrosurgical Generator and Smoke Evacuator in one unit
- Designed to work in small exam rooms
- Heavy-duty casters ensure easy mobility
- Convenient interior storage shelves
- Elegant design for modern medical facilities

Section 2 Unpacking and Assembly

Locate all three product boxes.

2.1 Unpacking the LEEP Cart (box #1)

When it comes to unpacking the carton that contains the Cart, great care must be taken. The box must be on its side (following the markings on the carton for which side is up) and the Cart rolled out. Orient the box as shown in Photo A.

CAUTION: DO NOT PICK THE CART UP OUT OF THE SHIPPING BOX BY ITS HANDLES. THE HANDLES ARE NOT DESIGNED TO ACCOMMODATE THE FULL WEIGHT OF THE CART.

1. After removing the Cart from the shipping box, lock the front two casters.
2. Install the Smoke Evacuator next.

Photo A


Figure 1

2.2 Unpacking and Installing the Smoke Evacuator (box #2)

Required Tools: Phillips Screwdriver

After unpacking the Smoke Evacuator box,

1. Remove the screws on the base of the Cart in preparation for attaching the Smoke Evacuator.
2. Then align the pivot hinges on the Smoke Evacuator with the screw holes on the Cart itself and fasten the Smoke Evacuator to the Cart. Make sure that the screws are tight, but DO NOT OVERTIGHTEN.
3. Remove the screws on the side of the Smoke Evacuator and then fasten down the bracket of the drop hinge to the Smoke Evacuator. See Figure 1.

NOTE: The Smoke Evacuator Controls are located at the top front of the LEEP Cart for easy access.

4. Install the Electro-surgical Generator next.

2.3 Unpacking and Installing the LEEP System 1000® Electro-surgical Generator (box #3)

After unpacking the Electro-surgical Generator box, align the pins on the bottom of the Electro-surgical Generator with the corresponding holes located on the top of the Cart. See Figure 1.

Face the Electro-surgical Generator toward the front of the Cart (where the door is) and drop the cord through the rectangular cutout located at the top of the Cart. Do not hook up this cord until you have installed the filter and tubing in the Smoke Evacuator.

2.4 Installing the Filters and Tubing to the Smoke Evacuator

2.4.1 ULPA Filter Installation

Tilt the Smoke Evacuator forward and insert the large ULPA Filter Cylinder with the air-flow arrow pointing down. (See Photo B)

Photo B

Photo C

2.4.2 Installing the Pre-Filter

Insert a clean, disposable Pre-Filter onto the ULPA Filter Cylinder. (See Photo C). Be sure this device is firmly seated.

2.4.3 Connecting the Tubing

There are two tubing hookup options:

- For procedures requiring close-proximity plume removal (i.e., Vaginal Speculum)
- For procedures requiring open-area plume removal (i.e., external lesions)

For procedures requiring close-proximity plume removal (i.e., Vaginal Speculum)

Assemble the $\frac{3}{8}$ " Reducer (REF 6083) onto the port on the disposable Pre-Filter (REF 6081) top with a slight twisting motion. Attach one end of an appropriate length of $\frac{3}{8}$ " ID Evacuation Tubing (REF 6084) to the Reducer connector and direct the other end to the patient and any appropriate device being used, such as a Vaginal Speculum equipped with a Smoke Evacuation Adapter. (See Photo D)

Photo D

For procedures requiring open-area plume removal (i.e., external lesions)


Assemble the sterile, disposable $1\frac{1}{4}$ " ID Evacuation Tubing (REF 6085) directly into the top of the Pre-Filter. Position the opposite end over the site to be treated.

2.5 Connecting the Smoke Evacuator and LEEP System 1000® Electrosurgical Generator together and to the Wall Outlet

1. Plug in the cord running from the Electrosurgical Generator to the receptacle at the back of the Smoke Evacuator [refer to (A) in Figure 1].
2. **Make sure both On/Off Switches on the Electrosurgical Generator panel are in the OFF position.** Refer to **1** and **12** in Figure 2. Next, attach the power cord to the receptacle at the back of the Electrosurgical Generator and then into a hospital-grade grounded wall outlet to obtain power [refer to (B) in Figure 1].

2.6 Installing the Foot Pedal Switch on the LEEP System 1000 Electrosurgical Generator

Connect the Foot Pedal Switch to the socket **2**, shown in Figure 2, without activating the Foot Pedal, and tighten the threaded plug. This is an air (pneumatic) operated control. There is no electric current, offering maximum safety.

2.7 Installing the Electrodes and Dispersive Patient Plate

2.7.1 Placement of the Electrode

- Connect the active Electrode Handpiece to the socket **3**, shown in Figure 2, and tighten the electrode of choice in the Handpiece.

2.7.2 Placement of the Patient Plate or Dispersive Electrode

When using an electrosurgical system, it is very important that all the current delivered to the patient returns correctly to the Electrosurgical Generator via the Dispersive Patient Plate only.

- Connect Dispersive Patient Plate to the socket **4**. Refer to Figure 2.
- The patient must be positioned correctly on the operating table. The patient and operator must not come in contact with any metal conductive surfaces.
- The Patient Plate must securely contact a vascular area close to the operating site. For a gynecology procedure the preferred sites are the patient's thigh (disposable adhesive pads) or under the patient's buttocks (reusable metal plate). The contact area must be clean, free of body lotions, shaved, and massaged for good circulation. The contact area of the Patient Plate must be maximized and frequently checked for uniform contact during the procedure, especially if the patient has moved or if liquids have contacted the Patient Plate. A CONDUCTIVE GEL IS RECOMMENDED. The Patient Plate MUST NEVER be placed so as to allow the patient's heart to be in the pathway from the active electrode.
- Power delivery to the operative site may be decreased appreciably if alternate pathways exist; for example, through the metal operating table, crossed Handpiece/Patient Plate Cables, etc.

Figures 3 through 5 show the proper and improper ways of hooking up and using the various electrodes and pads on the patient.

PROPER


Figure 3

IMPROPER


Figure 4

IMPROPER


Figure 5

Section 3 LEEP System 1000[®] Electrosurgical Generator Features

- Microprocessor-controlled for increased precision, accuracy, repeatability and safety
- Adequate power for all LEEP monopolar electrosurgical procedures
- Accurate selection of discrete power levels
- Digital display of output power levels
- Choice of radio-frequency waveforms including CUT, BLEND and COAG to accommodate subtle differences in technique and accessory performance
- Patient Plate continuity monitoring with audible alarm
- Distinct audible tones for CUT/BLEND modes and COAG mode with associated MODE light
- Fully regulated isolated output power
- Meets or exceeds IEC 601-2-2, second edition
- Non-electric pneumatic Foot Pedal to maximize safety
- Choice of reusable or disposable Patient Plate
- Choice of reusable or disposable Handpiece
- Choice of reusable or disposable electrodes
- Output power safety audible alarm with automatic power shut-off
- Class 1, type BF, protected for use with defibrillator
- Membrane switching to maximize cleanliness and ease of use

Section 4 Front Panel of the LEEP System 1000® Electrosurgical Generator

(Colored, numbered boxes are also located later in this manual.)



Figure 6

Workstation Controls

- | | |
|---|---|
| <ul style="list-style-type: none"> 1. Main Switch (On/Off) 2. Socket for pedal switch 3. Socket for active electrodes 4. Socket for neutral electrode 5. Warning light of neutral electrode alarm (red) 6. Coagulation light (blue) 7. Pure cut and blend light (yellow) | <ul style="list-style-type: none"> 8. Mode control: <div style="display: flex; align-items: center; margin-top: 10px;">  <ul style="list-style-type: none"> Pure Cut Blend Cut Coag </div> <ul style="list-style-type: none"> 9. Power control 10. Display |
|---|---|

Smoke Evacuator Controls

- 11. Suction Control Knob
- 12. Smoke Evacuator Power Switch (On/Off)

Symbols on the LEEP System 1000 Electrosurgical Generator

	Classification I Type BF protected against defibrillator effects
	Floating output circuit (Applied Part)
	Cautions – consult this manual for safety precautions
	Pedal connection
	Active handle connection
	Patient Plate connection
	High voltage
	Temperature limitation

IMPORTANT

The user of the LEEP System 1000[®] Electrosurgical Generator should be thoroughly trained in the techniques of Loop Electrosurgical Excision Procedures. This equipment has been designed for use with LEEP Electrosurgical Accessories. DO NOT use this equipment for any purpose other than that for which it has been designed. See warnings and caution statements throughout this manual.

Section 5 Professional Use Guide

This manual contains information about the proper procedures for inspecting and preparing the Electrosurgical Generator before its use, and for its care and storage after use.

This manual does not describe how an actual procedure is to be performed, nor is it meant to teach a beginner the proper technique or any of the medical considerations regarding the use of this equipment. CooperSurgical recommends that the prospective user obtain appropriate training before using this equipment, as improper use could be potentially hazardous to the patient and the user.

This device SHOULD NOT be used without proper training.

Training in the use of electrosurgical equipment should include:

- 1. A review of the published literature regarding the procedure of interest**
- 2. A review of the Loop Electrosurgical Excision Procedure (LEEP) Syllabus available from CooperSurgical**
- 3. Attendance at a course or courses offered by physicians experienced with the Loop Electrosurgical Excision Procedure**
- 4. Hands-on preceptor training from an experienced practitioner**

Read this entire manual carefully to become familiar with each of the controls and features before making any attempt to use the equipment clinically.

Instructions contained in the operating manuals of any equipment to be used in conjunction with this equipment must be followed to avoid any possible hazard from incompatibility.

Failure to thoroughly understand and follow the instructions given in this manual may result in serious injury to the patient and/or the operator. Failure to follow the instructions given in this manual may result in damage to or malfunction of this equipment.

No long-term follow-up studies with this device have been performed as to recurrence rates. The effects of Loop Electrosurgical Excision Procedure on pregnancy outcome are not known.

SAFETY PRECAUTIONS MUST ALWAYS BE EXERCISED WHEN USING ELECTRICAL EQUIPMENT TO PREVENT OPERATOR/PATIENT SHOCK, FIRE HAZARD AND EQUIPMENT DAMAGE.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. This device SHOULD NOT be used without proper training and preceptorship.

If any questions arise regarding the information contained in this manual, the operation or safety of the equipment or service, please contact CooperSurgical.

5.1 LEEP Cart Use

VERY IMPORTANT: The LEEP System 1000[®] Workstation must only be moved by grasping and holding firmly onto the Cart Handles to ensure the Cart does not tip over.

The Cart must only be pushed or pulled using the handles on the front or back of the Cart to ensure stability during transportation.

5.2 Indications

The LEEP procedure is indicated in the diagnosis and treatment of some Cervical Intraepithelial Neoplasia (CIN) in patients where there is:

- Cytological or colposcopic suspicion of CIN 2 or worse (including micro-invasion)
- Persistent CIN 1 (of more than 12 months duration)
- CIN 1 where the likelihood of follow-up is low or when the patient requests treatment
- A suspicion (cytological or colposcopic) of a glandular intraepithelial abnormality
- A disparity between the cytological and colposcopic diagnoses
- External anogenital lesion
- Large vaginal intraepithelial neoplastic (VAIN) lesions
- Cervical conization indications

5.3 Contraindications

The following are typical contraindications for performing the LEEP procedure. It is imperative that the physician carefully weigh the risks and benefits of treatment versus non-treatment in contraindicated patients:

- Pregnancy
- Gross invasive carcinoma of the cervix
- A bleeding disorder
- Acute or active inflammation of the cervix, endometrium, fallopian tube, ovary or peritoneum (cervicitis, endometritis, tubo-ovarian inflammatory disease or pelvic inflammatory disease)
- “Positive” endocervical curettage or a lesion in which the endocervical limit cannot be visualized colposcopically
- Less than three months postpartum
- Equivocal cervical abnormality

5.4 LEEP Procedure and Technique

It is recommended that the patient be provided with a brief description of the procedure and the equipment that will be used (ACOG, CooperSurgical and other professional organizations and equipment manufacturers have produced patient information brochures on the LEEP procedure that address many of the questions and concerns that your patients may have regarding the procedure).

5.5 Safety Precautions

1. This equipment should only be used by a thoroughly trained physician in an adequately equipped medical facility.
2. Replacement accessories and patient return pads should be kept on hand since defective active accessory or patient return pads can result in substandard performance of this equipment.
3. This equipment should only be connected to a properly grounded receptacle. NEVER use an adapter that defeats the ground of the built-in three (3) prong plug.
4. Care must be exercised when handling liquids around electrical equipment. DO NOT attempt to operate this equipment if liquids have spilled on the Electrosurgical Generator. DO NOT use flammable liquids around electrical equipment.
5. This equipment should never be used in conjunction with other equipment for which safety against leakage current has not been established.

6. When this equipment is operated:
 - a. A Patient Return Pad (dispersive pad) of adequate surface area MUST be properly attached to the patient or the risk of accidental burns will exist.
 - b. The Patient Return Pad (dispersive pad) should be placed as close as possible to the site of use of the active accessory – but MUST NEVER be placed so as to allow the patient's heart to be in the pathway from the active accessory to the return electrode!
7. The user should thoroughly understand the principles and use of radio frequency (RF) current before using this equipment. This understanding is essential to avoid the hazard of shocks or burns to the user and/or the patient.
8. The instructions for use described in this manual must be followed; otherwise, compromised safety, malfunction, injury to the operator and/or patient, or costly damage to the Electrosurgical Generator may occur.
9. There are no user-serviceable parts within the housing. Repairs to this equipment should only be performed by authorized CooperSurgical service personnel. For service information, please contact CooperSurgical (refer to Section 15).

5.6 Electrosurgical Procedures

This section provides only general information about the use of electrosurgical devices. Only the user can evaluate the clinical factors involved with each patient and determine if the use of this equipment is indicated. The user must then decide on the specific technique and procedure that will accomplish the desired clinical effect.

WARNING

Electrosurgical generators are designed to allow the controlled destruction of tissue and are inherently dangerous if operated improperly.

REPORTED PROBLEMS DUE TO IMPROPER OPERATION DURING ELECTROSURGICAL PROCEDURES HAVE INCLUDED:

- Inadvertent activation with resultant tissue damage at the wrong site and/or equipment damage
- Alternate current pathways resulting in burns where the patient or physician or assistant is in contact with exposed metal
- Explosions caused by electrosurgical sparking in a flammable gas mixture (i.e., explosive anesthetic gases and the inappropriate use of alcohol and other flammable liquids)
- Perforation and massive hemorrhage

A proper patient return pad pathway is extremely important during any monopolar electrosurgical procedure. Every effort must be made to ensure that, throughout the electrosurgical procedure, an adequate surface area is provided and remains in proper contact with the patient to reduce the current density below a level that might cause inadvertent tissue damage where the patient return pad has been applied.

5.7 Turning on the Smoke Evacuator and LEEP System 1000[®] Electrosurgical Generator

NOTE: Numbers in **blue boxes** refer to those component parts shown in Figure 6.

Push the POWER switch **12**, located under the Electrosurgical Generator (for the Smoke Evacuator), to the "On" position. This will start the Smoke Evacuator.

Turn ON the Electrosurgical Generator using the power switch **1**. The Electrosurgical Generator automatically performs a SELF-TEST that checks RAM memory, EPROM memory, supply voltage, signal modulation and the following displays: function selector green lights, digital display, cut and coagulant yellow light, coagulation blue light, and the audible signal. When the Electrosurgical Generator passes SELF-TEST, the display will show current software revision, i.e., r2A, r2B, for several seconds, then will go blank.

Section 6 Operating the System

NOTE: Numbers in blue boxes refer to those component parts shown in Figure 6.

- CUT MODE:** Control **8** sets the cut mode power output. The Electrosurgical Generator automatically powers up to cut mode when turned ON. The top LED lights on the Selector **8** and Control **9** sets the power.
- BLEND MODE:** Control **8** selects blend mode (center LED) and Control **9** sets the power.
- COAGULATION MODE:** Control **8** selects coagulation mode (bottom LED) and Control **9** sets the power.

Power settings are stored when the system is ON and will appear automatically on the Display **10** according to the selection of power mode by the Control **8** during the procedure.

Power may be changed at any time during the operation, except when the Electrosurgical Generator is activated by the Foot Pedal.

At the end of the procedure turn the Electrosurgical Generator OFF, and safely store the equipment and accessories. Power will reset to zero.

6.1 Setting the Controls and Output Mode

6.1.1 Using the Smoke Evacuator and Disposing of the Filters

Adjust the Suction Control Knob **11** to the desired level.

At the completion of each procedure activate the system to ensure safe particle containment. Using gloves and a mask, remove the Pre-Filter, the Reducer and the used section of Suction Tubing and discard into an infectious waste receptacle (See the **CAUTIONS** box). The Smoke Evacuator should be stored with a new Pre-Filter and Reducer in place on the ULPA Filter.

NOTE: The expected life of the ULPA Filter is three to six months, depending on usage. It must be discarded into an infectious waste receptacle if plume odor is detected or the suction is diminished.

CAUTIONS

This device produces a strong vacuum force, and care should be exercised to ensure that the suction control and the position of the inlet end of the Suction Tubing are adjusted properly to prevent injury to the patient or inadvertent damage to surgical materials.

The materials removed from the plume by this device are potentially hazardous. Handle according to 29 CFR 1910.1030 and OSHA 3127.1992 (Occupational Exposure to Blood Borne Pathogens) guidelines.

To prevent a fire or explosion hazard, do not use the system in the presence of flammable or potentially flammable materials.

Do not allow fluid to be pulled into the system.

To prevent premature failure of the ULPA Filter Cylinder, do not operate this device without a disposable Pre-Filter in place.

6.1.2 Electrosurgical Tissue Effect

Delivery of continuous sinusoid waveform currents through a small electrode at appropriate power levels can cause rapid heating of the intracellular fluids in the cells in close proximity to the electrode, turning these fluids to steam. The significant increase in volume (approximately five times) causes cellular structure to rupture, creating the clinical effect of "CUT", with little or no hemostatic effect along the margin of the divided tissue. Delivery of short-duration pulses of RF currents through a small electrode at appropriate power levels can cause heating of intracellular fluids at a more gradual pace. This allows evaporation of these fluids without rupturing the cellular structure, creating the clinical effect of desiccation, or "COAG", without the division of tissue.

By varying the pulse to an intermediate duration, it is possible to get a clinical effect that combines, or “blends”, the clinical characteristics of CUT and COAG, yielding the effect referred to as “BLEND”, where tissue is divided with a desirable amount of hemostatis along the margins of the divided tissue.

The electrosurgical effect may vary throughout the procedure, requiring the operator to adjust the relative power setting of the Electrosurgical Generator.

6.1.3 Select the Output Mode (i.e., “CUT”, “BLEND” or “COAG”, by Pushing the Corresponding Buttons

Output Mode	Waveform Description	General Effect
CUT	Continuous 450 kHz sinusoid with minimal modulation	Cutting without Hemostatis
BLEND	Interrupted 450 kHz sinusoid intermediate-duty cycle	Cutting with minimal Hemostatis
COAG	Bursts of 450 kHz sinusoid short-duty cycle	Coagulation without Cutting

6.1.4 Set the level of Output Power (confirmed on the digital display) by using the Output Power Selector Buttons as Desired

WARNING

The degree and speed of electrosurgical effect is largely dependent on Current Density at the point of contact of the active electrode. Loop electrosurgical excision procedure electrodes from other manufacturers may vary in the diameter thickness, size and configuration of the cutting wire. This may result in SIGNIFICANT changes in the electrosurgical effect at a given output power level setting. The use of CooperSurgical LEEP Electrodes is recommended.

6.2 Guidelines for Power Settings

The following guidelines for power settings may vary due to the technique, clinical circumstances, accessory style, cutting wire diameter, size, configuration and user preference.

Recommended Power Settings (Watts) for the CooperSurgical LEEP System 1000® Electrodes

Style	Loop Width (cm)			Ball Electrodes	Needle Electrodes	
	1.0	1.5	2.0		Macro	Micro
CUT*	—	—	—	—	—	14–24
BLEND	22–36	30–40	34–46	—	—	—
COAG	—	—	—	40–56 ⁺	30–36	—

NOTE

- * If the cut mode is desired, use the recommended settings for the blend mode.
- + The power setting can be increased beyond 56 to coagulate any bleeding points as required.

Wattage Suggestions for Fischer Cone Biopsy Excisor™

Part Number	Size	Wattage Rating	Part Number	Size	Wattage Rating
900-150	Small	30-35 Watts	900-152	Large	45-50 Watts
900-151	Medium	40-45 Watts	900-157	Small Wide	40-45 Watts
900-154	Medium Extended	40-45 Watts	900-158	Small Wide Extended	40-45 Watts
900-155	Large Shallow	40-45 Watts			

The above settings are for reference only and can be varied based on specific situations and the experience of the operator.

REMEMBER, THIS IS NOT AN ATTEMPT TO TEACH ELECTROSURGICAL TECHNIQUE.

The practitioner who lacks experience should not attempt the procedures described below based solely on this information; instead, the skills required should be acquired in the time-honored preceptor manner. Call CooperSurgical for information on courses that offer instruction on the proper use of electro-surgical generators and accessories.

NOTE: The best initial effect is accomplished with the cutting wire in only light contact with tissue. Tight pressure may cause desiccation of the tissue and will delay the start of the cutting effect.

If the use of other output modes is anticipated, repeat steps in Sections 6.1.3, 6.1.4 and 6.2 as desired. The output power level settings selected for each output mode will be retained as long as the Electro-surgical Generator remains ON.

IMPORTANT

The initial use of any electro-surgical generator always involves some degree of “trial and error”. This is true even when only changing from numbered dials to digital display models within the same manufacturer’s product line. As with any other therapeutic device, it is very helpful to experiment IN VITRO or on animal sample tissue before using any electro-surgical generator or methods which are not familiar.

The microprocessor control system of the Electro-surgical Generator was developed specifically to provide the best possible performance for loop electro-surgical excision procedures. By exhibiting patience and following the guidelines offered, the practitioner should easily become familiar with the performance characteristics of the LEEP System 1000[®] Electro-surgical Generator.

6.3 Thermal Effects on Tissue Treated with Loop Electrodes

Thermal effects on tissue specimens may include:

- 1) Thermal coagulation injury of the cervix, up to one-third the thickness of normal epithelium of the cervix;
- 2) Fragmentation of squamous epithelium of the cervix attributable to long exposure periods along the excision site that allows heat to dissipate laterally;
- 3) Partial coagulation of the endocervical epithelium because of lateral radiation of heat. Therefore, the Loop Electro-surgical Excision Procedure (LEEP) may produce thermal effects at the periphery of the excised tissue and may make histopathologic interpretation difficult or impossible and therefore, may not allow accurate diagnosis and may obscure the need for further treatment.

Section 7 Electro-surgical Precautions

The safety and effectiveness of electro-surgery is dependent, to a large degree, upon the skill of the user/operator. It is important that the user/operator read, understand and follow the operating instructions supplied with the CooperSurgical LEEP System 1000[®] Electro-surgical Generator and thoroughly understand the principles and use of radio-frequency (RF) electro-surgical systems.

WARNING: Electro-surgery uses radio-frequency energy to cut and coagulate tissue. Because of the sparking and heat associated with electro-surgery, do not use with flammable anesthetics or other flammable gases, near flammable fluids or objects, or with oxidizing agents.

- DO NOT use electro-surgery in the presence of flammable gases, flammable liquids or flammable objects in oxygen-enriched atmospheres, in nitrous oxide (N₂O) atmosphere or in the presence of other oxidizing agents.
- Prevent accumulation of oxygen, nitrous oxide (N₂O) and flammable gases under surgical drapes or within the area where electro-surgery is performed; moreover, avoid such accumulation in cases of thorax or head operations unless safely aspirated.
- Verify that all oxygen connections are leak-free before and during the use of electro-surgery.
- DO NOT use electro-surgery in the presence of naturally occurring flammable gases which may accumulate in body cavities such as the bowel.

- DO NOT use electrosurgery in the presence of flammable liquids, such as skin-prepping agents. Avoid pooling of flammable liquids near the electrosurgery site or in human body cavities such as the umbilicus or vagina.
- DO NOT place the electrosurgery active electrode near or in contact with flammable materials, such as cotton, wool or gauze. The active electrode is hot from use and can cause fire.
- It is possible that the radio frequency can interfere with the electronic circuitry in pacemakers. To reduce the risk, locate the patient-return electrode as close as possible to the treatment site and ensure that the current path between the surgical site and the patient-return electrode is as far removed from the heart as possible. For the gynecology procedures locate the patient-return electrode on the patient's upper thigh or under the buttocks. **Always monitor pacemaker patients during surgery.** In case of doubt ask the pacemaker manufacturer and/or cardiology department.
- In case of loss of power, turn the system off.
- The possibility exists that the radio frequency can interfere with other medical equipment when the electrosurgical system is operating. To reduce the interference, physically separate the device, utilize different electrical outlets that are hospital-grounded, do not allow cables to come in contact with each other, and utilize shielded devices where possible.
- The fixed output power level should be adjusted to the lowest power setting that will successfully complete the procedure. Refer to the following recommended power settings for the CooperSurgical LEEP System 1000® Electrodes.

Recommended Power Settings (Watts) for the CooperSurgical LEEP System 1000® Electrodes

Style	Loop Width (cm)			Ball Electrodes	Needle Electrodes	
	1.0	1.5	2.0		Macro	Micro
CUT*	—	—	—	—	—	14–24
BLEND	22–36	30–40	34–46	—	—	—
COAG	—	—	—	40–56 ⁺	30–36	—

NOTE
* If the cut mode is desired, use the recommended settings for the blend mode.
+ The power setting can be increased beyond 56 to coagulate any bleeding points as required.

Wattage Suggestions for Fischer Cone Biopsy Excisor™

Part Number	Size	Wattage Rating	Part Number	Size	Wattage Rating
900-150	Small	30-35 Watts	900-152	Large	45-50 Watts
900-151	Medium	40-45 Watts	900-157	Small Wide	40-45 Watts
900-154	Medium Extended	40-45 Watts	900-158	Small Wide Extended	40-45 Watts
900-155	Large Shallow	40-45 Watts			

The above settings are for reference only and can be varied based on specific situations and the experience of the operator.

- Skin to skin contact, for instance between the patient's arm and body, should be avoided by the placement of an appropriate separating device such as two to three inches of dry gauze. This will reduce the potential for alternate site burns.
- If monitoring, stimulation, imaging or similar devices are used simultaneously with electrosurgery, the monitoring electrodes must be placed as far as possible from the electrosurgery site and the patient-return electrode. Position the patient-return electrode close to the electrosurgery site, for example, on the thigh when treating the cervix. **NOTE:** Monitoring needle electrodes are not recommended.

- The electrical cord of the Electrosurgical Generator should be connected to a properly grounded receptacle. Extension cords and/or adapter plugs should not be used.
- The connecting cables to the electrosurgery electrodes should be placed so that they do not come in contact with the patient, or with other cables, or cross each other.
- Remove any metal items from the patient: for example, rings, chains, etc.
- Use accessories supplied by CooperSurgical; they are specifically designed for the LEEP System 1000 Workstation.
- Do not use old or worn accessories.
- Avoid High Frequency output settings where the Maximum Output Voltage may exceed the Rated Accessory Voltage.
- For surgical procedure where the High Frequency current could flow through parts of the body having a small cross-sectional area, the use of Bipolar techniques may be desirable in order to avoid unwanted tissue damage.
- Apparent low output or failure of the High Frequency surgical equipment to function correctly at normal operating settings may indicate faulty application of the Neutral Electrode or poor contact in its connections. In this case, the application of the Neutral Electrode and its connections should be checked before selecting a higher output power.
- Regularly inspect the accessories, including the electrode cables and any endoscopically used accessories, for possible damage.
- Failure of the High Frequency Surgical Equipment could result in an unintended increase of output power.

Accessory Precautions

- Rated Accessory Voltages should be determined as follows using the Maximum Output Voltage for each High Frequency Surgical Mode:
 - > For situations where the MAXIMUM OUTPUT VOLTAGE is less than or equal to 1600V, ASSOCIATED EQUIPMENT and ACTIVE ACCESSORIES should be selected that have a RATED ACCESSORY VOLTAGE equal to or greater than the MAXIMUM OUTPUT VOLTAGE.
 - > ASSOCIATED EQUIPMENT and ACTIVE accessories should be selected with RATED ACCESSORY VOLTAGE \geq MAXIMUM OUTPUT VOLTAGE when the smaller variable y [see below] or the number 6 is \leq CREST FACTOR for that High Frequency SURGICAL MODE.
 - > When MAXIMUM OUTPUT VOLTAGE (U_{max}) is $> 1600V$, and the CREST FACTOR is $<$ the variable y calculated below, indicating that any ASSOCIATED EQUIPMENT and ACTIVE ACCESSORIES used with such mode or setting must be rated to withstand the combination of actual voltage and CREST FACTOR.

$$y = \frac{U_{max} - 400 \text{ [volts]}}{600 \text{ [volts]}}$$

- Neutral Electrode Patient Return Pad requirements:
 - » Adult Single Dispersive Electrode Pad
 - » Single-Use only
 - » Foam pad with adhesive backing
 - » 19-square-inch gel surface area minimum
 - » 2-conductor wire connection for continuity monitoring
- ESU Pencil requirements:
 - » Single-Use
 - » Crimped wire Electrosurgical Generator Connector
 - » Footswitch Activated
 - » IEC 60601-2-2 Compliant

Section 8 KH1000A Accessories

NOTE: Use only genuine CooperSurgical KH1000A accessories (patient return pads, hand pencils, electrodes and smoke evacuation disposables) for optimal system performance and patient safety. The use of non-CooperSurgical authorized and genuine accessories is not recommended nor have been tested and verified to the level of safety and performance of a genuine CooperSurgical accessory. Contact CooperSurgical for a list of current genuine KH1000A accessories.

Section 9 Periodic Safety Checks and Maintenance

9.1 LEEP System 1000® Electrosurgical Generator

The following safety checks should be performed at least every 24 months by a qualified person who has adequate training, knowledge and practical experience to perform these tests.

- Inspect the equipment and accessories for mechanical and functional damage
- Inspect the relevant safety labels for legibility
- Inspect the fuse to verify compliance with rated current and breaking characteristics
- Inspect acoustical and visual alarms/displays
- Verify that the device functions properly as described in the Directions for Use
- Verify that the device shuts down the applied part circuit if the neutral electrode is disconnected
- Test the protection earth resistance according to IEC 601-1/1988: Limit 0.2 ohm
- Check that the output power is within tolerance versus the output control setting at a specified load resistance
- Check the power output at full and half setting of the output control over the range of load resistance as specified in the instructions for use (max. deviation is $\pm 20\%$)
- Test the enclosure leakage current according to IEC 601-1/1988: Limit 100 μ A
- Test the leakage current according to IEC 601-1/1988: Limit 100 μ A (BF)
- Test the patient leakage current under single fault condition with main voltage on the applied part according to IEC 601-1/1988: Limit: 5mA (BF)

The leakage current should never exceed the limit. The data should be recorded in an equipment log. If the device is not functioning properly or fails any of the above tests, the device must be repaired.

9.2 Smoke Evacuator

No maintenance is required during normal operation other than ensuring that ample space is maintained around the Smoke Evacuator to allow for free air flow and adequate cooling.

CAUTION

Make sure that the Smoke Evacuator unit is not operated for more than 15 minutes at a time. Allow the Smoke Evacuator to cool off a minimum of 15 minutes before turning back on.

Section 10 Safety Circuits

The Electrosurgical Generator is equipped with two safety circuits. The first one checks the Dispersive Patient Plate connection. The second one turns off the power in case of an internal failure. When activating the Electrosurgical Generator by the Foot Pedal, a power delivery higher than the one selected will stop power delivery and at the same time give an audible signal similar to the Patient Plate alarm, but at a higher frequency.

Section 11 Practical Suggestions

To optimize the performance when using an electrosurgical unit, the active electrode must be kept clean and use the lowest possible power setting required. Some sparks or superficial carbonization of the tissue may occur and the delivered power may decrease as a result of the electrical insulation caused by the tissue charring.

Too high of a power setting results in a shorter surgical procedure, but may cause discharges and/or superficial carbonization, sparking, arcs, etc.

Section 12 Cleaning the LEEP System 1000® Electrosurgical Generator and Smoke Evacuator

LEEP System 1000 Electrosurgical Generator: The Electrosurgical Generator may be cleaned with mild soap solution – but be sure that fluid does not enter the system. Wipe dry.

Smoke Evacuator: The outside of the Smoke Evacuator may be cleaned as necessary with a soft cloth dampened (not wet) with isopropyl alcohol. Just as for the Electrosurgical Generator, do not allow liquid to get inside of the unit.

Section 13 Troubleshooting

In the event of a failure during SELF-TEST, the display will show one of the error codes listed below:

ALARM SIGNALS		
PROBLEM	AUDIBLE TONE	DISPLAY
RAM Memory (during self-test phase)	1 kHz 100ms ON 250ms OFF	“Er0”
EPROM Memory (during self-test phase)	1 kHz 100ms ON 250ms OFF	“Er1”
Signal Modulation (during self-test phase)	1 kHz 200ms ON 250ms OFF	“Er2”
Signal Modulation (during activation phase)	1 kHz 40ms ON 60ms OFF	“Er2”
Supply Voltage (during self-test phase)	1 kHz 100ms ON 250ms OFF	“Er3”
Supply Voltage (during activation phase)	1 kHz 40ms ON 60ms OFF	“Er3”
Power Output (during activation phase)	1 kHz 100ms ON 250ms OFF	“Er5”
Foot Pedal Circuit (during self-test phase)	1 kHz 100ms ON 250ms OFF	“Er6”
Microcontroller Power Supply	1 kHz 100ms ON 250ms OFF	“Er7”
Dispersive Electrode	1 kHz 80ms ON 125ms OFF	“nP” & Red Light

- A. If upon completion of Self-Test the Electrosurgical Generator display reads Er0, Er1, Er2, Er3, Er5, Er6 or Er7, it must be returned to CooperSurgical for repair (refer to Section 15).
- B. If upon completion of Self-Test or before, if the Electrosurgical Generator stops functioning and the display reads nP and the red warning light is on, check the Dispersive Patient Plate to ensure that it is properly connected to the Electrosurgical Generator.
- C. If, following correct set-up of the system, it does not operate, works in an intermittent way, or stops working after a few seconds (without an audible signal), check for correct connection of pedal and its condition. The pedal is pneumatic, so even a slight leak can cause a performance problem.

Proceed as follows:

1. Tighten the threaded plug into the Foot Pedal Socket.
 2. Then, push the pedal hard repeatedly to detect possible breaks in the tubing or in the pedal.
- D. If the Electrosurgical Generator, correctly connected, appears to deliver a lower power output than usual, check:
1. The Dispersive Patient Plate for complete contact (refer to Dispersive Plate sections).
 2. The condition of active electrodes (refer to “Section 10, Practical Suggestions”).
 3. The condition of the handpiece (cable continuity, contact of the electrode in the handpiece and the connector) by moving the cable and the connector and electrode to detect possible breaks and poor contact in the socket of the system.
- E. The Electrosurgical Generator has a thermal protection circuit which will shut it off when internal operating temperatures exceed safe limits. Should the Electrosurgical Generator stop functioning with no alarm signal, ensure that the system has adequate ventilation and that you have not exceeded the recommended duty cycle of 10/30 seconds. Should the thermal protection circuit be activated under normal operating conditions, the Electrosurgical Generator must be returned to CooperSurgical for servicing (see Section 15).

Section 14 Accessories

NOTE: Use only genuine CooperSurgical KH1000A accessories (patient return pads, hand pencils, electrodes and smoke evacuation disposables) for optimal system performance and patient safety. The use of non-CooperSurgical authorized and genuine accessories is not recommended nor have been tested and verified to the level of safety and performance of a genuine CooperSurgical accessory. Contact CooperSurgical for a list of current genuine KH1000A accessories.

Section 15 Liability Statement

CooperSurgical guarantees the safety, reliability and performance of the LEEP System 1000® Workstation only if the installation, recalibrations and repairs are performed by personnel authorized by CooperSurgical and if it is used in compliance with given instructions in an area that meets all the applicable IEC requirements.

Section 16 Warranty

CooperSurgical, Inc., warrants that the LEEP System 1000 Workstation (the “Product”) will be free from defects in materials and workmanship for a period of one (1) year from the original date of purchase.

If the product should become inoperable due to a defect in material or workmanship during this one-year warranty period, CooperSurgical will, at its option, repair or replace the product. This limited warranty does not include replacement or service to repair damage resulting from improper installation, external electrical fault, accident, disaster, use for a purpose other than that for which originally designed or indicated in this manual, negligence, modification, service or repair by personnel not authorized by CooperSurgical or normal wear and tear, and also does not apply to disposable or single- or limited-use items or components. The sole and exclusive remedy under this limited warranty shall be repair or replacement as provided herein.

The foregoing limited warranty states the sole warranty made by CooperSurgical with respect to the product and all parts thereof, and is in lieu of any other warranty by CooperSurgical with respect to product. COOPERSURGICAL NEITHER MAKES NOR GRANTS ANY OTHER WARRANTY, EITHER EXPRESS OR IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL COOPERSURGICAL BE LIABLE FOR ANY DAMAGES ARISING OUT OF THE LOSS OF USE OF THE PRODUCT, OR ANY OTHER INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER OR NOT COOPERSURGICAL HAS ADVANCE KNOWLEDGE OF THE POSSIBILITY OF SAME.

No person, agent, distributor, dealer or company is authorized to change or modify the terms of this limited warranty.

Section 17 Service and Repair

LEEP System 1000® Electrosurgical Generator: There are no customer-serviceable parts

Smoke Evacuator: There are no user-serviceable parts (excluding intermittent filter changes)

LEEP Cart: There are no user-serviceable parts

Only CooperSurgical, Inc., is authorized to service or repair the Electrosurgical Generator or Smoke Evacuator. If repair is attempted outside the factory, the warranty is considered void. CooperSurgical is not responsible for any injury resulting from repairs made by other individuals or organizations not certified by CooperSurgical, Inc. If a repair is needed, the equipment *must be sanitized* before it is returned to the factory and carefully packaged in a protective carton.

On the note inserted in the box, please provide the following information:

- Customer and contact information, on the Repair Authorization Form (download from CooperSurgical Web site) or on company letterhead
- Nature of the problem
- Description of the item returned
- Serial Number (if applicable)

All shipments must be prepaid. COD packages will not be accepted. Return carton to:

CooperSurgical, Inc.

Attention: Repair Department
 95 Corporate Drive • Trumbull, CT 06611 USA
 Phone: (203) 601-5200 • (800) 444-8456
 Fax: (203) 601-4743

Section 18 LEEP System 1000 Electrosurgical Generator EMC Compliance Information

- MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The CooperSurgical Electrosurgical Generator is intended for use in the electromagnetic environment specified below. The customer or the end user of the Electrosurgical Generator should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	CooperSurgical Electrosurgical Generators use RF energy only for their internal function. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	CooperSurgical Electrosurgical Generators are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations / Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Electrosurgical Generator is intended for use in the electromagnetic environment specified below. The customer or the end user of the Electrosurgical Generator should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environmental - guidance
Electromagnetic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T^* (>95% dip in U_T^*) for 0.5 cycle 40% U_T^* (60% dip in U_T^*) for 5 cycles 70% U_T^* (30% dip in U_T^*) for 25 cycles <5% U_T^* (>95% dip in U_T^*) for 5 seconds	<5% U_T^* (>95% dip in U_T^*) for 0.5 cycle 40% U_T^* (60% dip in U_T^*) for 5 cycles 70% U_T^* (30% dip in U_T^*) for 25 cycles <5% U_T^* (>95% dip in U_T^*) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Electrosurgical Generator requires continued operation during power mains interruptions, it is recommended that the Electrosurgical Generator be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* U_T is the AC mains voltage prior to application of the test level.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environmental – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Electrosurgical Generator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_I} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_I} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_I} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LEEP System 1000[®] Workstation is used exceeds the applicable RF compliance level above, the Electrosurgical Generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the LEEP System 1000 Workstation.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distance between portable and mobile RF communications equipment and the LEEP System 1000® Electrosurgical Generator.

The Electrosurgical Generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Electrosurgical Generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Electrosurgical Generator as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <u>Watts</u>	Separation distance according to frequency of transmitter		
	<u>Meters</u>		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1167	0.1167	0.2333
0.1	0.3689	0.3689	0.7379
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3789
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

For the Electrosurgical Generator,
 $V_1 = 3$ Vrms
 $E_1 = 3$ V/m

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Section 19 Specifications

19.1 Electrical – LEEP System 1000® Workstation

- INPUT VOLTAGE 190–250 VAC, 50 Hz
- CURRENT 3.0 amps max.
- FUSES 8.0 amps, 250 V, T type
- WORKING FREQUENCY 450 kHz
- OUTPUT POWER 100 watts RMS (500 ohm load)
- DUTY CYCLE Intermittent 10/30 seconds
- LOW FREQUENCY LEAKAGE Less than 50 micro-amps
- POWER CORD Various
- WEIGHT Approximately 98 pounds (44 kg)

19.2 General Specifications – LEEP System 1000® Electrosurgical Generator

ENVIRONMENTAL CONDITIONS

(Usage)

ENVIRONMENTAL TEMPERATURE ... Between 50 °F and 113 °F (10 °C and 45 °C)
 RELATIVE HUMIDITY Between 30% and 75%

(Shipping and Storage)

ENVIRONMENTAL TEMPERATURE ... Between -40 °F and 158 °F (-40 °C and 70 °C)
 RELATIVE HUMIDITY Between 10% and 100%

PHYSICAL - Electrosurgical Generator

DIMENSIONS 13" wide x 14" deep x 7¼" high (330 mm x 356 mm x 184 mm)
 WEIGHT 30 pounds (14 kg)

ELECTROSURGICAL OUTPUT

RF OUTPUT FREQUENCY 450 kHz

RF OUTPUT POWER:		Volts p-p Max (open circuit)	Duty Cycle	Crest Factor
CUT	0-100 watts RMS*	830	---	1.4
BLEND	0-100 watts RMS*	1200	60%	2.0
COAG	0-80 watts RMS*	3800	10%	5.5

*Stable to >800 ohms (calibration at 500 ohms)

RF ISOLATION Less than 150 milliamps at 200 ohms

CLASSIFICATION I-Type BF

OUTPUT CIRCUIT Floating output. Protected against the effects of the defibrillator.

WORKING MODE Discontinuous maximum duty cycle: 10/30 sec.

COOLING Convection cooling without fan

CONTROL Foot pedal operated (pneumatic) with audible signals
 and mode lights

AUDIBLE SIGNALS AND LIGHTS FOR OPERATION AND ALARM:

MAIN Green light

ALARM, PATIENT PLATE CONTINUITY . Low pitch intermittent audible alarm - red light 5

ALARM, OUTPUT POWER Higher pitch intermittent audible alarm

CUT AND BLEND MODES Low pitch audible signal - yellow light 7

COAGULATION MODE High pitch audible signal - blue light 6

NOTE: Specifications subject to change.

OUTPUT POWER AT 500 OHMS:

Pure Cut

100 W RMS (open circuit 830 Vp-p; crest factor 1.43)

Waveform: sinusoidal at 450 kHz

Blend Cut

100 W RMS (open circuit 1200 Vp-p; crest factor 2)

Waveform: sinusoidal at 450 kHz

Duty cycle: 60%





Coagulation

80 W RMS (open circuit 3800 Vp-p; crest factor 5.5)

Waveform: sinusoidal at 450 kHz

Duty cycle: 10%

Output Power Diagrams according to PAR.6.8.3 (IEC 601-2-2)
 Tolerance: 20% according to PAR.50.2 (IEC 601-2-2)

19.3 Smoke Evacuator

INPUT VOLTAGE	93 – 120 VAC, 50 Hz
CURRENT	3.15 amps
FUSES	6.25 amps, 250V, T-Type
DIMENSIONS	8½" wide x 9" deep x 22" high (216 mm x 229 mm x 559 mm)
WEIGHT	24 pounds (11 kg)
AIR FLOW	35 CFM minimum at maximum setting
FILTER	Pre-Filter, charcoal filter for odor removal and hydrophobic ULPA filter rated at 99.999% efficiency level for 0.014 micron particles*.

*Based on Standard CNC Evaluations

19.4 LEEP Cart

DIMENSIONS	14" wide x 23" deep x 37" high (356 mm x 584 mm x 940 mm)
WEIGHT	44 pounds (20 kg)

Section 20 Explanation of Symbols

Reorder Number

This side up

Serial Number

Caution, risk of electric shock

ATTENTION:
See instructions for use.

Fragile

Medical equipment WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARD ONLY IN ACCORDANCE WITH UL60601-1 AND CAN/CSA C22.2 No.601.1

Keep Dry

DEFIBRILLATOR PROOF SYMBOL = indicates that the device will not be damaged if the defibrillator is active.

TYPE BF EQUIPMENT = Type of equipment is a Type B equipment with an F Type applied part.

Type B equipment is a piece of equipment providing a particular degree of protection against electric shock, particularly regarding:

- Allowable leakage current
- Reliability; of the protective earth connection.

Type F = Isolated applied part. Applied part isolated from all other parts of the equipment to such a degree that the patient leakage current allowable in single fault condition is not exceeded when a voltage equals to 1.1 times the highest rated mains voltage is applied between the applied part and earth.

APPLIED PART SYMBOL = Electrode is isolated from other parts of the equipment.

Rx ONLY

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

LEEP System 1000® is a registered trademark of CooperSurgical, Inc.
Fischer Cone Biopsy Excisor™ is a trademark of CooperSurgical, Inc.
CooperSurgical is a registered trademark of CooperSurgical, Inc.

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CooperSurgical

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: +1 (800) 243-2974
Fax: +1 (800) 262-0105
www.coopersurgical.com

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© 2013 CooperSurgical, Inc.
LEEP System 1000® é uma marca comercial registada da CooperSurgical, Inc.
Fischer Cone Biopsy Excisor™ é uma marca comercial da CooperSurgical, Inc.
CooperSurgical é uma marca comercial registada da CooperSurgical, Inc.