

ConMed Linvatec Shaver Handpieces



D4200/D9820



D4240/D9824

Proprietary Information

This manual contains information deemed proprietary to Linvatec Corporation. The information contained herein, including all of the designs and related materials, is the sole property of ConMed Linvatec and/or its licensors. ConMed Linvatec and/or its licensors reserve all patent, copyright and other proprietary rights to this document, including all design, manufacturing methodology and reproduction.

This document, and any related materials, is confidential and is protected by copyright laws and shall not be duplicated, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form or by any means, electronic, mechanical, magnetic, manual or otherwise, or disclosed to third parties, in whole or in part, without the prior express written consent of ConMed Linvatec.

ConMed Linvatec reserves the right to revise this publication and to make changes from time to time in the contents hereof without obligation to notify any person of such revision or changes, unless otherwise required by law.

Linvatec, Hall, Advantage, PowerPro, E9000 and Ergo are trademarks or registered trademarks of Linvatec Corporation.

© Linvatec Corporation 2014. All Rights Reserved.

1.0 INTRODUCTION

1.1	Operating Principle	1
1.2	Indications for Use.	1
1.3	Intended Use	1
1.4	Contraindications.	2
1.5	Warnings and Precautions	2
1.5.1	Warnings	2
1.5.2	Precautions	4
1.6	Environmental Directives	6
1.7	Product Photographs and Drawings	6
1.8	Symbol Definitions	7
1.8.1	Product Symbols	7
1.8.2	Warnings and Information Symbols	7

2.0 SYSTEM INSTALLATION AND OPERATION

2.1	Product Description	10
2.1.1	Handpiece Descriptions	10
2.1.2	Shaver Blades, Bur Description	11
2.2	Assembly/Installation Instructions	11
2.2.1	Shaver Blade Installation and Removal	11
2.3	Operating Instructions	12
2.4	Preoperative Functional Test	13

3.0 MAINTENANCE

3.1	Cleaning Information.	14
3.1.1	Warnings, Precautions and Notes	14
3.1.2	Manual Cleaning Instructions	15
3.2	Sterilization Information	16

Table of Contents

Page

3.2.1	Warnings, Precautions and Notes	16
3.2.2	Sterilization Instructions	17
3.3	Troubleshooting	18
3.4	Maintenance Schedule	19

4.0 TECHNICAL SPECIFICATIONS

4.1	Product Technical Specifications.	20
4.1.1	Handpieces	20
4.2	Product Environmental Requirements	22
4.2.1	Environmental Technical Specifications.	22
4.2.2	Electromagnetic Requirements	22
4.3	Accessories	22

5.0 CUSTOMER SERVICE

5.1	Assistance and Repair	23
-----	---------------------------------	----

1.0 INTRODUCTION

It is recommended that personnel study this manual before attempting to operate, clean, or sterilize this or associated equipment. The safe and effective use of this equipment requires the understanding of and compliance with all warnings, precautionary notices, and instructions marked on the product, and included in this manual.



This equipment is designed for use by medical professionals completely familiar with the required techniques and instructions for use of the equipment. Service intervals, as listed in Section “**3.4 Maintenance Schedule**”, are required to keep the equipment at its optimum operating performance.

1.1 Operating Principle

The shaver handpiece is electrically powered by the surgical console to provide high speed, rotary force to the accessory (shaver blade, shaver bur or attachment) for tissue resection. The shaver handpiece is controlled by operating buttons on the handpiece or using a footswitch. A suction source is connected to the shaver handpiece to draw tissue into the shaver blade and remove tissue from the surgical site.

Consult the associated ConMed Linvatec surgical console instruction manual prior to operating this equipment.

1.2 Indications for Use

The D9820/D9824 Shaver Handpieces, and their accessories perform cutting of soft tissue and bone. The fields of application include Arthroscopic, Foot, Hand, and Orthopedic surgical procedures.

The D4200/D4240 Shaver Handpieces, and their accessories perform cutting of soft tissue and bone. The fields of application include Arthroscopic, Foot, Hand, Orthopedic, and Otolaryngological surgical procedures.


1.3 Intended Use

Same as Indications for Use above.


1.4 Contraindications


None

1.5 Warnings and Precautions

Do not bypass this section. It contains warnings and precautions that must be thoroughly understood before operating any of the equipment. Lack of understanding or adherence to these warnings and precautions may result in injury or even death to the patient. 



The words **WARNING**, **PRECAUTION**, and **NOTE** carry special meanings and they must be read carefully.






WARNING: A warning contains critical information regarding serious adverse reactions and potential safety hazards that can occur in proper use or misuse of the equipment. Failure to observe the information or procedures presented in a Warning may result in injury or other serious adverse reactions to the patient and/or surgical staff. 

PRECAUTION: A precaution contains instructions for any special care to be exercised by the practitioner for the safe and effective use of the equipment. Failure to observe the information or procedures presented in a Precaution may result in damage to the equipment. 

NOTE: A note is added to provide additional focused information. This information has no critical effect on the patient or equipment.

1.5.1 Warnings

1. Eye protection is recommended when operating equipment. Eye injury may result. 
2. It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of the equipment and its associated accessories.
3. Do not use equipment if, upon receipt, package is opened, damaged, or shows any signs of tampering. 
4. Do not use equipment in the presence of flammable anesthetics, gases, disinfecting agents, cleaning solutions, or any material susceptible to ignition due to electrical sparking.

5. Do not use sterile equipment beyond the expiration date listed on the label. Sterility of the product cannot be assured beyond the expiration date.
6. Do not excessively bend or kink the handpiece cord. Always inspect cords for signs of excessive wear or damage, or bent, broken or missing pins within the connector(s). If any wear or damage is found, discontinue use and replace handpiece immediately. Using a damaged power cord could possibly cause injury.
7. Handpieces are supplied non-sterile. Clean and sterilize prior to each use. 
8. Do not contact the moving parts on the handpieces. Injury to the operator may occur.
9. Continually check handpiece for overheating. If overheating is sensed, immediately discontinue use and return equipment for service. Overheating of the blade or bur may cause damage to the blade or bur and may cause burns or thermal necrosis.
10. While handpiece is not in use do not place on patient/surgical drapes. Place handpiece on mayo stand.
11. Do not immerse the equipment in fluids. Immersion may render the device inoperable. 
12. Failure to follow the specified service interval could result in reduced instrument performance or overheating of the handpiece. Overheating can lead to possible burn injury to the patient or medical personnel. Rotation of handpiece usage per day will assist with proper performance. (Refer to section “**3.4 Maintenance Schedule**”). 
13. Do not attach, insert or remove accessories or attachments while the handpiece is operating. Injury to the operator and/or damage to the equipment may occur. Place the handpiece safety mechanism in the “safe” position prior to installation or removal of items.
14. Avoid contact with cutting tip of blade or bur when locking into handpiece. Tips are sharp and may cause injury.
15. After use blades, burs and tubing sets may be a potential biohazard and should be handled and disposed in accordance with acceptable medical practice and applicable local and national requirements. 
16. Do not use burs for plunge cutting. Injury or damage may occur.
17. Disposable blades and burs are supplied sterile and are for single-use only. Do not resterilize or reuse. The ability to effectively clean and re-sterilize these single use devices has not been established and subsequent re-use may adversely affect the performance, safety and/or sterility of the devices. 
18. Do not use shaver blade or bur if they show any signs of damage.

1.5.2 Precautions



1. United States Federal law restricts sale of this device to or on the order of a physician.
2. This device should only be used in compliance with its intended use.
3. Handle all equipment carefully. If any equipment is dropped or damaged in any way, return it immediately for service.
4. Use only associated ConMed Linvatec approved equipment and accessories. Using unapproved accessories may result in improper operation, and may result in non-compliance to medical standards.
5. The warranty becomes void and the manufacturer is not liable for direct or resulting damage if:
 - The device or the accessories are improperly used, prepared or maintained;
 - The instructions in the manual are not adhered to;
 - Non-authorized persons perform repairs, adjustments or alterations to the device or accessories.
6. There are no user-serviceable parts inside. No modification of this equipment is allowed.
7. Prior to each use, perform the following:
 - Ensure all accessories are correctly and completely attached. (Refer to section **“2.2 Assembly/Installation Instructions”**).
 - Perform the required Preoperative Functional Tests for the equipment and accessories. (Refer to section **“2.4 Preoperative Functional Test”**).
8. Clean and sterilize all the equipment and associated accessories according to instructions for use. (Refer to section **“3.1 Cleaning Information”** and section **“3.2 Sterilization Information”**).
9. Handpieces are factory sealed. Do not disassemble or lubricate, as this may void the warranty.
10. Always inspect for bent, dull or damaged blades or burs before each use. Do not attempt to straighten or sharpen. Do not use if damaged.
11. After each use, thoroughly clean the handpiece, attachments and accessories. (Refer to Section **“3.0 MAINTENANCE”**).
12. Do not stall handpieces, damage can occur.
13. Do not handle the handpiece by the cord. Do not pull on the cord to remove it from the surgical console. Damage can occur.



14. Running the shaver blade or bur without fluid flow (dry) may cause damage to the handpiece or may cause the bur hubs to melt due to excessive heat.
15. Running standard shaver blades above 6,000 rpm may generate excessive particulate. Only Ultra Series shaver blades are intended to run 12,000 rpm maximum speed.
16. Direct contact of the rotating cutting edge of shaver blades and burs with metallic surfaces and/or other hard surfaces such as arthroscopes, cannulas, or other instruments can cause damage to the devices. If contact does occur, shaver blades can break, seize, or shed metal particles. Shaver blade or bur should be examined for damage and replaced if necessary.
17. Do not apply excessive loading on the shaver blade or bur. Cutting performance is not increased with force. Excessive force or using shaver blades or burs as a lever can cause damage to the device including permanent deformation, shedding of metal (wear), motor seizure, and overheating.
18. If shaver blade clogs, use ConMed Linvatec Declogger (E9315) to declog shaver blade. Shaver blades or burs may also be declogged by removing the inner assembly of the shaver blade from the outer assembly and clearing any tissue from the window or aspiration hole. DO NOT disassemble Prebent Shaver Blades with catalog number ending in "M" or "MT" (i.e., C9299M and 42ULT-ZZ-MT).

1.6 Environmental Directives

WEEE Directive [2002/96/EC] on Waste Electrical and Electronic Equipment. This statement only applies to European countries with regard to the Waste Electrical and Electric Equipment (WEEE) European Directive.







The WEEE symbol on the product or its packaging indicates that this product must not be disposed of with other waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of Waste Electrical and Electronic Equipment. The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your medical equipment at the end of its useful life for recycling, please contact ConMed Linvatec.

1.7 Product Photographs and Drawings













The pictures in this manual are for reference only. Items shown may not represent the actual product. However, procedural steps are identical, unless otherwise specified. When necessary, the actual pictures will be represented.





























1.8 Symbol Definitions










1.8.1 Product Symbols

	Handpiece direction selection button (forward, reverse, oscillate). Will also convert to Speed Increase button (see step 4, “ 2.1.1 Handpiece Descriptions ” for more information).		Indicates blade collet unlock position
	Handpiece ON/OFF button. Will also convert to Speed Decrease button (see step 4, “ 2.1.1 Handpiece Descriptions ” for more information).		Suction ON/OFF

1.8.2 Warnings and Information Symbols

	Catalog Number		Serial Number
	Manufacturer		Date of Manufacture
	Consult Instructions for Use		Refer to Instruction Manual/Booklet (for critical safety instruction)
	Caution		DEHP Symbol
	Authorized Representative in the European Community		CE Mark of Conformity
	Prescription Only: Federal Law restricts this device to sale by or on the order of a physician		No User Service Recommended. Refer servicing to qualified ConMed Linvatec service personnel


	Non Sterile		Sterile
	Sterile - Sterilized Using EO		Sterile - Sterilized Using Irradiation
	Do Not Steam Sterilize		Do Not Sterilize
	Do Not Resterilize		Do Not Reuse (for Single Use Only)
	Do Not Use Oil		Do Not Use for Plunge Cutting
	Eye Protection Required		Biohazard Risk
	Do Not Immerse		Quantity
	Type B Applied Part		Type BF Applied Part
	UL Classification Mark		UL Recognized Components
	Rating Fuse		Fuse Location
	Alternating Current		Protective Earth Ground
	Equipotentiality (Equipment Potential)		Non-Ionizing Electromagnetic Radiation (RF Symbol)
	Temperature Limitation		Humidity Limitation
	Atmospheric Pressure Limitation		Use by Date

	Fragile		This Side Up
	Do Not Use if Package is Damaged		Keep Dry
	Warning: Corrosive Substance		Warning: Electrical Hazard/ High Voltage
	Waste Electrical and Electronic Equipment (WEEE) Symbol. Regarding European Union end-of-life of product, indicating separate collection for electrical and electronic equipment		
 	Recycle. Batteries contain materials which must be recycled or disposed of properly. The disposal of batteries as municipal waste is prohibited. Dispose or recycle in accordance with your local, state and governmental regulations. In the U.S. call 1-800-237-0169, or outside the U.S. contact your local ConMed Linvatec representative for additional information on battery disposal or recycling.		

2.0 SYSTEM INSTALLATION AND OPERATION

2.1 Product Description

2.1.1 Handpiece Descriptions




All these handpieces are equipped with sensors that automatically select the proper speed range for the attached blade or bur. For operating speeds, please refer to Section “**4.0 TECHNICAL SPECIFICATIONS**”. The 2-Button Shaver Handpieces (D4240, D9824) are also equipped with a **DIRECTION SELECT** button to change between forward, reverse, and oscillate, and an **ON/OFF** button to control the handpiece from within the sterile field.

1. **Blade Locking Collet** - To insert a blade or bur, align the key on the blade hub with one of the slots in the collet. Completely push the hub into the collet to lock it into place. Pull on the blade to ensure it is locked. To remove the blade, rotate and hold the collet in the “UNLOCK” direction and remove the blade.
2. **Suction Port Valve** - Manually controls suction flow. Move towards the front of the handpiece to open, or increase flow. Move towards the back of the handpiece to restrict or stop flow. To attain maximum performance while resecting tissue or bone, the suction port valve must be in the fully open position.
3. **Suction Port** - Connects to the suction line for aspiration.
4. **DIRECTION SELECT Button** - Press and release to select the direction of blade rotation between forward, reverse and oscillate.
5. **ON/OFF Button** - Turns handpiece on or off.



2.1.2 Shaver Blades, Bur Description



1. Cutting tip.
2. Shaver blade/bur hub- connects to the shaver handpiece.
3. Alignment key.

2.2 Assembly/Installation Instructions

2.2.1 Shaver Blade Installation and Removal

To insert a shaver blade or bur:

1. Align the key on the blade hub with one of the slots in the collet.
2. Completely push the hub into the collet to lock it into place.
3. Pull on the blade to ensure it is locked.



To remove the blade or bur:

1. Rotate and hold the collet in the “UNLOCK” direction.
2. Remove the blade.



2.3 Operating Instructions

To control the suction flow rate, move the suction port valve towards the front of the handpiece to increase flow and move the suction port valve towards the back of the handpiece to restrict or stop flow. To attain maximum performance while resecting tissue or bone, the suction port valve must be in the fully open position.

These shaver handpieces can be operated using a footswitch or the control buttons on the handpiece. For the footswitch operating instructions, please refer to the surgical console instruction manual.

For shaver handpiece operation using the buttons on the handpiece (D4240 and D9824):

Press and release to turn the shaver handpiece on and off.

Controlling window indexing using the buttons: The following operation applies only when the handpiece is NOT running. When the handpiece ON/OFF button is pressed and held, the console beeps once and the blade window starts to rotate (index) slowly at 6 rpm. When the window position of the blade is found (generally window completely open or window completely shut), release the ON/OFF button. The handpiece will return the blade window to this location when it is stopped.

Controlling speed using the buttons: The following operation applies only when the handpiece is NOT running. When the handpiece DIRECTION SELECT button is held for 1.5 seconds, the console beeps twice. Then the handpiece DIRECTION SELECT button changes to the Increase Speed button. The ON/OFF button changes to the Decrease Speed button. This enables the surgeon to change the handpiece operating speed from the sterile field. If no button is pressed for 1.5 seconds, the console beeps twice and returns to normal button functionality.

Turbo Mode Using Ultra Series Shaver Blades and D4000/D4000A Drive Console with console software revision V 2.0 or higher: The following applies to the D9824, and D4240 Shaver handpieces and only when it is actively running in oscillate mode. When the ON/OFF button is pressed and held for approximately 1 second, the handpiece will run in the preset forward mode for as long as the ON/OFF button is held down. Release the ON/OFF button and the handpiece will resume running in the oscillate mode. When the DIRECTION SELECT button is pressed and held for approximately 1 second, the handpiece will run in the preset reverse mode for as long as the DIRECTION SELECT button is held down. Release the DIRECTION SELECT button and the handpiece will resume running in the oscillate mode.

2.4 Preoperative Functional Test

Prior to each use, perform the following preoperative functional test:

1. Insert a bur or blade into the handpiece. Gently pull on the cutting accessory to ensure it is properly seated.
2. Verify that the console properly recognizes the handpiece and cutting accessory.
3. Press the direction select button or footswitch to ensure that it functions properly.
4. Run the handpiece for less than 5 seconds at a speed less than 1500 rpm to observe any abnormal noises, vibrations or heat rise.
5. If any operating difficulties occur, return the handpiece for service.

3.0 MAINTENANCE

3.1 Cleaning Information

3.1.1 Warnings, Precautions and Notes



1. Follow universal precautions for protective apparel when handling and cleaning contaminated instruments.
2. Clean instruments within 30 minutes after use to minimize the potential of blood and debris drying.
3. Never clean equipment in an ultrasonic cleaner.
4. Always disconnect the handpiece cord prior to cleaning.
5. Always detach accessories from equipment prior to cleaning.
6. Never clean handpieces with bleach, chlorine-based detergents, liquid or chemical disinfectants, or any products containing sodium hydroxide (such as, INSTRU-KLENZ or Buell Cleaner). These products degrade the anodized aluminum coating and may result in reduced handpiece reliability.
7. For aluminum surfaces, a neutral-pH agent should be used. To prevent corrosion, avoid contact with strong alkaline solutions (pH over 10.5) or agents containing iodine or chlorine. Follow universal precautions for protective apparel when handling and cleaning contaminated instruments.
8. Prior to using a washer/sanitizer, consult product labeling on all washer/sanitizer cleaning solutions for compatibility with aluminum. Use of a washer/sanitizer may decrease the life expectancy of the handpiece.

Refer also to section **“1.5 Warnings and Precautions”** and section.

3.1.2 Manual Cleaning Instructions

1. Thoroughly scrub the handpiece and attachments with a clean, soft brush dampened, not saturated, with a mild, pH-balanced detergent. Remove all traces of blood, debris, and stains.
2. Place the suction control valve in the fully open position. Clean the suction tube thoroughly with a cleaning brush until all debris is removed.
3. Using a low-velocity spray of water, spray into the collet. Ensure the collet operates smoothly and that no debris is binding the internal mechanism of the collet. Clean until all debris is removed.
4. To clean the cannulated section of the handpiece and attachment, feed the wire end of a cleaning brush through the cannulation of the handpiece or attachment. Pull the brush completely through and repeat until all debris is removed.
5. Manipulate all moving parts of the handpiece to ensure all debris is removed. If not, clean again until all debris is removed.
6. Keeping the nose of the handpiece pointed downward, rinse under running water to remove all traces of soap. Rinse all attachments in the same manner.
7. Flush the surfaces free of tap water with distilled water to prevent metal discoloration.
8. Gently shake the equipment free of water and wipe the surfaces with a clean, lint-free towel.
9. Visually inspect handpiece and accessories under good light conditions to check for visible soil and/or corrosion.
10. Perform functional checks according to section **“2.4 Preoperative Functional Test”**. Check mating accessories closely for proper assembly.
11. Remove and repair any damaged equipment. Repeat Manual Cleaning Instructions as necessary.

3.2 Sterilization Information

3.2.1 Warnings, Precautions and Notes



1. The use of disinfecting solutions for an exterior instrument wipe will not sterilize the equipment and is not recommended.
2. Do not sterilize equipment or accessories using Ethylene Oxide (EtO).
3. Do not sterilize equipment or accessories using a STERIS system or by comparable sterilization methods.
4. Do not sterilize equipment or accessories in cold sterilants like CIDEX.
5. Always detach accessories from equipment prior to sterilization.
6. Do not use equipment while warm. Allow adequate time for cooling prior to use. Cool by exposure to room temperature. Operation of equipment that is not completely cool or dry may decrease performance and/or reliability.
7. When sterilizing equipment and attachments in a system sterilization tray, use the recommended sterilization parameters in its instruction manual.
8. Do not “Peel Pack” handpieces for sterilization. Sterilization in a sealed pouch traps moisture which can cause damage.
9. The suction control valve, if one exists, must be in the fully open position during sterilization.
10. Attachments with collet mechanisms must be sterilized with the collet fully open.
11. An eight (8) minute minimum dry cycle must be run on all equipment and attachments every time the product is sterilized. Failure to use a dry cycle on the products may lead to reduced product performance or premature product failure.

Refer also to section “1.5 Warnings and Precautions”.

3.2.2 Sterilization Instructions

Steam sterilization is safe and effective and has no contraindications for sterilizing this equipment.

1. Individually wrap handpiece and accessories.
2. Follow the recommended minimum sterilization exposure times listed below.
3. Use the recommended tray parameters in its instruction manual when sterilizing handpieces and attachments in a system sterilization tray.

Handpieces and attachments may be processed in a pre-vacuum steam sterilizer (Steam Pre-vacuum) or in a gravity (downward) displacement sterilizer (Steam Gravity).

Minimum recommended sterilization exposure times are as follows:

Table 1: Sterilization Parameters

Method	Cycle	Minimum Temperature	Minimum Exposure	Minimum Dry Cycle
Steam (wrapped)	Pre-Vacuum	270°F (132°C)	4 minutes	8 minutes

3.3 Troubleshooting

Table 2: Troubleshooting Guide

Symptom	Possible Cause	Corrective Action
Handpiece does not operate.	<ul style="list-style-type: none">• Handpiece cord not connected securely.• Handpiece faulty.	<ul style="list-style-type: none">• Securely connect handpiece cord to controller.• Return for service.
Error Messages		
<p>“MOTOR OVER SPEED”</p> <p>“HP SERVICE REQUIRED”</p> <p>“CANNOT ID HANDPIECE”</p> messages display.	<ul style="list-style-type: none">• Handpiece and/or System faulty.	<ul style="list-style-type: none">• Replace handpiece. If symptom is corrected, return handpiece for service.• Else, return system for service.
CONNECT HANDPIECE message displays when the handpiece is already plugged into the console.	<ul style="list-style-type: none">• Handpiece not connected securely.• Handpiece and/or System faulty.	<ul style="list-style-type: none">• Securely connect handpiece to controller.• Replace handpiece. If symptom is corrected, return handpiece for service.• Else, return system for service.

3.4 Maintenance Schedule

Regular and proper maintenance of your equipment is the best way to protect your investment. It is essential that you have your equipment serviced as scheduled in order to retain its optimum performance and reliability, which will reward you with safer, less problematic product performance over time.

The equipment is not field repairable. Your ConMed Linvatec authorized service department is the most knowledgeable about this equipment and its accessories and will provide competent and efficient services. Service at ConMed Linvatec at the recommended service interval is mandatory to keep your product warranties in effect. Any services and/or repairs done by any unauthorized repair facility may result in reduced performance of the equipment or equipment failure. (Refer to section **“5.0 CUSTOMER SERVICE”**).

The Shaver Handpieces shall be returned every 12 months for servicing.

4.0 TECHNICAL SPECIFICATIONS

Medical electrical equipment complies with and was tested with respect to electric shock, fire, electromagnetic compatibility, mechanical and other specified hazards only, in accordance with UL60601-1, CAN/CSA C22.2 No. 601.1-M90, IEC60601-1:1988 +A1:1991 +A2:1995, ES60601-1:2005 +A1: 2009 +A2: 2010, CAN/CSA C22.2 No.60601-1-1-08 and IEC60601-1:2005 +C1:2006 +C2:2009.

Tested to IEC60601-1-2:2007 and Part 15 of the FCC Rules as follows: 1) The system may not cause harmful interference: 2) The system will accept interference, including interference that might cause undesired operation. If interference occurs, separate the instruments. For more information, contact customer service.

4.1 Product Technical Specifications

4.1.1 Handpieces

Advantage No-button (D9820) and 2-Button (D9824) Shaver Handpieces	
Classification:	Class I (intended for connection to Class I surgical console). Type rating of applied part is dependent on surgical console.
Protection Ingress of Fluids:	IPX0 (ordinary)
Mode of Operation:	Intermittent Loading
Operating Speed:*	500 - 10,000 rpm, Forward/Reverse 500 - 4,000 cpm, Single-and Multi-turn Oscillate
Torque (Nominal):	22 in. oz. (15.5 N•cm)
Length:	6.1 in (15.5 cm)
Diameter:	Approx. 1.125 in. tear drop (2.9 cm)
Weight:	15.8 oz. (450g)
Cord length:	10 ft. (3m) Non-detachable
Material:	Anodized Aluminum and Stainless Steel
Duty Cycle (once daily) with fluid flow:	30 seconds ON, 15 seconds OFF (3x) 3 minutes ON, 30 seconds OFF (4x)







*Operating speeds may vary depending on console, console software revision, and blade/bur selection.

Ergo No-button (D4200) and 2-Button (D4240) Shaver Handpieces	
Classification:	Class I (intended for connection to Class I surgical console). Type rating of applied part is dependent on surgical console.
Protection Ingress of Fluids:	IPX0 (ordinary)
Mode of Operation:	Intermittent Loading
Operating Speed:*	500 - 12,000 rpm, Forward/Reverse 500 - 4,000 cpm, Single-and Multi-turn Oscillate
Torque (Nominal):	35 in. oz. (24.7 N•cm)
Length:	6.1 in (15.5 cm)
Diameter:	Approx. 1.25 in. tear drop (3.2 cm)
Weight:	19.3 oz. (548g)
Cord length:	10 ft. (3m) Non-detachable
Material:	Anodized Aluminum and Stainless Steel
Duty Cycle (once daily) with fluid flow:	30 seconds ON, 15 seconds OFF (3x) 3 minutes ON, 30 seconds OFF (4x)

*Operating speeds may vary depending on console, console software revision, and blade/bur selection.

4.2 Product Environmental Requirements

4.2.1 Environmental Technical Specifications

Environmental Conditions	Operating	Storage and Transport
Temperature:		
Relative Humidity:	 30% — 75% Non-Condensing	 10% — 93% Non-Condensing
Atmospheric Pressure:	 700 hPa — 1060 hPa	 500 hPa — 1060 hPa

4.2.2 Electromagnetic Requirements

Refer to appropriate drive console instruction manual for electromagnetic requirements.

4.3 Accessories

Only use ConMed Linvatec shaver blades and burs.

Description

Ref Attachments

D9900 Thorne Chuck 5/32 in Attachment (for use with D4200 and D4240 only)

5.0 CUSTOMER SERVICE

5.1 Assistance and Repair

If you need technical assistance regarding the use or application of this product, or you encounter a problem that requires servicing or repair, contact ConMed Linvatec Customer Service at 1-800-237-0169 or your ConMed Linvatec Sales Representative. Outside the U.S. contact your local ConMed Linvatec Representative.

Report any events involving injuries or malfunctions to the ConMed Linvatec Regulatory Product Support.

Returning products for any reason requires an authorized Service Request (S.R.) number prominently displayed on the box and included on all paperwork. Refer to this number if making inquiries about the repair status. Please call ConMed Linvatec Customer Service and provide the following information to obtain an S.R. number prior to returning any product for repair:

- Product Number
- Serial Number/Lot Number
- Reason for Return
- Original Invoice Number
- Date of Purchase
- Detailed description of the problem

ConMed Linvatec
Attn.: Customer Service Dept.
11311 Concept Boulevard
Largo, Florida 33773-4908 USA

Customer Service

(within U.S.)	Phone:	1-800-237-0169
----------------------	---------------	-----------------------

	FAX:	(727)-399-5256
--	-------------	-----------------------

(outside U.S.)	Phone:	+1 (727)-392-6464
-----------------------	---------------	--------------------------

	FAX:	+1 (727)-397-4540
--	-------------	--------------------------

ConMed Linvatec Regulatory Product Support

(within U.S.)	Phone:	1-800-325-5900
----------------------	---------------	-----------------------

(outside U.S.)	Phone:	+1 (727)-392-6464
-----------------------	---------------	--------------------------



525 French Road
Utica, NY 13502-5994 USA
(727) 392-6464

EC REP

MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany

Customer Service: 1-800-237-0169

FAX: (727) 399-5256

International FAX: +1 (727) 397-4540

email: customer_service@conmed.com

www.conmed.com



Rx ONLY



All rights reserved.

W41-260-004 Rev. AC

10/2014