

# TROUBLESHOOTING REPAIR MANUAL:

Eppendorf 5417R Refrigerated Centrifuge

this guide brought to you by



***Begin Your Repair Now***





Fixing electronic or electromechanical devices on your own, without the proper training, knowledge, or tools, can pose significant risks not just to the device itself but also to your personal safety. These devices often operate under high voltage or store energy in capacitors that can cause severe electrical shocks. Additionally, tampering with internal components can lead to short circuits, fires, or even explosions if batteries are improperly handled. The complexity of these systems means that without a thorough understanding, it's easy to misdiagnose problems or create new issues that could render the device inoperable. Moreover, such unauthorized repairs can void warranties and violate safety regulations, leading to legal and financial repercussions. For these reasons, it's crucial to assess the risks and consider professional services for repairs, especially for sophisticated devices where precision and safety are paramount.

This educational document is provided for informational purposes only and is not intended to serve as a comprehensive guide for conducting repairs or other technical procedures. It does not suggest that repairs should be attempted without seeking further instruction or professional advice. The information contained herein is meant to supplement, not replace, the knowledge and expertise of qualified professionals in the field. Users are advised to exercise caution and consult with experts before undertaking any actions based on the information provided in this document.

The information presented in this document is derived from sources published by the original manufacturer and may not be entirely accurate or up to date. While we have made efforts to ensure the reliability of this information, we cannot guarantee its correctness or completeness. The details provided herein should be considered as a general guide and not as a definitive source of information. Users are encouraged to verify the currentness and accuracy of the data with the respective manufacturers or other authoritative sources. We accept no liability for any errors or omissions, or for any decisions made based on the information provided in this document.



The "Acme Revival" mark and brand are registered trademarks with the United States Patent and Trademark Office (USPTO). All rights associated with this trademark, including but not limited to the use of the mark in commerce, are exclusively owned by the registrant. Unauthorized use of the "Acme Revival" trademark, or any mark that is confusingly similar, is strictly prohibited and may result in legal action for trademark infringement.

This document may contain references to trademarks belonging to other entities. The inclusion of such trademarks is purely for informational purposes and aimed at fostering knowledge and education among our users. It is not our intention to create any confusion or mislead anyone regarding the ownership of these trademarks. Our goal is to provide valuable, educational content while respecting the rights of trademark owners..

## 1-1. Troubleshooting:

**No display:** Caused by no main connection or power failure. Solution involves checking the power supply cable and the main power fuse on the device and in the laboratory.

**CLEAn ro:** Triggered after 200 runs, indicating the rotor needs cleaning and lubrication of the axle.

**Error 3:** No rotor detected. Solution is to insert the rotor.

**Error 23:** Indicates an error in the drive or rotor recognition. The device should be switched off and then switched back on.

**Lid cannot be opened (No display):** Caused by a power failure. The solution is to bring the rotor to a standstill and activate the emergency lid release.

**Lid not closed completely (Display: LID):** Lid latch not engaged properly. Solution is to press down the lid again.

**Centrifuge shakes during acceleration and switches off (Display: InbAL):** Caused by the rotor not being loaded symmetrically or the rotor not being fastened properly. Check rotor loading and fasten the rotor correctly.

**Flashing temperature display and Warning signal:** Indicates a deviation from the nominal temperature value of more than 3°C or 5°C, suggesting the nominal value set is too low or refrigeration is defective. Solutions include adjusting the nominal speed, ensuring the rotor is recognized correctly, checking for an imbalance sensor defect, ensuring there is a rotor installed, and addressing errors in the rotational speed measuring system.

**Errors 1 through 9-25, including Error 18:** Various causes such as nominal speed too high, power failure during the run, rotor not recognized, imbalance sensor defective, no rotor, errors in rotational speed measuring system, unauthorized opening of lid or defective lid switch, converter overloaded, brake defective, rotor loose, overspeed, drive error, and electronics error. Solutions range from entering new nominal speed, repeating the run with the correct settings, ensuring the rotor is inserted properly, restarting the device after a wait, and tightening the rotor. For drive and electronics errors, contacting SERVICE is recommended.

### LEVEL 1 ERRORS

#### 0 System Error: NMI

**1 System Error: Trap**  
**2 System Error: ESC**  
**3 Memory Allocation Error**  
**LEVEL 2 ERRORS**  
**100 Fan Driver Failure**  
**101 TE Driver Failure**  
**102 Solenoid Driver Failure**  
**103 Block Probe Open?**  
**104 Block Probe Failure**  
**105 Sample Probe Open?**  
**106 Sample Probe Failure**  
**107 Cooling System Error**  
**108 Event Record Lost**  
**LEVEL 3 ERRORS**  
**200 Sample Pre-Freeze**  
**201 Sample Freeze Error**  
**202 Sample Did Not Freeze**  
**203 No Plateau**  
**204 Test Cancelled**  
**205 Calibration Cancelled**  
**206 Test Time-Out Error**  
**207 Recalibration Needed**  
**208 Standards Reversed?**  
**209 Calibration Out of Range**  
**210 Reset Probe Config.**  
**211 Parameter RAM Failed**  
**212 No Parameters in RAM**  
**213 Count Error: Parameter**  
**214 ROM Serial Number Error**  
**215 Baudrate Error**  
**216 Key Input Timeout**  
**217 Lift Sample Probe**  
**218 A/D Init Failure**  
**219 A/D High Filter Error**  
**ERRORS (continued)**  
**220 A/D Low Filter Error**  
**221 A/D Cal Mode Error**  
**222 Bar Code Cancelled**  
**223 Diagnostics Cancelled**  
**224 Low Battery**  
**225 Out of Range**

**Abrupt Loss of Power (No Display):**

Check power availability and fuses.

Suggested actions include restoring power, checking connections, and replacing blown fuses.

**Block Probe Open?:**

Check the block probe by running A/D tests and using an ohmmeter.

Actions include switching the instrument off then on and contacting the product service department if the problem persists.

**Fan Driver Failure:**

Indicates a problem with the fan under the sample handling assembly.

Check the fan wire connection.

**No Plateau Detected...:**

Check the sampler's calibration for proper volume.

Retry the test and verify sample osmolality within range limits.

**Parameter RAM Failed or No Parameters in RAM:**

Switch the instrument off, wait 10 seconds, then turn it on again.

Possible battery replacement required.

**Recalibration Needed:**

Calibrate the instrument.

**Results Not Repeatable (Too Scattered):**

Ensure cleanliness and calibration of the sample cell and sampler.

Replace and calibrate the plunger wire as needed.

**Sample Did Not Freeze...:**

Check sample configuration and block probe bin number.

Run the solenoid test for no impacts.

**Sample Pre-freeze...:**

Ensure the sample chamber is clean and the sample probe bin number is correct.

**Sample Probe Open?:**

Check the sample probe by running A/D tests and using an ohmmeter.

**Standards Reversed? Please Repeat...:**

Recalibrate with correct standards.

**T E Driver Failure:**

Check connections to the cooling assembly with an ohmmeter.

**Test Time-out Error:**

Verify the correct probe bin number and sampler tip positioning.

**Other Error Messages:**

Turn the instrument off and on again.

**Date and Time Not Maintained When Instrument Is Powered Off:**

Contact the product service department for assistance.

**"Pull Cradle Out":**

Adjust the cradle position towards the instrument front.

**"Cooling System Error":**

Use the A/D Test to check sample and block channels.

Replace probes if needed.

**Low Battery:**

Consider replacement of the memory backup battery or processor control board.

**Wrong Model Number/Operation:**

Verify and correct processor board option switch SW1 positions.