

# Vista® basic Infusion Pump

Easy, convenient, cost-effective infusion therapy for ambulatory and home-care settings



Designed to ensure simplicity during every phase of the infusion process — from programming and set installation to delivery of therapy.

## B. Flexible.

**B | BRAUN**  
SHARING EXPERTISE

# Vista® basic Infusion Pump

## Easy to Use and Convenient

- Saves time loading with straight-line gravity tubing
- Quick start-up saves time for continuous infusions—simply enter rate and volume, then press START
- Helpful prompts walk the user through the programming process
- Easy to read large display provides all information at a glance
- Ergonomic design is lightweight yet durable
- Bag hook for pump offers space and cost savings, eliminates separate IV pole

## Supports Safety

- Needle-free tubing — capless Ultrasite® needle-free valve
- Enhances patient safety with pump and set-based free-flow protection

## Flexible

### Continuous Mode

- Enter rate, volume and run to initiate

### Program Mode

- Flexible options allow up to nine steps in the program mode
- Sequential step delivery is ideal for circadian and intermittent therapies

### Ramp/Taper Mode

- Easy-to-program Ramp/Taper Mode uses step-by-step prompts to shorten learning curve
- Enhanced safety feature requires maximum plateau rate be approved by the clinician before initial infusion

## SPECIFICATIONS

Mechanism	Linear peristaltic
Width	5.5 in (14 cm)
Height	9.4 in (24 cm)
Depth	7.8 in (20 cm)
Weight	6.8 lbs (3.1 kg)
Power	100-120 VAC, 50-60 Hz, 13 Watts
Leakage Current	max 30 A
Battery	7.2V, 1.2Ah NiCd (rechargeable)
Battery Operating Time	Approx. 3.5 hours at 125 mL/hr.
Rate Range	0.1 to 800.0 mL/hr in 0.1 increments with Vista pump tubing
Air-in-Line	Adjustable 0.5-3.5 mL/hr.
Volume to be Infused Range	0.1 to 9999.9 mL
Occlusion Pressure	Variable settings: Low: 8 PSI and High: 17 PSI
KVO Rate	1 mL/hr or selected rate, whichever is lower
Flow Rate Accuracy	+/-5% measurement time over 8 hours at a temperature of 71.6 degrees F
Alarms	Air-in Line, Low Battery, Plug In Now, Bag Empty, Downstream Occlusion, Door Open, Pump Stopped, KVO, Unit Malfunction.
Memory	User option to retain last Continuous and Piggyback Mode data upon power up. Permanent data retention and Resume Feature for all Models after the pump is turned off.
Tubing	B. Braun Vista® basic pump administration sets. Needle-free sets available
Modes/Features	Continuous and Piggyback Modes, Program Mode for variable infusions, Ramp/Taper Mode for TPN infusions, and Delay A Mode for postponing infusion start and finish.



## Sharing Expertise™

It's not just a philosophy. It's a promise.

Sharing Expertise includes listening to the wisdom of colleagues, customers and clinicians. It's not just a philosophy. It's our promise to the healthcare community and ourselves to continuously seek knowledge. To expand it. To apply it. To share it. All to help improve and save lives.

This is Sharing Expertise. This is B. Braun.