LASAIR® II Models 310A, 310B, 510A

The LASAIR II from Particle Measuring Systems provides complete features in a compact, convenient particle counter. Its great portability, ease of use, and chemically-resistant, easy to clean cover have made it the de facto particle monitoring standard for pharmaceutical and high tech industries.

BENEFITS

Reduce Defects

- Use as a portable counter to localize source of particles
- Real-time measurement of defect-causing particles
- Patented sample-flow control gives more accurate measurements, eliminating errors associated with manifold sampling

Increase Productivity

- · Download via optional USB port
- Stores sampling recipes, reducing labor and errors
- Calculations and reports for ISO 14644-1, FS-209E, or EC GMP testing
- Comprehensive validation manual makes it easier for pharmaceutical to meet regulatory requirements
- Can operate in cleanroom from office PC via web browser: Set-up, sample, display, print, download data, update software
- Removable battery with optional external charger for continuous mobile use
- Quiet pump improves work environment
- Choice of 9 languages for display and printouts

Cost Effective

- Easy to clean/wipe down; designed to minimize particle traps
- \bullet Rugged, chemical-resistant, ESD-minimizing KYDEX $^{\textcircled{R}}$ casing
- Longer life pump eliminates annual vane replacement
- Diode laser for reduced maintenance



FEATURES

- Compact, convenient, lightweight
- Only 15.5 lb. (7.0 kg) with battery
- Sizing sensitivities from 0.3 25.0 μm
- Meets JIS accuracy standards
- Stores 3,000 samples
- Integrates up to 6 external analog sensors
- Provides range of data communications options
- · Controllable from PC web browser
- Built-in thermal printer
- Large, color LCD for higher visibility, more data displayed
- NiMH battery affords 4 7 hours of use

APPLICATIONS

- Cleanroom monitoring
- Facility certification
- Trend analysis
- Statistical process control
- Troubleshooting
- Manifold compatible
- Portable or dedicated use



Without measurement, there is no control.

Specifications			
LASAIR II	310A	310B	510A
Channal thresholds:	0.2.0.5.1.0.5.0.10.0.25.0.um	0 2 0 5 1 0 2 0 5 0 10 um	0.5, 1.0, 2.0, 5.0, 10.0, 25.0,
Flow rate	$28.3 \text{ I PM} (1 \text{ CFM}) \text{ with } \pm 5\% \text{ acc}$	0.5, 0.5, 1.0, 5.0, 5.0, 10 μm	0.5, 1.0, 2.0, 5.0, 10.0, 25.0 µm
Laser source:	Laser diode		
Calibration:	Materials traceable to US National Institute for Standards and Technology (NIST)		
Maximum concentration @ 5% coincidence loss:	> 375,000/ft ³	>375,000//ft ³	> 425.000/ft ³
Counting efficiency, sizing accuracy:	Manufactured units meet ISO 13323-2 standards		
Zero count level:	Meets JIS standard		
Data storage:	3,000 data sets (includes particle and environmental data, plus location, time, etc.). Cannot be altered		
Communication modes:	Ethernet or RS-232; optional 4-20 mA output board; optional wireless Ethernet; optional USB download		
Controlling software:	Facility Net, Pharmaceutical Net, Microsoft [®] Internet Explorer Rev. 5.0+		
Remote operation:	Use web browser to operate unit, save/load recipes; real-time download to Facility Net, Pharmaceutical Net		
Environmental sensors:	Temp/RH probe, plus integrates 4 external environmental sensors via 4-20 mA inputs		
Display:	1/4 VGA color LCD screen		
Languages:	English, French, German, Italian, Japanese (Kanji), Korean, Mandarin (Classical or Simplified), Spanish		
Printer:	Built-in thermal printer		
Reports:	Sample printouts; cleanroom certification reports for ISO, FS 209, EC GMP, or British Stds; averages of multiple runs		
	print on alarm; sampling configuration		
Key software features:	Recipes; pharm mode; scheduled sampling; historical data filtering, sorting and review; password protection		
External surface:	ESD-minimizing KYDEX chassis, polyester display screen, and zinc-plated steel bottom		
Cleaning materials:	Bleach, formaldehyde, ethyl/isopropyl alcohol, peroxide/quaternary ammonium solutions		
Sampling tubing:	Tubing ID: 3/8 in. Maximum tubing length: 10 ft. (3 M)		
Sample output filtering:	Internally filtered to HEPA standards (>99.97% @ 0.3 µm)		
Power:	85-264 V, 50-60 Hz		
Battery:	NiMH; expected operation 4-7 hrs. External recharge: 2-3 hrs. Internal recharge (@ 24° C): 80% by 4 hrs., 100% by 5-10 hrs.		
Dimensions (h, w, d):	11.5 x 10.7 x 9.0 in (29 x 27 x 23 cm)		
Weight:	12 lb. (5.4 kg) without battery; 15.5 lb. (7.0 kg) with optional battery		
Operating environment:	Temperature: 15-35° C, Humidity: 10-85% non-condensing		

Accessories

Included: Printer paper, Operations Manual, sample probe, probe adapters, sample tubing, zero count filter, power cord, spare fuse. Optional: USB data download, battery, external battery charger, external filter, HHIPA, hand-held/tripod/wall-mountable isokinetic sampling probe, tripod, temperature and relative humidity probe, air velocity probe, differential pressure probe, 4-20 mA data output board, TouchRAM wand and buttons.

LASAIR[®] is a registered trademark of Particle Measuring Systems, Inc. KYDEX[®] is a registered trademark of Kleerdex Company. Microsoft[®] is a registered trademark of Microsoft Corporation. U. S. patents: 6, 167,107. Foreign patents pending. Particle Measuring Systems, Inc. reserves the right to change specifications without notice.

AUTHORIZED REPRESENTATIVE



Registration applies to the Boulder, Colorado facility



Particle Measuring Systems Headquarters 5475 Airport Blvd., Boulder, CO 80301 (303)443-7100 1-800-238-1801 FAX: (303)546-7331 Instrument Service & Support: 1-800-557-6363 Customer Response Center: 1-877-475-3317 Particle Measuring Systems Europe Tel: +44 (0)1684-581000

Particle Measuring Systems Japan Tel: 813-5298-8175

Particle Measuring Systems Singapore Tel: (65) 68460-500

Particle Measuring Systems China Tel: 86-21-53855332-18