

Efficient Compliance

Highly Flexible

Advanced Diagnostics

Time Saving, Faster and Easier

Protects Your Investment



Specifications

9703+ Sampler

<i>Temperature Range</i>	5 to 40°C (40 to 104°F) 0 to 80% relative humidity, non-condensing
<i>Temperature Range of Sample</i>	5 to 40°C (40 to 104°F)
<i>Viscosity Limit</i>	<50 cP
<i>Power Requirement</i>	100/230 VAC, 50/60 Hz
<i>Maximum Power</i>	67 W at 115/230 VAC
<i>Dimensions</i>	343 W X 337 D X 482 H mm (13.5 X 13.3 X 19.0 inches)
<i>Weight</i>	10.7 kg (23.6 lbs)
<i>Sample Bottle Clearance</i>	153 mm (6.02 inches)
<i>Volume Accuracy</i>	> 95%
<i>Flow Rate Accuracy</i>	> 95%
<i>Sample Flow Rate</i>	10 to 100 ml/min; Actual flow rate for the system is determined by the sensor's flow rate with step-by-step instructions and Pass/Fail results.
<i>Tare Volumes for Probes</i>	0.091 ml for 81 L X 1.2 ID mm (3.2 X 0.047 inches) 0.172 ml for 154 L X 1.2 ID mm (6.1 X 0.047 inches) 1.57 ml for 140 L X 6.35 ID mm (5.5 X 0.25 inches) - large bore

PharmSpec 3 Software

<i>Processor</i>	Pentium® 1.0 GHz
<i>Operating System</i>	Windows® XP with Service Pack 3 or later Windows Vista® with Service Pack 2 or later Windows 7 platforms
<i>Memory</i>	512 MB of RAM
<i>Hard Disk Drives</i>	2 GB of available hard-disk space CD drive, DVD drive, USB port or network connection is required for installation
<i>Display</i>	Super VGA video adapter and monitor, 1024 x 768 or higher resolution
<i>Power Supply</i>	Uninterrupted power supply (UPS)
<i>Removable Drive</i>	Removable drive for backing up data, such as a recordable CD or a network server that can back up the computer hard drive Color printer capable of printing at 300 dpi
<i>Printer</i>	USB port to connect to the particle counting sensor
<i>COM Port Software</i>	Internet Explorer 6.0 or higher and Microsoft® Excel®, Word® or Adobe® Acrobat® (to access PharmSpec reports)

Typical System Configurations - other options available

<i>Part Number</i>	<i>Description</i>
2089105-01	9703+ Sampler
2089250-150	HRLD-150 Sensor
700060GP	PharmSpec 3 Software
DOC026.53.80099	Validation Documents - English (other languages available)
2089007-025	25 ml Calibration Flow Rate (other flow rates available)

<i>Standard Sensors</i>	<i>Range</i>
HRLD-150	1.3 µm to 150 µm
HRLD-150JA	1.3 µm to 150 µm (chemical compatibility)
HRLD-400	2 µm to 400 µm
HRLD-600JS	2 µm to 600 µm (chemical compatibility)

Lit. No. 5223
C10 Printed in USA
© Hach Company 2010. All rights reserved.
In the interest of improving and updating its equipment, Hach Company reserves the right to alter specifications to equipment at any time.

At Hach, it's about learning from our customers and providing the right answers.

Keep it pure.

Make it simple.

Be right.

For current price information, technical support, and ordering assistance, contact the Hach office or distributor serving your area.

In the United States, contact:

HACH COMPANY World Headquarters
PO Box 608
Loveland, Colorado 80539 USA
Telephone: 800-866-7889
Fax: 970-461-3914
E-mail: customersupport@hach.com
www.hach.com/hiac

U.S. exporters and customers in Canada, Latin America, sub-Saharan Africa, Asia, and Australia/ New Zealand, contact:

HACH COMPANY World Headquarters
PO Box 608
Loveland, Colorado 80539 USA
Telephone: 970-663-9760
Fax: 970-461-3914
E-mail: intl@hach.com
www.hach.com/hiac

In Europe, the Middle East, and Mediterranean Africa, contact:

HACH LANGE GmbH
Willstätterstraße 11
D-40549 Düsseldorf, GERMANY
Tel: +49 (0) 211 5288-0
Fax: +49 (0) 211 5288-143
E-mail: info@hach-lange.de
www.hach-lange.com



HIAC 9703+ Liquid Particle Counting System

Easily configurable for maximum flexibility



Be Right™



Be Right™



Be Right™

HIAC just got easier!

PharmSpec

Customized test routines and reporting

- Customized particle size, sample volume and run test recipes for your unique application needs
- User configurable reports

Added flexibility to configure security access

- Allows administration to define multi-level security access
- Meets requirements for 21 CFR Part 11 compliance

Save time utilizing pre-configured compendial test recipes

Eliminate manual calculations – Pass/Fail limits

Injectable Drugs

- USP <788>
- EP <2.9.19>
- JP <11>, <24>
- KP <52>

Ophthalmic

- USP <789>

Inhalants

- USP<788>/Levalbuterol Inhalation

Eliminate uncertainty from data anomalies

- Alarm notification when bubbles may have impacted particle counting results
- Operator alerted if the sensors have been contaminated
- Advanced notification when routine service or calibration is needed
- Operator notification when regulatory updates become available
- All alarms provide recommended actions

Eliminate anxiety of data loss

- Data loss recovery feature protects test data in the event of a power interruption or loss
- PharmSpec allows the operator to schedule data base back up to a secure network location

Use less sample, save valuable time – protect your investment

- Industry's smallest tare volume
- Small vial clamp ensures that your sample does not spill during testing

Enjoy more walk away time

- Automated flushing and cleaning routines



Manage all of your application needs

Designed to handle the full range from very small (1 ml) to very large (>1000 ml) volumes

- Small volume injectable (SVI)
- Large volume injectable (LVI)
- Research and Development
- Troubleshooting
- Stability testing

Easily configurable for maximum flexibility

- Interchangeable sensors from 0.5 µm to 600 µm
- Interchangeable sampling probes
- Interchangeable syringes – based on sampling size
- Configurable flow rates from 10 ml/min to 100 ml/min



Turn-key validation and after sale service

Comprehensive local and on-site HIAC service programs with multiple field service offices and authorized distributors are available throughout the world

Avoid costly down time

- On-site service significantly reduces down time – no shipping of your unit required
- Notification with a phone call when your unit is in need of routine service

Instrument Validation Services

- SOP assistance and IQ/OQ validation services
- Utilize our expertise to get your system out-of-the-box and ready-to-use
- Ensure the accuracy of your instrument with documented validation and ISO 21501-3 compliant calibrations