

# Trophon EPR

Distributed exclusively by GE Healthcare, the Trophon® EPR is a complete ultrasound transducer High-Level Disinfection system that's fast, easy to use, environmentally friendly and quality assured. It uses a unique platform technology to effectively disinfect the transducer, including the shaft and handle—in just seven minutes between patients.

When used as directed the Trophon EPR meets guidelines set forth by the FDA and CDC for High-Level Disinfection. For sale in the U.S. and Canada only.



Trophon EPR  
Catalog # E8350NA

## Accessories

The Trophon EPR can stand alone or for convenient point-of-care options, we offer both the Trophon Mobile Cart and Wall Mount solution.

## Consumables

We are your one-stop shop for Trophon EPR consumable supplies.



Trophon Mobile Cart  
Catalog # E8350NH  
34.5" H x 20" W x 20.5" D  
45 lbs. 21 kg.



Trophon Wall Mount  
Catalog # E8350NG  
8.6" H x 6.7" W x 2.0" D  
3 lbs. 1.35 kg.



Trophon Sonex-HL™  
Catalog # E8350NJ  
6 cartridges/box  
Approx. 240 cycles



Trophon Chemical Indicators  
Catalog # E8350NC  
300/pkg.  
1" 26 mm



## System specifications

System Specifications	
Cycle Time	7 minutes
Cycle Verification	Via User Screen, Chemical Indicator and Electronic Log.
Sonex-HL	Minimum Concentration 31.5% Minimum Disinfectant Dose 1.0 g Volume – 80 ml, approx. 40 cycles/cartridge Shelf Life – 2 years
Chemical Indicators	The Chemical Indicators are used exclusively for monitoring the High-Level Disinfection process when placed within the Trophon EPR chamber. Refer to the Chemical Indicator Color Assessment Chart for cycle verification information.
Probe Compatibility	Please check with your transducer manufacturer regarding compatibility with the Trophon EPR. You can also find a validated transducer list at <a href="http://www.trophon.com">www.trophon.com</a>
Electrical Specifications	
Electrical Requirements	Input Voltage: 120V AC Input Current: 5 Amp/50-60 Hz
Agency Approvals	UL61010-1:2004 R10.08, IEC61010-2-040:2005, CAN/CSA-C22.2 NO. 61010-1-04+G11 (R2009), CAN/CSA-C22.2 NO. 61010-2-040-07
EMC Compliance	EN61326-1:2006 FCC Rules Part 15 Subpart B Unintentional Radiators, Class B Digital Devices ANSI C63.4:2009
Mechanical Specifications	
Trophon EPR	19.3" H x 13.6" W x 13.6" D / 38 lbs. 490 mm x 345 mm x 345 mm / 17 kg

## Ordering is simple.

Order your Trophon EPR system, accessories, and consumable supplies with the method most convenient for you.

Call: 800-558-5102, opt. 2, opt. 1

Fax: 877-279-6990

Email: [gehcaccessorysales@ge.com](mailto:gehcaccessorysales@ge.com)

Contact: Your local GE representative

Visit [www.gehealthcare.com/trophon](http://www.gehealthcare.com/trophon)

GE Healthcare  
3000 North Grandview  
Waukesha, WI 53188  
U.S.A.

[www.gehealthcare.com](http://www.gehealthcare.com)



imagination at work

©2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

\*GE and GE Monogram are trademarks of General Electric Company.

Trophon EPR is a registered trademark of Nanosonics Limited.

Sonex-HL is a trademark of Nanosonics Limited.

GE Healthcare, a division of General Electric Company.

ACCD-0198-05.12-EN-US  
DOC0992431