

# High Flow

## 40-Liter Insufflator



# Superior Performance

## Unprecedented Safety & Reliability

### Stryker 40-Liter Insufflator

Stryker's 40-liter insufflator provides superior performance with unprecedented safety and reliability. This insufflator is designed to handle the needs of today's dynamic surgical environment. Its ability to maintain pneumoperitonium under the most extreme conditions, coupled with a multitude of safety features, place the Stryker 40-liter insufflator in a class by itself.



### The 40-liter insufflator provides a host of features including:

- The highest flow rate available on the market
- Neonatal mode offers increased sensitivity and flow rates as low as .1 liters per minute
- LCD display provides surgeons with ongoing status messages
- Sidne compatibility allows the surgeon to control the insufflator with voice commands
- Software driven to allow for future upgrades
- Built in gas heating technology
- Video on-screen display



High Flow Filtered Tube Set



Heated Tube Set

### Safety

- Surgivent-electronic venting system for added safety
- Visual and audible warning signals
- Safety Relief Valve
- Internal Leakage Detection

## System Specifications

40 Liter Insufflator	SYSTEM SPECIFICATIONS
<b>Power Supply Input:</b>	100-240 Vac; 3.5A
<b>Frequency:</b>	50-60 Hz
<b>Power Consumption:</b>	150VA
<b>Dimensions:</b>	5.83" (h) x 18.7" (d) x 12.52"(w)
<b>Weight:</b>	19.84 lbs
<b>Operation Conditions - temp:</b>	50 - 104 F
<b>Operation Conditions - humidity:</b>	30% - 75%

## Product Ordering information

ORDERING INFO.	PRODUCTS
<b>620-040-503</b>	40L Insufflator
<b>620-040-050</b>	40L Insufflator Set (hose and yoke)
<b>620-040-504</b>	Sidne Ready 40L Insufflator
<b>620-040-051</b>	Sidne Ready 40L Set (hose and yoke)
<b>620-030-301</b>	High-flow Tubing (disposable)
<b>620-030-407</b>	Heated Tubing (disposable)



Stryker Endoscopy  
5900 Optical Court  
San Jose, CA 95138  
t: 1 800 624 4422 f: 1 800 729 2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with <sup>TM</sup> designation are trademarks of Stryker.  
Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 1000900366 (Rev. B)  
VID/FWG 3k 07/04