



Hall® Powered Instruments System™


Product Catalog




MICROPOWER®
SAW



MICRO 100™
SURGAIRTOME TWO®



M-CLASS™
BLADES



HALL® LITHIUM
BATTERIES



HALL®
POWERED INSTRUMENTS

Ordering Information.....	2	Micro Sagittal Saw Blades	45 - 49
Product Support	3	Micro Oscillating Saw Blades	50 - 52
Product Support & Service Programs	4	Micro Reciprocating Saw Blades & Rasps	53
Return/Repair Information.....	5	Micro Burs	54 - 56
Hall 50"	6	Sternum Saws.....	57
Indications for Use: Large Bone & Small Bone Orthopaedics		Sternum Saw Blades	58 - 60
MPower® 2	7	Indications for Use: Primary & Repeat Sternotomy	
Indications for Use: Large Bone & Small Bone Orthopaedics		Ultrapower®	61
Batteries & Accessories	8, 9	Indications for Use: Spine & Neuro	
PowerPro® & PowerPro Electric II"	10	UltraPower Burs	62 - 66
Indications for Use: Large Bone & Small Bone Orthopaedics		Indications for Use: Otology (Round Burs Only), Spine & Neuro	
PowerPro® Pneumatic	11	E9000® System	67
Indications for Use: Large Bone & Small Bone Orthopaedics		Indications for Use: ENT, Spine & Neuro	
Pin and Wire Drivers.....	12	E9005 ENT Shaver Blades.....	68 - 71
Indications for Use: Large Bone & Small Bone Orthopaedics		Indications for Use: Sinus Surgery	
Attachments.....	13	E Burs.....	72, 73
Indications for Use: Large Bone & Small Bone Orthopaedics		Indications for Use: Otology, Spine & Neuro	
Sterilization Trays.....	14	Osteon™	74
Indications for Use: Large Bone Orthopaedics		Indications for use: ENT	
M-Class® Blades	15	Osteon Burs	75
Indication for Use: Total Knee Arthroplasty		Indications for use: Otology	
Hall® Classic Blades	16	Cebotome®	76
Indication for Use: Total Knee Arthroplasty		Indications for use: Large Bone Revision Surgery	
Intrex® Blade System.....	17	Cebotome® Burs & Metal Cutting Accessories.....	77
Indication for Use: Total Knee Arthroplasty		Indication for Use: Large Bone Revisions Surgery	
Stryker Replacement Blades	18	Wires & Staples.....	78
Indication for Use: Total Knee Arthroplasty		Indication for Use: Small Bone Trauma	
Large Bone Knee Blades.....	19 - 30	Miscellaneous Accessories	79
Indication for Use: Total Knee Arthroplasty		Electric Consoles & Accessories.....	80, 81
Large Bone Oscillating Blades	31 - 33	Pneumatic Regulators, Hoses & Accessories	82
ACL Blades		Kit Listings	83 - 86
Large Bone Oscillating Blades	34, 37	Numerical Index.....	88 - 107
Indications for Use: Total Hip Arthroplasty & Large Bone Trauma		Alphabetical Index	108 - 127
Large Bone Reciprocating Blades.....	38, 39		
Indications for Use: Total Joint Arthroplasty & Large Bone Trauma			
MicroPower™ +	40		
Indications for Use: Small Bone			
MicroPower®	41		
Indications for Use: Small Bone			
Micro 100"	42		
Indications for Use: Small Bone			
Attachments & Bur Guards	43		
For Use with Micro Handpieces			
Sterilization Trays.....	44		
For Use with Micro Handpieces			

Terms

ConMed standard terms are net thirty (30) days, F.O.B. shipping point, unless otherwise approved in writing by ConMed. All orders and transactions are subject to the laws of the state of Florida.

To Place Your Order

Please provide a purchase order number, billing and shipping instructions, product catalog number and description, valid ConMed quotation, and quantity and pricing to one of the following:

Your Sales Representative

Customer Service

Phone: 800-237-0169

Fax: 727-399-5256

E-mail: custserv1@conmed.com

ConMed

Attn: Customer Service

11311 Concept Boulevard

Largo, FL 33773-4908 USA

True Purchasing Flexibility

ConMed Financial Services offers a variety of flexible financing programs, from leases to Pre-Procedure Programs to delayed payment plans. Ask your Sales Representative for more information on how you can acquire the latest technology through a financing program that fits your unique budgetary requirements.

Rental Program for Instrument Sets & Capital Equipment

In every O.R., there are times when you have to schedule a case but do not have the equipment needed. Either the capital funding for the equipment has yet to be approved or you have been unable to justify the cost of the equipment due to the procedure's low volume.

For these situations, ConMed offers a rental program for procedure-specific instruments and capital equipment. Simply call ConMed at 888-587-9742 and everything that you need for the case is on its way.

To enroll in a rental program, contact your Sales Representative or call the ConMed field inventory team. Your representative will help you customize the kit(s) to fit your procedural requirements. After the case, simply return the instruments, cleaned and sterilized if appropriate, along with any disposables you did not use. The rental fee will be billed for the disposables not returned, and customer will be responsible for any expedited shipping cost.

Failure to return instrumentation/equipment at the end of the rental period will result in the customer being charged the standard rental fee automatically for each week that the instrumentation/equipment is not received by ConMed.

At ConMed, we want to offer more than a quality product. We take pride in being able to provide you with value added services to support our products. The information given in this catalog is designed to address your questions regarding product support.

Customer Support

If you encounter a problem or need technical assistance, contact Customer Service (Product Support) or your Sales Representative. Refer to the Customer Service telephone numbers listed.

Regulatory Compliance

Report any injuries or malfunctions to ConMed Customer Service. Refer to the Customer Service telephone numbers listed.

Original Warranty

ConMed (the "Company") warrants to the first purchaser or lessee ("You" or "Customer") that the instruments, attachments and parts manufactured by or for the Company (hereinafter collectively "Instruments") have been tested, inspected, and shipped in proper working order.

The Company warrants all new instruments to be free from defects in materials and workmanship for the following periods, measured from date of sale*:

Product Categories	Months
All Single Use Products.....	Upon Receipt
Powered Surgical Instruments.....	12
Power Attachments and Accessories.....	12
Handpiece Cords and Power Cords	6
Pneumatic Hoses	6
Bur Guards, Blade Guards, Attachments	6
Battery Packs and Accessories	3
Battery Chargers	12
Camera Consoles	12
Video Components	12
Video Cables and Light Guides.....	3
Scopes	12
Camera Heads	12
Instruments.....	12
Instrument Tray	3
Footswitches	12
Irrigation Systems (Pumps and Remotes).....	12
Reusable Procedure Specific Instruments.....	6

Following the Original Warranty, ConMed offers additional warranty and service programs that support our products throughout their lifetime. Service programs address the necessity to control cost and greatly reduce repair turnaround time by eliminating the need for repair quotes. Products not covered under a service program are subject to variable repair charges. To find out which service option best meets your needs, contact your local Sales Representative.

Service Agreements

Service Agreements provide additional service coverage on capital equipment. This additional warranty period covers all repair costs, periodic maintenance, free loaner usage during repair time and freight charges out from ConMed both on the loaner product and the repaired instrument.

Service Warranty

After expiration of the Original or Service Agreement period, instruments may be returned to ConMed for repair at the customer's expense. ConMed warrants that any service or repair work performed by ConMed shall be free from defects in workmanship and materials for a period of 90 days after the date the customer receives the Instrument(s). This limited warranty applies only to the actual service or repair work performed by ConMed or its authorized service centers. ConMed warrants that all parts and assemblies used in the repair or service of Instruments meet new part functional specifications, although some parts or assemblies may have been remanufactured.

ConMed components are manufactured to exacting tolerances and design specifications. As a standard practice, all repairs are subject to strict product performance and testing specifications prior to returning the Instrument to the customer. ConMed ensures that each Instrument is inspected according to the guidelines used for new manufactured Instruments. If an Instrument is determined unrepairable, ConMed will offer to sell to the customer a refurbished or brand new Instrument or comparable equivalent at a preferential price. All parts and assemblies replaced by ConMed shall become the property of ConMed.

Repair Quotes

Upon request Customer Service can estimate the repair charges for non-warranty repairs. If a firm repair quote is requested, please return your product for inspection and you will be contacted before the repair is initiated. If a product is returned and a firm quote is not requested, the repair will be processed and your account billed accordingly.

Minimum Service Charge

There is a minimum repair charge for the evaluation and testing of all returns (except for products covered under warranty).

Return Shipping

All repaired products will be returned to you via second day air delivery unless you request an alternate shipping method. The cost of return shipping is your responsibility unless the product is under warranty or service agreement.

* Effective February 1, 2014

Warranty Terms

If the customer discovers within the warranty period (Original or Service Agreement) that an Instrument has a defect in materials and/or workmanship, the customer must promptly notify the Company. Also, if the customer discovers within the service warranty period that any service or repair work performed by ConMed has a defect in workmanship and/or materials, the customer must promptly notify the Company. If it becomes necessary to return the Instrument to the Company, the customer must (a) acquire a "Service Request" (SR) authorization from the Company Customer Service Department, (b) pack the unit carefully, and (c) return it to the Company.

Within a reasonable time after receipt of the Instrument, the Company will investigate and correct any defect covered by warranty by providing, at its option, one of the following: service or repair of the Instrument, or replacement of the Instrument. All parts and assemblies replaced by ConMed shall become the property of ConMed. These remedies are the customer's exclusive remedies under the warranty.

The foregoing limited warranties do not apply to:

- » Instruments that have been tampered with, altered, neglected, or misused,
- » Instruments damaged through use with other than Company authorized accessories, attachments, burs or blades,
- » Instruments not manufactured by or for the Company,
- » Instruments used for purposes other than those for which they were designed and manufactured, including use in any way inconsistent with the instructions and warnings contained in the Company instruction manuals and package inserts, and
- » Instruments that were last serviced, refurbished, reprocessed or reconditioned by a nonauthorized service entity.

THE FOREGOING LIMITED WARRANTIES ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

In no case shall the Company be liable for any special, incidental, indirect or consequential damages based upon breach of warranty or any other legal theory. Some jurisdictions do not allow limits on warranties, or on remedies, and, in such jurisdictions, the limits in this and the preceding paragraphs may not apply.

The Company reserves the right (a) to make design changes to Instruments at any time without notice to customer and without incurring any obligation to incorporate those changes into Instruments previously purchased or leased, and (b) to make changes from time to time in the contents of any publication, instruction manual or package insert without any obligation to notify customers of such revisions or changes, and (c) to change price without notice.

U.S. Procedure for Returns/Repairs

1. All Instruments must be cleaned and sterilized, if appropriate, prior to return of any Instrument to ConMed. Instruments contaminated with biohazardous materials will be returned to you at your expense.
2. An authorized Service Request (SR) Number must be obtained prior to returning any Instrument for repair/return.
3. To return a product for repair, contact ConMed Customer Service. The following information is needed to process a repair/return:
 - » Catalog Number
 - » Serial/Lot Number
 - » Original Invoice Number
 - » Detailed description of the problem or reason for return
4. The SR number must be displayed on the outside of the box, included with all paperwork, and must be referred to when making an inquiry about the status of a repair/return.
5. ConMed maintains various repair centers around the world. When you are issued your SR number, you will be advised of the location where to send your product.

Return Terms

If it becomes necessary to return a product after purchase, the product must be returned within 45 days of its ship date. The product must be returned unused, in the original packaging and be restockable. Returns must be shipped prepaid freight. Returns that are not shipped prepaid freight will not be accepted. All product returns, except those attributable to a ConMed error, are subject to a restocking fee of 15% of the purchase price or \$25, whichever is greater.

Products Not Eligible for Credit:

- » Sterile products/single use devices returned with the outer packaging seals broken
- » Non-sterile products that have been sterilized and/or placed in the surgical field
- » Products damaged by customer misuse or negligence
- » Altered or customized products
- » Products damaged in return transit due to inadequate packaging
- » Less than full packages – all products must be in original packaging
- » Obsolete products

Packaging

Whenever it is required to return your product, be sure to package it in a protective carton. We recommend that you save the original shipping container for this purpose. In-transit damage is not covered by the warranty; therefore, it is best to always insure shipments.

Hall 50™

Indications for Use: Large Bone
& Small Bone Orthopaedics



Attachments Listed on Pages 12 and 13

Cat. No.

Hall 50 1-Trigger Modular HandpiecePR07100B

Hall 50 2-Trigger Modular HandpiecePR07200B

Hall 50 Oscillating Saw HandpiecePR07300B

Hall 50 Reciprocating Saw HandpiecePR07400B



PR07100B



PR07200B



PR07300B



PR07400B

All products may not be available in all markets.

Attachments Listed on Pages 12 and 13

	Cat. No.
Single Trigger Modular Handpiece	PRO6200M
2-Trigger Modular Handpiece	PRO6202M

Oscillating Saw

Reciprocating Saw



PRO6200M

PRO6202M



PRO6350M

PRO6400M

All products may not be available in all markets.



L3000



L3000SM

L3000LG



PRO3600



PRO3020



PRO3010



PRO3011

Cat. No.

Hall® Lithium Charger L3000

Hall Large Lithium Battery, 13.2 Volt L3000LG

Hall Small Lithium Battery, 13.2 Volt L3000SM

NiCad Battery Charger PRO3600

Small NiCad Battery, 12 Volt PRO3020

Large NiCad Battery, 9.6 Volt PRO3010

X-Large NiCad Battery, 12 Volt PRO3011

All products may not be available in all markets.

	Cat. No.
Small NiMH Sterile Transfer Battery, 12 Volt	PRO3520
Large NiMH Sterile Transfer Battery, 9.6 Volt.....	PRO3115

Battery Case for PRO3520	PRO3521
Battery Case for PRO3115.....	PRO3110

Transfer Shroud for PRO3521	PRO3531
Transfer Shroud for PRO3110	PRO3130

Charging Platform for PRO3520.....	PRO3535
Charging Platform for PRO3115	PRO3135



PRO3520



PRO3115



PRO3521



PRO3110



PRO3531



PRO3130



PRO3535



PRO3135

All products may not be available in all markets.

PowerPro® & PowerPro Electric II™

Indications for Use: Large Bone
& Small Bone Orthopaedics



PRO6100



PRO6125



PRO6135

Attachments Listed on Pages 12 and 13

Cat. No.

Single Trigger Modular HandpiecePR06100

Oscillating SawPR06125

Reciprocating SawPR06135

All products may not be available in all markets.

Attachments Listed on Pages 12 and 13

	Cat. No.
Single Trigger Modular Handpiece	PRO6150



PRO6150

Oscillating Saw	PRO6175
-----------------------	---------



PRO6175

Reciprocating Saw	PRO6185
-------------------------	---------



PRO6185

All products may not be available in all markets.

Pin and Wire Drivers

Indications for Use: Large Bone
& Small Bone Orthopaedics



PRO6128



PRO6140



PRO6228



PRO6240

Single Trigger Wire and Pin Drivers

Cat. No.

Wire Driver, Range: 0.7 – 1.6mm (0.028 – 0.062")..... PR06128

Pin Driver, Range: 1.8 – 4.0mm (0.071 – 0.157") PR06140

2-Trigger Wire and Pin Drivers

Wire Driver Range: 0.7 – 1.6mm (0.028 – 0.062")..... PR06228

Pin Driver Range: 1.8 – 4.0mm (0.071 – 0.157")..... PR06240

Drill-Ream 1:1/5:1 Combination Attachments



Drill-Ream Attachment, AO
Cat. No. **PRO1540**



Drill-Ream Attachment, 1/4" (6.35mm) Chuck
Cat. No. **PRO1541**



Drill-Ream Attachment, Zimmer/Hudson
Cat. No. **PRO1547**



Drill-Ream Attachment, DePuy/Hudson
Cat. No. **PRO1548**

1:1 Multipurpose Attachments



1/2" (4.0mm) Jacobs Chuck, with Chuck Key
Cat. No. **PRO2030**



1/2" (4.0mm) Replacement Chuck Key
Cat. No. **5044-999-53**



1/4" (6.35mm) Jacobs Chuck
Cat. No. **PRO2041**



1/4" (6.35mm) Replacement Chuck Key
Cat. No. **5044-999-52**



1/4" Keyless Chuck
Cat. No. **PRO6050**



1/4" Keyless Chuck
Cat. No. **PRO2250**

1:1 Quick Connect Drill Attachments



AO/Trinkle
Cat. No. **PRO2029**



Small AO—Low Profile
Cat. No. **PRO2038**



Small AO—One Handed Operation
Cat. No. **PRO2138**



Zimmer/Hudson
Cat. No. **PRO2046**



Large AO
Cat. No. **PRO2070**



Standard Hudson (DePuy)
Cat. No. **PRO2075**



Mini AO
Cat. No. **PRO2039**

3:1 Reamers



3:1 1/4" High Torque Jacobs Chuck
Cat. No. **PRO2042**



3:1 Zimmer/Hudson Reamer
Cat. No. **PRO2047**



3:1 Hudson Reamer (DePuy)
Cat. No. **PRO2065**



3:1 Large AO Reamer
Cat. No. **PRO2040**



3:1 Aesculap Reamer
Cat. No. **PRO2060**



3:1 DHS Reamer
Cat. No. **PRO2080**



3:1 Trinkle Reamer
Cat. No. **PRO2329**

5:1 Reamers



5:1 1/4" High Torque Jacobs Chuck
Cat. No. **PRO6042**



5:1 Zimmer/Hudson Reamer
Cat. No. **PRO6047**



5:1 Hudson Reamer(DePuy)
Cat. No. **PRO6048**



5:1 Large AO Reamer
Cat. No. **PRO6040**



5:1 Aesculap Reamer
Cat. No. **PRO6060**

Miscellaneous Attachments



Radiolucent Drive Adaptor
Cat. No. **PRO2514**



High Speed Drill Attachment
Cat. No. **PRO2033**



Sagittal Saw Attachment
Cat. No. **PRO2043**



Multipurpose High Speed
Cat. No. **PRO7020**



Reciprocating Saw Attachment for PRO6100, PRO6150, PRO6200M, & PRO6202M
Cat. No. **PRO6045**

Sterilization Trays

Indications for Use: Large Bone Orthopaedics



TR34R with PRO7000T




**TR34R with PRO7001T,
and PRO7002T**



TR12R with PRO7001T



TR12R with PRO7002T



PRO6001

Hall 50™ Sterilization Containers

Cat. No.

$\frac{3}{4}$ Size Rigid Sterilization Container TR34R

$\frac{1}{2}$ Size Rigid Sterilization Container TR12R


Hall 50 Sterilization Inner Trays

$\frac{3}{4}$ Size Inner Tray, Full Set PRO7000T

$\frac{1}{2}$ Size Inner Tray, Drill Set PRO7001T

$\frac{1}{2}$ Size Inner Tray, Saw Set PRO7002T

MPower® 2, Large Bone Tray with Battery Holders PRO6001




Blades images are 1:1 in length and width.

Body Width	19.0mm	19.0mm	13.0mm	19.0/13.0mm	25.0mm
Length	105.0mm	90.0mm	90.0mm	90.0mm	90.0mm
Thickness	Blade Catalog #				
0.89mm		TN190-089-90	TN130-089-90		TN250-089-90
1.00mm		TN190-100-90	TN130-100-90		TN250-100-90
1.07mm		TN190-107-90	TN130-107-90		TN250-107-90
1.19mm	TN190-119-05	TN190-119-90	TN130-119-90	TN193-119-90	TN250-119-90
1.27mm	TN190-127-05	TN190-127-90	TN130-127-90	TN193-127-90	TN250-127-90
1.37mm	TN190-137-05	TN190-137-90	TN130-137-90	TN193-137-90	TN250-137-90
1.47mm		TN190-147-90	TN130-147-90		TN250-147-90

Hall® Classic Blades


Indication for Use: Total Knee Arthroplasty



Blades images are 1:1 in length and width.

Body Width	19.0mm	19.0mm	13.0mm	19.0/13.0mm	25.0mm
Length	105.0mm	90.0mm	90.0mm	90.0mm	90.0mm
Thickness	Blade Catalog #				
0.89mm	V190-089-05	V190-089-90	V130-089-90		V250-089-90
1.00mm	V190-100-05	V190-100-90	V130-100-90		V250-100-90
1.07mm	V190-107-05	V190-107-90	V130-107-90		V250-107-90
1.19mm	V190-119-05	V190-119-90	V130-119-90	V193-119-90	V250-119-90
1.27mm	V190-127-05	V190-127-90	V130-127-90	V193-127-90	V250-127-90
1.37mm	V190-137-05	V190-137-90	V130-137-90	V193-137-90	V250-137-90
1.47mm	V190-147-05	V190-147-90	V130-147-90		V250-147-90

All products may not be available in all markets.




Blades images are 1:1 in length and width.

	Standard	MIS	Standard/ Wide	Standard/Long	Standard/ Narrow	Unicondylar
Head Width	19.5mm	25.0mm	25.0mm	25.0mm	12.5mm	13.0mm
Body Width	19.5mm	9.0mm	25.0mm	25.0mm	12.5mm	8.0mm
Length	90.0mm	90.0mm	90.0mm	100.0mm	90.0mm	70.0mm
Thickness	Blade Catalog #					
0.89mm	5071-300	5071-310	5071-320	5071-330	5071-340	5071-350
1.00mm	5071-301	5071-311	5071-321	5071-331	5071-341	5071-351
1.07mm	5071-302	5071-312	5071-322	5071-332	5071-342	5071-352
1.19mm	5071-303	5071-313	5071-323	5071-333	5071-343	5071-353
1.27mm	5071-304	5071-314	5071-324	5071-334	5071-344	5071-354
1.37mm	5071-305	5071-315	5071-325	5071-335	5071-345	5071-355
1.47mm	5071-306	5071-316	5071-326	5071-336	5071-346	5071-356

All products may not be available in all markets.

Stryker Replacement Blades

Indication for Use: Total Knee Arthroplasty









Blades images are 1:1 in length and width.

Body Width	19.0mm	19.0mm	13.0mm	19.0/13.0mm	25.0mm
Length	105.0mm	90.0mm	90.0mm	90.0mm	90.0mm
Thickness	Blade Catalog #				
0.89mm		STK190-089-90	STK130-089-90		STK193-137-90
1.00mm					STK250-100-90
1.07mm					
1.19mm		STK190-089-90	STK130-119-90		STK250-119-90
1.27mm	STK190-127-05	STK190-127-90	STK130-127-90		STK250-127-90
1.37mm		STK190-137-90		STK193-137-90	STK250-137-90
1.47mm					STK250-147-90

All products may not be available in all markets.

Hall® Blades

Thickness: 0.60mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-650	49.00	63.00	0.60	0.60	
Thickness: 0.64mm					
5071-594	29.46	58.70	0.64	0.64	
5271-502	29.00	74.00	0.64	0.64	
Thickness: 0.89mm					
5071-528	25.00	73.00	0.89	0.89	
5071-533	21.00	84.50	0.89	0.89	
5071-583	12.70	90.00	0.89	0.89	








Products shown at actual size. All products may not be available in all markets.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty










Hall® Blades

Thickness: 0.89mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-589	13.00	90.00	0.89	0.89	
Thickness: 0.90mm					
5071-532	12.00	80.00	0.90	0.90	
5071-185	19.50	80.00	0.90	0.90	
5071-195	19.50	80.00	0.90	0.90	
5071-146	25.40	90.00	0.90	0.90	
5071-196	25.40	90.00	0.90	0.90	
5071-120	19.50	95.00	0.90	0.90	

Grey shaded blades indicate coating.
Products shown at actual size.

Hall® Blades

Thickness: 0.90mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-121	25.40	95.00	0.90	0.90	
Thickness: 1.00mm					
5071-562	25.40	90.00	1.00	1.00	
5071-538	25.40	95.00	1.00	1.00	
Thickness: 1.19mm					
5071-548	12.70	86.00	1.19	1.19	
5071-547	25.40	86.00	1.19	1.19	
Thickness: 1.20mm					
5071-167	19.50	90.00	1.20	1.20	
5071-184	19.50	90.00	1.20	1.20	

Grey shaded blades indicate coating. Products shown at actual size. All products may not be available in all markets.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty










Hall® Blades

Thickness: 1.24mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-558	18.50	73.00	1.24	1.24	
5071-534	25.40	79.50	1.24	1.24	
5071-570-50	25.00	95.00	1.24	1.24	
5071-526	18.50	107.00	1.24	1.24	
Thickness: 1.26mm					
5071-573	18.50	90.00	1.26	1.26	
5071-575	31.00	90.00	1.26	1.26	

Products shown at actual size.

Hall® Blades

Thickness: 1.27mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-527	19.00	80.00	1.27	1.27	
5071-180	10.00	85.00	1.27	1.27	
5271-580	10.00	85.00	1.27	1.27	
5071-531	25.40	85.00	1.27	1.27	
5071-181	19.50	86.00	1.27	1.27	
5071-201	19.50	86.00	1.27	1.27	
5071-582	13.00	90.00	1.27	1.27	

Grey shaded blades indicate coating.
Products shown at actual size.

All products may not be available in all markets.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty






Hall® Blades

Thickness: 1.27mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-161	25.40	90.00	1.27	1.27	
5071-564	25.40	90.00	1.27	1.27	
5271-561	25.40	90.00	1.27	1.27	
5071-590	13.00	95.00	1.27	1.27	
5071-568	16.50	95.00	1.27	1.27	
5071-576	16.50	95.00	1.27	1.27	
5071-123	19.50	95.00	1.27	1.27	

Products shown at actual size.

Hall® Blades

Thickness: 1.27mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-519	21.00	95.00	1.27	1.27	
5071-124	25.40	95.00	1.27	1.27	
5071-129	19.50	105.00	1.27	1.27	
5071-529	19.50	105.00	1.27	1.27	

Products shown at actual size.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty








Hall® Blades

Thickness: 1.37mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-163	25.40	90.00	1.37	1.37	
5071-525	19.00	95.00	1.37	1.37	
5071-530	19.00	95.00	1.37	1.37	
5071-539	25.40	95.00	1.37	1.37	
Thickness: 1.47mm					
5071-565	25.40	90.00	1.47	1.47	
Thickness: 1.50mm					
5071-144	25.40	90.00	1.50	1.50	

Grey shaded blades indicate coating.
Products shown at actual size.

Stryker Blades

For use with: System 2000, System 4 & 5

Thickness: 0.64mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5072-145	9.12	25.22	0.64	0.64	
5072-140	29.46	58.70	0.64	0.64	
5072-100	25.00	73.79	0.64	0.64	
5072-150	12.50	73.80	0.64	0.64	
Thickness: 0.80mm					
5072-143	25.40	90.00	0.80	0.80	

Products shown at actual size.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty



Stryker Blades








For use with: System 2000, System 4 & 5

Thickness: 0.89mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5072-105	24.82	73.79	0.89	0.89	
5072-109	20.70	85.04	0.89	0.89	
5072-102	24.82	90.47	0.89	0.89	
Thickness: 0.90mm					
5072-195	19.50	80.00	0.90	0.90	
Thickness: 1.20mm					
5072-152	12.50	82.18	1.20	1.20	
5072-151	25.00	85.67	1.20	1.20	
5072-142	25.40	90.00	1.20	1.20	

Grey shaded blades indicate coating.
Products shown at actual size.

Stryker Blades

For use with: System 2000, System 4 & 5

Thickness: 1.24mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5072-218	25.00	98.50	1.24	1.24	
5072-185	18.50	104.00	1.24	1.24	
Thickness: 1.27mm					
5072-110	24.82	73.79	1.27	1.27	
5072-111	24.82	80.14	1.27	1.27	
5072-158	12.50	82.20	1.27	1.27	
5072-183	29.00	84.47	1.27	1.27	
5072-101	19.50	86.00	1.27	1.27	
5072-181	19.50	86.00	1.27	1.27	

Grey shaded blades indicate coating.
Products shown at actual size. All products may not be available in all markets.

Large Bone Knee Blades

Indication for Use: Total Knee Arthroplasty








Stryker Blades

For use with: System 2000, System 4 & 5

Thickness: 1.37mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5072-177	25.00	82.42	1.37	1.37	
Thickness: 1.50mm					
5072-176	25.00	80.14	1.50	1.50	

Products shown at actual size.

Graft Harvesting Blades

Cat. No.	Depth Stop (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)		
5071-105	20.00	16.50	0.40	5.00		5mm
5071-106	20.00	16.50	0.40	6.00		6mm
5071-107	20.00	16.50	0.40	7.00		7mm
5071-108	20.00	16.50	0.40	8.00		8mm
5071-109	20.00	16.50	0.40	9.00		9mm
5071-110	20.00	16.50	0.40	10.00		10mm
5071-111	20.00	16.50	0.40	11.00		11mm
5071-112	20.00	16.50	0.40	12.00		12mm
5071-113	20.00	16.50	0.40	13.00		13mm

Products shown at actual size.

Large Bone Oscillating Blades

ACL Blades





ACL Blades

Thickness: 0.60mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-133	9.50	25.50	0.40	0.60	
5071-134	9.00	41.00	0.40	0.60	
Thickness: 0.64mm					
5071-236	9.12	25.22	0.64	0.64	
Thickness: 0.80mm					
5071-127	8.50	38.00	0.40	0.80	
5071-246	5.50	25.50	0.60	0.80	

Products shown at actual size.

ACL Blades

Thickness: 0.84mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-571	9.50	30.00	0.64	0.84	
Thickness: 1.10mm					
5071-135	5.00	30.00	0.60	1.10	

Products shown at actual size.

Large Bone Oscillating Blades

Indications for Use: Total Hip Arthroplasty
& Large Bone Trauma










Oscillating Blades

Thickness: 0.80mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-132	19.50	41.00	0.40	0.80	
Thickness: 0.90mm					
5071-553	19.50	63.00	0.60	0.90	
5071-551	32.50	63.00	0.60	0.90	
5071-550	49.00	63.00	0.60	0.90	
5071-552	13.50	64.00	0.60	0.90	
5071-166	32.00	64.00	0.80	0.90	

All products may not be available in all markets.

Grey shaded blades indicate coating.
Products shown at actual size.

Oscillating Blades

Thickness: 0.90mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-186	32.00	64.00	0.80	0.90	
5071-556	33.50	64.00	0.80	0.90	
5071-177	19.50	71.00	0.80	0.90	
5071-187	19.50	71.00	0.80	0.90	
5071-572	19.50	71.00	0.80	0.90	
5071-178	19.50	80.00	0.80	0.90	
5071-188	19.50	80.00	0.80	0.90	

Grey shaded blades indicate coating.

Products shown at actual size.

All products may not be available in all markets.

Large Bone Oscillating Blades

Indications for Use: Total Hip Arthroplasty
& Large Bone Trauma







Oscillating Blades


Thickness: 0.90mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-174	19.50	86.00	0.80	0.90	
5071-143	25.40	90.00	0.80	0.90	
Thickness: 1.00mm					
5071-560	25.40	95.00	0.86	1.00	
Thickness: 1.06mm					
5071-557	12.50	73.80	0.63	1.06	
Thickness: 1.40mm					
5071-554	19.50	71.00	1.20	1.40	
5071-141	19.50	80.00	1.20	1.40	

Grey shaded blades indicate coating.
Products shown at actual size.

Oscillating Blades

Thickness: 1.40mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5071-140	25.40	80.00	1.20	1.40	
5071-142	25.40	90.00	1.20	1.40	
Thickness: 1.50mm					
5071-555	33.00	71.00	1.20	1.50	
5071-145	13.00	72.00	1.20	1.50	

Oscillating Abrader

Cat. No.	Width (mm)	Length (mm)	
5071-175	20.00	45.70	

Grey shaded blades indicate coating.
Products shown at actual size.

All products may not be available in all markets.

Large Bone Reciprocating Blades

Indications for Use: Total Joint Arthroplasty
& Large Bone Trauma







Hall® Reciprocating Blades

Thickness: 0.60mm					
Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5052-260	12.50	38.00	0.40	0.60	
Thickness: 0.90mm					
5052-179	12.50	76.00	0.90	0.90	
Thickness: 1.00mm					
5052-259	9.40	54.50	0.60	1.00	
5052-173	12.50	73.50	0.80	1.00	
5052-258	12.50	89.00	0.60	1.00	
Thickness: 1.27mm					
5052-291	12.50	76.00	1.27	1.27	
Thickness: 1.50mm					
5052-276	13.00	48.00	1.50	1.50	
5052-281	9.40	54.50	1.50	1.50	
Thickness: 1.60mm					
5052-261	12.50	89.00	1.50	1.60	

Grey shaded blades indicate coating.
Products shown at actual size.

Stryker Reciprocating Blades

For use with: System 2000, System 4 & 5

Cat. No.	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5052-334	12.50	73.50	0.60	0.60	
5052-332	12.00	76.00	0.80	0.80	
5052-135	12.50	76.00	0.90	0.90	
5052-333	12.00	76.00	1.30	1.30	


Products shown at actual size.



**Oral Max™ Elite
High Speed Drill**



**Orthopedic High
Speed Drill**




Medium Speed Drill



Sagittal Saw



Oscillating Saw



Reciprocating Saw

Attachments on page 43

Handpieces

Removable Lever

Oral Max Elite High Speed Drill.....PRO7000DE

Orthopedic High Speed DrillPRO7000SE

Medium Speed Drill PRO7100SE

Sagittal Saw.....PRO7200SE

Oscillating SawPRO7300SE

Reciprocating Saw PRO7400SE

All products may not be available in all markets.

Attachments on page 43

Handpieces

	Removable Lever	Permanently Attached Lever
Oral Max High Speed Drill	6020-025	6021-025
Orthopedic High Speed Drill	6020-026	6021-026
Medium Speed Drill	6020-021	6021-021
Sagittal Saw	6020-022	6021-022
Reciprocating Saw	6020-023	6021-023
Oscillating Saw	6020-024	6021-024



Oral Max™ High Speed Drill



Orthopedic High Speed Drill



Medium Speed Drill



Sagittal Saw



Reciprocating Saw




Oscillating Saw

All products may not be available in all markets.

Micro 100™

Indications for Use: Small Bone



Attachments on Page 43

Handpieces

Cat. No.

Surgairtome Two® High Speed Drill Kit (System listing on page 83)5058-000

Surgairtome Two High Speed Drill5058-001

Medium Speed Drill5053-009

Reciprocating Saw5053-010

Sagittal Saw5053-011

Oscillating Saw5053-012

XL Oscillating Saw5053-022

Wire Driver5053-013

Throttle Clip (for lever) Used during footswitch operation5053-074

Attachments & Accessories

Wire Guard for use with 5053-0135053-123

Wire Driver Adapter Chuck for use with 5053-0135053-019

Replacement Hex Wrench for Micro 100™ Sagittal & Oscillating Saw1365-005

Burguard Sleeve Surgairtome Two1375-112

Smartguard®5020-160

Smartguard Irrigator6020-260

Attachments

	Cat. No.
Medium Precision Bur Guard.....	1375-012P
Long Precision Bur Guard	1375-011P
X-Long Precision Bur Guard.....	1375-023P
Medium Bur Guard	1375-012
Long Bur Guard.....	1375-011
X-Long Bur Guard.....	1375-023
Laminectomy Guard.....	1375-020
Use 5092-103/203 long carbide burs, the laminectomy burs are the only burs that can be inserted prior to placing bur guard.	
Rhinoplasty Guard	1375-021
Tissue Retractor	1387-022
Use only with long burs (5092 series) which have diameter of 4.0mm or less. Recommended burs 5092-136.	
20 degree Angled Attachment	1375-032
Use only 5092 series X-long burs.	
20 degree X-Long Angled Attachment.....	1375-033
Use only 5093 series X-long burs.	
70 degree Contra Angle with Dental Head.....	1375-034
Use friction grip burs only – size .0625" or 1.58mm diameter. Not available from ConMed.	
70 degree Contra Angle with Surgical Head	1375-035
Use only 5089 ultra-short burs and 5090 series short burs.	
90 degree Angled Attachment	1375-036
Use only 5089 X-short burs and 5090 series short burs.	
Bur Changer for 70 and 90 degree Angled Attachments.....	1375-003
Irrigation Tubing Set.....	E9418



All products may not be available in all markets.

Sterilization Trays

For Use with Micro Handpieces




Hall 50™ Sterilization Containers

Cat. No.

½ Size Rigid Sterilization ContainerTR12R

Sterilization Inner Trays

½ Size Inner Tray, Small Bone Set.....PR07003T












PRO6000

Small Bone TrayPR06000

All products may not be available in all markets.

Sagittal Saw Blades








Thickness: 0.60mm						
Cat. No.	Tooth	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5023-144	Fine	5.50	18.50	0.40	0.60	
5023-118	Fine	9.50	25.50	0.40	0.60	
5023-130	Coarse	19.50	25.50	0.40	0.60	
5023-131	Coarse	19.50	41.00	0.40	0.60	
5023-132	Coarse	5.50	25.50	0.40	0.60	
5023-133	Coarse	14.00	41.00	0.40	0.60	
5023-137	Coarse	14.00	25.50	0.40	0.60	
5023-138	Fine	9.50	25.50	0.40	0.60	
5023-139	Fine	4.50	25.50	0.40	0.60	
5023-140	Tapered/Fine	9.50	41.00	0.40	0.60	

Products shown at actual size.
All products may not be available in all markets.

Micro Sagittal Saw Blades








Sagittal Saw Blades

Thickness: 0.60mm						
Cat. No.	Tooth	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5023-141	Fine	6.00	34.00	0.40	0.60	
5023-142	Fine	4.00	34.00	0.40	0.60	
5023-143	Coarse	9.50	25.50	0.40	0.60	
5023-244	Fine	4.00	6.00	0.40	0.60	
5044-145	Fine	10.00	18.50	0.40	0.60	
5023-247	Fine	9.50	25.50	0.60	0.60	
5023-249	Fine	9.00	31.00	0.40	0.60	

Products shown at actual size.

Angled Sagittal Saw Blades










Cat. No.	Tooth	Angle	Width (mm)	Length (mm)	Thickness: 0.60mm		
					Material Thickness (mm)	Cut Thickness (mm)	
5023-128	Fine	70 Degree	9.50	7.50	0.40	0.60	
5023-129	Fine	70 Degree	9.50	13.50	0.40	0.60	
5023-134	Coarse	90 Degree	5.50	9.00	0.40	0.60	
5023-135	Coarse	45 Degree	5.50	18.00	0.40	0.60	
5023-136	Coarse	90 Degree	14.00	15.50	0.40	0.60	

Products shown at actual size.

Micro Sagittal Saw Blades






ACL Blades





Thickness: 0.60mm								
Cat. No.	Tooth	Edge	Depth Stop (mm)	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
C8803	Fine	Straight	8.0 with guide/12 without guide	9.40	41.00	0.38	0.60	
C8804	Fine	Straight	8.0 Freehand	9.40	41.00	0.38	0.60	
C8806	Fine	Straight	10.0 Freehand	9.40	41.00	0.38	0.60	
8703	Fine	Straight	8.0 with guide/12 without guide	9.50	25.50	0.40	0.60	
5023-183	Fine	Straight	None	9.50	25.50	0.40	0.60	
5023-184	Tapered/Fine	Straight	None	9.50	41.00	0.40	0.60	
5023-185	Tapered/Coarse	Straight	None	6.50	30.00	0.40	0.60	
Thickness: 0.69mm								
5023-182	Coarse	Straight	None	19.00	41.00	0.40	0.69	
Thickness: 0.74mm								
5023-186	Coarse	Straight	None	8.50	38.00	0.40	0.74	

Products shown at actual size.

Dermabraders

Thickness: 3.80mm						
Cat. No.	Tooth	Edge	Width (mm)	Length (mm)	Material Thickness (mm)	
5023-055	Abrader	Straight	25.00	6.30	3.80	
5023-056	Abrader	25 Degree	25.00	6.30	3.80	
5023-057	Abrader	70 Degree	25.00	6.30	3.80	

Graft Harvesting Blades

Cat. No.	Depth Stop (mm)	Width (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5023-405	20.00	16.50	0.40	5.00	
5023-406	20.00	16.50	0.40	6.00	
5023-407	20.00	16.50	0.40	7.00	
5023-408	20.00	16.50	0.40	8.00	
5023-409	20.00	16.50	0.40	9.00	
5023-410	20.00	16.50	0.40	10.00	
5023-411	20.00	16.50	0.40	11.00	
5023-412	20.00	16.50	0.40	12.00	
5023-413	20.00	16.50	0.40	13.00	

Products shown at actual size.

All products may not be available in all markets.

Micro Oscillating Saw Blades



Oscillating Saw Blades





Thickness: 0.60mm							
Cat. No.	Tooth	Edge	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5023-160	Fine	Straight	5.50	12.00	0.40	0.60	
5023-161	Fine	Straight	10.00	12.00	0.40	0.60	
5023-162	Fine	Straight	5.50	18.50	0.40	0.60	
5023-163	Fine	Straight	10.00	18.50	0.40	0.60	
5023-164	Fine	Straight	5.50	28.50	0.40	0.60	
5023-165	Fine	Radial	10.00	28.50	0.40	0.60	
5023-166	Fine	Radial	18.50	28.50	0.40	0.60	
5023-169	Fine	Straight	10.00	35.00	0.40	0.60	

Angled Oscillating Saw Blades

Thickness: 0.60mm							
Cat. No.	Tooth	Edge	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	
5023-258	Fine	Radial/20 Degree	10.00	9.00	0.40	0.60	
5023-259	Fine	Radial/20 Degree	10.00	13.50	0.40	0.60	
5023-167	Fine	20 Degree	5.50	13.00	0.40	0.60	
5023-168	Fine	Radial/20 Degree	10.00	13.00	0.40	0.60	

Products shown at actual size.

Intra-Oral Blades *Not for use with Micro 100™





Thickness: 0.60mm							
Cat. No.	Tooth	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	Cutting Depth (mm)	
5023-101*	Coarse	11.50	7.00	0.40	0.60	7.00	
5023-102*	Coarse	11.50	12.00	0.40	0.60	12.00	
5023-103*	Coarse	9.00	12.00	0.40	0.60	12.0	
5023-104*	Coarse	5.00	12.00	0.40	0.60	12.00	

Products shown at actual size. All products may not be available in all markets.

Micro Oscillating Saw Blades








Crescentic Blades


Thickness: 0.60mm							
Cat. No.	Tooth	Width (mm)	Length (mm)	Material Thickness (mm)	Cut Thickness (mm)	Cutting Depth (mm)	
5023-170	Fine	9.0	12.00	0.40	0.60	12.00	
5023-171	Fine	13.5	25.50	0.40	0.60	25.50	
5023-172	Fine	18.0	32.00	0.40	0.60	32.00	
5023-176	Fine	13.5	32.00	0.40	0.60	32.00	

Products shown at actual size.





Reciprocating Rasp

Cat. No.	Shank	Width (mm)	Length (mm)	
5053-080	Round	5.50	9.50	
5053-082	Round	5.50	9.50	
5053-081	Round	6.50	12.50	
5053-083	Round	6.50	12.50	
5053-084	Round	5.30	24.50	
5053-085	Round	9.50	30.50	

Reciprocating Saw Blades

Cat. No.	Shank	Length (mm)	Material Thickness (mm)	
5053-122	Flat	25.50	0.25	
5053-220	Flat	25.50	0.40	
5053-221	Flat	38.00	0.40	
5053-054	Round	25.00	0.60	
5053-126	Round	28.00	0.60	

Depth Limiting Reciprocating Saw Blades

Thickness: 4.00mm					
Cat. No.	Shank	Length (mm)	Material Thickness (mm)	Cut Depth (mm)	
5053-051	Round	17.00	0.60	4.00	
Thickness: 5.00mm					
5053-050	Round	7.00	0.60	5.00	
5053-052	Round	23.00	0.60	5.00	
Thickness: 7.00mm					
5053-053	Round	25.00	0.60	7.00	

Products shown at actual size.
All products may not be available in all markets.

Micro Burs



	Round Fluted												
	Flutes (mm)	4.00	6.00	6.00	8.00	6.00	8.00	6.00	8.00	8.00	10.00	10.00	10.00
MicroPower & Micro 100 Bur Guards	Diameter (mm)	0.50	1.00	1.50	1.80	2.00	2.30	2.40	3.00	4.00	5.00	6.00	6.50
1375-012				5091-222	5091-221	5091-224	5091-223	5091-225	5091-226	5091-228	5091-230	M60R	5091-232
1375-011			5092-220	5092-222	L18R	5092-224	L23R		5092-226	5092-228	5092-230	L60R	5092-232
1375-023		5093-221	5093-220			5093-224			5093-226	5093-228	5093-230	X60R	5093-232
1387-022									5092-236				

	Fine Round Diamond						
MicroPower & Micro 100 Bur Guards	Diameter (mm)	1.00	2.00	3.00	4.00	5.00	6.00
1375-012			5091-260		5091-264	5091-268	M60D
1375-011		L10D	5092-260	5092-262	5092-264	5092-266	5092-268
1375-023		X10D	5093-260	5093-262	5093-264	5093-266	

	Wire Pass Burs			
	Flutes (mm)	2.00	2.00	2.00
	Diameter (mm)	1.00	1.50	2.00
MicroPower & Micro 100 Bur Guards	Length (mm)	19.00	19.00	19.00
1375-035		5090-247	5090-248	5090-249
1375-036		5090-247	5090-248	5090-249
1375-011		5092-247		
1375-032		5092-247		
1375-012		5091-247	5091-248	5091-249

	Swanson Reamers				
	Flutes (mm)	4.00	6.00	6.00	8.00
MicroPower & Micro 100 Bur Guards	Diameter (mm)	2.00	3.00	4.00	5.00
1375-012		5091-242	5091-243	5091-244	5091-245

	Dermabraders – Fine Diamond					N/A	N/A	N/A	N/A	
	Diameter (mm)	2.00	4.00	6.00	4.00	4.00	7.00	12.00	10.00	2.00
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)				8.00	15.00	15.00	15.00	15.00	19.00
1375-012		5091-260	5091-264	5091-268	5091-270	5091-080	5091-081	5091-082	5091-083	
1375-011					5092-270					

Bur Shank Diameter 2.34mm (0.82")


Products shown at actual size.

All products may not be available in all markets.

	Straight Side Cutting Burs														
		Flutes (mm)	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00	6.00
		Diameter (mm)	1.00	1.10	1.20	1.30	1.30	1.50	1.50	1.50	1.60	1.60	1.70	1.70	1.70
	MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	5.00	5.00	4.90	9.90	9.90	5.10	8.00	10.00	4.30	7.50	12.00	9.90	9.90
1375-012		5091-201	M10S	5091-216		M13S	5091-204		5091-202	5091-217	M15S			M17S	
1375-011					L13S		5092-204					5092-203	L17S		
1375-020												5092-203			
1375-035 & 1375-036								5089-208							

	Tapered Side-cutting						
		Diameter (mm)	0.80	1.20	1.20	1.20	1.60
		Bur Head Length (mm)	4.90	9.90	3.70	9.90	21.60
1375-012		5091-215	M12T	5091-219		M20T	5091-218
1375-011					L12T		

	Tapered Side-cutting					
		Diameter (mm)	1.60	2.00	2.00	2.00
		Bur Head Length (mm)	21.90	7.50	19.00	7.50
1375-012			5091-206	5091-214		
1375-011		L20T			5092-206	5092-214

	Barrel Burs		
		Flutes (mm)	6.00
		Diameter (mm)	2.00
MicroPower & Micro 100 Bur Guards		Bur Head Length (mm)	10.00
1375-012		5091-252	5091-254
1375-011		5092-252	5092-254


	Oval Burs					
		Flutes (mm)	8.00	8.00	8.00	8.00
		Diameter (mm)	4.00	5.50	4.00	5.50
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	8.00	10.0	7.9	10.0	
1375-012		5091-236	5091-238	M40B	M55B	
1375-011		5092-236	5092-238	L40B	L55B	
1375-023		5093-236	5093-238	X40B	X55B	
1387-022		5092-236				


Products shown at actual size.


All products may not be available in all markets.


Micro Burs






	Lindemann Burs	
	Flutes (mm)	4.00
	Diameter (mm)	1.60
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	31.00
1375-012		5091-212

	Cylindrical	
	Flutes (mm)	2.00
	Diameter (mm)	3.00
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	4.00
1375-011		5092-291

	Cone Burs		
	Flutes (mm)	8.00	12.00
	Diameter (mm)	4.50	4.50
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	12.00	12.00
1375-012		M45C	
1375-011			L45C

	Tapered Cross-cut Fissure Burs			
	Flutes (mm)	701 Type	702 Type	703 Type
	Diameter (mm)	1.20	1.60	2.10
MicroPower & Micro 100 Bur Guards	Bur Head Length (mm)	2.30	3.80	4.90
1375-012		5091-271	5091-272	5091-213
1375-011		5092-271	5092-272	5092-273

	Heavy Grit Diamond- Microbur. X-Long					
MicroPower & Micro 100 bur guards	Diameter (mm)	2.00	3.00	4.00	5.00	6.00
	Grit	60-70	60-70	40-45	40-45	40-45
1375-033 & 1375-023		5093-282	5093-283	5093-284	5093-285	
1375-032 & 1375-011			5092-283	5092-284	5092-285	5092-286

Products shown at actual size.

All products may not be available in all markets.

	Cat. No.
Hall 50™ Sternal Saw Handpiece.....	PRO7450B
MPOWER® 2 Battery Sternal Saw.....	PRO6450M
See pages 8, and 9 for battery accessories.	
Saw does not include sternum guard or blade.	

Pneumatic Sternal Saw.....5059-005

Sternal Saw Wrench, for PRO6450, PRO6450M & 5059-0055059-007

Sternal Guard, 2.25" (5.72mm), for PRO6450, PRO6450M
& 5059-005..... 5059-006

Sternal Guard, 3.12" (7.92mm), for PRO6450, PRO6450M
& 5059-005..... 5059-011

Repeat Sternotomy Saw, Pneumatic.....5057-040







All products may not be available in all markets.

Sternum Saw Blades

Indications for Use: Primary & Repeat Sternotomy





Sternum Saw Blades



Cat. No.	Width (mm)	Length (mm)	Non-sterile	
5059-235**		32.00		
5059-233*		34.00		
5059-532	10.00	35.00		
5059-032-01	10.00	35.00	50 Pack	
5059-531	10.00	54.00		
5059-031	10.00	54.00		
5059-031-01	10.00	54.00	50 Pack	

* For use with Sarns Sternum Saw

** For use with Stryker Sternum Saw

Repeat Sternotomy Blades




Cat. No.	Width (mm)	Length (mm)	
1384-175	52.00	19.00	
1384-176	63.00	19.00	




Cat. No.	Width (mm)	Length (mm)	
5023-187	50.00	24.00	
5023-188	63.00	24.00	

All products may not be available in all markets.

Products shown at actual size.

Repeat Sternotomy Blades

Cat. No.	Width (mm)	Length (mm)	
5023-178	44.60	40.00	
5023-179	54.60	40.00	
5023-180	65.80	40.00	

Cat. No.	Width (mm)	Length (mm)	
5059-278	44.60	40.00	
5059-279	54.60	40.00	
5059-280	65.80	40.00	

Products shown at actual size.

Sternum Saw Blades

Indications for Use: Primary & Repeat Sternotomy



Repeat Sternotomy Blades

Cat. No.	Width (mm)	Length (mm)	
5071-125	50.00	19.00	
5071-126	63.00	19.00	
5071-228	44.50	40.00	

Cat. No.	Width (mm)	Length (mm)	
5071-229	54.50	40.00	
5071-230	66.00	40.00	

Products shown at actual size.

All products may not be available in all markets.

Handpieces

	Cat. No.
High Speed Drill.....	7020-001
Low Speed Drill	7020-002
Perforator Drive.....	7020-003

Attachments & Accessories

Ultrapower Footswitch	7020-045
Short Bur Guard (like Midas-X) Red.....	7020-010
Medium Bur Guard (like Midas-O) Copper	7020-011
Long Bur Guard (like Midas-M/B) Mustard	7020-012
X-Long Bur Guard (like Midas-AM/A) Lime	7020-013
Short Revision Extension (like Midas TU) Teal.....	7020-014
Medium Revision Extension (like Midas-R) Purple	7020-015
Long Revision Extension (like Midas RX) Blue.....	7020-016
Short Bur Guard (like Midas B-5) Black.....	7020-017
Medium Bur Guard (like Midas B-1) Green.....	7020-018
Long Bur Guard (like Midas S-1) Brown	7020-019
90 degree Angle Attachment (like Midas AF/BF)..... Uses the 5089 and 5090 Series Burs Only	7020-022
20 degree Long Keyless Attachment	7020-023
Uses the 5092 Series Burs Only.	
20 degree Angle X-Long Keyless Attachment	7020-024
Uses the 5093 Series Burs Only.	



7020-001




7020-002



7020-003



7020-045



All products may not be available in all markets.

UltraPower® Burs

Indications for Use: Otology
(Round Burs Only), Spine & Neuro



UltraPower Bur Guard No.	Round Fluted									
	Flutes (mm)	4.00	6.00	6.00	6.00	8.00	8.00	8.00	8.00	8.00
UltraPower Bur Guard No.	Diameter (mm)	0.50	1.00	1.50	2.00	2.50	3.00	4.00	5.00	6.00
		7020-010			7021-024	7021-025	7021-026	7021-027	7021-028	7021-029
UltraPower Bur Guard No.	Diameter (mm)	7020-011		7021-122	7021-123	7021-124	7021-125	7021-126	7021-127	7021-128
UltraPower Bur Guard No.	Diameter (mm)	7020-012	7021-220	7021-222	7021-223	7021-224	7021-225	7021-226	7021-227	7021-228
UltraPower Bur Guard No.	Diameter (mm)	7020-013		7021-322	7021-323	7021-324	7021-325	7021-326	7021-327	7021-328
UltraPower Bur Guard No.	Diameter (mm)	7020-014					7021-426	7021-427	7021-428	7021-429
UltraPower Bur Guard No.	Diameter (mm)	7020-015					7021-526		7021-528	7021-529

UltraPower Bur Guard No.	Round Fluted					
	Flutes (mm)	10	10	10	10	10
UltraPower Bur Guard No.	Diameter (mm)	6.5	7.5	8.0	9.0	10.0
		7020-010	7021-030	7021-031	7021-032	7021-033
UltraPower Bur Guard No.	Diameter (mm)	7020-011	7021-130	7021-131	7021-132	7021-133
UltraPower Bur Guard No.	Diameter (mm)	7020-012	7021-230	7021-231	7021-232	7021-233
UltraPower Bur Guard No.	Diameter (mm)	7020-013		7021-331	7021-332	7021-333
UltraPower Bur Guard No.	Diameter (mm)	7020-014		7021-431	7021-432	7021-433
UltraPower Bur Guard No.	Diameter (mm)	7020-015		7021-531	7021-532	7021-533
UltraPower Bur Guard No.	Diameter (mm)	7020-016		7021-631	7021-632	7021-633

UltraPower Bur Guard No.	Fine Round Diamond							
	Diameter (mm)	1.0	2.0	3.0	4.0	5.0	6.0	7.0
UltraPower Bur Guard No.	Diameter (mm)	7020-010		7021-062	7021-063	7021-064	7021-065	7021-066
UltraPower Bur Guard No.	Diameter (mm)	7020-011	7021-161	7021-162	7021-163	7021-164	7021-165	7021-166
UltraPower Bur Guard No.	Diameter (mm)	7020-012	7021-261	7021-262	7021-263	7021-264	7021-265	7021-266
UltraPower Bur Guard No.	Diameter (mm)	7020-013	7021-361	7021-362	7021-363	7021-364	7021-365	7021-366

Products shown at actual size.

All products may not be available in all markets.

	Coarse Round Diamond							
UltraPower Bur Guard No.	Diameter (mm)	1.00	2.00	3.00	4.00	5.00	6.00	7.00
7020-010					7021-074	7021-075		7021-077
7020-011		7021-171	7021-172	7021-173	7021-174	7021-175	7021-176	7021-177
7020-012			7021-272	7021-273	7021-274	7021-275	7021-276	7021-277
7020-013		7021-371	7021-372	7021-373	7021-374	7021-375	7021-376	7021-377
7020-014					7021-474			





	Side Cutting, Burs					
	Flutes (mm)	8.00	8.00	8.00	8.00	8.00
	Diameter (mm)	5.00	5.00	7.50	7.50	7.50
UltraPower Bur Guard No.	Bur Head Length (mm)	8.00	17.00	5.00	9.00	28.00
7020-012		7021-202			7021-206	7021-207
7020-013		7021-302	7021-303	7021-305		
7020-014						7021-427


	Side Cutting Burs			
	Flutes (mm)	8.00	8.00	12.00
	Diameter (mm)	9.00	9.00	10.50
UltraPower Bur Guard No.	Bur Head Length (mm)	5.00	9.00	11.00
7020-012				7021-210
7020-013			7021-309	7021-310
7020-014		7021-408	7021-409	
7020-015			7021-509	


UltraPower® Burs

Indications for Use: Otology
(Round Burs Only), Spine & Neuro




Ultrapower Bur Guard No.	Side Cutting Burs									
	Flutes (mm)		1.00	1.00	1.00	2.00	2.00	2.00	2.00	2.00
	Diameter (mm)		1.50	1.8	2.00	1.20	2.20	2.20	3.00	3.00
Bur Head Length (mm)		8.00	11.00	16.0	6.00	16.00	23.00	20.00	25.00	
7020-010		7021-012	7021-013		7021-043	7021-015				
7020-011		7021-112	7021-113	7021-114	7021-143	7021-115	7021-116			
7020-012			7021-713	7021-214	7021-243	7021-815				7021-917
7020-013				7021-314		7021-315		7021-301		

Ultrapower Bur Guard No.	Side Cutting Burs			
	Flutes (mm)		2.00	6.00
	Diameter (mm)		3.00	1.50
Bur Head Length (mm)		36.00	4.90	
7020-011				7021-119
7020-012			7021-218	
7020-013			7021-318	


Ultrapower Bur Guard No.	Fluted Cylindrical Burs					
	Flutes (mm)		2.00	2.00	2.00	2.00
	Diameter (mm)		1.80	1.80	2.20	3.00
7020-010					7021-036	7021-037
7020-011					7021-136	7021-137
7020-012				7021-235	7021-236	7021-237
7020-013			7021-338	7021-335	7021-336	7021-337
7020-014						7021-437


Products shown at actual size.

All products may not be available in all markets.

	Fine Diamond Cylindrical Burs		
Ultrapower Bur Guard No.	Diameter (mm)	1.90	3.00
7020-012			7021-279
7020-013		7021-378	7021-379

	Fluted Acorn Burs				
	Flutes (mm)	8.00	8.00	8.00	8.00
Ultrapower Bur Guard No.	Diameter (mm)	5.00	6.00	7.50	9.00
7020-010			7021-081		7021-083
7020-011			7021-181	7021-182	7021-183
7020-012		7021-280	7021-281	7021-282	7021-283
7020-013		7021-380	7021-381	7021-382	7021-383
7020-014				7021-482	7021-483
7020-015			7021-581		7021-583
7020-016			7021-681		7021-683

	Trephine Burs				
Ultrapower Bur Guard No.	Diameter (mm)	3.20	5.30	10.70	13.50
7020-010			7021-051		
7020-012		7021-250		7021-252	7021-253

	Trephine Burs	
Ultrapower Bur Guard No.	Diameter (mm)	16.80
7020-012		7021-255

Products shown at actual size.

All products may not be available in all markets.

UltraPower® Burs

Indications for Use: Otology
(Round Burs Only), Spine & Neuro





	Wire Pass Burs			
	Flutes (mm)	2.00	2.00	2.00
UltraPower Bur Guard No.	Diameter (mm)	1.00	1.50	2.00
7020-010		7021-040	7021-041	7021-042
7020-011		7021-140	7021-141	7021-142


	Oval Fluted Burs		
	Flutes (mm)	8.00	8.00
	Diameter (mm)	4.00	5.50
UltraPower Bur Guard No.	Bur Head Length (mm)	8.00	10.00
7020-010		7021-095	7021-096
7020-011		7021-195	7021-196
7020-012		7021-295	7021-296
7020-013		7021-395	7021-396

	Fluted Frustrum Bur	
	Flutes (mm)	8.00
UltraPower Bur Guard No.	Diameter (mm)	6.50
7020-013		7021-391

	Carbide Metal Cutting Burs	
	Diameter (mm)	3.20
UltraPower Bur Guard No.	Bur Head Length (mm)	15.20
7020-010		7021-094
7021-012		7021-294

	Wire Pass Burs – Depth Limiting				
	Diameter (mm)	1.10	1.10	1.10	1.50
UltraPower Bur Guard No.	Depth (mm)	3.50	5.00	9.00	7.00
7020-011		7021-141-35	7021-141-50	7021-141-90	7021-145-70

	Carbide Metal Cutting Burs	
	Flutes (mm)	8.00
	Diameter (mm)	6.00
UltraPower Bur Guard No.	Bur Head Length (mm)	4.00
7020-012		7021-297

	Fine Diamond Wheels for Metal Cutting	
UltraPower Bur Guard No.	Diameter (mm)	25.40
7020-010		7021-092
7020-012		7021-292
7020-013		7021-392
7020-014		7021-492

Products shown at actual size.

All products may not be available in all markets.

Handpieces

CoolFlex® High Speed DrillE9010

High Speed Shaver E9005

Cranial Perforator Drive.....E9015

Attachments & Accessories

Medium Bur Guard for E9010 DrillE9011-1

Long Bur Guard for E9010 DrillE9011-2

X-Long Bur Guard for E9010 Drill.....E9011-3

Short Neuro Guard for E9010 Drill.....E9011-7

Medium Neuro Guard for E9010 Drill.....E9011-8

Long Neuro Guard for E9010 DrillE9011-9

Locking Collar for E9010 DrillE9012

Bur Rack, Holds up to 20 burs for E9010 DrillE9323

Drill Tubing for E9010 Drill without irrigation tipE9414

Drill Tubing for E9010 Drill with irrigation tip.....E9415A

4.2mm FESS Irrigation Tubing Set for E9005 ShaverE9416

3.7mm FESS Irrigation Tubing Set for E9005 Shaver (Not Shown)E9417

DeClogger for removing tissue and bone from shaver blades and bursE9315

Bur Cabinet for Neuro Burs.....7020-046



E9010



E9005



E9015



E9011-1



E9011-2



E9011-3



E9011-7



E9011-8



E9011-9



E9012



E9323



E9414



E9415A



E9416



E9315



7020-046

All products may not be available in all markets.

E9005 ENT Shaver Blades

Indications for Use: Sinus Surgery



Ideal for aggressive removal of ethmoid bone, cartilage and soft tissue



Ideal for aggressive removal of ethmoid bone, cartilage and soft tissue



Ideal for removal of soft tissue

Products not shown at actual size.

All products may not be available in all markets.

For use with:

E9005 High Speed Shaver

MicroXcisor® Blade

Cat. No.

4.2mm, Irrigated, 10.5cm Working Length, 6/Box N9330

Lynx® Blade

4.2mm, Irrigated, 10.5cm Working Length, 6/Box N9150

RhinoTec® Cuda® Blade

4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box..... E9350

3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box..... E9351

4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box..... A9350

4.2mm, Irrigated, 10.5cm Working Length, 6/Box N9210

3.5mm, 7.6cm Working Length, 6/Box..... C9952

2.9mm, 7.6cm Working Length, 6/Box..... C9951

2.0mm, 7.6cm Working Length, 6/Box..... C9950

Rhinotec® Gator® Blade

Cat. No.

4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box.....	E9360
3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box.....	E9361
4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box.....	A9360
4.2mm, Irrigated, 10.5cm Working Length, 6/Box.....	N9220
3.5mm, 7.6cm Working Length, 6/Box.....	C9962
2.9mm, 7.6cm Working Length, 6/Box.....	C9961
2.0mm, 7.6cm Working Length, 6/Box.....	C9960

Head & Neck Blade

4.2mm, 11.5cm Working Length, 6/Box.....	N9321
3.5mm, 7.0cm Working Length, 6/Box.....	N9320

Laryngeal Blade

3.7mm, Bendable to 45°, 20.0cm Working Length, 6/Box.....	E9363
---	-------



Ideal for aggressive removal of soft tissue, cartilage and thin bone



Ideal for the most gentle removal of soft tissue



Ideal for removal of soft tissue in laryngeal procedures

Products not shown at actual size.
All products may not be available in all markets.

E9005 ENT Shaver Blades


Indications for Use: Sinus Surgery



Ideal for aggressive removal of hard bone and cartilage



45° fixed curve (pre-bent)
Ideal for aggressive removal of hard bone



Ideal for less aggressive removal of hard bone Ideal for rhinoplasty

Products not shown at actual size.
All products may not be available in all markets.

Rhinotec® Vortex® Router

Cat. No.

- 3.5mm, Hooded, 7.6cm Working Length, 6/Box C9991
- 2.9mm, Hooded, 7.6cm Working Length, 6/Box C9990

Curved Vortex Router

- 4.5mm, 45° Pre-Bent, 11.5cm Working Length, 6/Box..... N9375

Otolaryngology Bur

- 4.0mm, Hooded, 10.8cm Working Length, 6/Box N9390

Rhinotec® Spherical Bur, Hooded

Cat. No.

4.5mm, Hooded, 10.8cm Working Length, 6/Box N9380

3.5mm, Hooded, 7.6cm Working Length, 6/Box C9912

2.9mm, Hooded, 7.6cm Working Length, 6/ Box C9911

Rhinotec® Spherical Bur, Unhooded

3.5mm, Unhooded, 7.6cm Working Length, 6/Box..... N9382

Rhinotec Spherical Diamond Bur, Hooded

3.5mm, Hooded, 10.8cm Working Length, 6/Box N9385



Ideal for aggressive removal of hard bone and cartilage



Ideal for less aggressive removal of hard bone and cartilage



Ideal for less aggressive removal of hard bone and cartilage


Products not shown at actual size.
All products may not be available in all markets.

E Burs

Indications for Use: Otology, Spine & Neuro



*These burs may be used with multiple bur guard options to change the working length of the bur.



	Round Fluted						
	Flutes (mm)	4.00	4.00	6.00	6.00	6.00	8.00
E9010 Bur Guard No.	Diameter (mm)	0.50	0.70	1.00	1.50	2.00	2.50
E9011-1 or E9011-2			E8007	E8010	E8015	E8020	E8025
E9011-1 or E9011-2		E8105		E8110	E8115	E8120	E8125
E9011-1 or E9011-2		E8205	E8207	E8210	E8215	E8220	E8225

	LSO Round Fluted				
	Flutes (mm)	8.00	8.00	8.00	10.00
Bur Guard No.	Diameter (mm)	3.00	4.00	5.00	6.50
E9011-1 or E9011-2		E8030LSO	E8040LSO	E8050LSO	E8065LSO
E9011-1 or E9011-2		E8130LSO	E8140LSO	E8150LSO	E8165LSO
E9011-1 or E9011-2		E8230LSO	E8240LSO		E8265LSO




	Fine Round Diamond										
E9010 Bur Guard No.	Diameter (mm)	0.50	0.70	1.00	1.50	2.00	2.50	3.00	4.00	5.00	6.00
E9011-1 or E9011-2									E6040		
E9011-1 or E9011-2		E6105		E6110	E6115	E6120	E6125	E6130	E6140	E6150	E6160
E9011-1 or E9011-2		E6205	E6207	E6210	E6215	E6220	E6225	E6230	E6240	E6250	


	Coarse Round Diamond									
E9010 Bur Guard No.	Diameter (mm)	1.00	1.50	2.00	2.50	3.00	4.00	5.00	6.00	7.00
E9011-1 or E9011-2				E5120		E5160	E5140	E5150	E5130	E5170
E9011-1 or E9011-2		E5210	E5215	E5220	E5225	E5230	E5240	E5050		


Products shown at actual size.


All products may not be available in all markets.



*These burs may be used with multiple bur guard options to change the working length of the bur.

	Wirepass burs			
	Diameter (mm)	1.00	1.50	2.00
E9010 Bur Guard No	Length (mm)	38.00	38.00	38.00
E9011-1 or E9011-2		E4110	E4115	E4120

	Neuro Burs			
	Flutes (mm)	1.00	2.00	2.00
	Diameter (mm)	1.80	2.20	3.00
E9010 Bur Guard No.	Bur Head Length (mm)	11.00	16.00	25.00
E9011-7 or 7020-017		7021-713		
E9011-8 or 7020-018			7021-815	
7020-019				7021-917

	Helical Bur for E9010 Drill	
E9010 Bur Guard No.	Bur Head Length (mm)	25.00
E9011-9, E9011-2 and 7020-012		7021-960

	E9010 Lindemann Bur	
	Diameter (mm)	1.60
E9010 Bur Guard No.	Bur Head Length (mm)	31.00
E9011-1		E4217

	E9010 Side Cutting Burs		
	Diameter (mm)	0.80	2.00
E9010 Bur Guard No.	Bur Head Length (mm)	4.90	7.50
E9011-1		E4208	E4220

5052-012



5038-009



5038-001



5038-002



5038-004



5038-005



5038-600



Attachments & Accessories

Cat. No.

Standard Sterilization Case	5052-012
Bur Rack, for Osteon™ and Ototom Burs.....	5038-009
Osteon/Ototome Bur Set, includes 16 burs (see System Listings Page 26)	5056-109

Handpieces

Osteon Handpiece, Straight.....	5038-001
---------------------------------	----------

Osteon Handpiece, Angled.....	5038-002
-------------------------------	----------

Attachments & Accessories

Osteon Footswitch with Hose	5038-004
Osteon Hose	5038-005
Osteon Tubing Set with irrigation tip	5038-600

	Fluted Round Burs											
		Flutes (mm)	4.00	4.00	6.00	6.00	6.00	6.00	8.00	8.00	8.00	10.00
Working Length (mm)	Diameter (mm)	0.50	0.70	1.00	1.50	2.00	2.50	2.50	3.00	4.00	5.00	6.50
30			5056-173	5056-174	5056-177	5056-178	5056-179		5056-116	5056-119	5056-118	5056-117
38									5056-128	5056-127	5056-126	5056-125
48		5056-166	5056-165	5056-164	5056-140	5056-142		5056-149	5056-146	5056-145	5056-144	

	Fine Diamond Round Burs										
	Diameter (mm)		0.50	0.70	1.00	1.50	2.00	3.00	4.00	5.00	6.00
Working Length (mm)	Diameter (mm)										
30									5056-168	5056-147	5056-167
38									5056-148	5056-150	
48			5056-161	5056-160	5056-158	5056-163	5056-162	5056-159	5056-152		

	Coarse Diamond Round Burs							
	Diameter (mm)		1.00	2.00	3.00	4.00	5.00	6.00
Working Length (mm)	Diameter (mm)							
48			5056-958	5056-962	5056-959	5056-952	5056-950	5056-967

	Deep Access Burs						
	Flutes (mm)	Diameter (mm)	10.00	12.00	6.00		
Working Length (mm)	Diamond					Fine	Fine
	Diameter (mm)		3.00	4.00	5.50	3.00	4.00
90			5056-170	5056-171		5056-175	5056-176
95					5056-172		

	Wire Pass Burs				
	Diameter (mm)		1.00	1.50	2.00
Working Length (mm)	Bur Head Length (mm)		19.00	19.00	19.00
38			5056-184	5056-185	5056-186

	80 Grit Dermabrasion Burs					
	Diameter (mm)		4.00	7.00	12.00	10.00
Working Length (mm)	Bur Head Length (mm)		15.00	15.00	15.00	15.00
30			5056-080	5056-081	5056-082	5056-083

Products shown at actual size.
All products may not be available in all markets.



Handpieces

Cat. No.

Cebotome® Revision Drill Handpiece 5052-041

Attachments & Accessories

Cebotome Bur Wrench.....5052-046

Cebotome Metal Cutting Bur Wrench 5052-143

Standard Bur Extension, Short, use with standard Cebotome Burs5052-042



Standard Bur Extension, Medium, use with standard Cebotome Burs5052-048

Standard Bur Extension, Long, use with standard Cebotome Burs.....5052-049



Metal Cutting Bur Extension, use with metal
cutting wheels and diamond metal cutting burs..... 5052-142



Universal Hose, 10' (3.05m)..... 5052-010

All products may not be available in all markets.

Cebotome® Burs								
Flutes (mm)	2.00	2.00	2.00	10.00	10.00	10.00	6.00	6.00
Diameter (mm)	5.00	6.00	7.00	7.00	8.00	9.00	8.50	10.00
	5052-181	5052-182	5052-183	5052-184	5052-185	5052-186	5052-187	5052-188

Cebotome Burs				
Flutes (mm)	6.00	10.00	6.00	8.00
Diameter (mm)	5.00	10.00	7.00	8.50
	5052-189	5052-190	5052-191	5052-192

Cebotome Metal Cutting Wheels		
Grit (mm)	150.00	150.00
Diameter (mm)	25.40	25.40
	5052-125	5052-126

Cebotome Diamond Metal Cutting Burs		
Diameter (mm)	3.2	3.2
Diamond	Flat	Round
	5052-144	5052-145

Wires & Staples


Indication for Use: Small Bone Trauma



C-Wires

Cat. No.	Tip	Diameter	Length	Quantity
5050-041	Spade	.028" (.70mm)	5" (13.00cm)	50 per Pkg.
5050-042	Spade	.035" (.90mm)	5" (13.00cm)	50 per Pkg.
5050-043	Spade	.045" (1.10mm)	5" (13.00cm)	50 per Pkg.
5050-044	Spade	.062" (1.60mm)	5" (13.00cm)	50 per Pkg.
5050-046	Trocar	.028" (.70mm)	5" (13.00cm)	50 per Pkg.
5050-047	Trocar	.035" (.90mm)	5" (13.00cm)	50 per Pkg.
5050-048	Trocar	.045" (1.10mm)	5" (13.00cm)	50 per Pkg.
5050-049	Trocar	.062" (1.60mm)	5" (13.00cm)	50 per Pkg.


Note: To order C-Wires in a box of 5, instead of 50 per package, please order using catalog number 5050-2XX. For example, 5051-241 instead of 5050-041.



Staplizer® Cartridges

Cat. No.	Width (mm)	Length (mm)
7707-03	7.00	7.00
7710-03	7.00	10.00
7007-03	10.00	7.00
7010-03	10.00	10.00
7015-03	10.00	15.00
7307-03	13.00	7.00
7310-03	13.00	10.00
7315-03	13.00	15.00
7320-03	13.00	20.00
7610-03	16.00	10.00
7615-03	16.00	15.00
7620-03	16.00	20.00
7625-03	16.00	25.00

Note: 3 cartridges per Pkg.
5 Staples per cartridge.



Products shown at actual size.

All products may not be available in all markets.

Sterilization Containers

	Cat. No.
Standard Sterilization Case	5052-012



5052-012

Holding Clamps for use with 5030-010 & 5030-020

Small	5030-028
Medium.....	5030-029



5030-029

5030-028

Bur Cleaning Brush.....	1375-015
-------------------------	----------



1375-015

Cleaning Brush for Cannulated Handpieces	5053-124
--	----------



5053-124

Electric Consoles & Accessories



D4000



D3000



E9000



PRO2000

Cat. No.

D4000 Console D4000

D4000 Console with Irrigation D4000A

Advantage® Console D3000

E9000® Console E9000

PowerPro® Console PRO2000

2-Pedal Footswitch.....Cat. No. 5020-053

3-Pedal Footswitch.....C9863

Handpiece Cord.....MC5057

AC Power Cord.....5039-003

Fluid Bag Support Rod for D3000 & E9000 5040-180

Console Stand.....5020-054



5020-053



C9863



MC5057



5039-003



5040-180



5020-054

All products may not be available in all markets.

Pneumatic Regulators, Hoses & Accessories



1375-028



Universal Hose



Extension Hose



Hose for Staplizer®



M207



M208



M322



5144-002



1375-040



M105A



M317



1365-016



1375-037

Hall® offers nitrogen specific products. Please contact customer service for information.

Cat. No.

Micro 100™ Footswitch..... 1375-028

Universal Hose, 10' (3.05m) for Hall Handpieces 5052-010

Extension Hose, 10' (3.05m), Male Schrader to Female Schrader..... 5052-018

Extension Hose, 10' (3.05m), Male Diss and Female Schrader 5052-029

10' (3.05m) Hose for Staplizer® Handpiece A201

Regulator, Single Connection.....M207

Regulator, Dual Connection.....M208

DISS to Schrader Adapter.....M322

Schrader Coupler..... 5144-002

Y Connector 1375-040

Blitz® Instrument Cleaner (Use with Staplizer Handpiece)M105A

Instrument Lubricant for Staplizer.....M317

Hall Lubricant for Angled Attachments..... 1365-016

Hall Lubricant for Surgical Angled Attachments 1375-037

All products may not be available in all markets.

Surgairtome Two® System

Components Ordered Separately

	Cat. No.
Surgairtome Two Handpiece	5058-001
Medium Bur Guard	1375-012
Standard Sterilization Case	5052-012
Universal Hose, 10' (3.05m).....	5052-010
Bur Cleaning Brush.....	1375-015

Ultrapower® Complete Drill System Product Listing

High Speed Drill.....	7020-001
Low Speed Drill	7020-002
Perforator Drive	7020-003
Universal Hose, 10' (3.05m).....	5052-010
Short Bur Guard like Midas-X, Red.....	7020-010
Medium Bur Guard like Midas-O, Copper	7020-011
Long Bur Guard like Midas-M/B, Mustard	7020-012
X-Long Bur Guard like Midas-AM/A, Lime	7020-013
Short Revision Extension like Midas TU, Teal	7020-014
Medium Revision Extension like Midas-R, Purple	7020-015
Long Revision Extension like Midas RX, Blue	7020-016
Short Bur Guard like Midas B-5, Black	7020-017
Medium Bur Guard like Midas B-1, Green.....	7020-018
Long Bur Guard like Midas S-1, Brown	7020-019
90 degree Angle Attachment like Midas AF/BF.....	7020-022
20 degree Long Keyless Attachment	7020-023
20 degree Angle X-Long Keyless Attachment	7020-024
Key Wrench, 2 per box	7020-040
Footswitch.....	7020-045
Storage Cabinet for Neuro Burs.....	7020-046
Attachment Spray	7020-041



All products may not be available in all markets.



Ultrapower® Neuro/Spine Drill System Product Listing

Components Ordered Separately

	Cat. No.
High Speed Drill.....	7020-001
Low Speed Drill	7020-002
Perforator Drive	7020-003
Universal Hose, 10' (3.05m).....	5052-010
Short Bur Guard like Midas-X, Red.....	7020-010
Medium Bur Guard like Midas-O, Copper	7020-011
Long Bur Guard like Midas-M/B, Mustard	7020-012
X-Long Bur Guard like Midas-AM/A, Lime	7020-013
Short Bur Guard like Midas B-5, Black.....	7020-017
Medium Bur Guard like Midas B-1, Green.....	7020-018
Long Bur Guard like Midas S-1, Brown	7020-019
90 degree Angle Attachment like Midas AF/BF.....	7020-022
20 degree Long Keyless Attachment	7020-023
20 degree Angle X-Long Keyless Attachment	7020-024
Key Wrench, 2 per box	7020-040
Footswitch.....	7020-045
Storage Cabinet for Neuro Burs.....	7020-046
Attachment Spray	7020-041

Ultrapower Basic Drill System Product Listing

Components Ordered Separately

High Speed Drill.....	7020-001
Universal Hose, 10' (3.05m).....	5052-010
Short Bur Guard like Midas-X, Red.....	7020-010
Medium Bur Guard like Midas-O, Copper	7020-011
Long Bur Guard like Midas-M/B, Mustard	7020-012
X-Long Bur Guard like Midas-AM/A, Lime	7020-013
Short Revision Extension like Midas TU, Teal.....	7020-014
Medium Revision Extension like Midas-R, Purple	7020-015
Long Revision Extension like Midas RX, Blue.....	7020-016
Short Bur Guard like Midas B-5, Black.....	7020-017
Medium Bur Guard like Midas B-1, Green.....	7020-018
Long Bur Guard like Midas S-1, Brown	7020-019
90 degree Angle Attachment like Midas AF/BF.....	7020-022
20 degree Long Keyless Attachment	7020-023
20 degree Angle X-Long Keyless Attachment	7020-024
Key Wrench, 2 per box	7020-040
Footswitch.....	7020-045
Storage Cabinet for Neuro Burs.....	7020-046
Attachment Spray	7020-041

Osteon™/Ototome Burs

Components Ordered Separately

	Cat. No.
Round Carbide Bur, 6.5 x 30 mm, 2 each	5056-117
Round Carbide Bur, 5.0 x 30mm.....	5056-118
Round Carbide Bur, 4.0 x 30mm.....	5056-119
Round Carbide Bur, 5.0 x 38mm.....	5056-126
Round Steel Bur, 1.5 x 48mm.....	5056-140
Round Steel Bur, 2.0 x 48mm.....	5056-142
Round Carbide Bur, 4.0 x 48mm.....	5056-145
Round Carbide Bur, 3.0 x 48mm.....	5056-146
Round Diamond Bur, 4.0 x 48mm	5056-152
Round Diamond Bur, 1.0 x 48mm	5056-158
Round Diamond Bur, 3.0 x 48mm	5056-159
Round Diamond Bur, 0.5 x 48mm	5056-161
Round Diamond Bur, 2.0 x 48mm	5056-162
Round Steel Bur, 1.0 x 48mm	5056-164
Round Steel Bur, 0.5 x 48mm.....	5056-166

Osteon Drill System with Burs

Components Ordered Separately

Bur Cleaning Brush.....	1375-015
Osteon Footswitch with Universal Hose 5052-010	5038-004
Osteon Hose	5038-005
Osteon/Ototome Bur Rack	5038-009
Osteon Handpiece, Straight.....	5038-001
Osteon Handpiece, Angled.....	5038-002



All products may not be available in all markets.



Osteon™ Drill System without Bur System

Components Ordered Separately

	Cat. No.
Bur Cleaning Brush.....	1375-015
Osteon Footswitch with Universal Hose 5052-010.....	5038-004
Osteon Hose	5038-005
Osteon/Ototome Bur Rack	5038-009
Osteon Handpiece, Straight.....	5038-001
Osteon Handpiece, Angled.....	5038-002

Revision Drill System

Components Ordered Separately

Universal Hose, 10' (3.05m).....	5052-010
Cebotome Revision Drill Handpiece	5052-041
Cebotome Bur Extension, Short.....	5052-042
Cebotome Wrench.....	5052-046
Cebotome Bur Extension, Medium.....	5052-048

Cebotome Metal Cutting Wheel System

Components Ordered Separately

Cebotome Metal Cutting Wheel.....	5052-125
Cebotome Metal Trimming Wheel.....	5052-126
Cebotome Metal Cutting Wheel Wrench	5052-129

Cebotome Metal Cutting Bur System

Components Ordered Separately

Cebotome Metal Cutting Bur Extension	5052-142
Cebotome Metal Cutting Bur Wrench	5052-143
ST Cebotome Metal Cutting Bur, Side, Flat.....	5052-144
ST Cebotome Metal Cutting Bur, Side, Round.....	5052-145

HALL®
**POWERED
INSTRUMENTS
INDEX**

All products may not be available in all markets.

Cat. No.	Description	Page No.
1365-005	Replacement Hex Wrench for Micro 100 Sagittal & Oscillating Saw	42
1365-016	Hall Lubricant for Angled Attachments	82
1375-003	Bur Changer for 70 and 90 degree Angled Attachments	43
1375-011	Long Bur Guard	43, 54 - 56
1375-011P	Long Precision Bur Guard	43
1375-012	Medium Bur Guard	43, 54 - 56, 83
1375-012P	Medium Precision Bur Guard	43
1375-015	Bur Cleaning Brush	79, 83, 85, 86
1375-020	Laminectomy Guard	43, 55
1375-021	Rhinoplasty Guard	43
1375-023	X-Long Bur Guard	43, 54, 55, 56
1375-023P	X-Long Precision Bur Guard	43
1375-028	Micro 100 Footswitch	82
1375-032	20 degree Angled Attachment	43, 54, 56
1375-033	20 degree X-Long Angled Attachment	43, 56
1375-034	70 degree Contra Angle with Dental Head	43
1375-035	70 degree Contra Angle with Surgical Head	43, 54, 55
1375-036	90 degree Angled Attachment	43, 54, 55
1375-037	Hall Lubricant for Surgical Angled Attachments	82
1375-040	Y Connector	82
1375-112	Burguard Sleeve Surgairtome Two	42
1384-175	Repeat Sternotomy Blade	58
1384-176	Repeat Sternotomy Blade	58
1387-022	Tissue Retractor	43, 54, 55
5020-053	2-Pedal Footswitch	81
5020-054	Console Stand	81
5020-160	Smartguard	42
5023-055	Micro Dermabrader Sagittal Saw Blade	49
5023-056	Micro Dermabrader Sagittal Saw Blade	49
5023-057	Micro Dermabrader Sagittal Saw Blade	49
5023-101	Micro Intra-Oral Oscillating Saw Blade	51
5023-102	Micro Intra-Oral Oscillating Saw Blade	51
5023-103	Micro Intra-Oral Oscillating Saw Blade	51
5023-104	Micro Intra-Oral Oscillating Saw Blade	51
5023-118	Micro Sagittal Saw Blade	45
5023-128	Micro Angled Sagittal Saw Blade	47
5023-129	Micro Angled Sagittal Saw Blade	47
5023-130	Micro Sagittal Saw Blade	45
5023-131	Micro Sagittal Saw Blade	45
5023-132	Micro Sagittal Saw Blade	45
5023-133	Micro Sagittal Saw Blade	45
5023-134	Micro Angled Sagittal Saw Blade	47
5023-135	Micro Angled Sagittal Saw Blade	47
5023-136	Micro Angled Sagittal Saw Blade	47
5023-137	Micro Sagittal Saw Blade	45
5023-138	Micro Sagittal Saw Blade	45
5023-139	Micro Sagittal Saw Blade	45
5023-140	Micro Sagittal Saw Blade	45
5023-141	Micro Sagittal Saw Blade	46
5023-142	Micro Sagittal Saw Blade	46
5023-143	Micro Sagittal Saw Blade	46
5023-144	Micro Sagittal Saw Blade	45
5023-160	Micro Oscillating Saw Blade	50
5023-161	Micro Oscillating Saw Blade	50
5023-162	Micro Oscillating Saw Blade	50

Cat. No.	Description	Page No.
5023-163	Micro Oscillating Saw Blade	50
5023-164	Micro Oscillating Saw Blade	50
5023-165	Micro Oscillating Saw Blade	50
5023-166	Micro Oscillating Saw Blade	50
5023-167	Micro Angled Oscillating Saw Blade	50
5023-168	Micro Angled Oscillating Saw Blade	50
5023-169	Micro Oscillating Saw Blade	50
5023-170	Micro Crescentic Oscillating Saw Blade	52
5023-171	Micro Crescentic Oscillating Saw Blade	52
5023-172	Micro Crescentic Oscillating Saw Blade	52
5023-176	Micro Crescentic Oscillating Saw Blade	52
5023-178	Repeat Sternotomy Blade	59
5023-179	Repeat Sternotomy Blade	59
5023-180	Repeat Sternotomy Blade	59
5023-182	Micro ACL Sagittal Saw Blade	48
5023-183	Micro ACL Sagittal Saw Blade	48
5023-184	Micro ACL Sagittal Saw Blade	48
5023-185	Micro ACL Sagittal Saw Blade	48
5023-186	Micro ACL Sagittal Saw Blade	48
5023-187	Repeat Sternotomy Blade	58
5023-188	Repeat Sternotomy Blade	58
5023-244	Micro Sagittal Saw Blade	46
5023-247	Micro Sagittal Saw Blade	46
5023-249	Micro Sagittal Saw Blade	46
5023-258	Micro Angled Oscillating Saw Blade	50
5023-259	Micro Angled Oscillating Saw Blade	50
5023-405	Micro Graft Harvesting Sagittal Saw Blade	49
5023-406	Micro Graft Harvesting Sagittal Saw Blade	49
5023-407	Micro Graft Harvesting Sagittal Saw Blade	49
5023-408	Micro Graft Harvesting Sagittal Saw Blade	49
5023-409	Micro Graft Harvesting Sagittal Saw Blade	49
5023-410	Micro Graft Harvesting Sagittal Saw Blade	49
5023-411	Micro Graft Harvesting Sagittal Saw Blade	49
5023-412	Micro Graft Harvesting Sagittal Saw Blade	49
5023-413	Micro Graft Harvesting Sagittal Saw Blade	49
5030-028	Small Holding Clamps	79
5030-029	Medium Holding Clamps	79
5038-001	Osteon Handpiece, Straight	74, 85, 86
5038-002	Osteon Handpiece, Angled	74, 85, 86
5038-004	Osteon Footswitch with Hose	74, 85, 86
5038-005	Osteon Hose	74, 85, 86
5038-009	Bur Rack, for Osteon and Ototom Burs	74, 85, 86
5038-600	Osteon Tubing Set with irrigation tip	74
5039-003	AC Power Cord	81
5040-180	Fluid Bag Support Rod for D3000 & E9000	81
5044-145	Micro Sagittal Saw Blade	46
5044-999-52	1/4" (6.35mm) Replacement Chuck Key	13
5044-999-53	5/32" (4.0mm) Replacement Chuck Key	13
5050-041	C-Wire	78
5050-042	C-Wire	78
5050-043	C-Wire	78
5050-044	C-Wire	78
5050-046	C-Wire	78
5050-047	C-Wire	78
5050-048	C-Wire	78

Numerical Index



Cat. No.	Description	Page No.
5050-049	C-Wire	78
5052-010	Universal Hose, 10' (3.05m)	76, 82 - 84, 86
5052-012	Standard Sterilization Case	74, 79, 83
5052-018	Extension Hose, 10' (3.05m), Male Schrader to Female Schrader	82
5052-029	Extension Hose, 10' (3.05m), Male Diss and Female Schrader	82
5052-041	Cebotome Revision Drill Handpiece	76, 86
5052-042	Standard Bur Extension, Short, use with standard Cebotome Burs	76, 86
5052-046	Cebotome Bur Wrench	76, 86
5052-048	Standard Bur Extension, Medium, use with standard Cebotome Burs	76, 86
5052-049	Standard Bur Extension, Long, use with standard Cebotome Burs	76
5052-125	Cebotome Metal Cutting Wheels	77, 86
5052-126	Cebotome Metal Cutting Wheels	77, 86
5052-129	Cebotome Metal Cutting Wheel Wrench	86
5052-135	Stryker Large Bone Reciprocating Blade	39
5052-142	Metal Cutting Bur Extension	76, 86
5052-143	Cebotome Metal Cutting Bur Wrench	76, 86
5052-144	Cebotome Diamond Metal Cutting Bur	77, 86
5052-145	Cebotome Diamond Metal Cutting Bur	77, 86
5052-173	Hall Large Bone Reciprocating Blade	38
5052-179	Hall Large Bone Reciprocating Blade	38
5052-181	Cebotome Bur	77
5052-182	Cebotome Bur	77
5052-183	Cebotome Bur	77
5052-184	Cebotome Bur	77
5052-185	Cebotome Bur	77
5052-186	Cebotome Bur	77
5052-187	Cebotome Bur	77
5052-188	Cebotome Bur	77
5052-189	Cebotome Bur	77
5052-190	Cebotome Bur	77
5052-191	Cebotome Bur	77
5052-192	Cebotome Bur	77
5052-258	Hall Large Bone Reciprocating Blade	38
5052-259	Hall Large Bone Reciprocating Blade	38
5052-260	Hall Large Bone Reciprocating Blade	38
5052-261	Hall Large Bone Reciprocating Blade	38
5052-276	Hall Large Bone Reciprocating Blade	38
5052-281	Hall Large Bone Reciprocating Blade	38
5052-291	Hall Large Bone Reciprocating Blade	38
5052-332	Stryker Large Bone Reciprocating Blade	39
5052-333	Stryker Large Bone Reciprocating Blade	39
5052-334	Stryker Large Bone Reciprocating Blade	39
5053-009	Micro 100 Medium Speed Drill	42
5053-010	Micro 100 Reciprocating Saw	42
5053-011	Micro 100 Sagittal Saw	42
5053-012	Micro 100 Oscillating Saw	42
5053-013	Micro 100 Wire Driver	42
5053-019	Wire Driver Adapter Chuck for use with 5053-013	42
5053-022	Micro 100 XL Oscillating Saw	42
5053-050	Micro Depth Limiting Reciprocating Saw Blade	53
5053-051	Micro Depth Limiting Reciprocating Saw Blade	53
5053-052	Micro Depth Limiting Reciprocating Saw Blade	53
5053-053	Micro Depth Limiting Reciprocating Saw Blade	53
5053-054	Micro Reciprocating Saw Blade	53
5053-074	Micro 100 Throttle Clip (for lever) Used during footswitch operation	42

Cat. No.	Description	Page No.
5053-080	Micro Reciprocating Rasp	53
5053-081	Micro Reciprocating Rasp	53
5053-082	Micro Reciprocating Rasp	53
5053-083	Micro Reciprocating Rasp	53
5053-084	Micro Reciprocating Rasp	53
5053-085	Micro Reciprocating Rasp	53
5053-122	Micro Reciprocating Saw Blade	53
5053-123	Wire Guard for use with 5053-013	42
5053-124	Cleaning Brush for Cannulated Handpieces	79
5053-126	Micro Reciprocating Saw Blade	53
5053-220	Micro Reciprocating Saw Blade	53
5053-221	Micro Reciprocating Saw Blade	53
5056-080	80 Grit Dermabrasion Osteon Bur	75
5056-081	80 Grit Dermabrasion Osteon Bur	75
5056-082	80 Grit Dermabrasion Osteon Bur	75
5056-083	80 Grit Dermabrasion Osteon Bur	75
5056-109	Osteon/Ototome Bur Set, includes 16 burs (see System Listings Page 26)	74
5056-116	Fluted Round Osteon Bur	75
5056-117	Fluted Round Osteon Bur	75, 85
5056-118	Fluted Round Osteon Bur	75, 85
5056-119	Fluted Round Osteon Bur	75, 85
5056-125	Fluted Round Osteon Bur	75
5056-126	Fluted Round Osteon Bur	75, 85
5056-127	Fluted Round Osteon Bur	75
5056-128	Fluted Round Osteon Bur	75
5056-140	Fluted Round Osteon Bur	75, 85
5056-142	Fluted Round Osteon Bur	75, 85
5056-144	Fluted Round Osteon Bur	75
5056-145	Fluted Round Osteon Bur	75, 85
5056-146	Fluted Round Osteon Bur	75, 85
5056-147	Fine Diamond Round Osteon Bur	75
5056-148	Fine Diamond Round Osteon Bur	75
5056-149	Fluted Round Osteon Bur	75
5056-150	Fine Diamond Round Osteon Bur	75
5056-152	Fine Diamond Round Osteon Bur	75, 85
5056-158	Fine Diamond Round Osteon Bur	75, 85
5056-159	Fine Diamond Round Osteon Bur	75, 85
5056-160	Fine Diamond Round Osteon Bur	75
5056-161	Fine Diamond Round Osteon Bur	75, 85
5056-162	Fine Diamond Round Osteon Bur	75, 85
5056-163	Fine Diamond Round Osteon Bur	75
5056-164	Fluted Round Osteon Bur	75, 85
5056-165	Fluted Round Osteon Bur	75
5056-166	Fluted Round Osteon Bur	75, 85
5056-167	Fine Diamond Round Osteon Bur	75
5056-168	Fine Diamond Round Osteon Bur	75
5056-170	Deep Access Osteon Bur	75
5056-171	Deep Access Osteon Bur	75
5056-172	Deep Access Osteon Bur	75
5056-173	Fluted Round Osteon Bur	75
5056-174	Fluted Round Osteon Bur	75
5056-175	Deep Access Osteon Bur	75
5056-176	Deep Access Osteon Bur	75
5056-177	Fluted Round Osteon Bur	75
5056-178	Fluted Round Osteon Bur	75

Cat. No.	Description	Page No.
5056-179	Fluted Round Osteon Bur	75
5056-184	Wire Pass Osteon Bur	75
5056-185	Wire Pass Osteon Bur	75
5056-186	Wire Pass Osteon Bur	75
5056-950	Coarse Diamond Round Osteon Bur	75
5056-952	Coarse Diamond Round Osteon Bur	75
5056-958	Coarse Diamond Round Osteon Bur	75
5056-959	Coarse Diamond Round Osteon Bur	75
5056-962	Coarse Diamond Round Osteon Bur	75
5056-967	Coarse Diamond Round Osteon Bur	75
5057-040	Repeat Sternotomy Saw, Pneumatic	57
5058-000	Surgairtome Two® High Speed Drill Kit (System listing on page xx)	42
5058-001	Surgairtome Two High Speed Drill	42
5058-001	Surgairtome Two Handpiece	83
5059-005	Pneumatic Sternum Saw	57
5059-006	Sternum Guard, 2.25" (5.72mm), for PRO6450, PRO6450M	57
5059-007	Sternum Saw Wrench, for PRO6450, PRO6450M & 5059-005	57
5059-011	Sternum Guard, 3.12" (7.92mm), for PRO6450, PRO6450M	57
5059-031	Sternum Saw Blade	58
5059-031-01	Sternum Saw Blade	58
5059-032-01	Sternum Saw Blade	58
5059-233	Sternum Saw Blade	58
5059-235	Sternum Saw Blade	58
5059-278	Repeat Sternotomy Blade	59
5059-279	Repeat Sternotomy Blade	59
5059-280	Repeat Sternotomy Blade	59
5059-531	Sternum Saw Blade	58
5059-532	Sternum Saw Blade	58
5071-105	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-106	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-107	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-108	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-109	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-110	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-111	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-112	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-113	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-120	Hall Large Bone Knee Blade	20
5071-121	Hall Large Bone Knee Blade	21
5071-123	Hall Large Bone Knee Blade	24
5071-124	Hall Large Bone Knee Blade	25
5071-125	Repeat Sternotomy Blade	60
5071-126	Repeat Sternotomy Blade	60
5071-127	ACL Large Bone Oscillating Blade	32
5071-129	Hall Large Bone Knee Blade	25
5071-132	Large Bone Oscillating Blade	34
5071-133	ACL Large Bone Oscillating Blade	32
5071-134	ACL Large Bone Oscillating Blade	32
5071-135	ACL Large Bone Oscillating Blade	33
5071-140	Large Bone Oscillating Blade	37
5071-141	Large Bone Oscillating Blade	36
5071-142	Large Bone Oscillating Blade	37
5071-143	Large Bone Oscillating Blade	36
5071-144	Hall Large Bone Knee Blade	26
5071-145	Large Bone Oscillating Blade	37

Cat. No.	Description	Page No.
5071-146	Hall Large Bone Knee Blade	20
5071-161	Hall Large Bone Knee Blade	24
5071-163	Hall Large Bone Knee Blade	26
5071-166	Large Bone Oscillating Blade	34
5071-167	Hall Large Bone Knee Blade	21
5071-174	Large Bone Oscillating Blade	36
5071-175	Large Bone Oscillating Abrader	37
5071-177	Large Bone Oscillating Blade	35
5071-178	Large Bone Oscillating Blade	35
5071-180	Hall Large Bone Knee Blade	23
5071-181	Hall Large Bone Knee Blade	23
5071-184	Hall Large Bone Knee Blade	21
5071-185	Hall Large Bone Knee Blade	20
5071-186	Large Bone Oscillating Blade	35
5071-187	Large Bone Oscillating Blade	35
5071-188	Large Bone Oscillating Blade	35
5071-195	Hall Large Bone Knee Blade	20
5071-196	Hall Large Bone Knee Blade	20
5071-201	Hall Large Bone Knee Blade	23
5071-228	Repeat Sternotomy Blade	60
5071-229	Repeat Sternotomy Blade	60
5071-230	Repeat Sternotomy Blade	60
5071-236	ACL Large Bone Oscillating Blade	32
5071-246	ACL Large Bone Oscillating Blade	32
5071-300	Intrex Blade	17
5071-301	Intrex Blade	17
5071-302	Intrex Blade	17
5071-303	Intrex Blade	17
5071-304	Intrex Blade	17
5071-305	Intrex Blade	17
5071-306	Intrex Blade	17
5071-310	Intrex Blade	17
5071-311	Intrex Blade	17
5071-312	Intrex Blade	17
5071-313	Intrex Blade	17
5071-314	Intrex Blade	17
5071-315	Intrex Blade	17
5071-316	Intrex Blade	17
5071-320	Intrex Blade	17
5071-321	Intrex Blade	17
5071-322	Intrex Blade	17
5071-323	Intrex Blade	17
5071-324	Intrex Blade	17
5071-325	Intrex Blade	17
5071-326	Intrex Blade	17
5071-330	Intrex Blade	17
5071-331	Intrex Blade	17
5071-332	Intrex Blade	17
5071-333	Intrex Blade	17
5071-334	Intrex Blade	17
5071-335	Intrex Blade	17
5071-336	Intrex Blade	17
5071-340	Intrex Blade	17
5071-341	Intrex Blade	17
5071-342	Intrex Blade	17

Cat. No.	Description	Page No.
5071-343	Intrex Blade	17
5071-344	Intrex Blade	17
5071-345	Intrex Blade	17
5071-346	Intrex Blade	17
5071-350	Intrex Blade	17
5071-351	Intrex Blade	17
5071-352	Intrex Blade	17
5071-353	Intrex Blade	17
5071-354	Intrex Blade	17
5071-355	Intrex Blade	17
5071-356	Intrex Blade	17
5071-519	Hall Large Bone Knee Blade	25
5071-525	Hall Large Bone Knee Blade	26
5071-526	Hall Large Bone Knee Blade	22
5071-527	Hall Large Bone Knee Blade	23
5071-528	Hall Large Bone Knee Blade	19
5071-529	Hall Large Bone Knee Blade	25
5071-530	Hall Large Bone Knee Blade	26
5071-531	Hall Large Bone Knee Blade	23
5071-532	Hall Large Bone Knee Blade	20
5071-533	Hall Large Bone Knee Blade	19
5071-534	Hall Large Bone Knee Blade	22
5071-538	Hall Large Bone Knee Blade	21
5071-539	Hall Large Bone Knee Blade	26
5071-547	Hall Large Bone Knee Blade	21
5071-548	Hall Large Bone Knee Blade	21
5071-550	Large Bone Oscillating Blade	34
5071-551	Large Bone Oscillating Blade	34
5071-552	Large Bone Oscillating Blade	34
5071-553	Large Bone Oscillating Blade	34
5071-554	Large Bone Oscillating Blade	36
5071-555	Large Bone Oscillating Blade	37
5071-556	Large Bone Oscillating Blade	35
5071-557	Large Bone Oscillating Blade	36
5071-558	Hall Large Bone Knee Blade	22
5071-560	Large Bone Oscillating Blade	36
5071-562	Hall Large Bone Knee Blade	21
5071-564	Hall Large Bone Knee Blade	24
5071-565	Hall Large Bone Knee Blade	26
5071-568	Hall Large Bone Knee Blade	24
5071-570-50	Hall Large Bone Knee Blade	22
5071-571	ACL Large Bone Oscillating Blade	33
5071-572	Large Bone Oscillating Blade	35
5071-573	Hall Large Bone Knee Blade	22
5071-575	Hall Large Bone Knee Blade	22
5071-576	Hall Large Bone Knee Blade	24
5071-582	Hall Large Bone Knee Blade	23
5071-583	Hall Large Bone Knee Blade	19
5071-589	Hall Large Bone Knee Blade	20
5071-590	Hall Large Bone Knee Blade	24
5071-594	Hall Large Bone Knee Blade	19
5071-650	Hall Large Bone Knee Blade	19
5072-100	Stryker Large Bone Knee Blade	27
5072-101	Stryker Large Bone Knee Blade	29
5072-102	Stryker Large Bone Knee Blade	27, 28

Cat. No.	Description	Page No.
5072-105	Stryker Large Bone Knee Blade	28
5072-109	Stryker Large Bone Knee Blade	28
5072-110	Stryker Large Bone Knee Blade	29
5072-111	Stryker Large Bone Knee Blade	29
5072-140	Stryker Large Bone Knee Blade	27
5072-142	Stryker Large Bone Knee Blade	28
5072-143	Stryker Large Bone Knee Blade	27
5072-145	Stryker Large Bone Knee Blade	27
5072-150	Stryker Large Bone Knee Blade	27
5072-151	Stryker Large Bone Knee Blade	28
5072-152	Stryker Large Bone Knee Blade	28
5072-158	Stryker Large Bone Knee Blade	29
5072-176	Stryker Large Bone Knee Blade	30
5072-177	Stryker Large Bone Knee Blade	30
5072-181	Stryker Large Bone Knee Blade	29
5072-183	Stryker Large Bone Knee Blade	29
5072-185	Stryker Large Bone Knee Blade	29
5072-195	Stryker Large Bone Knee Blade	28
5072-218	Stryker Large Bone Knee Blade	29
5089-208	Micro Straight Side – Cutting Bur	55
5090-247	Micro Wire Pass Bur	54
5090-248	Micro Wire Pass Bur	54
5090-249	Micro Wire Pass Bur	54
5091-080	Micro Dermabraders – Fine Diamond Bur	54
5091-081	Micro Dermabraders – Fine Diamond Bur	54
5091-082	Micro Dermabraders – Fine Diamond Bur	54
5091-083	Micro Dermabraders – Fine Diamond Bur	54
5091-201	Micro Straight Side – Cutting Bur	55
5091-202	Micro Straight Side – Cutting Bur	55
5091-204	Micro Straight Side – Cutting Bur	55
5091-206	Micro Tapered Side – Cutting Bur	55
5091-212	Micro Lindemann Bur	56
5091-213	Micro Tapered Cross-cut Fissure Bur	56
5091-214	Micro Tapered Side – Cutting Bur	55
5091-215	Micro Tapered Side – Cutting Bur	55
5091-216	Micro Straight Side – Cutting Bur	55
5091-217	Micro Straight Side – Cutting Bur	55
5091-218	Micro Tapered Side – Cutting Bur	55
5091-219	Micro Tapered Side – Cutting Bur	55
5091-221	Micro Round Fluted Bur	54
5091-222	Micro Round Fluted Bur	54
5091-223	Micro Round Fluted Bur	54
5091-224	Micro Round Fluted Bur	54
5091-225	Micro Round Fluted Bur	54
5091-226	Micro Round Fluted Bur	54
5091-228	Micro Round Fluted Bur	54
5091-230	Micro Round Fluted Bur	54
5091-232	Micro Round Fluted Bur	54
5091-236	Micro Oval Bur	55
5091-238	Micro Oval Bur	55
5091-242	Micro Swanson Reamer Bur	54
5091-243	Micro Swanson Reamer Bur	54
5091-244	Micro Swanson Reamer Bur	54
5091-245	Micro Swanson Reamer Bur	54
5091-247	Micro Wire Pass Bur	54

Cat. No.	Description	Page No.
5091-248	Micro Wire Pass Bur	54
5091-249	Micro Wire Pass Bur	54
5091-252	Micro Barrel Bur	55
5091-254	Micro Barrel Bur	55
5091-260	Micro Fine Round Diamond Bur	54
5091-264	Micro Fine Round Diamond Bur	54
5091-268	Micro Fine Round Diamond Bur	54
5091-270	Micro Dermabraders – Fine Diamond Bur	54
5091-271	Micro Tapered Cross-cut Fissure Bur	56
5091-272	Micro Tapered Cross-cut Fissure Bur	56
5092-203	Micro Straight Side – Cutting Bur	55
5092-204	Micro Straight Side – Cutting Bur	55
5092-206	Micro Tapered Side – Cutting Bur	55
5092-214	Micro Tapered Side – Cutting Bur	55
5092-220	Micro Round Fluted Bur	54
5092-222	Micro Round Fluted Bur	54
5092-224	Micro Round Fluted Bur	54
5092-226	Micro Round Fluted Bur	54
5092-228	Micro Round Fluted Bur	54
5092-230	Micro Round Fluted Bur	54
5092-232	Micro Round Fluted Bur	54
5092-236	Micro Round Fluted Bur	54, 55
5092-238	Micro Oval Bur	55
5092-247	Micro Wire Pass Bur	54
5092-252	Micro Barrel Bur	55
5092-254	Micro Barrel Bur	55
5092-260	Micro Fine Round Diamond Bur	54
5092-262	Micro Fine Round Diamond Bur	54
5092-264	Micro Fine Round Diamond Bur	54
5092-266	Micro Fine Round Diamond Bur	54
5092-268	Micro Fine Round Diamond Bur	54
5092-270	Micro Dermabraders – Fine Diamond Bur	54
5092-271	Micro Tapered Cross-cut Fissure Bur	56
5092-272	Micro Tapered Cross-cut Fissure Bur	56
5092-273	Micro Tapered Cross-cut Fissure Bur	56
5092-283	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-284	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-285	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-286	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-291	Micro Cylindrical Bur	56
5093-220	Micro Round Fluted Bur	54
5093-221	Micro Round Fluted Bur	54
5093-224	Micro Round Fluted Bur	54
5093-226	Micro Round Fluted Bur	54
5093-228	Micro Round Fluted Bur	54
5093-230	Micro Round Fluted Bur	54
5093-232	Micro Round Fluted Bur	54
5093-236	Micro Oval Bur	55
5093-238	Micro Oval Bur	55
5093-260	Micro Fine Round Diamond Bur	54
5093-262	Micro Fine Round Diamond Bur	54
5093-264	Micro Fine Round Diamond Bur	54
5093-266	Micro Fine Round Diamond Bur	54
5093-282	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-283	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56

Cat. No.	Description	Page No.
5093-284	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-285	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5144-002	Schrader Coupler	82
5271-502	Hall Large Bone Knee Blade	19
5271-561	Hall Large Bone Knee Blade	24
5271-580	Hall Large Bone Knee Blade	23
6020-021	MicroPower Medium Speed Drill, Removable Lever	41
6020-022	MicroPower Sagittal Saw, Removable Lever	41
6020-023	MicroPower Reciprocating Saw, Removable Lever	41
6020-024	MicroPower Oscillating Saw, Permanently Attached Lever	41
6020-025	Oral Max High Speed Drill, Removable Lever	41
6020-026	MicroPower Orthopedic High Speed Drill, Removable Lever	41
6020-260	Smartguard Irrigator	42
6021-021	MicroPower Medium Speed Drill, Permanently Attached Lever	41
6021-022	MicroPower Sagittal Saw, Permanently Attached Lever	41
6021-023	MicroPower Reciprocating Saw, Permanently Attached Lever	41
6021-024	MicroPower Oscillating Saw, Permanently Attached Lever	41
6021-025	Oral Max High Speed Drill, Permanently Attached Lever	41
6021-026	MicroPower Orthopedic High Speed Drill, Permanently Attached Lever	41
7007-03	Staplizer Cartridge	78
7010-03	Staplizer Cartridge	78
7015-03	Staplizer Cartridge	78
7020-001	Ultrapower High Speed Drill	61, 83, 84
7020-002	Ultrapower Low Speed Drill	61, 83, 84
7020-003	Ultrapower Perforator Drive	61, 83, 84
7020-010	Short Bur Guard (like Midas-X) Red	61 - 66, 83, 84
7020-011	Medium Bur Guard (like Midas-O) Copper	61 - 66, 83, 84
7020-012	Long Bur Guard (like Midas-M/B) Mustard	61 - 66, 73, 83, 84
7020-013	X-Long Bur Guard (like Midas-AM/A) Lime	61 - 66, 83, 84
7020-014	Short Revision Extension (like Midas TU) Teal	61 - 66, 83, 84
7020-015	Medium Revision Extension (like Midas-R) Purple	61 - 63, 65, 83, 84
7020-016	Long Revision Extension (like Midas RX) Blue	61, 62, 65, 83, 84
7020-017	Short Bur Guard (like Midas B-5) Black	61, 67, 83, 84
7020-018	Medium Bur Guard (like Midas B-1) Green	61, 67, 83, 84
7020-019	Long Bur Guard (like Midas S-1) Brown	61, 73, 83, 84
7020-022	90 degree Angle Attachment (like Midas AF/BF)	61, 83, 84
7020-023	20 degree Long Keyless Attachment	61, 83, 84
7020-024	20 degree Angle X-Long Keyless Attachment	83, 84
7020-040	Key Wrench, 2 per box	83, 84
7020-041	Attachment Spray	61, 83, 84
7020-045	Ultrapower Footswitch	67, 83, 84
7020-046	Bur Cabinet for Neuro Burs	64, 66, 84
7021-012	UltraPower Side Cutting Bur	62, 64
7021-013	UltraPower Side Cutting Bur	62
7021-015	UltraPower Side Cutting Bur	62
7021-024	UltraPower Round Fluted Bur	62
7021-025	UltraPower Round Fluted Bur	62
7021-026	UltraPower Round Fluted Bur	62
7021-027	UltraPower Round Fluted Bur	62
7021-028	UltraPower Round Fluted Bur	62
7021-029	UltraPower Round Fluted Bur	62
7021-030	UltraPower Round Fluted Bur	62
7021-031	UltraPower Round Fluted Bur	62
7021-032	UltraPower Round Fluted Bur	64
7021-033	UltraPower Round Fluted Bur	64

Cat. No.	Description	Page No.
7021-034	UltraPower Round Fluted Bur	66
7021-036	UltraPower Fluted Cylindrical Bur	66
7021-037	UltraPower Fluted Cylindrical Bur	66
7021-040	UltraPower Wire Pass Bur	64
7021-041	UltraPower Wire Pass Bur	65
7021-042	UltraPower Wire Pass Bur	62
7021-043	UltraPower Side Cutting Bur	62
7021-051	UltraPower Trephine Bur	62
7021-062	UltraPower Fine Round Diamond Bur	62
7021-063	UltraPower Fine Round Diamond Bur	62
7021-064	UltraPower Fine Round Diamond Bur	63
7021-065	UltraPower Fine Round Diamond Bur	63
7021-066	UltraPower Fine Round Diamond Bur	63
7021-074	UltraPower Coarse Round Diamond Bur	65
7021-075	UltraPower Coarse Round Diamond Bur	65
7021-077	UltraPower Coarse Round Diamond Bur	66
7021-081	UltraPower Fluted Acorn Bur	66
7021-083	UltraPower Fluted Acorn Bur	66
7021-092	UltraPower Fine Diamond Wheel for Metal Cutting Bur	66
7021-094	UltraPower Carbide Metal Cutting Bur	64
7021-095	UltraPower Oval Fluted Bur	64
7021-096	UltraPower Oval Fluted Bur	64
7021-112	UltraPower Side Cutting Bur	64
7021-113	UltraPower Side Cutting Bur	64
7021-114	UltraPower Side Cutting Bur	64
7021-115	UltraPower Side Cutting Bur	62
7021-116	UltraPower Side Cutting Bur	62
7021-119	UltraPower Side Cutting Bur	62
7021-122	UltraPower Round Fluted Bur	62
7021-123	UltraPower Round Fluted Bur	62
7021-124	UltraPower Round Fluted Bur	62
7021-125	UltraPower Round Fluted Bur	62
7021-126	UltraPower Round Fluted Bur	62
7021-127	UltraPower Round Fluted Bur	62
7021-128	UltraPower Round Fluted Bur	62
7021-129	UltraPower Round Fluted Bur	62
7021-130	UltraPower Round Fluted Bur	62
7021-131	UltraPower Round Fluted Bur	62
7021-132	UltraPower Round Fluted Bur	64
7021-133	UltraPower Round Fluted Bur	64
7021-134	UltraPower Round Fluted Bur	66
7021-136	UltraPower Fluted Cylindrical Bur	66
7021-137	UltraPower Fluted Cylindrical Bur	66
7021-140	UltraPower Wire Pass Bur	66
7021-141	UltraPower Wire Pass Bur	66
7021-141-35	UltraPower Fluted Frustrum Bur – Depth Limiting	66
7021-141-50	UltraPower Fluted Frustrum Bur – Depth Limiting	64
7021-141-90	UltraPower Fluted Frustrum Bur – Depth Limiting	66
7021-142	UltraPower Wire Pass Bur	62
7021-143	UltraPower Side Cutting Bur	62
7021-145-70	UltraPower Fluted Frustrum Bur – Depth Limiting	62
7021-161	UltraPower Fine Round Diamond Bur	62
7021-162	UltraPower Fine Round Diamond Bur	62
7021-163	UltraPower Fine Round Diamond Bur	62
7021-164	UltraPower Fine Round Diamond Bur	62

Cat. No.	Description	Page No.
7021-165	UltraPower Fine Round Diamond Bur	63
7021-166	UltraPower Fine Round Diamond Bur	63
7021-167	UltraPower Fine Round Diamond Bur	63
7021-171	UltraPower Coarse Round Diamond Bur	63
7021-172	UltraPower Coarse Round Diamond Bur	63
7021-173	UltraPower Coarse Round Diamond Bur	63
7021-174	UltraPower Coarse Round Diamond Bur	63
7021-175	UltraPower Coarse Round Diamond Bur	65
7021-176	UltraPower Coarse Round Diamond Bur	65
7021-177	UltraPower Coarse Round Diamond Bur	65
7021-181	UltraPower Fluted Acorn Bur	66
7021-182	UltraPower Fluted Acorn Bur	66
7021-183	UltraPower Fluted Acorn Bur	63
7021-195	UltraPower Oval Fluted Bur	63
7021-196	UltraPower Oval Fluted Bur	63
7021-202	UltraPower Side Cutting Bur	63
7021-206	UltraPower Side Cutting Bur	64
7021-207	UltraPower Side Cutting Bur	64
7021-210	UltraPower Side Cutting Bur	62
7021-214	UltraPower Side Cutting Bur	62
7021-218	UltraPower Side Cutting Bur	62
7021-220	UltraPower Round Fluted Bur	62
7021-222	UltraPower Round Fluted Bur	62
7021-223	UltraPower Round Fluted Bur	62
7021-224	UltraPower Round Fluted Bur	62
7021-225	UltraPower Round Fluted Bur	62
7021-226	UltraPower Round Fluted Bur	62
7021-227	UltraPower Round Fluted Bur	62
7021-228	UltraPower Round Fluted Bur	62
7021-229	UltraPower Round Fluted Bur	62
7021-230	UltraPower Round Fluted Bur	62
7021-231	UltraPower Round Fluted Bur	64
7021-232	UltraPower Round Fluted Bur	64
7021-233	UltraPower Round Fluted Bur	64
7021-235	UltraPower Fluted Cylindrical Bur	64
7021-236	UltraPower Fluted Cylindrical Bur	65
7021-237	UltraPower Fluted Cylindrical Bur	65
7021-243	UltraPower Side Cutting Bur	65
7021-250	UltraPower Trephine Bur	65
7021-252	UltraPower Trephine Bur	62
7021-253	UltraPower Trephine Bur	62
7021-255	UltraPower Trephine Bur	62
7021-261	UltraPower Fine Round Diamond Bur	62
7021-262	UltraPower Fine Round Diamond Bur	62
7021-263	UltraPower Fine Round Diamond Bur	62
7021-264	UltraPower Fine Round Diamond Bur	62
7021-265	UltraPower Fine Round Diamond Bur	63
7021-266	UltraPower Fine Round Diamond Bur	63
7021-267	UltraPower Fine Round Diamond Bur	63
7021-272	UltraPower Coarse Round Diamond Bur	63
7021-273	UltraPower Coarse Round Diamond Bur	63
7021-274	UltraPower Coarse Round Diamond Bur	63
7021-275	UltraPower Coarse Round Diamond Bur	65
7021-276	UltraPower Coarse Round Diamond Bur	65
7021-277	UltraPower Coarse Round Diamond Bur	65

Cat. No.	Description	Page No.
7021-279	UltraPower Fine Diamond Cylindrical Bur	65
7021-280	UltraPower Fluted Acorn Bur	65
7021-281	UltraPower Fluted Acorn Bur	66
7021-282	UltraPower Fluted Acorn Bur	66
7021-283	UltraPower Fluted Acorn Bur	66
7021-292	UltraPower Fine Diamond Wheel for Metal Cutting Bur	66
7021-294	UltraPower Carbide Metal Cutting Bur	66
7021-295	UltraPower Oval Fluted Bur	64
7021-296	UltraPower Oval Fluted Bur	63
7021-297	UltraPower Carbide Metal Cutting Bur	63
7021-301	UltraPower Side Cutting Bur	63
7021-302	UltraPower Side Cutting Bur	63
7021-303	UltraPower Side Cutting Bur	63
7021-305	UltraPower Side Cutting Bur	64
7021-309	UltraPower Side Cutting Bur	64
7021-310	UltraPower Side Cutting Bur	64
7021-314	UltraPower Side Cutting Bur	62
7021-315	UltraPower Side Cutting Bur	62
7021-318	UltraPower Side Cutting Bur	62
7021-322	UltraPower Round Fluted Bur	62
7021-323	UltraPower Round Fluted Bur	62
7021-324	UltraPower Round Fluted Bur	62
7021-325	UltraPower Round Fluted Bur	62
7021-326	UltraPower Round Fluted Bur	62
7021-327	UltraPower Round Fluted Bur	62
7021-328	UltraPower Round Fluted Bur	62
7021-329	UltraPower Round Fluted Bur	62
7021-331	UltraPower Round Fluted Bur	64
7021-332	UltraPower Round Fluted Bur	64
7021-333	UltraPower Round Fluted Bur	64
7021-335	UltraPower Fluted Cylindrical Bur	64
7021-336	UltraPower Fluted Cylindrical Bur	62
7021-337	UltraPower Fluted Cylindrical Bur	62
7021-338	UltraPower Fluted Cylindrical Bur	62
7021-361	UltraPower Fine Round Diamond Bur	62
7021-362	UltraPower Fine Round Diamond Bur	62
7021-363	UltraPower Fine Round Diamond Bur	62
7021-364	UltraPower Fine Round Diamond Bur	62
7021-365	UltraPower Fine Round Diamond Bur	63
7021-366	UltraPower Fine Round Diamond Bur	63
7021-367	UltraPower Fine Round Diamond Bur	63
7021-371	UltraPower Coarse Round Diamond Bur	63
7021-372	UltraPower Coarse Round Diamond Bur	63
7021-373	UltraPower Coarse Round Diamond Bur	63
7021-374	UltraPower Coarse Round Diamond Bur	63
7021-375	UltraPower Coarse Round Diamond Bur	65
7021-376	UltraPower Coarse Round Diamond Bur	65
7021-377	UltraPower Coarse Round Diamond Bur	65
7021-378	UltraPower Fine Diamond Cylindrical Bur	65
7021-379	UltraPower Fine Diamond Cylindrical Bur	65
7021-380	UltraPower Fluted Acorn Bur	65
7021-381	UltraPower Fluted Acorn Bur	66
7021-382	UltraPower Fluted Acorn Bur	66
7021-383	UltraPower Fluted Acorn Bur	66
7021-391	UltraPower Fluted Frustrum Bur	66

Cat. No.	Description	Page No.
7021-392	UltraPower Fine Diamond Wheel for Metal Cutting Bur	63
7021-395	UltraPower Oval Fluted Bur	63
7021-396	UltraPower Oval Fluted Bur	62
7021-408	UltraPower Side Cutting Bur	62
7021-409	UltraPower Side Cutting Bur	63
7021-426	UltraPower Round Fluted Bur	62
7021-427	UltraPower Round Fluted Bur	62
7021-427	UltraPower Side Cutting Bur	62
7021-428	UltraPower Round Fluted Bur	62
7021-429	UltraPower Round Fluted Bur	62
7021-431	UltraPower Round Fluted Bur	64
7021-432	UltraPower Round Fluted Bur	63
7021-433	UltraPower Round Fluted Bur	65
7021-437	UltraPower Fluted Cylindrical Bur	65
7021-474	UltraPower Coarse Round Diamond Bur	66
7021-482	UltraPower Fluted Acorn Bur	63
7021-483	UltraPower Fluted Acorn Bur	62
7021-492	UltraPower Fine Diamond Wheel for Metal Cutting Bur	62
7021-509	UltraPower Side Cutting Bur	62
7021-526	UltraPower Round Fluted Bur	62
7021-528	UltraPower Round Fluted Bur	62
7021-529	UltraPower Round Fluted Bur	62
7021-531	UltraPower Round Fluted Bur	65
7021-532	UltraPower Round Fluted Bur	65
7021-533	UltraPower Round Fluted Bur	62
7021-581	UltraPower Fluted Acorn Bur	62
7021-583	UltraPower Fluted Acorn Bur	62
7021-631	UltraPower Round Fluted Bur	65
7021-632	UltraPower Round Fluted Bur	65
7021-633	UltraPower Round Fluted Bur	64
7021-681	UltraPower Fluted Acorn Bur	73
7021-683	UltraPower Fluted Acorn Bur	64
7021-713	UltraPower Side Cutting Bur	64, 73
7021-815	UltraPower Side Cutting Bur	73
7021-917	UltraPower Side Cutting Bur	78
7021-960	Helical Bur for E9010 Drill	78
7307-03	Staplizer Cartridge	78
7310-03	Staplizer Cartridge	78
7315-03	Staplizer Cartridge	78
7320-03	Staplizer Cartridge	78
7610-03	Staplizer Cartridge	78
7615-03	Staplizer Cartridge	78
7620-03	Staplizer Cartridge	78
7625-03	Staplizer Cartridge	82
7707-03	Staplizer Cartridge	68
7710-03	Staplizer Cartridge	69
A201	10' (3.05m) Hose for Staplizer® Handpiece	48
A9350	Rhinotec Cuda Blade 4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box	48
A9360	Rhinotec Gator Blade 4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box	48
C8803	Micro ACL Sagittal Saw Blade	81
C8804	Micro ACL Sagittal Saw Blade	71
C8806	Micro ACL Sagittal Saw Blade	71
C9863	3-Pedal Footswitch	68
C9911	Rhinotec Spherical Bur, 2.9mm, Hooded, 7.6cm Working Length, 6/ Box	68
C9912	Rhinotec Spherical Bur, 3.5mm, Hooded, 7.6cm Working Length, 6/Box	68

Cat. No.	Description	Page No.
C9950	Rhinotec Cuda Blade 2.0mm, 7.6cm Working Length, 6/Box	69
C9951	Rhinotec Cuda Blade 2.9mm, 7.6cm Working Length, 6/Box	69
C9952	Rhinotec Cuda Blade 3.5mm, 7.6cm Working Length, 6/Box	69
C9960	Rhinotec Gator Blade 2.0mm, 7.6cm Working Length, 6/Box	70
C9961	Rhinotec Gator Blade 2.9mm, 7.6cm Working Length, 6/Box	70
C9962	Rhinotec Gator Blade 3.5mm, 7.6cm Working Length, 6/Box	80
C9990	Rhinotec Vortex Router 2.9mm, Hooded, 7.6cm Working Length, 6/Box	80
C9991	Rhinotec Vortex Router 3.5mm, Hooded, 7.6cm Working Length, 6/Box	80
D3000	Advantage Console	73
D4000	D4000 Console	73
D4000A	D4000 Console with Irrigation	73
E4110	Wirepass E Bur	73
E4115	Wirepass E Bur	73
E4120	Wirepass E Bur	73
E4208	E9010 Side Cutting Bur	72
E4217	E9010 Lindemann Bur	72
E4220	E9010 Side Cutting Bur	72
E5050	Coarse Round Diamond E Bur	72
E5120	Coarse Round Diamond E Bur	72
E5130	Coarse Round Diamond E Bur	72
E5140	Coarse Round Diamond E Bur	72
E5150	Coarse Round Diamond E Bur	72
E5160	Coarse Round Diamond E Bur	72
E5170	Coarse Round Diamond E Bur	72
E5210	Coarse Round Diamond E Bur	72
E5215	Coarse Round Diamond E Bur	72
E5220	Coarse Round Diamond E Bur	72
E5225	Coarse Round Diamond E Bur	72
E5230	Coarse Round Diamond E Bur	72
E5240	Coarse Round Diamond E Bur	72
E6040	Fine Round Diamond E Bur	72
E6105	Fine Round Diamond E Bur	72
E6110	Fine Round Diamond E Bur	72
E6115	Fine Round Diamond E Bur	72
E6120	Fine Round Diamond E Bur	72
E6125	Fine Round Diamond E Bur	72
E6130	Fine Round Diamond E Bur	72
E6140	Fine Round Diamond E Bur	72
E6150	Fine Round Diamond E Bur	72
E6160	Fine Round Diamond E Bur	72
E6205	Fine Round Diamond E Bur	72
E6207	Fine Round Diamond E Bur	72
E6210	Fine Round Diamond E Bur	72
E6215	Fine Round Diamond E Bur	72
E6220	Fine Round Diamond E Bur	72
E6225	Fine Round Diamond E Bur	72
E6230	Fine Round Diamond E Bur	72
E6240	Fine Round Diamond E Bur	72
E6250	Fine Round Diamond E Bur	72
E8007	Round Fluted E Bur	72
E8010	Round Fluted E Bur	72
E8015	Round Fluted E Bur	72
E8020	Round Fluted E Bur	72
E8025	Round Fluted E Bur	72
E8030LSO	LSO Round Fluted E Bur	72

Cat. No.	Description	Page No.
E8040LS0	LS0 Round Fluted E Bur	72
E8050LS0	LS0 Round Fluted E Bur	72
E8065LS0	LS0 Round Fluted E Bur	72
E8105	Round Fluted E Bur	72
E8110	Round Fluted E Bur	72
E8115	Round Fluted E Bur	72
E8120	Round Fluted E Bur	72
E8125	Round Fluted E Bur	72
E8130LS0	LS0 Round Fluted E Bur	72
E8140LS0	LS0 Round Fluted E Bur	72
E8150LS0	LS0 Round Fluted E Bur	72
E8165LS0	LS0 Round Fluted E Bur	72
E8205	Round Fluted E Bur	72
E8207	Round Fluted E Bur	72
E8210	Round Fluted E Bur	72
E8215	Round Fluted E Bur	72
E8220	Round Fluted E Bur	72
E8225	Round Fluted E Bur	72
E8230LS0	LS0 Round Fluted E Bur	80
E8240LS0	LS0 Round Fluted E Bur	67
E8265LS0	LS0 Round Fluted E Bur	67
E9000	E9000 Console	67
E9005	High Speed Shaver	73
E9010	CoolFlex® High Speed Drill	73
E9011-1	Medium Bur Guard for E9010 Drill	67, 72, 73
E9011-2	Long Bur Guard for E9010 Drill	67, 72, 73
E9011-3	X-Long Bur Guard for E9010 Drill	67
E9011-7	Short Neuro Guard for E9010 Drill	67, 73
E9011-8	Medium Neuro Guard for E9010 Drill	67
E9011-9	Long Neuro Guard for E9010 Drill	67
E9012	Locking Collar for E9010 Drill	67
E9015	Cranial Perforator Drive	68
E9315	DeClogger for removing tissue and bone from shaver blades and burs	68
E9323	Bur Rack, Holds up to 20 burs for E9010 Drill	69
E9350	Rhinotec Cuda Blade 4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	69
E9351	Rhinotec Cuda Blade 3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	69
E9360	Rhinotec Gator Blade 4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	67
E9361	Rhinotec Gator Blade 3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	67
E9363	Laryngeal Blade 3.7mm, Bendable to 45°, 20.0cm Working Length, 6/Box	67
E9414	Drill Tubing for E9010 Drill without irrigation tip	67
E9415A	Drill Tubing for E9010 Drill with irrigation tip	43
E9416	4.2mm FESS Irrigation Tubing Set for E9005 Shaver	54
E9417	3.7mm FESS Irrigation Tubing Set for E9005 Shaver (Not Shown)	55
E9418	Irrigation Tubing Set	55
L10D	Micro Fine Round Diamond Bur	55
L12T	Micro Tapered Side – Cutting Bur	54
L13S	Micro Straight Side – Cutting Bur	55
L17S	Micro Straight Side – Cutting Bur	54
L18R	Micro Round Fluted Bur	8
L20T	Micro Tapered Side – Cutting Bur	8
L23R	Micro Round Fluted Bur	8
L3000	Hall Lithium Charger	55
L3000LG	Hall Large Lithium Battery, 13.2 Volt	56
L300SM	Hall Small Lithium Battery, 13.2 Volt	55
L40B	Micro Oval Bur	54

Cat. No.	Description	Page No.
L45C	Micro Cone Bur	82
L55B	Micro Oval Bur	55
L60R	Micro Round Fluted Bur	55
M105A	Blitz Instrument Cleaner (Use with Staplizer Handpiece)	55
M10S	Micro Straight Side – Cutting Bur	55
M12T	Micro Tapered Side – Cutting Bur	55
M13S	Micro Straight Side – Cutting Bur	82
M15S	Micro Straight Side – Cutting Bur	82
M17S	Micro Straight Side – Cutting Bur	55
M207	Regulator, Single Connection	82
M208	Regulator, Dual Connection	82
M20T	Micro Tapered Side – Cutting Bur	55
M317	Instrument Lubricant for Staplizer	56
M322	DISS to Schrader Adapter	55
M40B	Micro Oval Bur	54
M45C	Micro Cone Bur	54
M55B	Micro Oval Bur	81
M60D	Micro Fine Round Diamond Bur	68
M60R	Micro Round Fluted Bur	68
MC5057	Handpiece Cord	69
N9150	Lynx Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	69
N9210	Rhinotec Cuda Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	69
N9220	Rhinotec Gator Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	68
N9320	Head & Neck Blade 3.5mm, 7.0cm Working Length, 6/Box	70
N9321	Head & Neck Blade 4.2mm, 11.5cm Working Length, 6/Box	71
N9330	MicroXcisor Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	71
N9375	Curved Vortex Router 4.5mm, 45° Pre-Bent, 11.5cm Working Length, 6/Box	71
N9380	Rhinotec Spherical Bur, 4.5mm, Hooded, 10.8cm Working Length, 6/Box	70
N9382	Rhinotec Spherical Bur, 3.5mm, Unhooded, 7.6cm Working Length, 6/Box	13
N9385	Rhinotec Spherical Diamond Bur, 3.5mm, Hooded, 10.8cm Working Length, 6/Box	13
N9390	Otolaryngology Bur 4.0mm, Hooded, 10.8cm Working Length, 6/Box	13
PR01540	Drill-Ream Attachment, AO	13
PR01541	Drill-Ream Attachment, 1/4" (6.35mm) Chuck	80
PR01547	Drill-Ream Attachment, Zimmer/Hudson	13
PR01548	Drill-Ream Attachment, DePuy/Hudson	13
PR02000	PowerPro Console	13
PR02029	AO/Trinkle	13
PR02030	5/32" (4.0mm) Jacobs Chuck, with Chuck Key	13
PR02033	High Speed Drill Attachment	13
PR02038	Small AO—Low Profile	13
PR02039	Mini AO	13
PR02040	3:1 Large AO Reamer	13
PR02041	1/4" (6.35mm) Jacobs Chuck	13
PR02042	3:1 1/4" High Torque Jacobs Chuck	13
PR02043	Sagittal Saw Attachment	13
PR02046	Zimmer/Hudson	13
PR02047	3:1 Zimmer/Hudson Reamer	13
PR02060	3:1 Aesculap Reamer	13
PR02065	3:1 Hudson Reamer (DePuy)	13
PR02070	Large AO	13
PR02075	Standard Hudson (DePuy)	13
PR02080	3:1 DHS Reamer	13
PR02138	Small AO—One Handed Operation	13
PR02250	1/4" Keyless Chuck	8
PR02329	3:1 Trinkle Reamer	8

Cat. No.	Description	Page No.
PR02514	Radiolucent Drive Adaptor	8
PR03010	Large NiCad Battery, 9.6 Volt	9
PR03011	X-Large NiCad Battery, 12 Volt	9
PR03020	Small NiCad Battery, 12 Volt	9
PR03110	Battery Case for PR03115	9
PR03115	Large NiMH Sterile Transfer Battery, 9.6 Volt	9
PR03130	Transfer Shroud for PR03110	9
PR03135	Charging Platform for PR03115	9
PR03520	Small NiMH Sterile Transfer Battery, 12 Volt	9
PR03521	Battery Case for PR03520	8
PR03531	Transfer Shroud for PR03521	44
PR03535	Charging Platform for PR03520	13
PR03600	NiCad Battery Charger	13
PR06000	Small Bone Tray	13
PR06040	5:1 Large AO Reamer	13
PR06042	5:1 1/4" High Torque Jacobs Chuck	13
PR06045	Reciprocating Saw Attachment for PR06100, PR06150, PR06200M, & PR06202M	13
PR06047	5:1 Zimmer/Hudson Reamer	13
PR06048	5:1 Hudson Reamer(DePuy)	10
PR06050	1/8" Keyless Chuck	10
PR06060	5:1 Aesculap Reamer	12
PR06100	Single Trigger Modular Handpiece	10
PR06125	Oscillating Saw	12
PR06128	Wire Driver, Range: 0.7 – 1.6mm (0.028 – 0.062")	11
PR06135	Reciprocating Saw	11
PR06140	Pin Driver, Range: 1.8 – 4.0mm (0.071 – 0.157")	11
PR06150	Single Trigger Modular Handpiece	7
PR06175	Oscillating Saw	7
PR06185	Reciprocating Saw	12
PR06200M	Single Trigger Modular Handpiece	12
PR06202M	2-Trigger Modular Handpiece	7
PR06228	Wire Driver Range: 0.7 – 1.6mm (0.028 – 0.062")	7
PR06240	Pin Driver Range: 1.8 – 4.0mm (0.071 – 0.157")	57
PR06350M	Oscillating Saw	40
PR06400M	Reciprocating Saw	40
PR06450M	MPower 2 Battery Sternum Saw	14
PR07000DE	Oral Max Elite High Speed Drill	14
PR07000SE	MicroPower+ Orthopedic High Speed Drill	14
PR07000T	3/4 Size Inner Tray, Full Set	44
PR07001T	1/2 Size Inner Tray, Drill Set	13
PR07002T	1/2 Size Inner Tray, Saw Set	6
PR07003T	1/2 Size Inner Tray, Small Bone Set	40
PR07020	Multipurpose High Speed	6
PR07100B	Hall 50 1-Trigger Modular Handpiece	40
PR07100SE	MicroPower+ Medium Speed Drill	6
PR07200B	Hall 50 2-Trigger Modular Handpiece	40
PR07200SE	MicroPower+ Sagittal Saw	6
PR07300B	Hall 50 Oscillating Saw Handpiece	40
PR07300SE	MicroPower+ Oscillating Saw	57
PR07400B	Hall 50 Reciprocating Saw Handpiece	18
PR07400SE	MicroPower+ Reciprocating Saw	18
PR07450B	Hall 50 Sternum Saw Handpiece	18
STK130-089-90	Stryker Replacement Blade	18
STK130-119-90	Stryker Replacement Blade	18
STK130-127-90	Stryker Replacement Blade	18

Cat. No.	Description	Page No.
STK190-089-90	Stryker Replacement Blade	18
STK190-089-90	Stryker Replacement Blade	18
STK190-127-05	Stryker Replacement Blade	18
STK190-127-90	Stryker Replacement Blade	18
STK190-137-90	Stryker Replacement Blade	18
STK193-137-90	Stryker Replacement Blade	18
STK193-137-90	Stryker Replacement Blade	18
STK250-100-90	Stryker Replacement Blade	18
STK250-119-90	Stryker Replacement Blade	18
STK250-127-90	Stryker Replacement Blade	15
STK250-137-90	Stryker Replacement Blade	15
STK250-147-90	Stryker Replacement Blade	15
TN130-089-90	M-Class Blade	15
TN130-100-90	M-Class Blade	15
TN130-107-90	M-Class Blade	15
TN130-119-90	M-Class Blade	15
TN130-127-90	M-Class Blade	15
TN130-137-90	M-Class Blade	15
TN130-147-90	M-Class Blade	15
TN190-089-90	M-Class Blade	15
TN190-100-90	M-Class Blade	15
TN190-107-90	M-Class Blade	15
TN190-119-05	M-Class Blade	15
TN190-119-90	M-Class Blade	15
TN190-127-05	M-Class Blade	15
TN190-127-90	M-Class Blade	15
TN190-137-05	M-Class Blade	15
TN190-137-90	M-Class Blade	15
TN190-147-90	M-Class Blade	15
TN193-119-90	M-Class Blade	15
TN193-127-90	M-Class Blade	15
TN193-137-90	M-Class Blade	15
TN250-089-90	M-Class Blade	15
TN250-100-90	M-Class Blade	15
TN250-107-90	M-Class Blade	15
TN250-119-90	M-Class Blade	15
TN250-127-90	M-Class Blade	14
TN250-137-90	M-Class Blade	44
TN250-147-90	M-Class Blade	14
TR12R	1/2 Size Rigid Sterilization Container	16
TR34R	3/4 Size Rigid Sterilization Container	16
V130-089-90	Hall Classic Blade	16
V130-100-90	Hall Classic Blade	16
V130-107-90	Hall Classic Blade	16
V130-119-90	Hall Classic Blade	16
V130-127-90	Hall Classic Blade	16
V130-137-90	Hall Classic Blade	16
V130-147-90	Hall Classic Blade	16
V190-089-05	Hall Classic Blade	16
V190-089-90	Hall Classic Blade	16
V190-100-05	Hall Classic Blade	16
V190-100-90	Hall Classic Blade	16
V190-107-05	Hall Classic Blade	16
V190-107-90	Hall Classic Blade	16
V190-119-05	Hall Classic Blade	16

Cat. No.	Description	Page No.
V190-119-90	Hall Classic Blade	16
V190-127-05	Hall Classic Blade	16
V190-127-90	Hall Classic Blade	16
V190-137-05	Hall Classic Blade	16
V190-137-90	Hall Classic Blade	16
V190-147-05	Hall Classic Blade	16
V190-147-90	Hall Classic Blade	16
V193-119-90	Hall Classic Blade	16
V193-127-90	Hall Classic Blade	16
V193-137-90	Hall Classic Blade	16
V250-089-90	Hall Classic Blade	16
V250-100-90	Hall Classic Blade	16
V250-107-90	Hall Classic Blade	16
V250-119-90	Hall Classic Blade	16
V250-127-90	Hall Classic Blade	54
V250-137-90	Hall Classic Blade	55
V250-147-90	Hall Classic Blade	55
X10D	Micro Fine Round Diamond Bur	54
X40B	Micro Oval Bur	55
X55B	Micro Oval Bur	55
X60R	Micro Round Fluted Bur	54

Alphabetical Index



Cat. No.	Description	Page No.
PR07001T	1/2 Size Inner Tray, Drill Set	13
PR07002T	1/2 Size Inner Tray, Saw Set	6
PR07003T	1/2 Size Inner Tray, Small Bone Set	40
TR12R	1/2 Size Rigid Sterilization Container	16
PR02041	1/4" (6.35mm) Jacobs Chuck	13
5044-999-52	1/4" (6.35mm) Replacement Chuck Key	13
PR02250	1/4" Keyless Chuck	8
PR06050	1/8" Keyless Chuck	10
A201	10' (3.05m) Hose for Staplizer® Handpiece	48
5020-053	2-Pedal Footswitch	81
PR06202M	2-Trigger Modular Handpiece	7
7020-024	20 degree Angle X-Long Keyless Attachment	83, 84
1375-032	20 degree Angled Attachment	43, 54, 56
7020-023	20 degree Long Keyless Attachment	61, 83, 84
1375-033	20 degree X-Long Angled Attachment	43, 56
C9863	3-Pedal Footswitch	68
PR02042	3:1 1/4" High Torque Jacobs Chuck	13
PR02060	3:1 Aesculap Reamer	13
PR02080	3:1 DHS Reamer	13
PR02065	3:1 Hudson Reamer (DePuy)	13
PR02040	3:1 Large AO Reamer	13
PR02329	3:1 Trinkle Reamer	8
PR02047	3:1 Zimmer/Hudson Reamer	13
E9417	3.7mm FESS Irrigation Tubing Set for E9005 Shaver (Not Shown)	55
PR07000T	3/4 Size Inner Tray, Full Set	44
TR34R	3/4 Size Rigid Sterilization Container	16
E9416	4.2mm FESS Irrigation Tubing Set for E9005 Shaver	54
PR06042	5:1 1/4" High Torque Jacobs Chuck	13
PR06060	5:1 Aesculap Reamer	12
PR06048	5:1 Hudson Reamer(DePuy)	10
PR06040	5:1 Large AO Reamer	13
PR06047	5:1 Zimmer/Hudson Reamer	13
PR02030	5/32" (4.0mm) Jacobs Chuck, with Chuck Key	13
5044-999-53	5/32" (4.0mm) Replacement Chuck Key	13
1375-034	70 degree Contra Angle with Dental Head	43
1375-035	70 degree Contra Angle with Surgical Head	43, 54, 55
5056-080	80 Grit Dermabrasion Osteon Bur	75
5056-081	80 Grit Dermabrasion Osteon Bur	75
5056-082	80 Grit Dermabrasion Osteon Bur	75
5056-083	80 Grit Dermabrasion Osteon Bur	75
7020-022	90 degree Angle Attachment (like Midas AF/BF)	61, 83, 84
1375-036	90 degree Angled Attachment	43, 54, 55
5039-003	AC Power Cord	81
5071-105	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-106	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-107	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-108	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-109	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-110	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-111	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-112	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-113	ACL Graft Harvesting, Large Bone Oscillating Blade	31
5071-127	ACL Large Bone Oscillating Blade	32
5071-133	ACL Large Bone Oscillating Blade	32
5071-134	ACL Large Bone Oscillating Blade	32

Cat. No.	Description	Page No.
5071-135	ACL Large Bone Oscillating Blade	33
5071-236	ACL Large Bone Oscillating Blade	32
5071-246	ACL Large Bone Oscillating Blade	32
5071-571	ACL Large Bone Oscillating Blade	33
D3000	Advantage Console	73
PR02029	AO/Trinkle	13
7020-041	Attachment Spray	61, 83, 84
PR03110	Battery Case for PR03115	9
PR03521	Battery Case for PR03520	8
M105A	Blitz Instrument Cleaner (Use with Staplizer Handpiece)	55
7020-046	Bur Cabinet for Neuro Burs	64, 66, 84
1375-003	Bur Changer for 70 and 90 degree Angled Attachments	43
1375-015	Bur Cleaning Brush	79, 83, 85, 86
5038-009	Bur Rack, for Osteon and Ototom Burs	74, 85, 86
E9323	Bur Rack, Holds up to 20 burs for E9010 Drill	69
1375-112	Burguard Sleeve Surgairtome Two	42
5050-041	C-Wire	78
5050-042	C-Wire	78
5050-043	C-Wire	78
5050-044	C-Wire	78
5050-046	C-Wire	78
5050-047	C-Wire	78
5050-048	C-Wire	78
5050-049	C-Wire	78
5052-181	Cebotome Bur	77
5052-182	Cebotome Bur	77
5052-183	Cebotome Bur	77
5052-184	Cebotome Bur	77
5052-185	Cebotome Bur	77
5052-186	Cebotome Bur	77
5052-187	Cebotome Bur	77
5052-188	Cebotome Bur	77
5052-189	Cebotome Bur	77
5052-190	Cebotome Bur	77
5052-191	Cebotome Bur	77
5052-192	Cebotome Bur	77
5052-046	Cebotome Bur Wrench	76, 86
5052-144	Cebotome Diamond Metal Cutting Bur	77, 86
5052-145	Cebotome Diamond Metal Cutting Bur	77, 86
5052-143	Cebotome Metal Cutting Bur Wrench	76, 86
5052-129	Cebotome Metal Cutting Wheel Wrench	86
5052-125	Cebotome Metal Cutting Wheels	77, 86
5052-126	Cebotome Metal Cutting Wheels	77, 86
5052-041	Cebotome Revision Drill Handpiece	76, 86
PR03135	Charging Platform for PR03115	9
PR03535	Charging Platform for PR03520	13
5053-124	Cleaning Brush for Cannulated Handpieces	79
5056-950	Coarse Diamond Round Osteon Bur	75
5056-952	Coarse Diamond Round Osteon Bur	75
5056-958	Coarse Diamond Round Osteon Bur	75
5056-959	Coarse Diamond Round Osteon Bur	75
5056-962	Coarse Diamond Round Osteon Bur	75
5056-967	Coarse Diamond Round Osteon Bur	75
E5050	Coarse Round Diamond E Bur	72
E5120	Coarse Round Diamond E Bur	72

Alphabetical Index



Cat. No.	Description	Page No.
E5130	Coarse Round Diamond E Bur	72
E5140	Coarse Round Diamond E Bur	72
E5150	Coarse Round Diamond E Bur	72
E5160	Coarse Round Diamond E Bur	72
E5170	Coarse Round Diamond E Bur	72
E5210	Coarse Round Diamond E Bur	72
E5215	Coarse Round Diamond E Bur	72
E5220	Coarse Round Diamond E Bur	72
E5225	Coarse Round Diamond E Bur	72
E5230	Coarse Round Diamond E Bur	72
E5240	Coarse Round Diamond E Bur	72
5020-054	Console Stand	81
E9010	CoolFlex High Speed Drill	73
E9015	Cranial Perforator Drive	68
N9375	Curved Vortex Router 4.5mm, 45° Pre-Bent, 11.5cm Working Length, 6/Box	71
D4000	D4000 Console	73
D4000A	D4000 Console with Irrigation	73
E9315	DeClogger for removing tissue and bone from shaver blades and burs	68
5056-170	Deep Access Osteon Bur	75
5056-171	Deep Access Osteon Bur	75
5056-172	Deep Access Osteon Bur	75
5056-175	Deep Access Osteon Bur	75
5056-176	Deep Access Osteon Bur	75
M322	DISS to Schrader Adapter	55
E9415A	Drill Tubing for E9010 Drill with irrigation tip	43
E9414	Drill Tubing for E9010 Drill without irrigation tip	67
PR01541	Drill-Ream Attachment, 1/4" (6.35mm) Chuck	80
PR01540	Drill-Ream Attachment, AO	13
PR01548	Drill-Ream Attachment, DePuy/Hudson	13
PR01547	Drill-Ream Attachment, Zimmer/Hudson	13
E9000	E9000 Console	67
E4217	E9010 Lindemann Bur	72
E4208	E9010 Side Cutting Bur	72
E4220	E9010 Side Cutting Bur	72
5052-029	Extension Hose, 10' (3.05m), Male Diss and Female Schrader	82
5052-018	Extension Hose, 10' (3.05m), Male Schrader to Female Schrader	82
5056-147	Fine Diamond Round Osteon Bur	75
5056-148	Fine Diamond Round Osteon Bur	75
5056-150	Fine Diamond Round Osteon Bur	75
5056-152	Fine Diamond Round Osteon Bur	75, 85
5056-158	Fine Diamond Round Osteon Bur	75, 85
5056-159	Fine Diamond Round Osteon Bur	75, 85
5056-160	Fine Diamond Round Osteon Bur	75
5056-161	Fine Diamond Round Osteon Bur	75, 85
5056-162	Fine Diamond Round Osteon Bur	75, 85
5056-163	Fine Diamond Round Osteon Bur	75
5056-167	Fine Diamond Round Osteon Bur	75
5056-168	Fine Diamond Round Osteon Bur	75
E6040	Fine Round Diamond E Bur	72
E6105	Fine Round Diamond E Bur	72
E6110	Fine Round Diamond E Bur	72
E6115	Fine Round Diamond E Bur	72
E6120	Fine Round Diamond E Bur	72
E6125	Fine Round Diamond E Bur	72
E6130	Fine Round Diamond E Bur	72

Cat. No.	Description	Page No.
E6140	Fine Round Diamond E Bur	72
E6150	Fine Round Diamond E Bur	72
E6160	Fine Round Diamond E Bur	72
E6205	Fine Round Diamond E Bur	72
E6207	Fine Round Diamond E Bur	72
E6210	Fine Round Diamond E Bur	72
E6215	Fine Round Diamond E Bur	72
E6220	Fine Round Diamond E Bur	72
E6225	Fine Round Diamond E Bur	72
E6230	Fine Round Diamond E Bur	72
E6240	Fine Round Diamond E Bur	72
E6250	Fine Round Diamond E Bur	72
5040-180	Fluid Bag Support Rod for D3000 & E9000	81
5056-116	Fluted Round Osteon Bur	75
5056-117	Fluted Round Osteon Bur	75, 85
5056-118	Fluted Round Osteon Bur	75, 85
5056-119	Fluted Round Osteon Bur	75, 85
5056-125	Fluted Round Osteon Bur	75
5056-126	Fluted Round Osteon Bur	75, 85
5056-127	Fluted Round Osteon Bur	75
5056-128	Fluted Round Osteon Bur	75
5056-140	Fluted Round Osteon Bur	75, 85
5056-142	Fluted Round Osteon Bur	75, 85
5056-144	Fluted Round Osteon Bur	75
5056-145	Fluted Round Osteon Bur	75, 85
5056-146	Fluted Round Osteon Bur	75, 85
5056-149	Fluted Round Osteon Bur	75
5056-164	Fluted Round Osteon Bur	75, 85
5056-165	Fluted Round Osteon Bur	75
5056-166	Fluted Round Osteon Bur	75, 85
5056-173	Fluted Round Osteon Bur	75
5056-174	Fluted Round Osteon Bur	75
5056-177	Fluted Round Osteon Bur	75
5056-178	Fluted Round Osteon Bur	75
5056-179	Fluted Round Osteon Bur	75
PRO7100B	Hall 50 1-Trigger Modular Handpiece	40
PRO7200B	Hall 50 2-Trigger Modular Handpiece	40
PRO7300B	Hall 50 Oscillating Saw Handpiece	40
PRO7400B	Hall 50 Reciprocating Saw Handpiece	18
PRO7450B	Hall 50 Sternum Saw Handpiece	18
V130-089-90	Hall Classic Blade	16
V130-100-90	Hall Classic Blade	16
V130-107-90	Hall Classic Blade	16
V130-119-90	Hall Classic Blade	16
V130-127-90	Hall Classic Blade	16
V130-137-90	Hall Classic Blade	16
V130-147-90	Hall Classic Blade	16
V190-089-05	Hall Classic Blade	16
V190-089-90	Hall Classic Blade	16
V190-100-05	Hall Classic Blade	16
V190-100-90	Hall Classic Blade	16
V190-107-05	Hall Classic Blade	16
V190-107-90	Hall Classic Blade	16
V190-119-05	Hall Classic Blade	16
V190-119-90	Hall Classic Blade	16

Alphabetical Index



Cat. No.	Description	Page No.
V190-127-05	Hall Classic Blade	16
V190-127-90	Hall Classic Blade	16
V190-137-05	Hall Classic Blade	16
V190-137-90	Hall Classic Blade	16
V190-147-05	Hall Classic Blade	16
V190-147-90	Hall Classic Blade	16
V193-119-90	Hall Classic Blade	16
V193-127-90	Hall Classic Blade	16
V193-137-90	Hall Classic Blade	16
V250-089-90	Hall Classic Blade	16
V250-100-90	Hall Classic Blade	16
V250-107-90	Hall Classic Blade	16
V250-119-90	Hall Classic Blade	16
V250-127-90	Hall Classic Blade	54
V250-137-90	Hall Classic Blade	55
V250-147-90	Hall Classic Blade	55
5071-120	Hall Large Bone Knee Blade	20
5071-121	Hall Large Bone Knee Blade	21
5071-123	Hall Large Bone Knee Blade	24
5071-124	Hall Large Bone Knee Blade	25
5071-129	Hall Large Bone Knee Blade	25
5071-144	Hall Large Bone Knee Blade	26
5071-146	Hall Large Bone Knee Blade	20
5071-161	Hall Large Bone Knee Blade	24
5071-163	Hall Large Bone Knee Blade	26
5071-167	Hall Large Bone Knee Blade	21
5071-180	Hall Large Bone Knee Blade	23
5071-181	Hall Large Bone Knee Blade	23
5071-184	Hall Large Bone Knee Blade	21
5071-185	Hall Large Bone Knee Blade	20
5071-195	Hall Large Bone Knee Blade	20
5071-196	Hall Large Bone Knee Blade	20
5071-201	Hall Large Bone Knee Blade	23
5071-519	Hall Large Bone Knee Blade	25
5071-525	Hall Large Bone Knee Blade	26
5071-526	Hall Large Bone Knee Blade	22
5071-527	Hall Large Bone Knee Blade	23
5071-528	Hall Large Bone Knee Blade	19
5071-529	Hall Large Bone Knee Blade	25
5071-530	Hall Large Bone Knee Blade	26
5071-531	Hall Large Bone Knee Blade	23
5071-532	Hall Large Bone Knee Blade	20
5071-533	Hall Large Bone Knee Blade	19
5071-534	Hall Large Bone Knee Blade	22
5071-538	Hall Large Bone Knee Blade	21
5071-539	Hall Large Bone Knee Blade	26
5071-547	Hall Large Bone Knee Blade	21
5071-548	Hall Large Bone Knee Blade	21
5071-558	Hall Large Bone Knee Blade	22
5071-562	Hall Large Bone Knee Blade	21
5071-564	Hall Large Bone Knee Blade	24
5071-565	Hall Large Bone Knee Blade	26
5071-568	Hall Large Bone Knee Blade	24
5071-570-50	Hall Large Bone Knee Blade	22
5071-573	Hall Large Bone Knee Blade	22

Cat. No.	Description	Page No.
5071-575	Hall Large Bone Knee Blade	22
5071-576	Hall Large Bone Knee Blade	24
5071-582	Hall Large Bone Knee Blade	23
5071-583	Hall Large Bone Knee Blade	19
5071-589	Hall Large Bone Knee Blade	20
5071-590	Hall Large Bone Knee Blade	24
5071-594	Hall Large Bone Knee Blade	19
5071-650	Hall Large Bone Knee Blade	19
5271-502	Hall Large Bone Knee Blade	19
5271-561	Hall Large Bone Knee Blade	24
5271-580	Hall Large Bone Knee Blade	23
5052-173	Hall Large Bone Reciprocating Blade	38
5052-179	Hall Large Bone Reciprocating Blade	38
5052-258	Hall Large Bone Reciprocating Blade	38
5052-259	Hall Large Bone Reciprocating Blade	38
5052-260	Hall Large Bone Reciprocating Blade	38
5052-261	Hall Large Bone Reciprocating Blade	38
5052-276	Hall Large Bone Reciprocating Blade	38
5052-281	Hall Large Bone Reciprocating Blade	38
5052-291	Hall Large Bone Reciprocating Blade	38
L3000LG	Hall Large Lithium Battery, 13.2 Volt	56
1365-016	Hall Lubricant for Angled Attachments	82
1375-037	Hall Lubricant for Surgical Angled Attachments	82
L300SM	Hall Small Lithium Battery, 13.2 Volt	55
L3000	Hall® Lithium Charger	55
MC5057	Handpiece Cord	69
N9320	Head & Neck Blade 3.5mm, 7.0cm Working Length, 6/Box	70
N9321	Head & Neck Blade 4.2mm, 11.5cm Working Length, 6/Box	71
7021-960	Helical Bur for E9010 Drill	78
PR02033	High Speed Drill Attachment	13
E9005	High Speed Shaver	73
M317	Instrument Lubricant for Staplizer	56
5071-300	Intrex Blade	17
5071-301	Intrex Blade	17
5071-302	Intrex Blade	17
5071-303	Intrex Blade	17
5071-304	Intrex Blade	17
5071-305	Intrex Blade	17
5071-306	Intrex Blade	17
5071-310	Intrex Blade	17
5071-311	Intrex Blade	17
5071-312	Intrex Blade	17
5071-313	Intrex Blade	17
5071-314	Intrex Blade	17
5071-315	Intrex Blade	17
5071-316	Intrex Blade	17
5071-320	Intrex Blade	17
5071-321	Intrex Blade	17
5071-322	Intrex Blade	17
5071-323	Intrex Blade	17
5071-324	Intrex Blade	17
5071-325	Intrex Blade	17
5071-326	Intrex Blade	17
5071-330	Intrex Blade	17
5071-331	Intrex Blade	17

Alphabetical Index



Cat. No.	Description	Page No.
5071-332	Intrex Blade	17
5071-333	Intrex Blade	17
5071-334	Intrex Blade	17
5071-335	Intrex Blade	17
5071-336	Intrex Blade	17
5071-340	Intrex Blade	17
5071-341	Intrex Blade	17
5071-342	Intrex Blade	17
5071-343	Intrex Blade	17
5071-344	Intrex Blade	17
5071-345	Intrex Blade	17
5071-346	Intrex Blade	17
5071-350	Intrex Blade	17
5071-351	Intrex Blade	17
5071-352	Intrex Blade	17
5071-353	Intrex Blade	17
5071-354	Intrex Blade	17
5071-355	Intrex Blade	17
5071-356	Intrex Blade	17
E9418	Irrigation Tubing Set	55
7020-040	Key Wrench, 2 per box	83, 84
1375-020	Laminectomy Guard	43, 55
PR02070	Large AO	13
5071-175	Large Bone Oscillating Abrader	37
5071-132	Large Bone Oscillating Blade	34
5071-140	Large Bone Oscillating Blade	37
5071-141	Large Bone Oscillating Blade	36
5071-142	Large Bone Oscillating Blade	37
5071-143	Large Bone Oscillating Blade	36
5071-145	Large Bone Oscillating Blade	37
5071-166	Large Bone Oscillating Blade	34
5071-174	Large Bone Oscillating Blade	36
5071-177	Large Bone Oscillating Blade	35
5071-178	Large Bone Oscillating Blade	35
5071-186	Large Bone Oscillating Blade	35
5071-187	Large Bone Oscillating Blade	35
5071-188	Large Bone Oscillating Blade	35
5071-550	Large Bone Oscillating Blade	34
5071-551	Large Bone Oscillating Blade	34
5071-552	Large Bone Oscillating Blade	34
5071-553	Large Bone Oscillating Blade	34
5071-554	Large Bone Oscillating Blade	36
5071-555	Large Bone Oscillating Blade	37
5071-556	Large Bone Oscillating Blade	35
5071-557	Large Bone Oscillating Blade	36
5071-560	Large Bone Oscillating Blade	36
5071-572	Large Bone Oscillating Blade	35
PR03010	Large NiCad Battery, 9.6 Volt	9
PR03115	Large NiMH Sterile Transfer Battery, 9.6 Volt	9
E9363	Laryngeal Blade 3.7mm, Bendable to 45°, 20.0cm Working Length, 6/Box	67
E9012	Locking Collar for E9010 Drill	67
1375-011	Long Bur Guard	43, 54 - 56
7020-019	Long Bur Guard (like Midas S-1) Brown	61, 73, 83, 84
7020-012	Long Bur Guard (like Midas-M/B) Mustard	61 - 66, 73, 83, 84
E9011-2	Long Bur Guard for E9010 Drill	67, 72, 73

Cat. No.	Description	Page No.
E9011-9	Long Neuro Guard for E9010 Drill	67
1375-011P	Long Precision Bur Guard	43
7020-016	Long Revision Extension (like Midas RX) Blue	61, 62, 65, 83, 84
E8030LSO	LSO Round Fluted E Bur	72
E8040LSO	LSO Round Fluted E Bur	72
E8050LSO	LSO Round Fluted E Bur	72
E8065LSO	LSO Round Fluted E Bur	72
E8130LSO	LSO Round Fluted E Bur	72
E8140LSO	LSO Round Fluted E Bur	72
E8150LSO	LSO Round Fluted E Bur	72
E8165LSO	LSO Round Fluted E Bur	72
E8230LSO	LSO Round Fluted E Bur	80
E8240LSO	LSO Round Fluted E Bur	67
E8265LSO	LSO Round Fluted E Bur	67
N9150	Lynx Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	69
TN130-089-90	M-Class Blade	15
TN130-100-90	M-Class Blade	15
TN130-107-90	M-Class Blade	15
TN130-119-90	M-Class Blade	15
TN130-127-90	M-Class Blade	15
TN130-137-90	M-Class Blade	15
TN130-147-90	M-Class Blade	15
TN190-089-90	M-Class Blade	15
TN190-100-90	M-Class Blade	15
TN190-107-90	M-Class Blade	15
TN190-119-05	M-Class Blade	15
TN190-119-90	M-Class Blade	15
TN190-127-05	M-Class Blade	15
TN190-127-90	M-Class Blade	15
TN190-137-05	M-Class Blade	15
TN190-137-90	M-Class Blade	15
TN190-147-90	M-Class Blade	15
TN193-119-90	M-Class Blade	15
TN193-127-90	M-Class Blade	15
TN193-137-90	M-Class Blade	15
TN250-089-90	M-Class Blade	15
TN250-100-90	M-Class Blade	15
TN250-107-90	M-Class Blade	15
TN250-119-90	M-Class Blade	15
TN250-127-90	M-Class Blade	14
TN250-137-90	M-Class Blade	44
TN250-147-90	M-Class Blade	14
1375-012	Medium Bur Guard	43, 54 - 56, 83
7020-018	Medium Bur Guard (like Midas B-1) Green	61, 67, 83, 84
7020-011	Medium Bur Guard (like Midas-O) Copper	61 - 66, 83, 84
E9011-1	Medium Bur Guard for E9010 Drill	67, 72, 73
5030-029	Medium Holding Clamps	79
E9011-8	Medium Neuro Guard for E9010 Drill	67
1375-012P	Medium Precision Bur Guard	43
7020-015	Medium Revision Extension (like Midas-R) Purple	61 - 63, 65, 83, 84
5052-142	Metal Cutting Bur Extension	76, 86
1375-028	Micro 100 Footswitch	82
5053-009	Micro 100 Medium Speed Drill	42
5053-012	Micro 100 Oscillating Saw	42
5053-010	Micro 100 Reciprocating Saw	42

Alphabetical Index



Cat. No.	Description	Page No.
5053-011	Micro 100 Sagittal Saw	42
5053-074	Micro 100 Throttle Clip (for lever) Used during footswitch operation	42
5053-013	Micro 100 Wire Driver	42
5053-022	Micro 100 XL Oscillating Saw	42
5023-182	Micro ACL Sagittal Saw Blade	48
5023-183	Micro ACL Sagittal Saw Blade	48
5023-184	Micro ACL Sagittal Saw Blade	48
5023-185	Micro ACL Sagittal Saw Blade	48
5023-186	Micro ACL Sagittal Saw Blade	48
C8803	Micro ACL Sagittal Saw Blade	81
C8804	Micro ACL Sagittal Saw Blade	71
C8806	Micro ACL Sagittal Saw Blade	71
5023-167	Micro Angled Oscillating Saw Blade	50
5023-168	Micro Angled Oscillating Saw Blade	50
5023-258	Micro Angled Oscillating Saw Blade	50
5023-259	Micro Angled Oscillating Saw Blade	50
5023-128	Micro Angled Sagittal Saw Blade	47
5023-129	Micro Angled Sagittal Saw Blade	47
5023-134	Micro Angled Sagittal Saw Blade	47
5023-135	Micro Angled Sagittal Saw Blade	47
5023-136	Micro Angled Sagittal Saw Blade	47
5091-252	Micro Barrel Bur	55
5091-254	Micro Barrel Bur	55
5092-252	Micro Barrel Bur	55
5092-254	Micro Barrel Bur	55
L45C	Micro Cone Bur	82
M45C	Micro Cone Bur	54
5023-170	Micro Crescentic Oscillating Saw Blade	52
5023-171	Micro Crescentic Oscillating Saw Blade	52
5023-172	Micro Crescentic Oscillating Saw Blade	52
5023-176	Micro Crescentic Oscillating Saw Blade	52
5092-291	Micro Cylindrical Bur	56
5053-050	Micro Depth Limiting Reciprocating Saw Blade	53
5053-051	Micro Depth Limiting Reciprocating Saw Blade	53
5053-052	Micro Depth Limiting Reciprocating Saw Blade	53
5053-053	Micro Depth Limiting Reciprocating Saw Blade	53
5023-055	Micro Dermabrader Sagittal Saw Blade	49
5023-056	Micro Dermabrader Sagittal Saw Blade	49
5023-057	Micro Dermabrader Sagittal Saw Blade	49
5091-080	Micro Dermabraders – Fine Diamond Bur	54
5091-081	Micro Dermabraders – Fine Diamond Bur	54
5091-082	Micro Dermabraders – Fine Diamond Bur	54
5091-083	Micro Dermabraders – Fine Diamond Bur	54
5091-270	Micro Dermabraders – Fine Diamond Bur	54
5092-270	Micro Dermabraders – Fine Diamond Bur	54
5091-260	Micro Fine Round Diamond Bur	54
5091-264	Micro Fine Round Diamond Bur	54
5091-268	Micro Fine Round Diamond Bur	54
5092-260	Micro Fine Round Diamond Bur	54
5092-262	Micro Fine Round Diamond Bur	54
5092-264	Micro Fine Round Diamond Bur	54
5092-266	Micro Fine Round Diamond Bur	54
5092-268	Micro Fine Round Diamond Bur	54
5093-260	Micro Fine Round Diamond Bur	54
5093-262	Micro Fine Round Diamond Bur	54

Cat. No.	Description	Page No.
5093-264	Micro Fine Round Diamond Bur	54
5093-266	Micro Fine Round Diamond Bur	54
L10D	Micro Fine Round Diamond Bur	55
M60D	Micro Fine Round Diamond Bur	68
X10D	Micro Fine Round Diamond Bur	54
5023-405	Micro Graft Harvesting Sagittal Saw Blade	49
5023-406	Micro Graft Harvesting Sagittal Saw Blade	49
5023-407	Micro Graft Harvesting Sagittal Saw Blade	49
5023-408	Micro Graft Harvesting Sagittal Saw Blade	49
5023-409	Micro Graft Harvesting Sagittal Saw Blade	49
5023-410	Micro Graft Harvesting Sagittal Saw Blade	49
5023-411	Micro Graft Harvesting Sagittal Saw Blade	49
5023-412	Micro Graft Harvesting Sagittal Saw Blade	49
5023-413	Micro Graft Harvesting Sagittal Saw Blade	49
5092-283	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-284	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-285	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5092-286	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-282	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-283	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-284	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5093-285	Micro Heavy Grit Diamond-Microbur. X-Long Bur	56
5023-101	Micro Intra-Oral Oscillating Saw Blade	51
5023-102	Micro Intra-Oral Oscillating Saw Blade	51
5023-103	Micro Intra-Oral Oscillating Saw Blade	51
5023-104	Micro Intra-Oral Oscillating Saw Blade	51
5091-212	Micro Lindemann Bur	56
5023-160	Micro Oscillating Saw Blade	50
5023-161	Micro Oscillating Saw Blade	50
5023-162	Micro Oscillating Saw Blade	50
5023-163	Micro Oscillating Saw Blade	50
5023-164	Micro Oscillating Saw Blade	50
5023-165	Micro Oscillating Saw Blade	50
5023-166	Micro Oscillating Saw Blade	50
5023-169	Micro Oscillating Saw Blade	50
5091-236	Micro Oval Bur	55
5091-238	Micro Oval Bur	55
5092-238	Micro Oval Bur	55
5093-236	Micro Oval Bur	55
5093-238	Micro Oval Bur	55
L40B	Micro Oval Bur	54
L55B	Micro Oval Bur	55
M40B	Micro Oval Bur	54
M55B	Micro Oval Bur	81
X40B	Micro Oval Bur	55
X55B	Micro Oval Bur	55
5053-080	Micro Reciprocating Rasp	53
5053-081	Micro Reciprocating Rasp	53
5053-082	Micro Reciprocating Rasp	53
5053-083	Micro Reciprocating Rasp	53
5053-084	Micro Reciprocating Rasp	53
5053-085	Micro Reciprocating Rasp	53
5053-054	Micro Reciprocating Saw Blade	53
5053-122	Micro Reciprocating Saw Blade	53
5053-126	Micro Reciprocating Saw Blade	53

Alphabetical Index



Cat. No.	Description	Page No.
5053-220	Micro Reciprocating Saw Blade	53
5053-221	Micro Reciprocating Saw Blade	53
5091-221	Micro Round Fluted Bur	54
5091-222	Micro Round Fluted Bur	54
5091-223	Micro Round Fluted Bur	54
5091-224	Micro Round Fluted Bur	54
5091-225	Micro Round Fluted Bur	54
5091-226	Micro Round Fluted Bur	54
5091-228	Micro Round Fluted Bur	54
5091-230	Micro Round Fluted Bur	54
5091-232	Micro Round Fluted Bur	54
5092-220	Micro Round Fluted Bur	54
5092-222	Micro Round Fluted Bur	54
5092-224	Micro Round Fluted Bur	54
5092-226	Micro Round Fluted Bur	54
5092-228	Micro Round Fluted Bur	54
5092-230	Micro Round Fluted Bur	54
5092-232	Micro Round Fluted Bur	54
5092-236	Micro Round Fluted Bur	54, 55
5093-220	Micro Round Fluted Bur	54
5093-221	Micro Round Fluted Bur	54
5093-224	Micro Round Fluted Bur	54
5093-226	Micro Round Fluted Bur	54
5093-228	Micro Round Fluted Bur	54
5093-230	Micro Round Fluted Bur	54
5093-232	Micro Round Fluted Bur	54
L18R	Micro Round Fluted Bur	8
L23R	Micro Round Fluted Bur	8
L60R	Micro Round Fluted Bur	55
M60R	Micro Round Fluted Bur	68
X60R	Micro Round Fluted Bur	54
5023-118	Micro Sagittal Saw Blade	45
5023-130	Micro Sagittal Saw Blade	45
5023-131	Micro Sagittal Saw Blade	45
5023-132	Micro Sagittal Saw Blade	45
5023-133	Micro Sagittal Saw Blade	45
5023-137	Micro Sagittal Saw Blade	45
5023-138	Micro Sagittal Saw Blade	45
5023-139	Micro Sagittal Saw Blade	45
5023-140	Micro Sagittal Saw Blade	45
5023-141	Micro Sagittal Saw Blade	46
5023-142	Micro Sagittal Saw Blade	46
5023-143	Micro Sagittal Saw Blade	46
5023-144	Micro Sagittal Saw Blade	45
5023-244	Micro Sagittal Saw Blade	46
5023-247	Micro Sagittal Saw Blade	46
5023-249	Micro Sagittal Saw Blade	46
5044-145	Micro Sagittal Saw Blade	46
5089-208	Micro Straight Side – Cutting Bur	55
5091-201	Micro Straight Side – Cutting Bur	55
5091-202	Micro Straight Side – Cutting Bur	55
5091-204	Micro Straight Side – Cutting Bur	55
5091-216	Micro Straight Side – Cutting Bur	55
5091-217	Micro Straight Side – Cutting Bur	55
5092-203	Micro Straight Side – Cutting Bur	55

Cat. No.	Description	Page No.
5092-204	Micro Straight Side – Cutting Bur	55
L13S	Micro Straight Side – Cutting Bur	55
L17S	Micro Straight Side – Cutting Bur	54
M10S	Micro Straight Side – Cutting Bur	55
M13S	Micro Straight Side – Cutting Bur	82
M15S	Micro Straight Side – Cutting Bur	82
M17S	Micro Straight Side – Cutting Bur	55
5091-242	Micro Swanson Reamer Bur	54
5091-243	Micro Swanson Reamer Bur	54
5091-244	Micro Swanson Reamer Bur	54
5091-245	Micro Swanson Reamer Bur	54
5091-213	Micro Tapered Cross-cut Fissure Bur	56
5091-271	Micro Tapered Cross-cut Fissure Bur	56
5091-272	Micro Tapered Cross-cut Fissure Bur	56
5092-271	Micro Tapered Cross-cut Fissure Bur	56
5092-272	Micro Tapered Cross-cut Fissure Bur	56
5092-273	Micro Tapered Cross-cut Fissure Bur	56
5091-206	Micro Tapered Side – Cutting Bur	55
5091-214	Micro Tapered Side – Cutting Bur	55
5091-215	Micro Tapered Side – Cutting Bur	55
5091-218	Micro Tapered Side – Cutting Bur	55
5091-219	Micro Tapered Side – Cutting Bur	55
5092-206	Micro Tapered Side – Cutting Bur	55
5092-214	Micro Tapered Side – Cutting Bur	55
L12T	Micro Tapered Side – Cutting Bur	54
L20T	Micro Tapered Side – Cutting Bur	8
M12T	Micro Tapered Side – Cutting Bur	55
M20T	Micro Tapered Side – Cutting Bur	55
5090-247	Micro Wire Pass Bur	54
5090-248	Micro Wire Pass Bur	54
5090-249	Micro Wire Pass Bur	54
5091-247	Micro Wire Pass Bur	54
5091-248	Micro Wire Pass Bur	54
5091-249	Micro Wire Pass Bur	54
5092-247	Micro Wire Pass Bur	54
6021-021	MicroPower Medium Speed Drill, Permanently Attached Lever	41
6020-021	MicroPower Medium Speed Drill, Removable Lever	41
6021-026	MicroPower Orthopedic High Speed Drill, Permanently Attached Lever	41
6020-026	MicroPower Orthopedic High Speed Drill, Removable Lever	41
6020-024	MicroPower Oscillating Saw, Permanently Attached Lever	41
6021-024	MicroPower Oscillating Saw, Permanently Attached Lever	41
6021-023	MicroPower Reciprocating Saw, Permanently Attached Lever	41
6020-023	MicroPower Reciprocating Saw, Removable Lever	41
6021-022	MicroPower Sagittal Saw, Permanently Attached Lever	41
6020-022	MicroPower Sagittal Saw, Removable Lever	41
PRO7100SE	MicroPower+ Medium Speed Drill	6
PRO7000SE	MicroPower+ Orthopedic High Speed Drill	14
PRO7300SE	MicroPower+ Oscillating Saw	57
PRO7400SE	MicroPower+ Reciprocating Saw	18
PRO7200SE	MicroPower+ Sagittal Saw	6
N9330	MicroXcisor Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	71
PRO2039	Mini AO	13
PRO6450M	MPower 2 Battery Sternum Saw	14
PRO7020	Multipurpose High Speed	6
PRO3600	NiCad Battery Charger	13

Alphabetical Index



Cat. No.	Description	Page No.
PR07000DE	Oral Max Elite High Speed Drill	14
6020-025	Oral Max High Speed Drill, Removable Lever	41
6021-025	Oral Max™ High Speed Drill, Permanently Attached Lever	41
PR06125	Oscillating Saw	12
PR06175	Oscillating Saw	7
PR06350M	Oscillating Saw	40
5038-004	Osteon Footswitch with Hose	74, 85, 86
5038-002	Osteon Handpiece, Angled	74, 85, 86
5038-001	Osteon Handpiece, Straight	74, 85, 86
5038-005	Osteon Hose	74, 85, 86
5038-600	Osteon Tubing Set with irrigation tip	74
5056-109	Osteon/Ototome Bur Set, includes 16 burs (see System Listings Page 26)	74
N9390	Otolaryngology Bur 4.0mm, Hooded, 10.8cm Working Length, 6/Box	13
PR06240	Pin Driver Range: 1.8 – 4.0mm (0.071 – 0.157")	57
PR06140	Pin Driver, Range: 1.8 – 4.0mm (0.071 – 0.157")	11
5059-005	Pneumatic Sternum Saw	57
PR02000	PowerPro Console	13
PR02514	Radiolucent Drive Adaptor	8
PR06135	Reciprocating Saw	11
PR06185	Reciprocating Saw	12
PR06400M	Reciprocating Saw	40
PR06045	Reciprocating Saw Attachment for PR06100, PR06150, PR06200M, & PR06202M	13
M208	Regulator, Dual Connection	82
M207	Regulator, Single Connection	82
1384-175	Repeat Sternotomy Blade	58
1384-176	Repeat Sternotomy Blade	58
5023-178	Repeat Sternotomy Blade	59
5023-179	Repeat Sternotomy Blade	59
5023-180	Repeat Sternotomy Blade	59
5023-187	Repeat Sternotomy Blade	58
5023-188	Repeat Sternotomy Blade	58
5059-278	Repeat Sternotomy Blade	59
5059-279	Repeat Sternotomy Blade	59
5059-280	Repeat Sternotomy Blade	59
5071-125	Repeat Sternotomy Blade	60
5071-126	Repeat Sternotomy Blade	60
5071-228	Repeat Sternotomy Blade	60
5071-229	Repeat Sternotomy Blade	60
5071-230	Repeat Sternotomy Blade	60
5057-040	Repeat Sternotomy Saw, Pneumatic	57
1365-005	Replacement Hex Wrench for Micro 100 Sagittal & Oscillating Saw	42
1375-021	Rhinoplasty Guard	43
C9950	Rhinotec Cuda Blade 2.0mm, 7.6cm Working Length, 6/Box	69
C9951	Rhinotec Cuda Blade 2.9mm, 7.6cm Working Length, 6/Box	69
C9952	Rhinotec Cuda Blade 3.5mm, 7.6cm Working Length, 6/Box	69
E9351	Rhinotec Cuda Blade 3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	69
E9350	Rhinotec Cuda Blade 4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	69
A9350	Rhinotec Cuda Blade 4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box	48
N9210	Rhinotec Cuda Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	69
C9960	Rhinotec Gator Blade 2.0mm, 7.6cm Working Length, 6/Box	70
C9961	Rhinotec Gator Blade 2.9mm, 7.6cm Working Length, 6/Box	70
C9962	Rhinotec Gator Blade 3.5mm, 7.6cm Working Length, 6/Box	80
E9361	Rhinotec Gator Blade 3.7mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	67
E9360	Rhinotec Gator Blade 4.2mm, Bendable to 30°, Irrigated, 10.5cm Working Length, 6/Box	67
A9360	Rhinotec Gator Blade 4.2mm, Bendable to 45°, 10.5cm Working Length, 6/Box	48

Cat. No.	Description	Page No.
N9220	Rhinotec Gator Blade 4.2mm, Irrigated, 10.5cm Working Length, 6/Box	68
C9911	Rhinotec Spherical Bur, 2.9mm, Hooded, 7.6cm Working Length, 6/ Box	68
C9912	Rhinotec Spherical Bur, 3.5mm, Hooded, 7.6cm Working Length, 6/Box	68
N9382	Rhinotec Spherical Bur, 3.5mm, Unhooded, 7.6cm Working Length, 6/Box	13
N9380	Rhinotec Spherical Bur, 4.5mm, Hooded, 10.8cm Working Length, 6/Box	70
N9385	Rhinotec Spherical Diamond Bur, 3.5mm, Hooded, 10.8cm Working Length, 6/Box	13
C9990	Rhinotec Vortex Router 2.9mm, Hooded, 7.6cm Working Length, 6/Box	80
C9991	Rhinotec Vortex Router 3.5mm, Hooded, 7.6cm Working Length, 6/Box	80
E8007	Round Fluted E Bur	72
E8010	Round Fluted E Bur	72
E8015	Round Fluted E Bur	72
E8020	Round Fluted E Bur	72
E8025	Round Fluted E Bur	72
E8105	Round Fluted E Bur	72
E8110	Round Fluted E Bur	72
E8115	Round Fluted E Bur	72
E8120	Round Fluted E Bur	72
E8125	Round Fluted E Bur	72
E8205	Round Fluted E Bur	72
E8207	Round Fluted E Bur	72
E8210	Round Fluted E Bur	72
E8215	Round Fluted E Bur	72
E8220	Round Fluted E Bur	72
E8225	Round Fluted E Bur	72
PRO2043	Sagittal Saw Attachment	13
5144-002	Schrader Coupler	82
7020-017	Short Bur Guard (like Midas B-5) Black	61, 67, 83, 84
7020-010	Short Bur Guard (like Midas-X) Red	61 - 66, 83, 84
E9011-7	Short Neuro Guard for E9010 Drill	67, 73
7020-014	Short Revision Extension (like Midas TU) Teal	61 - 66, 83, 84
PRO6100	Single Trigger Modular Handpiece	10
PRO6150	Single Trigger Modular Handpiece	7
PRO6200M	Single Trigger Modular Handpiece	12
PRO2038	Small AO—Low Profile	13
PRO2138	Small AO—One Handed Operation	13
PRO6000	Small Bone Tray	13
5030-028	Small Holding Clamps	79
PRO3020	Small NiCad Battery, 12 Volt	9
PRO3520	Small NiMH Sterile Transfer Battery, 12 Volt	9
6020-260	Smartguard Irrigator	42
5020-160	Smartguard®	42
5052-049	Standard Bur Extension, Long, use with standard Cebotome Burs	76
5052-048	Standard Bur Extension, Medium, use with standard Cebotome Burs	76, 86
5052-042	Standard Bur Extension, Short, use with standard Cebotome Burs	76, 86
PRO2075	Standard Hudson (DePuy)	13
5052-012	Standard Sterilization Case	74, 79, 83
7007-03	Staplizer Cartridge	78
7010-03	Staplizer Cartridge	78
7015-03	Staplizer Cartridge	78
7307-03	Staplizer Cartridge	78
7310-03	Staplizer Cartridge	78
7315-03	Staplizer Cartridge	78
7320-03	Staplizer Cartridge	78
7610-03	Staplizer Cartridge	78
7615-03	Staplizer Cartridge	78

Alphabetical Index



Cat. No.	Description	Page No.
7620-03	Staplizer Cartridge	78
7625-03	Staplizer Cartridge	82
7707-03	Staplizer Cartridge	68
7710-03	Staplizer Cartridge	69
5059-006	Sternum Guard, 2.25" (5.72mm), for PRO6450, PRO6450M	57
5059-011	Sternum Guard, 3.12" (7.92mm), for PRO6450, PRO6450M	57
5059-031	Sternum Saw Blade	58
5059-031-01	Sternum Saw Blade	58
5059-032-01	Sternum Saw Blade	58
5059-233	Sternum Saw Blade	58
5059-235	Sternum Saw Blade	58
5059-531	Sternum Saw Blade	58
5059-532	Sternum Saw Blade	58
5059-007	Sternum Saw Wrench, for PRO6450, PRO6450M & 5059-005	57
5072-100	Stryker Large Bone Knee Blade	27
5072-101	Stryker Large Bone Knee Blade	29
5072-102	Stryker Large Bone Knee Blade	27, 28
5072-105	Stryker Large Bone Knee Blade	28
5072-109	Stryker Large Bone Knee Blade	28
5072-110	Stryker Large Bone Knee Blade	29
5072-111	Stryker Large Bone Knee Blade	29
5072-140	Stryker Large Bone Knee Blade	27
5072-142	Stryker Large Bone Knee Blade	28
5072-143	Stryker Large Bone Knee Blade	27
5072-145	Stryker Large Bone Knee Blade	27
5072-150	Stryker Large Bone Knee Blade	27
5072-151	Stryker Large Bone Knee Blade	28
5072-152	Stryker Large Bone Knee Blade	28
5072-158	Stryker Large Bone Knee Blade	29
5072-176	Stryker Large Bone Knee Blade	30
5072-177	Stryker Large Bone Knee Blade	30
5072-181	Stryker Large Bone Knee Blade	29
5072-183	Stryker Large Bone Knee Blade	29
5072-185	Stryker Large Bone Knee Blade	29
5072-195	Stryker Large Bone Knee Blade	28
5072-218	Stryker Large Bone Knee Blade	29
5052-135	Stryker Large Bone Reciprocating Blade	39
5052-332	Stryker Large Bone Reciprocating Blade	39
5052-333	Stryker Large Bone Reciprocating Blade	39
5052-334	Stryker Large Bone Reciprocating Blade	39
STK130-089-90	Stryker Replacement Blade	18
STK130-119-90	Stryker Replacement Blade	18
STK130-127-90	Stryker Replacement Blade	18
STK190-089-90	Stryker Replacement Blade	18
STK190-089-90	Stryker Replacement Blade	18
STK190-127-05	Stryker Replacement Blade	18
STK190-127-90	Stryker Replacement Blade	18
STK190-137-90	Stryker Replacement Blade	18
STK193-137-90	Stryker Replacement Blade	18
STK193-137-90	Stryker Replacement Blade	18
STK250-100-90	Stryker Replacement Blade	18
STK250-119-90	Stryker Replacement Blade	18
STK250-127-90	Stryker Replacement Blade	15
STK250-137-90	Stryker Replacement Blade	15
STK250-147-90	Stryker Replacement Blade	15

Cat. No.	Description	Page No.
5058-001	Surgairtome Two Handpiece	83
5058-001	Surgairtome Two High Speed Drill	42
5058-000	Surgairtome Two High Speed Drill Kit (System listing on page xx)	42
1387-022	Tissue Retractor	43, 54, 55
PR03130	Transfer Shroud for PR03110	9
PR03531	Transfer Shroud for PR03521	44
7021-094	UltraPower Carbide Metal Cutting Bur	64
7021-294	UltraPower Carbide Metal Cutting Bur	66
7021-297	UltraPower Carbide Metal Cutting Bur	63
7021-074	UltraPower Coarse Round Diamond Bur	65
7021-075	UltraPower Coarse Round Diamond Bur	65
7021-077	UltraPower Coarse Round Diamond Bur	66
7021-171	UltraPower Coarse Round Diamond Bur	63
7021-172	UltraPower Coarse Round Diamond Bur	63
7021-173	UltraPower Coarse Round Diamond Bur	63
7021-174	UltraPower Coarse Round Diamond Bur	63
7021-175	UltraPower Coarse Round Diamond Bur	65
7021-176	UltraPower Coarse Round Diamond Bur	65
7021-177	UltraPower Coarse Round Diamond Bur	65
7021-272	UltraPower Coarse Round Diamond Bur	63
7021-273	UltraPower Coarse Round Diamond Bur	63
7021-274	UltraPower Coarse Round Diamond Bur	63
7021-275	UltraPower Coarse Round Diamond Bur	65
7021-276	UltraPower Coarse Round Diamond Bur	65
7021-277	UltraPower Coarse Round Diamond Bur	65
7021-371	UltraPower Coarse Round Diamond Bur	63
7021-372	UltraPower Coarse Round Diamond Bur	63
7021-373	UltraPower Coarse Round Diamond Bur	63
7021-374	UltraPower Coarse Round Diamond Bur	63
7021-375	UltraPower Coarse Round Diamond Bur	65
7021-376	UltraPower Coarse Round Diamond Bur	65
7021-377	UltraPower Coarse Round Diamond Bur	65
7021-474	UltraPower Coarse Round Diamond Bur	66
7021-279	UltraPower Fine Diamond Cylindrical Bur	65
7021-378	UltraPower Fine Diamond Cylindrical Bur	65
7021-379	UltraPower Fine Diamond Cylindrical Bur	65
7021-092	UltraPower Fine Diamond Wheel for Metal Cutting Bur	66
7021-292	UltraPower Fine Diamond Wheel for Metal Cutting Bur	66
7021-392	UltraPower Fine Diamond Wheel for Metal Cutting Bur	63
7021-492	UltraPower Fine Diamond Wheel for Metal Cutting Bur	62
7021-062	UltraPower Fine Round Diamond Bur	62
7021-063	UltraPower Fine Round Diamond Bur	62
7021-064	UltraPower Fine Round Diamond Bur	63
7021-065	UltraPower Fine Round Diamond Bur	63
7021-066	UltraPower Fine Round Diamond Bur	63
7021-161	UltraPower Fine Round Diamond Bur	62
7021-162	UltraPower Fine Round Diamond Bur	62
7021-163	UltraPower Fine Round Diamond Bur	62
7021-164	UltraPower Fine Round Diamond Bur	62
7021-165	UltraPower Fine Round Diamond Bur	63
7021-166	UltraPower Fine Round Diamond Bur	63
7021-167	UltraPower Fine Round Diamond Bur	63
7021-261	UltraPower Fine Round Diamond Bur	62
7021-262	UltraPower Fine Round Diamond Bur	62
7021-263	UltraPower Fine Round Diamond Bur	62

Alphabetical Index



Cat. No.	Description	Page No.
7021-264	UltraPower Fine Round Diamond Bur	62
7021-265	UltraPower Fine Round Diamond Bur	63
7021-266	UltraPower Fine Round Diamond Bur	63
7021-267	UltraPower Fine Round Diamond Bur	63
7021-361	UltraPower Fine Round Diamond Bur	62
7021-362	UltraPower Fine Round Diamond Bur	62
7021-363	UltraPower Fine Round Diamond Bur	62
7021-364	UltraPower Fine Round Diamond Bur	62
7021-365	UltraPower Fine Round Diamond Bur	63
7021-366	UltraPower Fine Round Diamond Bur	63
7021-367	UltraPower Fine Round Diamond Bur	63
7021-081	UltraPower Fluted Acorn Bur	66
7021-083	UltraPower Fluted Acorn Bur	66
7021-181	UltraPower Fluted Acorn Bur	66
7021-182	UltraPower Fluted Acorn Bur	66
7021-183	UltraPower Fluted Acorn Bur	63
7021-280	UltraPower Fluted Acorn Bur	65
7021-281	UltraPower Fluted Acorn Bur	66
7021-282	UltraPower Fluted Acorn Bur	66
7021-283	UltraPower Fluted Acorn Bur	66
7021-380	UltraPower Fluted Acorn Bur	65
7021-381	UltraPower Fluted Acorn Bur	66
7021-382	UltraPower Fluted Acorn Bur	66
7021-383	UltraPower Fluted Acorn Bur	66
7021-482	UltraPower Fluted Acorn Bur	63
7021-483	UltraPower Fluted Acorn Bur	62
7021-581	UltraPower Fluted Acorn Bur	62
7021-583	UltraPower Fluted Acorn Bur	62
7021-681	UltraPower Fluted Acorn Bur	73
7021-683	UltraPower Fluted Acorn Bur	64
7021-036	UltraPower Fluted Cylindrical Bur	66
7021-037	UltraPower Fluted Cylindrical Bur	66
7021-136	UltraPower Fluted Cylindrical Bur	66
7021-137	UltraPower Fluted Cylindrical Bur	66
7021-235	UltraPower Fluted Cylindrical Bur	64
7021-236	UltraPower Fluted Cylindrical Bur	65
7021-237	UltraPower Fluted Cylindrical Bur	65
7021-335	UltraPower Fluted Cylindrical Bur	64
7021-336	UltraPower Fluted Cylindrical Bur	62
7021-337	UltraPower Fluted Cylindrical Bur	62
7021-338	UltraPower Fluted Cylindrical Bur	62
7021-437	UltraPower Fluted Cylindrical Bur	65
7021-391	UltraPower Fluted Frustrum Bur	66
7021-141-35	UltraPower Fluted Frustrum Bur – Depth Limiting	66
7021-141-50	UltraPower Fluted Frustrum Bur – Depth Limiting	64
7021-141-90	UltraPower Fluted Frustrum Bur – Depth Limiting	66
7021-145-70	UltraPower Fluted Frustrum Bur – Depth Limiting	62
7020-045	Ultrapower Footswitch	67, 83, 84
7020-001	Ultrapower High Speed Drill	61, 83, 84
7020-002	Ultrapower Low Speed Drill	61, 83, 84
7021-095	UltraPower Oval Fluted Bur	64
7021-096	UltraPower Oval Fluted Bur	64
7021-195	UltraPower Oval Fluted Bur	63
7021-196	UltraPower Oval Fluted Bur	63
7021-295	UltraPower Oval Fluted Bur	64

Cat. No.	Description	Page No.
7021-296	UltraPower Oval Fluted Bur	63
7021-395	UltraPower Oval Fluted Bur	63
7021-396	UltraPower Oval Fluted Bur	62
7020-003	Ultrapower Perforator Drive	61, 83, 84
7021-024	UltraPower Round Fluted Bur	62
7021-025	UltraPower Round Fluted Bur	62
7021-026	UltraPower Round Fluted Bur	62
7021-027	UltraPower Round Fluted Bur	62
7021-028	UltraPower Round Fluted Bur	62
7021-029	UltraPower Round Fluted Bur	62
7021-030	UltraPower Round Fluted Bur	62
7021-031	UltraPower Round Fluted Bur	62
7021-032	UltraPower Round Fluted Bur	64
7021-033	UltraPower Round Fluted Bur	64
7021-034	UltraPower Round Fluted Bur	66
7021-122	UltraPower Round Fluted Bur	62
7021-123	UltraPower Round Fluted Bur	62
7021-124	UltraPower Round Fluted Bur	62
7021-125	UltraPower Round Fluted Bur	62
7021-126	UltraPower Round Fluted Bur	62
7021-127	UltraPower Round Fluted Bur	62
7021-128	UltraPower Round Fluted Bur	62
7021-129	UltraPower Round Fluted Bur	62
7021-130	UltraPower Round Fluted Bur	62
7021-131	UltraPower Round Fluted Bur	62
7021-132	UltraPower Round Fluted Bur	64
7021-133	UltraPower Round Fluted Bur	64
7021-134	UltraPower Round Fluted Bur	66
7021-220	UltraPower Round Fluted Bur	62
7021-222	UltraPower Round Fluted Bur	62
7021-223	UltraPower Round Fluted Bur	62
7021-224	UltraPower Round Fluted Bur	62
7021-225	UltraPower Round Fluted Bur	62
7021-226	UltraPower Round Fluted Bur	62
7021-227	UltraPower Round Fluted Bur	62
7021-228	UltraPower Round Fluted Bur	62
7021-229	UltraPower Round Fluted Bur	62
7021-230	UltraPower Round Fluted Bur	62
7021-231	UltraPower Round Fluted Bur	64
7021-232	UltraPower Round Fluted Bur	64
7021-233	UltraPower Round Fluted Bur	64
7021-322	UltraPower Round Fluted Bur	62
7021-323	UltraPower Round Fluted Bur	62
7021-324	UltraPower Round Fluted Bur	62
7021-325	UltraPower Round Fluted Bur	62
7021-326	UltraPower Round Fluted Bur	62
7021-327	UltraPower Round Fluted Bur	62
7021-328	UltraPower Round Fluted Bur	62
7021-329	UltraPower Round Fluted Bur	62
7021-331	UltraPower Round Fluted Bur	64
7021-332	UltraPower Round Fluted Bur	64
7021-333	UltraPower Round Fluted Bur	64
7021-426	UltraPower Round Fluted Bur	62
7021-427	UltraPower Round Fluted Bur	62
7021-428	UltraPower Round Fluted Bur	62

Alphabetical Index



Cat. No.	Description	Page No.
7021-429	UltraPower Round Fluted Bur	62
7021-431	UltraPower Round Fluted Bur	64
7021-432	UltraPower Round Fluted Bur	63
7021-433	UltraPower Round Fluted Bur	65
7021-526	UltraPower Round Fluted Bur	62
7021-528	UltraPower Round Fluted Bur	62
7021-529	UltraPower Round Fluted Bur	62
7021-531	UltraPower Round Fluted Bur	65
7021-532	UltraPower Round Fluted Bur	65
7021-533	UltraPower Round Fluted Bur	62
7021-631	UltraPower Round Fluted Bur	65
7021-632	UltraPower Round Fluted Bur	65
7021-633	UltraPower Round Fluted Bur	64
7021-012	UltraPower Side Cutting Bur	62, 64
7021-013	UltraPower Side Cutting Bur	62
7021-015	UltraPower Side Cutting Bur	62
7021-043	UltraPower Side Cutting Bur	62
7021-112	UltraPower Side Cutting Bur	64
7021-113	UltraPower Side Cutting Bur	64
7021-114	UltraPower Side Cutting Bur	64
7021-115	UltraPower Side Cutting Bur	62
7021-116	UltraPower Side Cutting Bur	62
7021-119	UltraPower Side Cutting Bur	62
7021-143	UltraPower Side Cutting Bur	62
7021-202	UltraPower Side Cutting Bur	63
7021-206	UltraPower Side Cutting Bur	64
7021-207	UltraPower Side Cutting Bur	64
7021-210	UltraPower Side Cutting Bur	62
7021-214	UltraPower Side Cutting Bur	62
7021-218	UltraPower Side Cutting Bur	62
7021-243	UltraPower Side Cutting Bur	65
7021-301	UltraPower Side Cutting Bur	63
7021-302	UltraPower Side Cutting Bur	63
7021-303	UltraPower Side Cutting Bur	63
7021-305	UltraPower Side Cutting Bur	64
7021-309	UltraPower Side Cutting Bur	64
7021-310	UltraPower Side Cutting Bur	64
7021-314	UltraPower Side Cutting Bur	62
7021-315	UltraPower Side Cutting Bur	62
7021-318	UltraPower Side Cutting Bur	62
7021-408	UltraPower Side Cutting Bur	62
7021-409	UltraPower Side Cutting Bur	63
7021-427	UltraPower Side Cutting Bur	62
7021-509	UltraPower Side Cutting Bur	62
7021-713	UltraPower Side Cutting Bur	64, 73
7021-815	UltraPower Side Cutting Bur	73
7021-917	UltraPower Side Cutting Bur	78
7021-051	UltraPower Trephine Bur	62
7021-250	UltraPower Trephine Bur	65
7021-252	UltraPower Trephine Bur	62
7021-253	UltraPower Trephine Bur	62
7021-255	UltraPower Trephine Bur	62
7021-040	UltraPower Wire Pass Bur	64
7021-041	UltraPower Wire Pass Bur	65
7021-042	UltraPower Wire Pass Bur	62

Cat. No.	Description	Page No.
7021-140	UltraPower Wire Pass Bur	66
7021-141	UltraPower Wire Pass Bur	66
7021-142	UltraPower Wire Pass Bur	62
5052-010	Universal Hose, 10' (3.05m)	76, 82 - 84, 86
5053-019	Wire Driver Adapter Chuck for use with 5053-013	42
PRO6228	Wire Driver Range: 0.7 – 1.6mm (0.028 – 0.062")	7
PRO6128	Wire Driver, Range: 0.7 – 1.6mm (0.028 – 0.062")	11
5053-123	Wire Guard for use with 5053-013	42
5056-184	Wire Pass Osteon Bur	75
5056-185	Wire Pass Osteon Bur	75
5056-186	Wire Pass Osteon Bur	75
E4110	Wirepass E Bur	73
E4115	Wirepass E Bur	73
E4120	Wirepass E Bur	73
PRO3011	X-Large NiCad Battery, 12 Volt	9
1375-023	X-Long Bur Guard	43, 54, 55, 56
7020-013	X-Long Bur Guard (like Midas-AM/A) Lime	61 - 66, 83, 84
E9011-3	X-Long Bur Guard for E9010 Drill	67
1375-023P	X-Long Precision Bur Guard	43
1375-040	Y Connector	82
PRO2046	Zimmer/Hudson	13



Hall® Powered Instruments System™

Product Catalog

525 French Road
Utica, New York 13502

Local 727-392-6464
Toll Free 1-800-237-0169

ConMed.com
customer_service@conmed.com

HALL®
POWERED INSTRUMENTS