

Combining clinical performance with ergonomic design



The Tec® 7 Vaporizer from Datex-Ohmeda delivers isoflurane, sevoflurane, enflurane and halothane effectively and efficiently. Each agent specific Tec 7 incorporates significant performance, convenience and ergonomic benefits with important safety features to help improve agent delivery and reliability while reducing overall operating costs.

- Provides consistent output throughout the clinical flow range
- Easy to turn dial and fine graduations help control agent delivery
- Enhanced design helps improve reliability, operating excellence and product life
- Easy-fil™ filler is designed to simplify agent filling and help minimize agent leaks while filling. Tec 7 is also available in funnel fill and Quik-Fil™ (sevoflurane only) variants
- Improved ergonomics and contemporary design complements Datex-Ohmeda's anesthesia systems
- Attaches to Datex-Ohmeda interlocking Selectatec® manifold
- Planned factory service free
- Three year warranty period



Tec 7 agent specific vaporizers

Clinical performance

- Tec 7 Vaporizers are designed to provide consistent output throughout the clinical flow range from 200 mL/min to 15 L/min. Models are available for isoflurane, sevoflurane, enflurane and halothane.
- Tec 7 Vaporizers are equipped with a large diameter control dial which incorporates fine graduations of 0.2% between 0 and 1%, and 0.5% from 1% to 8%. The dial for sevoflurane is marked in steps of 0.2% up to 1% v/v, and in steps of 1% between 1% and 8%.
- The easy turning dials and small graduations on the Tec 7 Vaporizer help you fine tune anesthetic delivery over the full range of dial settings and flow rates.

Enhanced ergonomics

- The ergonomic agent dial release on the Tec 7 Vaporizer allows either left or right hand operation.
- A prismatic site glass on the front panel of the Tec 7 Vaporizer provides clear indication of liquid agent level.

Reliability, convenience and operating excellence

- Easy-fill filler is designed to simplify agent filling and help minimize agent leaks while filling. The Tec 7 Vaporizer is also available with funnel fill and Quik-Fil (sevoflurane only) filling mechanisms.
- To minimize filling frequency, the Tec 7 Vaporizer accommodates 225 mL of anesthetic agent.
- The Tec 7 Vaporizer can be mounted on Datex-Ohmeda anesthesia systems equipped with the Selectatec manifold without tools and without taking the anesthesia system out of service.
- The Tec 7 Vaporizer also interlocks in series with Datex-Ohmeda Tec 4, Tec 5, Tec 6 and Tec 6 Plus Vaporizers.
- An innovative non-spill system in the Tec 7 Vaporizer limits movement of liquid agent if the vaporizer is tilted or inverted, protecting internal components and helping maintain output within clinically acceptable limits.



Easy-fill helps simplify agent filling




Agent color coded Easy-fill bottle adapters

Lower overall ownership costs

- The Tec 7 Vaporizer is planned factory service free, which helps lower ownership costs and eliminate the logistical challenges associated with the return of vaporizers to the factory. For complete user maintenance requirements, refer to the User Reference Manual.
- The Tec 7 Vaporizer comes with a three year warranty.

For more information on how the Tec 7 Vaporizer can benefit your medical practice, contact your local Datex-Ohmeda representative.



Datex-Ohmeda Division • Instrumentarium Corp.
P.O. Box 900, FIN-00031 Datex-Ohmeda, Finland
Tel. +358 10 39411 • Fax +358 9 146 3310

www.datex-ohmeda.com

Datex-Ohmeda, Inc.
P.O. Box 7550, Madison, WI 53707-7550, USA
Tel. 800 345 2700 • Fax 608 221 4384

www.us.datex-ohmeda.com

