



**SLA-600TC**

SUPER-LITE™ Tissue Culture  
 Aluminum Rotor  
 12-Place, 50 ml  
 Tube Dimensions: 30 × 152 mm  
 Fixed Angle: 34°  
 Maximum Speed: 13,000 rpm  
 K Factor: 913 Rotor Mass: 9.8 kg  
 Autoclavable: 121°C

	RCF (g)	R (cm)
Max.	27,750	14.70
Min.	15,083	7.99

		RC-28S	RC-26 Plus	RC-24 <sup>1</sup>	RC-5C Plus RC-5C <sup>1</sup>	RC-5B Plus RC-5B <sup>1</sup> RC-5 <sup>1</sup>	RC-2B <sup>1</sup> RC-2 <sup>1</sup>	SS-3 <sup>1</sup> SS-4 <sup>1</sup>	RC-3B <sup>1</sup> RC-3B Plus RC-3C <sup>1</sup> RC-3C Plus
<b>SLA-600TC</b>	<b>Max Speed(rpm)</b>	n/a	13,000	13,000	12,500	12,500	12,500	n/a	n/a
	<b>Max RCF (g)</b>	n/a	27,750	27,750	25,656	25,656	25,656	n/a	n/a

**Catalog No. 74220**

Includes:  
 SLA-600TC Rotor  
 Rotor Cover Assembly **74221**  
 Vacuum Grease **65937**  
 Lubricant **61556**  
 Rotor Extractor **76818**

Instruction Manual **11463**  
 Rotor Cover O-ring (6 supplied) **92233**  
 Retaining Ring Pliers **65614**  
 Locking Knob O-ring (6 supplied) **91035**  
 Center Hole O-ring (6 supplied) **92231**

**Accessories**

Software Update Kit for RC-24<sup>2</sup> **12939**  
 Software Update Kit for RC-5C<sup>2</sup> **12814**  
 Software Update Kit for RC-5C Plus<sup>2</sup> **12996**  
 Software Update Kit for RC-26 Plus<sup>2</sup> **71366**

**TUBES AND BOTTLES**

NOMINAL CAPACITY (ml)	DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	ADAPTER CAT. NO. (Qty/Pkg)	NO. OF PLACES
50	Conical Tissue Culture Tube	*	30 × 115	—	—
50	Conical Tissue Culture Tube (Freestanding)	*	30 × 115	—	—
35	Nalgene® Conical Oak Ridge Tube	*	29 × 105	<b>74234</b> (6)	—
20	Filtron Macrosep Ultrafiltration Device	*	29 × 109	—	—
15	Conical Tissue Culture Tube	*	17 × 120	<b>74232</b> (6)	1
15	Amicon Centriplus Ultrafiltration Device	*	29 × (variable length)	—	—
15	Amicon Centriprep Ultrafiltration Device	*	28 × 135	—	—
3.5	Filtron Microsep Ultrafiltration Device	*	17 × 99	<b>74233</b> (6)	1
2	Amicon Centricon Ultrafiltration Device	*	16 × (variable length)	<b>74232</b> (6)	1
2	Millipore Ultrafree-CL Ultrafiltration Device	*	17 × 65	<b>74232</b> or <b>74233</b> (6)	1
1	Amicon Centrifree Ultrafiltration Device	*	17 × 95	<b>74233</b> (6)	1

<sup>(1)</sup>Instrument no longer available.

<sup>(2)</sup>Contact your SORVALL® representative to determine your need for a Software Update Kit.

\*Maximum g-force ratings of tubes varies by manufacturer. Verify correct maximum speed/g-force before use.

All trademarks are the property of their respective owners.