

## MA 40

### Air/Bone Conduction Audiometer

Portable air/bone conduction audiometer built for tough, rugged durability.

#### Features at a Glance

- Portable air and bone conduction with Masking (125 – 8000 Hz)
- AC adapter operated
- Standard headphones included, insert phone optional
- Large level and frequency dials
- Continuous, pulse, warble tones
- Sound room compatible
- Hard or soft carrying case included
- NOAH compatible

#### LED Display

- Large LED numeric displays shows the selected frequency and HL intensity.
- Bright individual LEDs are used to show the selected function.

#### Instrument Reliability

MA 40 has a rugged, integral carrying case or optional soft carrying case with ample storage of all accessories. It takes less than one minute to set up the device and start testing. The MA 40 is fast and easy to operate. Sensible keyboard layout with LEDs for each function ensure error free operations. The bright LED-display shows frequency and intensity level settings.



## Technical Data

Frequency Accuracy	+/- 3% maximum
Test Frequencies	125 Hz to 8000 Hz
Crosstalk	-100 dB minimum
Distortion	.5% typical, 3% maximum
On/Off Ratio	Greater than 90 dB
Rise/Fall Time	35 msec (typical)
Frequency Mod. Rate	+/-5% frequency modulating rate at 5 Hz
Dimensions	12.5" W x 6.25" H x 15.5" D (32cm W x 16cm H x 40cm D)
Weight	16.5 lb./ 7.5 kg Hard case: 10.5 kg Soft case: 4.76 kg
Power Supply	117/234 volts AC, 50/60 Hz, 25W

### STANDARD ACCESSORIES

- DD45 Headset
- B-71 bone conductor
- Bond headband
- Audiogram pad

### OPTIONAL ACCESSORIES

- Patient response switch
- Audio cup headset
- Insert phones
- Sound room kit
- Patch cord

### CPT CODES:

- 92551: Screening, pure tone
- 92552: Threshold, air
- 92553: Air and bone

Specifications subject to change without notice.

#### MAICO Diagnostics

10393 West 70th Street, Eden Prairie, MN 55344

Tel: (888) 941-4201 | Fax: (952) 903-4100

e-mail: [info@maico-diagnostics.com](mailto:info@maico-diagnostics.com) | Web: [www.maico-diagnostics.com](http://www.maico-diagnostics.com)

**ONE YEAR LIMITED WARRANTY:** This warranty is extended to the original purchaser of the instrument, by MAICO, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of invoice of the instrument to the original purchaser. Accessories which are purchased from MAICO at the same time as the instrument are also under warranty for one year from the date of invoice. For additional information contact your MAICO Distributor.