

Cast Removal Products

AMERICAN ORTHOPAEDIC™ VACUUM

- Lightweight and mobile; facilitates easy movement from room to room or patient to patient.
- Unique lift-out disposable filter cartridge keeps the system running at peak efficiency and makes waste disposal quicker, easier, safer and cleaner; no more messy bags.
- Quiet operation eases patient's fears through reduction of noise.
- Remote on/off switch allows the saw and vacuum to be turned on simultaneously.
- The ten (10) foot hose increases maneuverability.

*Note: This vacuum may also be adapted to the following cast cutters:

M-PACT #4183-106, Stryker #840, DePuy #3300-00, Zimmer #8909-01



Specifications

Product Measurements	Length	23" (58.4 cm)
	Width	15" (38.1 cm)
	Height	23.5" (59.7 cm)
	Weight	25 lbs. (11.3 kg)

The American Orthopaedic™ Cast Dust Vacuum is supplied with a pre-installed Disposable Filter Cartridge, 10' x 1.25" I.D. hose, 6 hose/cord straps, vacuum nozzle adapter, owner's manual and registration card.

Product #	Description and Electrical Requirements
0295-400*	115 VAC, 50-60 Hz, 8.5 amps
<i>for Europe</i>	
0295-401	230 VAC, 50-60 Hz, 4.8 amps, Europlug



ACCESSORIES

Product #	Description	Unit
TCFC1	Disposable Filter Cartridge for TCCVS1, TCCVS2 or 0295-400 vacuums	ea
0295-175	Portable filter cartridge for 0295-300 vacuum	ea
0295-283	Straight vacuum nozzle assembly connects 4183-106 or 0295-200 cutters to 0295-300 or 0295-400 vacuums	ea
4183-122	Vacuum nozzle assembly connects 4183-106 or 0295-200 cutters to Stryker vacuum w/cord inside the hose	ea
4183-123	Vacuum nozzle adapter connects the 0295-300 or 0295-400 vacuums to Stryker cutter	ea
4183-124	Right angle vacuum nozzle assembly connects 4183-106 or 0295-200 cutters to 0295-300, 0295-400, 0295-250 vacs	ea
0295-405	Vacuum nozzle adapter connects Stryker 940 cutter to 0295-300 or 0295-400 vacuums	ea